

# Chromatin Condensation/Dead Cell Apoptosis Kit with Hoechst 33342 and PI for Flow Cytometry

Catalog no. V13244

**Table 1.** Contents and storage information.

Material	Amount	Concentration	Storage	Stability
Hoechst 33342 (Component A)	200 µL	5.0 mg/mL solution in water	<ul style="list-style-type: none"> <li>• 2–6°C</li> <li>• Protect from light</li> </ul>	When stored as directed, kit components should be stable for at least 6 months.
Propidium iodide (PI, Component B)	200 µL	1.0 mg/mL solution in water		

\* Hoechst 33342 and propidium iodide may be handled in normal room light, but avoid prolonged exposure to light.


**Number of assays:** The kit provides sufficient reagents for 200 flow cytometry assays, each having  $2 \times 10^5$  to  $1 \times 10^6$  cells in a 1 mL volume.

**Approximate fluorescence excitation/emission maxima, when bound to DNA:** ~350/461 nm Hoechst dye; ~535/617 nm propidium iodide

## Introduction

Apoptosis is a carefully regulated process of cell death that occurs as a normal part of development. Inappropriately regulated apoptosis is implicated in disease states, such as Alzheimer's disease and cancer. Apoptosis is distinguished from necrosis, or accidental cell death, by characteristic morphological and biochemical changes, including compaction and fragmentation of the nuclear chromatin, shrinkage of the cytoplasm, and loss of membrane asymmetry.<sup>1-5</sup>

The Chromatin Condensation/Dead Cell Apoptosis Kit with Hoechst 33342 and PI provides a rapid and convenient assay for apoptosis based upon fluorescence detection of the compacted state of the chromatin in apoptotic cells. The kit contains ready-to-use solutions of the blue-fluorescent Hoechst 33342 dye (excitation/emission maxima ~350/461 nm, when bound to DNA), which stains the condensed chromatin of apoptotic cells more brightly than the chromatin of normal cells, and the red-fluorescent propidium iodide dye (excitation/emission maxima ~535/617 nm, when bound to DNA), which is permeant only to dead cells. The staining pattern resulting from the simultaneous use of these dyes makes it possible to distinguish normal, apoptotic, and dead cell populations by flow cytometry (Figure 1).<sup>2,6,7</sup> The 351 nm line of an argon-ion laser or other suitable UV source is required for excitation of the Hoechst 33342 dye, whereas propidium iodide may be excited with the 488 nm line of an argon-ion laser. We have optimized this assay using Jurkat cells, a human T-cell leukemia clone, treated with camptothecin to induce apoptosis. Some modifications may be required for use with other cell types.



**Figure 1.** Jurkat cells (T-cell leukemia, human) treated with 10  $\mu$ M camptothecin for four hours (bottom panel) or untreated (as control, top panel). Cells were then treated with the reagents in the kit and analyzed by flow cytometry using UV/488 nm dual excitation. Note that the camptothecin-treated cells (bottom panel) have a higher percentage of apoptotic cells (indicated by an "A") than the basal level of apoptosis seen in the control cells (top panel). L = live cells, D = dead cells.

## Before You Begin

---

### Materials Required but Not Provided

- Phosphate-buffered saline (PBS)
- Camptothecin or other agent to induce apoptosis

### Caution

Propidium iodide and Hoechst 33342 are known or suspected mutagens and should be used with appropriate precautions.

## Experimental Protocols

---

We have optimized this assay using Jurkat cells treated with camptothecin to induce apoptosis. Some modifications may be required for use with other cell types.

1. Induce apoptosis in cells using the desired method. Prepare a negative control by incubating cells in the absence of inducing agent.
2. Harvest the cells after the incubation period, wash in cold phosphate-buffered saline (PBS) and adjust the cell density to  $\sim 1 \times 10^6$  cells/mL in PBS. For each assay, a 1 mL volume will be used.
3. Add 1  $\mu$ L of the Hoechst 33342 stock solution (Component A) and 1  $\mu$ L of the PI stock solution (Component B) to each 1 mL of cell suspension.
4. Incubate the cells on ice for 20–30 minutes.
5. As soon as possible after the incubation period, analyze the stained cells by flow cytometry, using UV/488 nm dual excitation and measuring the fluorescence emission at  $\sim 460$  nm and  $>575$  nm. The population should separate into three groups: live cells will show only a low level of fluorescence, apoptotic cells will show a higher level of blue fluorescence, and dead cells will show both blue and red fluorescence (see Figure 1). Confirm the flow cytometry results by viewing the cells under a fluorescence microscope, using filters appropriate for Hoechst or DAPI and rhodamine (TRITC) or Texas Red<sup>®</sup> dye.

## References

---

1. Immunol Cell Biol 76, 1 (1998); 2. Cytometry 27, 1 (1997); 3. J Pharmacol Toxicol Methods 37, 215 (1997); 4. FASEB J 9, 1277 (1995); 5. Am J Pathol 146, 3 (1995); 6. Cytometry 17, 59 (1994); 7. Exp Cell Res 211, 322 (1994).

## Product List

Current prices may be obtained from our website or from our Customer Service Department.

---

Cat. no.	Product Name	Unit Size
V13244	Chromatin Condensation/Dead Cell Apoptosis Kit with Hoechst 33342 and PI *for flow cytometry* *200 assays* .....	1 kit

## Contact Information

---

### **Molecular Probes, Inc.**

29851 Willow Creek Road  
Eugene, OR 97402  
Phone: (541) 465-8300  
Fax: (541) 335-0504

### **Customer Service:**

6:00 am to 4:30 pm (Pacific Time)  
Phone: (541) 335-0338  
Fax: (541) 335-0305  
probesorder@invitrogen.com

### **Toll-Free Ordering for USA:**

Order Phone: (800) 438-2209  
Order Fax: (800) 438-0228

### **Technical Service:**

8:00 am to 4:00 pm (Pacific Time)  
Phone: (541) 335-0353  
Toll-Free (800) 438-2209  
Fax: (541) 335-0238  
probestech@invitrogen.com

### **Invitrogen European Headquarters**

Invitrogen, Ltd.  
3 Fountain Drive  
Inchinnan Business Park  
Paisley PA4 9RF, UK  
Phone: +44 (0) 141 814 6100  
Fax: +44 (0) 141 814 6260  
Email: euroinfo@invitrogen.com  
Technical Services: eurotech@invitrogen.com

**For country-specific contact information,  
visit [www.invitrogen.com](http://www.invitrogen.com).**

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

### **Limited Use Label License No. 223: Labeling and Detection Technology**

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2009, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.