

ThiolTracker™ Violet (Glutathione Detection Reagent)

Catalog nos. T10095, T10096

Table 1. Contents and storage information.

Material	Amount	Concentration	Storage*	Stability
ThiolTracker™ Violet	180 assays (T10095) 5 microplates (T10096)	Solid	<ul style="list-style-type: none"> • ≤ -20°C • Desiccate • Protect from light 	When stored as directed the product is stable for 1 year

Number of assays: Sufficient material is supplied for T10095 for staining suspension cells in 180 tubes (0.5 mL/tube) for flow cytometric analysis. Sufficient material is supplied for T10096 for staining cells in 5 × 96-well microplates for imaging, or cells on approximately 120 of 5 mm round coverslips for imaging.

Approximate fluorescence excitation/emission maxima: 404/526 in nm

Introduction

ThiolTracker™ Violet dye is ideal for cellular labeling due to the following properties of the dye: (1) reacts actively with reduced thiols in intact cells, (2) can be efficiently excited with 405 nm laser or traditional xenon or mercury arc lamps for use in cellular analysis using flow cytometry and fluorescence microscopy, (3) excitation and emission property is well suitable for imaging with fluorescence microscopes equipped with filter set for Hoechst 33342 dye.

Glutathione (GSH) plays a central role in protecting mammalian cells against damage incurred by free radicals, oxidants, and electrophiles. Since reduced glutathione represents the majority of intracellular free thiols in cells, ThiolTracker™ Violet dye can be used in estimating the cellular level of reduced glutathione.

ThiolTracker™ Violet dye is applied directly to live cells in thiol-free buffer to stain the cells. After labeling, live cells can be directly imaged or cells can be fixed with aldehyde before imaging. Preservation of ThiolTracker™ Violet dye signal upon Triton® X-100 permeabilization step is excellent allowing ThiolTracker™ Violet dye labeled cells to be probed with antibodies or other probes.

Before You Begin

Materials Required but Not Provided

- Dimethyl sulfoxide (DMSO), high quality DMSO (>99.7 %, such as Sigma D-2650, anhydrous) is recommended.

- Dulbecco's phosphate buffered saline, containing Ca⁺⁺ and Mg⁺⁺, glucose, and sodium pyruvate. (D-PBS with C/M, Invitrogen Cat. no. 14287-080)
- Dulbecco's phosphate buffered saline, without Ca⁺⁺ or Mg⁺⁺, (D-PBS, Invitrogen Cat. no. 14190-144)
- Formaldehyde 37% aqueous solution

Preparing Cells

Seed cells the day before addition of the test compound. For adherent cells, optimize the cell number and plate coating requirements for the chosen cell model and time span of test compound treatment before labeling with ThiolTracker™ Violet dye. While many cell lines tested do not require special coating of surface, it is highly recommended to seed HepG2 cells on poly-lysine coated surface for staining with ThiolTracker™ Violet dye. Similar coating requirements might be applicable to other cell lines. Detachment of cells upon labeling with ThiolTracker™ Violet dye may indicate the need for optimization of dye concentration, labeling time, and substratum requirements.

ThiolTracker™ Violet

Follow these guidelines for preparing and storing ThiolTracker™ Violet dye:

- Prepare the ThiolTracker™ Violet dye stock solution in DMSO as described in the protocol. Alternatively, you can use 1-methyl-2-pyrrolidine to prepare the dye.
- It may take up to ten minutes to completely dissolve ThiolTracker™ Violet dye in DMSO. Pipeting up and down facilitates dye solubilization into DMSO and appears to work better than vortexing.
- In the protocols below, the ThiolTracker™ Violet dye stock solution concentration in DMSO is 20 mM. The working concentration is 20 μM for imaging and 10 μM for flow cytometry. Lower concentrations can be chosen for the stock solution to suit the optimized working concentrations of ThiolTracker™ Violet dye for a cell line of interest.
- Once reconstituted in DMSO, use the reconstituted ThiolTracker™ Violet dye **immediately**. Store any remaining unused stock solution of ThiolTracker™ Violet dye at –20°C and use the stock solution within 2 weeks of preparation. Avoid repeated freeze-thaw steps.

Labeling Guidelines

Since ThiolTracker™ Violet is a thiol reactive dye, use **thiol-free buffer** for labeling live cells. Our studies indicate that labeling cells with 20 μM ThiolTracker™ Violet dye for 30 minutes at 37°C gives excellent results in reporting the increase or decrease of stain intensity after treatment of cells with agents that modulate the intracellular reduced glutathione levels. However, some optimization of the working dye concentration and labeling time for your cell line of interest may be necessary. Many cell lines tested were found to be labeled well with ThiolTracker™ Violet dye ranging from 10 μM to 20 μM with the protocols outlined below.

Experimental Protocol for Labeling Cells in 96-well Plates for Imaging

Preparing Dye Stock Solution

Prepare ThiolTracker™ Violet dye in DMSO as follows to obtain a final concentration of 20 mM:

- For Cat. no. T10095, add 15 μL DMSO to each vial
- For Cat. no. T10096, add 75 μL DMSO to each vial

Preparing Dye Working Solution

Prepare ThiolTracker™ Violet dye working solution **just before use** by diluting the ThiolTracker™ Violet dye stock solution into D-PBS with C/M to the working concentration.

For example, you need 100 µL of working dye solution for each well of a 96-well plate. Prepare 12 mL of **20 µM** ThiolTracker™ Violet dye working solution by diluting 12 µL of 20 mM ThiolTracker™ Violet stock solution into 12 mL of D-PBS with C/M for one 96-well plate.

Protocol for Labeling Cells in 96-well Plates

- 1.1 After test compound or drug treatment of cells (if required), remove the incubation medium from the wells of the 96-well plate.
- 1.2 Rinse cells twice with 100 µL D-PBS C/M per well. Remove the D-PBS C/M.
- 1.3 Add 100 µL prewarmed ThiolTracker™ Violet dye working solution in D-PBS C/M. Incubate the plate in a 37°C cell culture incubator for 30 minutes.
- 1.4 Replace the ThiolTracker™ Violet dye working solution with a suitable buffer or medium.
- 1.5 The cells are now ready for imaging using a fluorescence microscope or other suitable fluorescence imaging instrument. The approximate fluorescence excitation and emission wavelength for ThiolTracker™ Violet dye is 404/526 in nm.
- 1.6 **Optional:** Perform optional formaldehyde fixation as follows:
 - a. Prepare 3–4% formaldehyde solution in D-PBS. You need 100 µL/well for a 96-well plate.
 - b. Remove the ThiolTracker™ Violet dye working solution and add 100 µL 3–4% formaldehyde in D-PBS to each well.
 - c. Incubate the plate in a fume hood at room temperature for 30 minutes.
 - d. Remove the formaldehyde fixative and rinse cells twice with D-PBS. The cells are now ready for imaging.
- 1.7 **Optional:** Perform optional counterstaining for fixed samples. For example, counterstain for nucleus using nucleic acid dyes for cell identification purposes often employed in automated image analysis (high-content analysis).

Note: We recommend imaging the cells immediately after labeling and processing. If the plate cannot be imaged immediately, keep the processed plate at 4°C.

Experimental Protocol for Labeling Cells on Coverslips for imaging

Preparing Dye Stock Solution

Prepare ThiolTracker™ Violet dye in DMSO as follows to obtain a final concentration of 20 mM:

- For Cat. no. T10095, add 15 μL DMSO to each vial
- For Cat. no. T10096, add 75 μL DMSO to each vial

Preparing Dye Working Solution

Prepare ThiolTracker™ Violet dye working solution **just before use** by diluting the ThiolTracker™ Violet dye stock solution into D-PBS with C/M to the working concentration.

Prepare enough dye working solution for the type of cell culture dish holding the coverslips. For example, 1 mL of ThiolTracker™ Violet dye working solution in D-PBS with C/M is sufficient for coverslips in a well of a 6-well plate or a 3-cm cell culture dish. Required volume can be adjusted according to the size of coverslips and the type of cell culture vessel used.

Protocol for Labeling Cells on Coverslips

- 2.1 After test compound or drug treatment of cells (if required), remove the incubation medium.
- 2.2 Rinse cells twice with D-PBS C/M. Remove the D-PBS C/M.
- 2.3 Add the appropriate volume of prewarmed ThiolTracker™ Violet dye working solution in D-PBS C/M depending on the cell culture dish. Incubate cells in a 37°C cell culture incubator for 30 minutes.
- 2.4 Replace the ThiolTracker™ Violet dye working solution with a suitable buffer or medium.
- 2.5 The cells are now ready for imaging using a fluorescence microscope or other suitable fluorescence imaging instrument. The approximate fluorescence excitation and emission wavelength for ThiolTracker™ Violet dye is 404/526 nm.
- 2.6 **Optional:** Perform optional formaldehyde fixation as follows:
 - a. Prepare 3–4% formaldehyde solution in D-PBS.
 - b. Remove the ThiolTracker™ Violet dye working solution and add 3–4% formaldehyde in D-PBS to each coverslip.
 - c. Incubate cells in a fume hood at room temperature for 30 minutes.
 - d. Remove the formaldehyde fixative and rinse cells twice with D-PBS. The cells are now ready for imaging.

Note: We recommend imaging the cells immediately after labeling and processing. If the plate cannot be imaged immediately, keep the processed plate at 4°C.
- 2.7 **Optional:** Perform optional counterstaining for fixed samples. For example, counterstain for nucleus using nucleic acid dyes for cell identification purposes often employed in automated image analysis (high-content analysis).

Note: We recommend imaging the cells immediately after labeling and processing. If the plate cannot be imaged immediately, keep the processed plate at 4°C.

Experimental Protocol for Labeling Suspension Cells for Flow Cytometry

Preparing Dye Stock Solution Prepare ThiolTracker™ Violet dye in DMSO as follows to obtain a final concentration of 10 mM:

- For Cat. no. T10095, add 30 µL DMSO to each vial
- For Cat. no. T10096, add 150 µL DMSO to each vial

Preparing Dye Working Solution

Prepare ThiolTracker™ Violet dye working solution **just before use** by diluting the ThiolTracker™ Violet dye stock solution into D-PBS with C/M to 2X working concentration.

Use a working concentration of 10 µM as a starting concentration for optimization. Required volume of working dye solution is 250 µL per tube for a final volume of 500 µL.

Protocol for Labeling Suspension Cells

- 3.1 After test compound or drug treatment of cells (if required), remove the incubation medium by centrifugation.
- 3.2 Add the appropriate volume of D-PBS with C/M to wash cells. Pellet cells by centrifugation and discard the supernatant.
- 3.3 Resuspend the cells in an appropriate volume of D-PBS with C/M. Count cells and then adjust the cell density to 2X the desired final cell density. Dispense 250 µL cell suspension in a sterile test tube.
- 3.4 Add 250 µL 2X ThiolTracker™ Violet dye working solution in D-PBS with C/M to each tube containing 250 µL of cells in suspension. Gently mix. Incubate cells in a 37°C cell incubator for 30 minutes, protected from light.
- 3.5 Wash cells once with D-PBS or other buffer. Pellet cells by centrifugation and discard the supernatant. Resuspend cells in D-PBS with C/M or other buffer. Acquire data using a flow cytometer at 405 nm excitation and 525BP or similar emission filter.

Product List Current prices may be obtained from our website or from our Customer Service Department.

Cat. no.	Product Name	Unit Size
T10095	ThiolTracker™ Violet (Glutathione Detection Reagent) *for 180 assays*	1 set of 3 vials
T10096	ThiolTracker™ Violet (Glutathione Detection Reagent) *for 5 microplates*	each

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probestech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2008, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.