

SiteClick™ Antibody Labeling Kits

Table 1 Contents and storage

Material	Cap color	Amount	Storage*
Antibody preparation buffer (Component A)	Yellow	1.8 mL	<ul style="list-style-type: none">• 2–8°C• DO NOT FREEZE• Protect from light.
Antibody concentrator (small) (Component B)	N/A	each	
Collection tube (Component C)	N/A	each	
β-Galactosidase (Component D)	Green	12 μL	
UDP-GalNAz (Component E)	Blue	220 μg	
20X Tris pH 7.0 (Component F)	Red	1.8 mL	
Buffer additive (Component G)	Purple	30 μL	
β-1,4-galactosyltransferase (GalT) (Component H)	Orange	88 μL	
Antibody concentrator (large) (Component I)	N/A	each	
DIBO-modified label (Component J)†	Dark Orange	55 μL (Qdot®) or 80 μL (R-PE)	
Purification concentrator (Component K)‡	N/A	each	
* When stored as directed, this kit is stable for at least 3 months. † DIBO-modified Qdot® nanocrystal or DIBO-modified R-Phycoerythrin (R-PE). ‡ Only available with Cat. no. S10469. N/A = not applicable.			
Approximate fluorescence excitation and emission maxima: See Table 2, page 4.			

Introduction

The SiteClick™ Antibody Labeling Kits allow you to conjugate your own antibodies to DIBO-modified Qdot® nanocrystals (525, 565, 585, 605, 625, 655, 705, and 800 nm emission) or DIBO-modified R-Phycoerythrin (R-PE). The SiteClick™ conjugation workflow consists of three steps (antibody carbohydrate domain modification, azide attachment to the antibody, and conjugation with the DIBO-modified label) and relies on copper-free click chemistry to covalently link the label containing the DIBO moiety with the azide-modified antibody without reducing the protein. The antibody concentrators provided in the kits are used to purify and concentrate the antibody at each step of the SiteClick™ antibody labeling workflow (Figure 1, page 2).

In the first step of SiteClick™ conjugation, terminal galactose residues on the N-linked sugars in the Fc region of the antibody are removed by β -Galactosidase. The azide-containing sugar, GalNAz, is then added to the modified carbohydrate domain of the antibody via the β -1,4-galactosyltransferase (Gal-T)-catalyzed reaction targeting the terminal GlcNAc residues. This specific targeting maintains the integrity of the antigen binding site on the antibody. Finally, the antibody (now containing an azide moiety) is conjugated to the DIBO-modified label (Qdot® nanocrystals or R-PE) in a copper-free click reaction with simple overnight incubation (Figure 2, below).

Each SiteClick™ Antibody Conjugation Kit contains sufficient reagents to perform one conjugation reaction of Qdot® nanocrystals or R-PE to a primary IgG antibody sample. The protocol in this manual describes a conjugation reaction starting with 100–125 μ g of whole IgG from any host species.

Figure 1 SiteClick™ antibody labeling workflow




Figure 2 SiteClick™ conjugation reaction




Table 2 Optimal fluorescence excitation and emission maxima of DIBO-modified labels

Label	Ex *	Em *	Extinction coefficient (ϵ)	Measured At	Cat. no.
Qdot® 525	<525 nm	525 nm	200,000 M ⁻¹ cm ⁻¹	between 504 nm and 512nm	S10449
Qdot® 565	<565 nm	565 nm	300,000 M ⁻¹ cm ⁻¹	between 548 nm and 556 nm	S10450
Qdot® 585	<585 nm	585 nm	250,000 M ⁻¹ cm ⁻¹	between 572 nm and 580 nm	S10451
Qdot® 605	<605 nm	605 nm	400,000 M ⁻¹ cm ⁻¹	between 592 nm and 600 nm	S10469
Qdot® 625	<625 nm	625 nm	500,000 M ⁻¹ cm ⁻¹	between 605 nm and 612 nm	S10452
Qdot® 655	<655 nm	655 nm	1,700,000 M ⁻¹ cm ⁻¹	550 nm	S10453
Qdot® 705	<705 nm	705 nm	1,700,000 M ⁻¹ cm ⁻¹	550 nm	S10454
Qdot® 800	<800 nm	800 nm	1,700,000 M ⁻¹ cm ⁻¹	550 nm	S10455
R-Phycoerythrin	496, 546, 565 nm †	578 nm	1,960,000 M ⁻¹ cm ⁻¹	578 nm	S10467

* Qdot® nanocrystals are excitable (Ex, in nm) at any wavelength below their emission maxima (Em, in nm). For most practical applications, they should be excited at wavelengths below 405 nm.
† Multiple absorbance peaks.

Before Starting

Equipment Required

- Centrifuge with fixed angle rotor that can accommodate 1.5-mL centrifuge tubes.
- Centrifuge with swinging bucket or fixed angle rotor with 17 mm × 100 mm inserts

Materials Required but Not Provided

- 100 to 125 µg of whole IgG antibody, preferably at a concentration of 2 to 4 mg/mL in a Tris-based buffer, free of carrier proteins and/or azide.
- Centrifuge tubes: 1.5-mL and 15-mL
- Distilled water (dH₂O)

Caution

- β-Galactosidase (Component D) may cause an allergic skin reaction, and it may cause allergy or asthma symptoms or breathing difficulties, if inhaled. Read the Safety Data Sheet (SDS), available at www.lifetechnologies.com, before handling this reagent.
- The Qdot® conjugate contains cadmium and selenium in an inorganic crystalline form.
- Dispose of the reagents in compliance with all pertaining local regulations. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Always wear suitable laboratory protective clothing and gloves when handling these reagents.

Step 1. Antibody Concentration and/or Buffer Exchange (Optional)

Time Required: 1 hour

This antibody concentration and buffer exchange step is required if:

- Your antibody concentration is less than 2 mg/mL, and/or
- Your antibody is in a phosphate-based buffer (e.g. PBS), and/or
- Your antibody is in a buffer containing azide.

Before you begin, briefly centrifuge Components A, C, D, E, F, G, H, and J to ensure all material is at the bottom of the tubes.

Wash Antibody Concentrator

- 1.1 Add 450 μL of dH_2O to the small antibody concentrator (Component B) and cap the device as shown in Figure 3, below.
- 1.2 Centrifuge for 6 minutes at $5000 \times g$, ensuring that the cap strap and one membrane panel of the concentrator faces the center of the rotor.
- 1.3 Discard the flow through.


Concentrate Antibody and Exchange Buffer

- 1.4 Add a sufficient volume of antibody solution to contain 100–125 μg of antibody to the small antibody concentrator. For example, if the antibody concentration is 1 mg/mL , add 125 μL .
- 1.5 Dilute the added antibody to 500 μL using antibody preparation buffer (Component A).
- 1.6 Centrifuge for 6 minutes at $5000 \times g$, ensuring that the cap strap and one membrane panel of the concentrator faces the center of the rotor.
- 1.7 Discard the flow through.
- 1.8 Add 450 μL of antibody preparation buffer (Component A) to the small antibody concentrator (Component B) and centrifuge for 6 minutes at $5000 \times g$, ensuring that the cap strap and one membrane panel of the concentrator faces the center of the rotor.

Note: If antibody volume in concentrator is greater than 50 μL following Step 1.8, centrifuge for an additional 3 minutes at $5000 \times g$ or until the appropriate volume is achieved.

- 1.9 Invert the small antibody concentrator (Component B) into the collection tube (Component C) as shown in Figure 3, below.
- 1.10 Centrifuge for 3 minutes at $1000 \times g$ to collect the concentrated antibody. Following collection, you should have approximately 50 μL of concentrated antibody in the collection tube.

Figure 3 Antibody concentration and/or buffer exchange



Step 2. Modification of Antibody Carbohydrate Domain

Time Required: 4 hours, hands-off

Add β -galactosidase

- 2.1 Add 10 μ L of β -galactosidase (Component D) to the antibody collected in Step 1.10, as shown in Figure 4, below.
- 2.2 Wrap the tube cap with Parafilm[®] laboratory film or similar and incubate for 4 hours at 37°C.


Step 3. Azide Attachment

Time Required: 5 minutes hands-on, then overnight incubation

Add GalT enzyme

- 3.1 Prepare the azide modification solution by adding the following components to the tube containing UDP-GalNAz (Component E), as shown in Figure 4, below:
 - 75 μ L of dH₂O
 - 10 μ L of 20 \times Tris buffer, pH 7.0 (Component F)
 - 25 μ L of buffer additive (Component G)
 - 80 μ L of GalT enzyme (Component H)
- 3.2 Vortex the reaction components and then add the modified antibody from Step 2.2 to the tube.
- 3.3 Briefly centrifuge the tube, wrap the tube cap with Parafilm[®] laboratory film or similar, and incubate overnight at 30°C.

Figure 4 Modification of antibody carbohydrate domain and azide attachment



Step 4. Purification and Concentration of Azide-Modified Antibody

Time Required: 1 hour

Wash Antibody Concentrator

- 4.1. Prepare 10 mL of 1× Tris, pH 7.0 by adding 500 µL of 20× Tris, pH 7.0 (Component F) to 9.5 mL of dH₂O in a 15-mL conical tube. Vortex briefly to mix.


Note: TBS may also be used for the purification and collection of the modified antibody (Steps 4.2 – 4.12). 20× Tris, pH 7.0 is provided for convenience.

- 4.2 Remove the conical collection tube from the large antibody concentrator (Component I) as shown in Figure 5, below.
- 4.3 Add 1 mL of 1× Tris, pH 7.0 to the large antibody concentrator (Component I) and centrifuge for 10 minutes at 1200 × g, ensuring that one membrane panel of the concentrator faces the center of the rotor.
- 4.4 Discard flow through.

Purify the Antibody

- 4.5 Add 1.75 mL of 1× Tris, pH 7.0 and 250 µL of the azide modified antibody from Step 3.3 to the large antibody concentrator (Component I) as shown in Figure 5, below.
- 4.6. Centrifuge for 6 minutes at 1200 × g, ensuring one membrane panel of the concentrator faces the center of the rotor.
- 4.7 Discard flow through.
- 4.8 Add 1.8 mL of 1× Tris, pH 7.0 to the large antibody concentrator (Component I) and centrifuge for 10 minutes at 1200 × g, ensuring that one membrane panel of the concentrator faces the center of the rotor.
- 4.9 Discard flow through and repeat Step 4.8 once more.

Figure 5 Purification and concentration of azide-modified antibody



Concentrate and Collect Antibody

- 4.10** Add 1.8 mL of 1× Tris, pH 7.0 to the large antibody concentrator (Component I) and centrifuge for 10 minutes at $1400 \times g$. Discard flow through. The final volume in the concentrator should be approximately 80–120 μL

Note: If antibody volume in concentrator is greater than 100 μL , centrifuge for an additional 5 minutes at $1400 \times g$ or until the appropriate volume is achieved.


- 4.11** Invert the antibody concentrator into the conical collection tube as shown in Figure 6, below.

- 4.12** Centrifuge for 3 minutes at $1000 \times g$ to collect the concentrated antibody.

- 4.13** Transfer the antibody from the conical collection tube to a 1.5 mL centrifuge tube. If the final collected volume is less than 100 μL , dilute antibody to 100 μL with 20× Tris, pH 7.0.

Note: At this stage, the antibody can be stored at 2–8°C for attachment of the DIBO-modified label at a later time.

Figure 6 Collection of purified and concentrated azide-modified antibody



Step 5. Conjugation with DIBO-modified Label

Time Required: 5 minutes hands-on, then overnight incubation

Add DIBO-modified Label

- 5.1** Add the DIBO-modified label (Component J) to the azide-modified antibody in the 1.5-mL centrifuge tube:
- If using DIBO-modified Qdot[®] nanocrystal, add 50 μL of Qdot[®] DIBO (Component J).
 - If using DIBO-modified R-PE, add 75 μL of R-PE DIBO (Component J).
- 5.2** Vortex the reaction mixture, briefly centrifuge, and incubate overnight at 25°C.
- 5.3** The antibody conjugate can now be stored at 2–8°C, protected from light (see **Antibody Conjugate Storage**, page 10) or optionally purified of excess antibody (Step 6, page 9).

Step 6. Purification and Concentration of Antibody Conjugate (Optional)

Time Required: 1 hour

The purification step removes any excess antibody that has not been conjugated with the Qdot® DIBO-modified label as well as unconjugated Qdot® nanocrystals. Typical yield from this step is approximately 80%.

Note: The flow-through from the final concentration and purification step may be colored, and this color may vary. This is normal and will not impact the quality of the final purified conjugate.


Note: For many applications, purification of the antibody conjugate is not necessary because the efficiency of the SiteClick™ conjugation reaction yields a minimal quantity of free antibody.

Note: The purification concentrator (Component K) is only available in SiteClick™ antibody purification kits containing DIBO-modified Qdot® nanocrystals with emission maxima of 605, 625, 655, 705, and 800 nm. Antibody conjugates labeled with R-PE or with Qdot® nanocrystals that emit at 525, 565, or 585 nm are too small to be retained by the membrane in the purification concentrator.

Purify and Concentrate Antibody

- 6.1 Add 500 µL of dH₂O to the purification concentrator (Component K) and centrifuge for 5 minutes at 1500 × g. Discard flow through.
- 6.2 Add 150 µL of the labeled antibody conjugate from Step 5.2 to the purification concentrator (Component K) as shown in Figure 7, below.
- 6.3 Add buffer (TBS or PBS) to the purification concentrator to bring the volume to 500 µL and centrifuge for 10 minutes at 1500 × g. Discard flow through.
- 6.4 Add 500 µL of buffer (TBS or PBS) to the purification concentrator and centrifuge for 10 minutes at 1500 × g.
- 6.5 Discard flow through and repeat Step 6.4 once more.
- 6.6 Collect the purified labeled antibody conjugate from the top of the concentrator and transfer to a new 1.5 mL centrifuge tube. The antibody conjugate can now be stored at 2–8°C, protected from light (see **Antibody Conjugate Storage**, page 10).

Figure 7 Optional purification of conjugated antibody



Storing and Using Antibody Conjugates

Antibody Conjugate Storage	Store the antibody conjugate labeled with the Qdot® nanocrystal or R-PE at 2–8°C, in the dark, until needed. DO NOT FREEZE. Sodium azide or thimerosal can be added at this stage to a final concentration of 0.02% (w/v) for long term storage, if preferred.
Antibody Conjugate Use	<p>The concentration of the label in the final preparation can be determined by measuring the optical density of the label at the specified wavelength and then using the formula $A = \epsilon cL$, where A is the absorbance, ϵ is the molar extinction coefficient (Table 2, page 4), c is the molar concentration, and L is the pathlength.</p> <p>For example, for the Qdot® 605 nanocrystal, if the material eluting from the final column has A = 0.80 measured at peak between 592 nm and 600 nm in a cuvette with 1 cm pathlength , then</p> $c = A/\epsilon = 0.80/400,000 = 2 \mu\text{M of label, based on nanocrystal absorbance.}$ <p>The optimal concentration of any labeled antibody should be determined empirically for a particular application or experiment. Determine optimal working concentrations by performing a titration series for the application of interest.</p>

Product List

Current prices may be obtained from our website or from our Customer Service Department.

Cat. no.	Product Name	Unit Size
S10449	SiteClick™ Qdot® 525 Antibody Labeling Kit	1 kit
S10450	SiteClick™ Qdot® 565 Antibody Labeling Kit	1 kit
S10451	SiteClick™ Qdot® 585 Antibody Labeling Kit	1 kit
S10469	SiteClick™ Qdot® 605 Antibody Labeling Kit	1 kit
S10452	SiteClick™ Qdot® 625 Antibody Labeling Kit	1 kit
S10453	SiteClick™ Qdot® 655 Antibody Labeling Kit	1 kit
S10454	SiteClick™ Qdot® 705 Antibody Labeling Kit	1 kit
S10455	SiteClick™ Qdot® 800 Antibody Labeling Kit	1 kit
S10467	SiteClick™ R-PE Antibody Labeling Kit	1 kit

Purchaser Notification

Corporate Headquarters

5791 Van Allen Way
Carlsbad, CA 92008
USA
Phone: +1 760 603 7200
Fax: +1 760 602 6500
Email: techsupport@lifetech.com

European Headquarters

Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF
UK
Phone: +44 141 814 6100
Toll-Free Phone: 0800 269 210
Toll-Free Tech: 0800 838 380
Fax: +44 141 814 6260
Tech Fax: +44 141 814 6117
Email: euroinfo@invitrogen.com
Email Tech: eurotech@invitrogen.com

Japanese Headquarters

LOOP-X Bldg. 6F
3-9-15, Kaigan
Minato-ku, Tokyo 108-0022
Japan
Phone: +81 3 5730 6509
Fax: +81 3 5730 6519
Email: jpinfo@invitrogen.com

Additional international offices are listed at
www.lifetechnologies.com

These high-quality reagents and materials must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Read the Safety Data Sheet provided for each product; other regulatory considerations may apply.

Obtaining Support

For the latest services and support information for all locations, go to www.lifetechnologies.com.

At the website, you can:

- Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities
- Search through frequently asked questions (FAQs)
- Submit a question directly to Technical Support (techsupport@lifetech.com)
- Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents
- Obtain information about customer training
- Download software updates and patches

SDS

Safety Data Sheets (SDSs) are available at www.lifetechnologies.com/sds.

Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to www.lifetechnologies.com/support and search for the Certificate of Analysis by product lot number, which is printed on the product packaging (tube, pouch, or box).

Disclaimer

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial services of any kind, including, without limitation, reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008.

The trademarks mentioned herein are the property of Life Technologies Corporation and/or its affiliate(s) or their respective owners.

Parafilm is a registered trademark of Bemis Company, Inc.

©2013 Life Technologies Corporation. All rights reserved.

