

Rhodamine 110–Based Proteinase Substrates

Table 1. Contents and storage information.

Material	Amount	Storage	Stability
Bisamide substrate of rhodamine 110, lyophilized powder	See product label *	<ul style="list-style-type: none"> • $\leq -20^{\circ}\text{C}$ • Desiccate • Protect solutions from light 	When stored as directed, product is stable for at least 6 months.

* The molecular weight is indicated on the product label.

Approximate Fluorescence Excitation and Emission, in nm: 498/521

Introduction

Molecular Probes' bisamide derivatives of rhodamine 110 are sensitive and selective substrates for assaying proteinases in solution or inside living cells. Originally developed by Walter F. Mangel and colleagues, these fluorogenic substrates contain an amino acid or peptide covalently linked to each of rhodamine 110's amino groups (Figure 1).^{1,2} Upon enzymatic cleavage, the nonfluorescent bisamide substrate is converted first to the fluorescent monoamide and then to rhodamine 110, with a further increase in fluorescence (Figure 2). Both enzymatic products exhibit spectral properties similar to those of fluorescein, with peak excitation and emission wavelengths of 498 nm and 521 nm, respectively. Thus, proteinase assays that employ these substrates are compatible with flow cytometers and other argon laser–based instrumentation. Moreover, unlike fluorescein, the fluorescence intensity of the monoamide and rhodamine 110 is constant from pH 3–9.

The bis-(benzyloxycarbonyl-L-arginine amide) derivative of rhodamine 110 (bis-(CBZ-Arg)-R110, R6501) is a general substrate for serine proteinases. With several enzymes, this substrate has been shown to be 50- to 300-fold more sensitive than the analogous coumarin-based substrate.² The increased sensitivity of rhodamine 110–based assays can be attributed both to the greater fluorescence of the enzymatic product and to the enhanced reactivity of the cleavage site. The tripeptide derivative bis-(CBZ-Ile-Pro-Arg)-rhodamine 110 (known as BZiPAR, R6505) has also proven to be an excellent substrate for serine proteinases. Because this substrate allows direct and continuous monitoring of enzyme turnover, it was used by Leytus and co-workers to determine individual kinetic constants of fast-acting, irreversible trypsin inhibitors.³




Figure 1. Chemical structure of the bis-(CBZ-L-arginine amide) derivative of rhodamine 110 (R6501).

This strategy has led to the development of several peptide derivatives of rhodamine 110 (Table 2) that can be used to detect specific proteinase activity *in vitro* and *in vivo*. Reyftmann and colleagues have shown that BZiPAR enters intact cells and acts as a substrate for lysosomal proteinases that are released in response to porphyrin photosensitization.⁴ The bis-(CBZ-Phe-Arg) derivative of rhodamine 110 (R6502) has been employed for flow cytometric measurement of the cysteine proteinases cathepsin B and L in human monocytes and rat macrophages.⁵⁻⁷ In similar experiments, the bis-(CBZ-Ala-Ala) derivative of rhodamine 110 (R6504) was used to measure the activity of the lysosomal serine proteinase elastase.⁶

Rhodamine 110 can be derivatized with other peptides to create any number of proteinase-specific substrates. For example, rhodamine 110 has been derivatized with the consensus sequence for a human adenovirus proteinase and then employed to identify cofactors required for virion proteinase activity.⁸ Molecular Probes also offers rhodamine 110 derivatives which serve as effective substrates for various caspases.⁹ These products are described in more detail in our product information sheet entitled *Rhodamine 110-Based Caspase Substrates*, mp22120.

Table 2 lists the rhodamine 110-based substrates that we currently have available, along with the proteinases for which they were designed. Our Custom and Bulk Sales Department will be pleased to respond to your inquiries about other rhodamine 110-based proteinase substrates to meet your specific research needs.

Our quality analysis includes HPLC and for many substrates, spectrophotometric quantitation of the rhodamine 110 released by complete enzymatic hydrolysis. To ensure consistency between different lots of these products, we have designed our packaging protocol so that each vial contains the specified amount of anhydrous bis-amide-rhodamine 110. In addition, the product may also contain a variable amount of water of hydration and a small amount of inorganic salts that do not affect its use. For several substrates, the lot-specific weight purity is indicated on the product label. For these products, the total amount of material in the container is:

$$\text{Total material (mg)} = \left[\frac{(\text{substrate (mg)} \times 100)}{\text{weight purity (\%)}} \right]$$




Figure 2. Upon proteinase cleavage, the nonfluorescent bisamide derivative of rhodamine 110 is converted first to the fluorescent monoamide and then to rhodamine 110, with a further increase in fluorescence.

Table 2. Rhodamine 110–based bis-peptide substrates.

Cat #	Proteinase Substrates *	Enzymes
R6504	(CBZ-Ala-Ala) ₂ -R110	elastase ²
R6506	(CBZ-Ala-Ala-Ala-Ala) ₂ -R110	elastase ³
R6508	(CBZ-Ala-Arg) ₂ -R110	trypsin ⁴
R6501	(CBZ-Arg) ₂ -R110	trypsin ^{4,5}
R6505	(CBZ-Ile-Pro-Arg) ₂ -R110	trypsin ^{6,7}
R6502	(CBZ-Phe-Arg) ₂ -R110	plasmin, cathepsin L ^{1,8,9}
R22124	(<i>p</i> -tosyl-Gly-Pro-Arg) ₂ -R110	thrombin ¹²

* SC = succinoyl; CBZ = benzyloxycarbonyl; *p*-tosyl = toluenesulfonyl; R110 = rhodamine 110; **1.** Biochemistry 38, 13906 (1999); **2.** Biol Chem Hoppe-Seyler 373, 547 (1992); **3.** Anal Chem 65, 2352 (1993); **4.** Biochem J 215, 253 (1983); **5.** Biochem J 209, 299 (1983); **6.** Biochim Biophys Acta 788, 74 (1984); **7.** Photochem Photobiol 44, 461 (1986); **8.** Glia 7, 183 (1993); **9.** Biol Chem Hoppe-Seyler 373, 433 (1992); **10.** J Biol Chem 274, 17484 (1999); **11.** Exp Cell Res 250, 203 (1999); **12.** Biomed Instrum Technol 30, 245 (1996).

Guidelines for Use

Materials Required but Not Provided

- DMSO or DMF
- Buffer such as Tris or HEPES, pH 7.5
- EDTA
- Ethanol or detergents for cell lysis

Preparing the Rhodamine 110 Stock Solutions

A 5–10 mM stock solution may be prepared in high-quality, anhydrous dimethylsulfoxide (DMSO) or dimethylformamide (DMF). This solution can be stored desiccated at 2–6°C or below for at least six months.

In Vitro Assay of Proteinase Activity

In vitro proteinase assays with the bisamide derivatives of rhodamine 110 are typically performed at 22°C in 10 mM Tris or HEPES buffer, pH 7.5; 15% (v/v) ethanol^{1,2} or detergents used to promote cell lysis are sometimes present in the assay buffer. The substrate stock solution in DMSO or DMF is diluted into the assay buffer just prior to initiating the reaction by addition of the enzyme preparation or cell lysate. If reaction conditions are chosen such that less than 15% of the substrate is hydrolyzed, then the increase in fluorescence is due solely to the production of the monoamide derivative of rhodamine 110.² Thus, under these conditions, the interpretation of kinetic data is not complicated by the bifunctionality of the bisamide substrate.

Intracellular Proteinase Assays

In preparation for intracellular proteinase assays with the rhodamine 110–based substrates, cells should be sedimented, resuspended in HEPES-buffered saline (HBS; 5 mM HEPES, 0.15 M NaCl, pH 7.35) containing 2 mM EDTA (HBS-EDTA) to a density of approximately 0.5–1 × 10⁷ cells/mL and stored at 4°C for no longer than two hours.^{5,6} Just prior to adding substrate, the cell suspension is diluted 100-fold in HBS-EDTA and incubated for approximately 20 minutes in the presence of 10 μM substrate. If desired, dead cells can then be identified by briefly counterstaining the sample with propidium iodide or other cell-impermeant nucleic acid dye before analysis by flow cytometry or fluorescence microscopy.

Spectral Properties

The peak excitation and emission wavelengths of rhodamine 110 are 498 nm and 521 nm, respectively. Rhodamine 110 is reported to have an extinction coefficient of $81,000 \text{ cm}^{-1}\text{M}^{-1}$ at 498 nm, with a quantum yield of 0.91.² The monoamide, CBZ-(L-arginine amide)-rhodamine 110, has an extinction coefficient of about $23,500 \text{ cm}^{-1}\text{M}^{-1}$ at 492 nm, with a quantum yield of 0.29.²

References

1. Biochem J 215, 253 (1983); 2. Biochem J 209, 299 (1983); 3. Biochim Biophys Acta 788, 74 (1984); 4. Photochem Photobiol 44, 461 (1986); 5. Glia 7, 183 (1993); 6. Biol Chem Hoppe-Seyler 373, 547 (1992); 7. Biol Chem Hoppe-Seyler 373, 433 (1992); 8. Nature 361, 274 (1993); 9. Biochemistry 38, 13906 (1999).

Product List Current prices may be obtained from our website or from our Customer Service Department.

Cat #	Product Name	Unit Size
R22124	rhodamine 110, bis-(p-tosyl-L-glycyl-L-prolyl-L-arginine amide).....	2 mg
R6479	rhodamine 110 (R110) *reference standard*.....	25 mg
R6501	rhodamine 110, bis-(CBZ-L-arginine amide), dihydrochloride (BZAR).....	5 mg
R6502	rhodamine 110, bis-(CBZ-L-phenylalanyl-L-arginine amide), dihydrochloride.....	5 mg
R6504	rhodamine 110, bis-(CBZ-L-alanyl-L-alanine amide).....	5 mg
R6505	rhodamine 110, bis-(CBZ-L-isoleucyl-L-prolyl-L-arginine amide), dihydrochloride (BZIPAR).....	5 mg
R6506	rhodamine 110, bis-(CBZ-L-alanyl-L-alanyl-L-alanyl-L-alanine amide).....	5 mg
R6507	rhodamine 110, bis-(CBZ-L-prolyl-L-arginine amide), dihydrochloride.....	5 mg
R6508	rhodamine 110, bis-(CBZ-L-alanyl-L-arginine amide), dihydrochloride.....	5 mg
R6513	rhodamine 110, bis-(t-BOC-L-leucyl-L-methionine amide).....	5 mg
R6560	rhodamine 110, 4-(chloromethyl)benzoyl amide, CBZ-L-arginine amide, hydrochloride (CMB-R110, CBZ-Arg).....	1 mg
R6577	rhodamine 110, bis-(L-phenylalanine amide), di(trifluoroacetic acid) salt.....	5 mg

Contact Information

Molecular Probes, Inc.
29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:
8:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:
Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:
8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probetech@invitrogen.com

Invitrogen European Headquarters
Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Assistance Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2006, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.