

# BacMam GFP Transduction Control \*BacMam 2.0\*

Catalog no. B10383

**Table 1.** Contents and storage information.

Material	Amount	Concentration	Storage	Stability
BacMam GFP Transduction Control *BacMam 2.0*	1 mL	$\sim 1 \times 10^8$ particles/mL	<ul style="list-style-type: none"> <li>• 2–6°C</li> <li>• Protect from light</li> <li>• DO NOT FREEZE</li> </ul>	When stored as directed, this kit is stable for 6 months.
<b>Approximate fluorescence excitation/emission maxima:</b> GFP: 488/512 in nm.				
<b>Number of assays:</b> Sufficient material is supplied to transduce approximately $5 \times 10^6$ cells based on the protocol below.				


## Introduction

The BacMam GFP Transduction Control lets you experience the power and convenience of the enhanced BacMam 2.0 technology, with a convenient readout of emerald GFP (emGFP). The extremely bright and photostable *Aequorea victoria* emGFP lights up the entire cell in this non-targeted form. Evaluate BacMam 2.0 in your cell model before advancing to the more specialized BacMam reagents.

BacMam technology is based on the use of an insect cell virus (baculovirus) to efficiently deliver and express genes in mammalian cells.<sup>1,2</sup> Transgenes under mammalian promoter elements are expressed, while baculoviral genes and their promoters are not recognized. As mammalian cells do not support replication of baculoviruses, transduction is extremely well tolerated and generally lacking in cytopathic effects, even at high virus levels. The inability of baculoviruses to replicate renders them safe as research reagents.

BacMam 2.0 greatly expands the efficiency and utility of this popular gene delivery platform.<sup>3–5</sup> Cell types previously not compatible with the technology (primary neurons), or cells that were poorly transduced with version 1.0 (some stem cells, CHO) can now be transduced quantitatively in a simple, one step process. The improved performance is due to inclusion of elements that greatly enhance transduction efficiency and expression levels: a pseudotyped capsid protein for more efficient cell entry and genetic elements (enhanced CMV promoter and Woodchuck Post-transcriptional Regulatory Element) that boost expression levels.

For more information on a range of BacMam-based reagents, including CellLight™ cell labeling reagents, Premo™ Biosensors, ion channel drug targets, and pathway analysis kits that facilitate the study of live cells, visit [www.invitrogen.com/bacmam](http://www.invitrogen.com/bacmam).




**Figure 1.** BacMam 2.0 workflow. Using the BacMam 2.0 to express genes is as simple as it is efficient:

1. Add the BacMam reagent directly to the cells.
2. Analyze the transduction efficiency the next day or freeze the cells for future use.

**Table 2.** Generalized comparison of the transduction efficiencies of BacMam 1.0 and BacMam 2.0 across different cell types. Gene expression can be increased by using the BacMam enhancer (Cat. no. B10107).

	<b>BacMam 1.0 No enhancer</b>	<b>BacMam 1.0 With enhancer</b>	<b>BacMam 2.0 No enhancer</b>	<b>BacMam 2.0 With enhancer</b>
<b>Cell lines</b>	+/-	++	++	++
<b>Primary cells</b>	+	++	++	++
<b>Stem cells</b>	-	++	++	++
<b>Neurons</b>	-	-	++	++
<b>Immortalized T-cells</b>	-	-	+/-	++
<b>Primary T-cells and B-cells</b>	-	-	-	-

- = not transduced  
 +/- = transduction efficiency ~10%  
 + = transduction efficiency <50%  
 ++ = transduction efficiency >50%



**Figure 2.** Comparison of the transduction efficiency of (1) BacMam 1.0 and (2) BacMam 2.0 GFP in T84 adenocarcinoma (A) and adipocytes-derived stem cells (ADSC) (B), respectively. Transduction conditions including virus titer, particles per cell, treatment time, and cell density were identical.

## Before Starting

---

### Materials Required but Not Provided

- Phosphate-buffered saline (PBS) without  $\text{Ca}^{2+}/\text{Mg}^{2+}$  (Cat. no. 10010)
- Hemacytometer and Trypan blue, or the Countess® Automated Cell Counter
- *Optional:* TrypLE™ Express dissociation enzyme (Cat. no. 12604-013)
- *Optional:* Freezing medium such as Recovery™ Cell Culture Freezing Medium (Cat. no. 12648-010)

### Guidelines for Working with BacMam Reagents

- The standard protocol is based on a 1 mL labeling volume for a 35-mm dish or 1 well of a 6-well culture plate, with cells about 70% confluent at time of transduction, and 30 BacMam particles per cell.
- For applications that require a larger number of cells, such as high-content screening (HCS), we recommend transducing the cells in a 10-cm dish or a T-75 flask and increasing the labeling volume to 10 mL with a proportionate increase in the volume of the virus.
- Following an overnight incubation under normal growth conditions, trypsinize and count the cells for distribution to appropriate plates at the desired cell number.
- If the transduction efficiency needs to be optimized, we suggest adjusting the following variables: MOI (from 10 to 100), cell density (80,000 to 200,000 cells/mL), temperature (room temperature for 1 or 2 hours before moving to the incubator), transduction volume, and incubation time. For some very sensitive cell types a relatively high MOI for one or two hours followed by medium removal has also been reported to be effective.
- The BacMam Enhancer (Cat. no. B10107) is generally not required for BacMam 2.0 reagents. However, its use has been shown to boost expression in some challenging cell types such as Jurkat.
- We recommend transducing the cells at a confluence of about 70% for best results.
- For first time users of BacMam reagents, we recommend exceptionally well-transduced cells like U-2 OS (ATCC® Number: HTB-96™).

## Experimental Protocols

---

The BacMam GFP Transduction Control is provided as a  $1 \times 10^8$  particles/mL solution. BacMam 2.0 reagents work best at a multiplicity of infection (MOI) between 10 and 50 in most cell types.

At a concentration of  $1 \times 10^8$  particles/mL, for every 10,000 cells:

- 1  $\mu$ L BacMam GFP Transduction Control is required for an MOI of 10
- 2  $\mu$ L BacMam GFP Transduction Control is required for an MOI of 20
- 5  $\mu$ L BacMam GFP Transduction Control is required for an MOI of 50

**Note:** Some cell types, such as neurons, may need a higher MOI.

### Protocol for Adherent Cells

#### Day 1

- 1.0 Harvest cells for subculture in complete medium at the desired density, e.g. 200,000 cells/mL
- 1.1 Calculate the appropriate volume of BacMam GFP Transduction Control for the number of cells to be plated.

$$\text{Volume of BacMam Reagent (mL)} = \frac{\text{number of cells} \times \text{desired MOI}}{1 \times 10^8 \text{ particles/mL}}$$

- 1.2 Mix the BacMam GFP Transduction Control reagent several times by inversion to ensure a homogenous solution. When using more than one BacMam reagent, pre-mix the reagents ahead of addition to the cells.
- 1.3 Add the appropriate volume of BacMam GFP Transduction Control reagent directly to the cells in complete cell medium and mix gently.
- 1.4 Incubate the cells at room temperature for 10 minutes.
- 1.5 Return the cells to the culture incubator and incubate overnight ( $\geq 16$  hours) for GFP expression.

**Note:** Alternatively, add the BacMam GFP Transduction Control reagent to established cultures in complete medium at 1% to 10% (vol/vol), mix gently, and return to the incubator.

#### Day 2

- 1.6 Image the cells using the appropriate filters for GFP/FITC.
- 1.7 Replace medium if the cells are to be maintained in culture.

#### Notes:

- The BacMam GFP Transduction Control reagent was developed for use in live-cell studies. Should you prefer fixed cell analysis, the fluorescence from GFP has been shown to be compatible with fixation with 4% formaldehyde and permeabilization with 0.1% Triton<sup>®</sup> X-100.
- Increased transduction efficiencies are observed if the virus is added directly after plating, or while the cells are still in suspension.
- In most cases, the BacMam particles do not need to be removed, although a medium exchange the day following transduction may be desirable.
- Avoid exposing the BacMam reagent to light for more than 10 minutes.

## Frequently Asked Questions

---

**Q:** Will BacMam 2.0 transduce my cells?

**A:** The first generation BacMam reagents were shown to efficiently transduce over 90 cell types, including stable cell lines and primary cells. For the most up to date list of cells and transduction efficiencies, refer to [www.invitrogen.com/BacMamCompatible](http://www.invitrogen.com/BacMamCompatible). BacMam 2.0 GFP is based on a modified and enhanced BacMam vector with superior transduction properties. For instance, it is now possible to efficiently transduce primary neurons and stem cells.

**Q:** How long does expression last?

**A:** The duration of transgene expression depends on many factors, including transduction levels, cell division rates, mRNA and protein stability. In most transformed cell lines such as HeLa and CHO expression lasts about 5 days. In cells that divide more slowly or show contact inhibition, such as some stem cells, primary cells, and neurons, we have observed bright staining and transgene expression for more than two weeks. For non-dividing, terminally differentiated cells we have observed expression for 2 to 4 weeks.

**Q:** Can I transduce with more than one BacMam reagent at a time?

**A:** Yes, this is one of the advantages of the system. For instance the Premo™ FUCCI Cell Cycle Sensor and the BacMam Kv7.2/7.3 Potassium Ion Channel reagent are based on optimized ratios of two BacMam constructs that give rise to a two-color cell cycle sensor and a functional heterotetrameric K channel, respectively.

**Q:** Will BacMam transduction hurt my cells?

**A:** BacMam transduction is generally exceptionally well tolerated, even at very high number of viral particles to cell ratios (>1,000). However, we have occasionally observed apparent cytotoxic effects by some BacMam reagents at very high virus levels; this may be due to the nature of the transgene. For this reason, we recommend using no more virus than is needed.

**Q:** If I freeze my cells after transduction, how long can I store them without reducing expression levels?

**A:** Our data show that transduced cells can be stored at  $-80^{\circ}\text{C}$  for several months without reducing the level of transgene expression.

**Q:** Can transduction be optimized if my cells are difficult to transduce?

**A:** Yes. Try varying virus-to-cell ratio (MOI), incubation temperature and duration, and cell density (if adherent cells are transduced). For adherent cells, we recommend a confluence of about 70%. Media have also been shown to affect transduction efficiency; if your cells do not tolerate the recommended PBS without  $\text{Mg}^{2+}$  and  $\text{Ca}^{2+}$ , we suggest using RPMI1640.

**Q:** Can a cell be transduced more than once?

**A:** Yes. Because transduction is so well tolerated, you can readily add more BacMam reagent after a few days if expression levels need to be boosted or if a different BacMam-based assay is needed.

**Q:** It's a virus—is it safe to use?

**A:** Yes. Baculoviruses are insect viruses that do not replicate in mammalian cells and are generally used under the safety precautions common for standard cell-based reagents.

## References

---

1. Nature Biotechnol 23, 567 (2005); 2. Expert Opin Drug Discov 2, 1669 (2007); 3. Biochem Biophys Res Comm 349, 1220 (2006); 4. J Biotechnol 131, 1 (2007); 5. Mol Ther 17, 1585 (2009).

## Product List Current prices may be obtained from our website or from our Customer Service Department.

---

Cat. no.	Product Name	Unit Size
B10383	BacMam GFP Transduction Control *BacMam 2.0*	1 mL
<b>Related Products</b>		
A6455	anti-green fluorescent protein, rabbit serum (anti-GFP, serum)	100 µL
A11122	anti-green fluorescent protein, rabbit IgG fraction (anti-GFP, IgG) *2 mg/mL*	100 µL
B10019	BacMam-hERG *for 10 microplates* *BacMam 1.0*	1 kit
B10107	BacMam Enhancer Kit	1 kit
B10146	BacMam Kir2.1 *for 10 microplates* *BacMam 1.0*	1 kit
B10147	BacMam Kv7.2 and Kv7.3 *for 10 microplates* *BacMam 1.0*	1 kit
C10106	Cellular Lights™ Tubulin-GFP *BacMam 1.0*	1 kit
C10112	Cellular Lights™ Tubulin-RFP *BacMam 1.0*	1 kit
C10126	Cellular Lights™ Actin-GFP *BacMam 1.0*	1 kit
C10127	Cellular Lights™ Actin-RFP *BacMam 1.0*	1 kit
O10100	Organelle Lights™ Lysosomes-RFP *BacMam 1.0*	1 kit
O10104	Organelle Lights™ Endosomes-GFP *BacMam 1.0*	1 kit
O36210	Organelle Lights™ Mito-GFP *BacMam 1.0*	1 kit
O36231	Organelle Lights™ Endosomes-RFP *BacMam 1.0*	1 kit
P36232	Premo™ FUCCI Cell Cycle Sensor *BacMam 1.0*	1 kit
P36235	Premo™ Autophagy Sensor LC3B-GFP *BacMam 2.0*	1 kit
P36236	Premo™ Autophagy Sensor LC3B-RFP *BacMam 2.0*	1 kit

## Contact Information

---

### **Molecular Probes, Inc.**

29851 Willow Creek Road  
Eugene, OR 97402  
Phone: (541) 465-8300  
Fax: (541) 335-0504

### **Customer Service:**

6:00 am to 4:30 pm (Pacific Time)  
Phone: (541) 335-0338  
Fax: (541) 335-0305  
probestorder@invitrogen.com

### **Toll-Free Ordering for USA:**

Order Phone: (800) 438-2209  
Order Fax: (800) 438-0228

### **Technical Service:**

8:00 am to 4:00 pm (Pacific Time)  
Phone: (541) 335-0353  
Toll-Free (800) 438-2209  
Fax: (541) 335-0238  
probestech@invitrogen.com

### **Invitrogen European Headquarters**

Invitrogen, Ltd.  
3 Fountain Drive  
Inchinnan Business Park  
Paisley PA4 9RF, UK  
Phone: +44 (0) 141 814 6100  
Fax: +44 (0) 141 814 6260  
Email: euroinfo@invitrogen.com  
Technical Services: eurotech@invitrogen.com

**For country-specific contact information,  
visit [www.invitrogen.com](http://www.invitrogen.com).**

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

### **Limited Use Label License No. 223: Labeling and Detection Technology**

The manufacture, use, sale or import of this product may be subject to one or more patents or pending applications owned or licensed by Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity) in a manner consistent with the accompanying product literature. The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the most restrictive terms apply. Life Technologies Corporation will not assert a claim against the buyer of infringement of patents that are owned or controlled by Life Technologies Corporation and/or Molecular Probes, Inc. which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

### **Limited Use Label License No. 21: Bac-to-Bac® and Bac-to-Bac® HT**

This product is the subject of U.S. Patent No. 5,348,886. This product is sold under patent license from Monsanto for research purposes only and no license for commercial use is included. Requests for licenses for commercial manufacture or use should be directed to Director, Monsanto Corporate Research, 800 N. Lindbergh, St. Louis, Missouri 63167.

### **Limited Use Label License No: 127 GFP with Heterologous Promoter**

This product and its use is the subject of one or more of U.S. Patent Nos. 5,491,084 and 6,146,826, and foreign equivalents. This product is sold under license from Columbia University. Rights to use this product are limited to research use only, and expressly exclude the right to manufacture, use, sell or lease this product for use for measuring the level of toxicity for chemical agents and environmental samples in cells and transgenic animals. No other rights are conveyed. Not for human use or use in diagnostic or therapeutic procedures. Inquiry into the availability of a license to broader rights or the use of this product for commercial purposes should be directed to Columbia Innovation Enterprise, Columbia University, Engineering Terrace-Suite 363, New York, New York 10027.

### **Limited Use Label License No: 198 Fluorescent Protein Products**

This product and its use is the subject of one or more of U.S. Patent Nos. 5,777,079, 6,066,476, and 6,319,669 and foreign equivalents. Any use of this product by a commercial entity requires a separate license from either GE Healthcare or Life Technologies Corporation. For information on obtaining a commercial license to use this product, please refer to the contact information located at the bottom of this statement. No rights are conveyed to modify or clone the gene encoding GFP contained in this product. For information on licensing, contact Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008 ; Phone (760) 603-7200 or e-mail: outlicensing@lifetech.com.

### **Limited Use Label License No: 267 Mutant GFP Products**

This product and its use is the subject of one or more of U.S. Patent Nos. 6,090,919, 5,804,387, 5,994,077, and foreign equivalents.

### **Limited Use Label License No: 272 Humanized GFP**

This product is the subject of one or more of U.S. Patent Numbers 5,786,464, 5,795,737, 5,874,304, and 6,114,148 and foreign equivalents licensed by Life Technologies Corporation. This product is sold for research use only. Not for therapeutic or diagnostic use in humans.

### **Limited Use Label License No: 306 Baculovirus Vectors**

Certain methods that utilize the product associated with this limited use label license are covered by U.S. Patent Nos. 5,731,182 and 5,871,986. This product is for research use only by those researchers in laboratories of academic, government, industrial and/or clinical institutions engaged in the investigation of biological or biochemical processes, or research and development of biological products. This product is not to be used in the manufacture, use or sale of human or animal diagnostic, therapeutic or prophylactic products.

**Limited Use Label License No: 308 WPRE Element**

This product contains the Woodchuck Post-transcriptional Regulatory Element ("WPRE") which is the subject of intellectual property owned by The Salk Institute for Biological Studies, and licensed to Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; and/or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. In addition, any use of WPRE outside of this product or the product's authorized use requires a separate license from the Salk Institute. Life Technologies will not assert a claim against the buyer of infringement of patents owned by Life Technologies and claiming this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product or for a Commercial Purpose. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008, Phone (760) 603-7200. Fax (760) 602-6500, or The Salk Institute for Biological Studies, 10010 North Torrey Pines Road, La Jolla, CA 92037, Attn.: Office of Technology Management, Phone: (858) 453-4100 extension 1275, Fax: (858) 546-8093.

**Limited Use Label License No: 332 BacMam Virus Use**

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer solely in accordance with the accompanying product literature or manual. Purchase of this product does not convey a license to expand, amplify, or otherwise propagate the provided viral particles or to otherwise modify or alter the virus by any means.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Triton® is trademark of Union Carbide.

Copyright 2010, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.

**For research use only. Not intended for any animal or human therapeutic or diagnostic use.**