

CellEvent™ Caspase-3/7 Green Detection Reagent

Catalog no. C10423

Table 1 Contents and storage

Material	Amount	Concentration	Storage	Stability
CellEvent™ Caspase-3/7 Green Detection Reagent	100 µL	2.0 mM solution in DMSO	<ul style="list-style-type: none"> • ≤-20°C • Protect from light 	When stored as directed, the product is stable for 3 months from the date of receipt.
Approximate fluorescence excitation and emission maxima: 502/530 nm.				

Introduction

CellEvent™ Caspase-3/7 Green Detection Reagent is a novel fluorogenic substrate for activated caspases 3 and 7. The reagent consists of a four amino acid peptide (DEVD) conjugated to a nucleic acid binding dye. This cell-permeant substrate is intrinsically non-fluorescent, because the DEVD peptide inhibits the ability of the dye to bind to DNA. After activation of caspase-3 or caspase-7 in apoptotic cells, the DEVD peptide is cleaved, enabling the dye to bind to DNA and produce a bright, fluorogenic response with an absorption/emission maxima of ~502/530 nm (Figure 1).

There are several benefits of the CellEvent™ Caspase-3/7 Green Detection Reagent compared to other methods of detecting activated caspase 3/7. One important advantage of this assay is that no wash steps are required; this helps preserve the fragile apoptotic cells typically lost during wash steps. To use the CellEvent™ Caspase-3/7 Green Detection Reagent, add the substrate to your cells in complete growth medium or buffer, incubate for 30 minutes, and image (Figure 2). Apoptotic cells with activated caspase-3/7 show bright green nuclei, while cells without activated caspase 3/7 exhibit minimal fluorescence signal (Figure 3, page 2).

Figure 1 Fluorescence excitation and emission spectra of CellEvent™ Caspase-3/7 Green Detection Reagent after reaction with activated caspase-3 or 7, bound to DNA





Figure 2 Workflow for the CellEvent™ Caspase-3/7 Green Detection Reagent



This robust assay is highly specific for caspase-3/7 activation (Figure 4) and can be used to monitor caspase-3 or -7 activation with live-cell fluorescence imaging. Because the cleaved reagent labels nuclei of caspase 3/7–positive cells, this stain can also provide information on nuclear morphology, including condensed nuclei typical of late-stage apoptosis. Additionally, the fluorescent signal from CellEvent™ Caspase-3/7 Detection Reagent survives formaldehyde fixation and detergent permeabilization. This provides flexibility in assay workflow and extends the multiplexability of this probe for the detection other proteins of interest using immunocytochemistry.

In addition to traditional fluorescence microscopy, the CellEvent™ Caspase-3/7 Green Detection Reagent has been validated for high-content imaging and analysis. The drastic change in fluorescence between normal and apoptotic cells within a population provides an excellent assay window, and the z-factor value indicates the reagent is robust enough for use in high-content imaging assays (Figure 5).

Figure 3 HeLa cells were loaded with 7.5 μM CellEvent™ Caspase-3/7 Green Detection Reagent then treated with 0.5 μM staurosporine (right) or vehicle control (left) for 4 hours. Staurosporine-induced apoptosis was detected with CellEvent™ Caspase-3/7 Green Detection Reagent. Apoptotic cells (right) fluoresce bright green, while non-apoptotic control cells do not show any signal (left).




Figure 4 To determine the specificity of CellEvent™ Caspase-3/7 Green Detection Reagent, HeLa cells were treated with 0.5 μM staurosporine in the presence or absence of the Caspase 3/7 Inhibitor 1 (EMD Chemicals) at 0–30 μM for 4 hours. Cells were then labeled with 5 μM CellEvent™ Caspase-3/7 Green Detection Reagent, followed by Hoechst 33342 for 15 minutes in complete medium. Images were acquired and analyzed on a Thermo Fisher Cellomics ArrayScan® VTI. Quantitative analysis revealed a decrease in the percent of cells positive for active caspase 3/7 with increasing concentrations of inhibitor.





Figure 5 To generate a dose response curve, U-2 OS cells were plated onto a 96-well plate and treated with 0–0.75 μM staurosporine or vehicle control for 24 hours. Cells were then labeled with 7.5 μM CellEvent™ Caspase-3/7 Green Detection Reagent for 30 minutes at 37°C followed by Hoechst 33342 stain. Cells were analyzed on a Thermo Fisher Cellomics ArrayScan® VTI. The percent of cells positive for active caspase 3/7 was determined, and the EC50 was calculated.



Before Starting

Materials required but not provided

- Cells
- Complete medium or buffer suitable for the cell type used
- *Optional:* Fixative (i.e., 3.7% formaldehyde in PBS)

Caution

DMSO is known to facilitate the entry of organic molecules into tissues. Handle reagents containing DMSO using equipment and practices appropriate for the hazards posed by such materials. Dispose of the reagents in compliance with all pertaining local regulations. In case of contact with eyes, rinse immediately and thoroughly with water and seek medical advice. Always wear suitable laboratory protective clothing and gloves when handling this reagent.

Experimental Protocols

The following protocol was developed with HeLa and U-2 OS cells with an optimized concentration of 5 μM , but it can be adapted for any cell type. Growth medium, cell density, cell type variations, and other factors may influence labeling. In initial experiments, we recommend testing a concentration range of CellEvent™ Caspase-3/7 Green Detection Reagent to determine the optimal conditions for your model.

Staining protocol

1. Add CellEvent™ Caspase-3/7 Green Detection Reagent at a final concentration of between 2–8 μM to the sample and/or appropriately induced cells and incubate for 30 minutes at 37°C. We recommend initial testing between 2–10 μM of CellEvent™ Caspase-3/7 Green Detection Reagent; however, optimal concentration may be more or less depending on model.

Note: It is best to prepare an intermediate dilution of CellEvent™ Caspase-3/7 Detection Reagent in complete medium, mix it well by pipetting, and then add the diluted solution to the cells so that the final concentration of the reagent on the cells is 5 μM . Incubation can also be performed successfully at room temperature (25°C).

2. *Optional:* Cells can be preserved with a formaldehyde-based fixative at this stage. Fixation with 3.7% formaldehyde for 15 minutes is recommended, but this can be altered based on the cell type.
3. *Optional:* Cells can be stained with a nuclear or other counterstain at this step.
4. *Optional:* To stabilize and prolong the signal, ProLong® Gold antifade reagent (Cat. no. P36934) can be used for ultimate overnight mounting. For quick mounting, SlowFade® Gold antifade reagent (Cat. no. S36937) can be used.
6. Image the cells using the appropriate instrument filter sets such as those used for FITC and the Alexa Fluor® 488 dye. The excitation/emission maxima for the CellEvent™ Caspase-3/7 Green Detection Reagent is 502/530 (Figure 1, page 1).

Product List Current prices may be obtained at www.invitrogen.com or from our Customer Service Department.

Catalog no.	Product name	Unit Size
C10423	CellEvent™ Caspase-3/7 Green Detection Reagent *2 mM solution in DMSO*	100 µL
C10422	CellROX™ Deep Red Reagent *for oxidative stress detection*	5 × 50 µL
Related Products		
C10422	CellROX™ Deep Red Reagent *for oxidative stress detection*	5 × 50 µL
H3570	Hoechst 33342, trihydrochloride, trihydrate *10 mg/mL solution in water*	10 mL
T668	tetramethylrhodamine, methyl ester, perchlorate (TMRM)	25 mg
L10382	LC3B Antibody Kit for Autophagy *rabbit polyclonal LC3B* *includes autophagosome inducer*	1 kit

Contact Information

Corporate Headquarters

5791 Van Allen Way
Carlsbad, CA 92008
USA
Phone: +1 760 603 7200
Fax: +1 760 602 6500
Email: techsupport@lifetech.com

European Headquarters

Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF
UK
Phone: +44 141 814 6100
Toll-Free Phone: 0800 269 210
Toll-Free Tech: 0800 838 380
Fax: +44 141 814 6260
Tech Fax: +44 141 814 6117
Email: euroinfo@invitrogen.com
Email Tech: euotech@invitrogen.com

Japanese Headquarters

LOOP-X Bldg. 6F
3-9-15, Kaigan
Minato-ku, Tokyo 108-0022
Japan
Phone: +81 3 5730 6509
Fax: +81 3 5730 6519
Email: jpinfo@invitrogen.com

Additional international offices are listed at
www.invitrogen.com

These high-quality reagents and materials must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Read the Safety Data Sheet provided for each product; other regulatory considerations may apply.

Web Resources

Visit the Invitrogen website at www.invitrogen.com for:

- Technical resources, including manuals, vector maps and sequences, application notes, Meds, FAQs, formulations, citations, handbooks, etc.
- Complete technical support contact information
- Access to the Invitrogen Online Catalog
- Additional product information and special offers

SDS

Safety Data Sheets (SDSs) are available at www.invitrogen.com/sds.

Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to www.invitrogen.com/support and search for the Certificate of Analysis by product lot number, which is printed on the product packaging (tube, pouch, or box).

Limited Warranty

Invitrogen (a part of Life Technologies Corporation) is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Support Representatives.

All Invitrogen products are warranted to perform according to specifications stated on the certificate of analysis. The Company will replace, free of charge, any product that does not meet those specifications. This warranty limits the Company's liability to only the price of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Support Representatives.

Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Limited Use Label License No. 358: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial services of any kind, including, without limitation, reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

©2011 Life Technologies Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.