

Concanavalin A Conjugates

Table 1. Contents and storage information.

Product	Amount	Storage	Stability
Texas Red®, fluorescein, and tetramethylrhodamine conjugates of concanavalin A	10 mg, lyophilized	<ul style="list-style-type: none"> • ≤-20°C • Desiccate • Protect from air 	When stored as directed, products are stable for at least 1 year.
Oregon Green®, Marina Blue®, and Alexa Fluor® conjugates of concanavalin A	5 mg, lyophilized		
Alexa Fluor® conjugate of succinylated concanavalin A			
Approximate fluorescence excitation/emission maxima: See Table 2			

Introduction

Lectins are versatile probes for detecting glycoconjugates in histochemical and flow cytometric applications and for localizing glycoproteins in gels. Concanavalin A selectively binds to α -mannopyranosyl and α -glucopyranosyl residues. In neutral and alkaline solutions, concanavalin A exists as a tetramer with a molecular weight of approximately 104,000 daltons. In acidic solutions (pH below 5.0), concanavalin A exists as a dimer.

When concanavalin A is succinylated with succinic anhydride, it is irreversibly converted to a dimer that retains the same sugar-binding specificity as the parent lectin. However, succinyl-concanavalin A has a profile of biological activities different from the tetrameric form. In contrast to the tetramer, succinyl-concanavalin A does not induce capping of cell surface glycoprotein receptors, inhibit capping of cell surface immunoglobulin receptors, or strongly agglutinate erythrocytes or spleen cells. The mitogenic effect of succinyl-concanavalin A is similar to that of the native lectin, although it is mitogenic over a significantly wider range of concentrations than the tetramer.^{1,2}

Invitrogen offers a broad selection of concanavalin A conjugates. Table 2 provides a summary of these products, including their peak excitation and emission wavelengths (where applicable).

Table 2. Fluorescent concanavalin A conjugates and their spectral characteristics.

Cat #	Label	Ex *	Em *
C11254	Alexa Fluor® 350	346	442
C827	Fluorescein	494	518
C11252, C21401 †	Alexa Fluor® 488	495	519
C6741	Oregon Green® 488	496	524
C860	Tetramethylrhodamine	555	580
C11253	Alexa Fluor® 594	590	617
C825	Texas Red®	595	615
C21402	Alexa Fluor® 633	632	647
C21421	Alexa Fluor® 647	650	668

* Approximate fluorescence excitation (Ex) and emission (Em) maxima, in nm. Complete spectra for these dyes are available at our website (probes.invitrogen.com). † The concanavalin A of this conjugate is succinylated.

Before You Begin

Preparing the Stock Solutions

Stock solutions can be made at 1–5 mg/ml in 0.1 M sodium bicarbonate (approximate pH 8.3). In most cases a small percentage of the conjugate will remain as a visible aggregate in solution. Store solutions at 2–6°C with the addition of 2 mM sodium azide; for longer storage, divide the solution into aliquots and freeze at ≤–20°C. PROTECT FROM LIGHT, AVOID REPEATED FREEZING AND THAWING.

Guideline for Use

It is a good practice to centrifuge the protein conjugate solution briefly in a microcentrifuge before use; only the supernatant should then be added to the experiment. This step will eliminate any protein aggregates that may have formed during storage, thereby reducing nonspecific background staining.

Because staining protocols vary with application, the appropriate dilution of concanavalin A conjugate should be determined empirically. Typically, a final concentration of 50–200 µg/mL is satisfactory for immunohistochemical applications.

References

1. *Biochemistry* 32, 5116 (1993); 2. *Proc Nat Acad Sci USA* 70, 1012 (1973).

Product List Current prices may be obtained from our website or from our Customer Service Department.

Cat #	Product Name	Unit Size
C11254	concanavalin A, Alexa Fluor® 350 conjugate	5 mg
C11252	concanavalin A, Alexa Fluor® 488 conjugate	5 mg
C21401	concanavalin A, succinylated, Alexa Fluor® 488 conjugate	5 mg
C11253	concanavalin A, Alexa Fluor® 594 conjugate	5 mg
C21402	concanavalin A, Alexa Fluor® 633 conjugate	5 mg
C21421	concanavalin A, Alexa Fluor® 647 conjugate	5 mg
C827	concanavalin A, fluorescein conjugate	10 mg
C6741	concanavalin A, Oregon Green® 488 conjugate	5 mg
C860	concanavalin A, tetramethylrhodamine conjugate	10 mg
C825	concanavalin A, Texas Red® conjugate	10 mg

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probestech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2006, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.