

Fluo-4 Calcium Imaging Kit

Catalog no. F10489

Table 1 Contents and storage

Material	Amount	Storage	Stability
Fluo-4, AM, 1000X in DMSO (Component A)	50 µL	<ul style="list-style-type: none"> • 2–8°C • Dessicate • Protect from light • DO NOT FREEZE 	When stored as directed the product is stable for at least 1 year.
PowerLoad™ Concentrate, 100X (Component B)	500 µL	<ul style="list-style-type: none"> • 2–8°C • DO NOT FREEZE 	
Neuro Background Suppressor, 10X (Component C)	5 mL	<ul style="list-style-type: none"> • 2–8°C • Protect from light • DO NOT FREEZE 	
Probenecid, water soluble (Component D)	77 mg	<ul style="list-style-type: none"> • ≤25°C • Desiccate 	

Number of assays: Sufficient material is supplied for 25 assays based on the protocol below.

Approximate fluorescence excitation and emission maxima: 494/506 nm (see Figure 1, page 2).

K_d for Ca²⁺: 335 nM

Introduction

Changes in calcium gradients have been demonstrated to be important for several biological processes. As the intracellular and extracellular calcium gradient is vast (10–50 nM and 1–3 mM, respectively), small changes in intracellular calcium [Ca²⁺] can result in large cellular modifications.¹ Visible light-excitable calcium indicators have been established as important tools for signal transduction and cell-based pharmacological testing. These small chemical entities (e.g., Fura-2, Indo, Fluo-3, and Fluo-4) display high sensitivity and a large fluorescent increase upon calcium binding.

The Fluo-4 Calcium Imaging Kit has been formulated, optimized, and contains all the necessary components for the detection of calcium flux by imaging applications. In addition to Fluo-4, AM, the kit contains PowerLoad™ Concentrate, a 10X Neuro Background Suppressor solution, and Probenecid.

For easy cell loading, the Fluo-4 Calcium Kit contains PowerLoad™ Concentrate. Due to its unique nature, PowerLoad™ solution can be used in the presence of complete culture medium, thus reducing the negative effects of replacing medium or loading in serum-free medium.


For Research Use Only. Not for use in diagnostic procedures.

Baseline autofluorescence caused by components within the growth medium can be greatly reduced by the addition of the Neuro Background Suppressor solution. The Neuro Background Suppressor solution has been specifically formulated for use with neuronal cells and will not cause osmotic shock. Additionally, the Neuro Background Suppressor solution has been used successfully with many different cell types to efficiently suppress background fluorescence without sacrificing the specific cellular fluorescence generated in the assay.

A proprietary, water-soluble Probenecid (which is commonly used to inhibit cellular transport and thus reduce the baseline signal) is supplied with the Fluo-4 Calcium Imaging Kit.

- Properties**
- The emission from the calcium-bound Fluo-4 dye (excitation 494 nm/emission 506 nm) can be detected with standard FITC filters.
 - Upon binding calcium the fluorescence intensity increases >100 fold.
 - K_d for Ca^{2+} in buffer: ~335 nM

Figure 1 Fluorescence excitation and emission spectra of Ca^{2+} -saturated Fluo-4, AM in pH 7.2 buffer.



Before You Begin

Materials Required but Not Provided

- Cell line and culture medium of choice
- For optimal performance, we recommend Live Cell Imaging Solution (LCIS) (Cat. no. A14291DJ), a physiological buffered saline equivalent to Ringer's Solution.
- A sterile-filtered, 2 M Glucose Stock Solution will also be required for LCIS and many other solution formats to support cell health in longer term (hours) experiments, as well as studies with primary or differentiated neural cell types.
- Buffered and pH adjusted physiological saline solution for dye loading and imaging. Phosphate buffered saline (PBS), Hank's balanced salt solution (HBSS), Ringer's solution, or Krebs' solution are acceptable depending on cell type, and should be pH or osmotically adjusted for signaling or long term studies.

**Prepare 20 mM Glucose +
Live Cell Imaging Solution**

Dilute 2 M Glucose Stock Solution 1:100 into LCIS for a final glucose concentration of 20 mM. Keep this solution clean and free of contaminants to prevent bacterial, fungal, or yeast growth once glucose has been added.

Prepare Probenecid

Dissolve the contents of the Probenecid (Component D) vial into 1 mL of LCIS to prepare the 100X Probenecid stock solution. Use the solution the same day or store at $\leq -20^{\circ}\text{C}$ for up to 6 months.

Experimental Protocols

The protocol below provides instructions for performing the calcium flux detection assay using adherent cells grown in a 35-mm dish with 2 mL of culture medium.

- 1.1** To a 15-mL tube, add the following reagents in the order listed below to prepare fresh Fluo-4, AM Loading Solution:

100X PowerLoad™ concentrate (Component B)	100 μL
Fluo-4, AM, 1000X (Component A)	10 μL

Vortex to mix.

Physiological buffer of choice or 20 mM Glucose Stock + LCIS	10 mL
--	-------

Invert the tube to mix.

Optional: Add 100 μL of 100X Probenecid stock solution to prevent extrusion of cytosolic dye by anion pumps, which can decrease loading efficiency on some cell types.

- 1.2** Remove medium from adherent cells and wash cells once in physiological buffer of choice or LCIS.
- 1.3** Add 2 mL of Fluo-4, AM Loading Solution (from step 1.1) to cells, and incubate cells at 37°C for 15–30 minutes, followed by 15–30 minutes at room temperature.
- Note:** Cells may be loaded at room temperature for as long as 60 minutes.
- 1.4** Remove Fluo-4, AM Loading Solution, and wash cells once in physiological buffer of choice or LCIS.
- 1.5** Add 2 mL of physiological buffer of choice or 20 mM Glucose Stock in LCIS. Cells are now ready for live-cell imaging.

Optional: To suppress background fluorescence, add 1:10 diluted Neuro Backdrop Background Suppressor solution (Component C).

Imaging

Standard FITC settings may be used to visualize the cytosolic staining of the Fluo-4, AM dye. For a positive control, add Ionomycin (a calcium ionophore, Cat. no. I24222) to a final concentration of 10 μM for a large increase in cytosolic calcium concentration.

References

1. Nat Rev Mol Cell Biol 7, 517 (2003); 2. Anal Biochem 291, 175 (2001); 3. Biochem J 356, 345 (2001).

Product List

Current prices may be obtained from our website or from our Customer Service Department.

Cat. no.	Product Name	Unit Size
F10489	Fluo-4 Calcium Imaging Kit	1 kit
Related Products		
A14291DJ	Live Cell Imaging Solution	500 mL
I24222	Ionomycin, calcium salt	1 mg
14025-092	Hanks' Balanced Salt Solution (HBSS) (1X), liquid	500 mL
15630-106	HEPES Buffer Solution (1 M)	20 mL

Purchaser Notification

Corporate Headquarters
5791 Van Allen Way
Carlsbad, CA 92008
USA
Phone: +1 760 603 7200
Fax: +1 760 602 6500
Email: techsupport@lifetech.com

European Headquarters
Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF
UK
Phone: +44 141 814 6100
Toll-Free Phone: 0800 269 210
Toll-Free Tech: 0800 838 380
Fax: +44 141 814 6260
Tech Fax: +44 141 814 6117
Email: euroinfo@invitrogen.com
Email Tech: eurotech@invitrogen.com

Japanese Headquarters
LOOP-X Bldg. 6F
3-9-15, Kaigan
Minato-ku, Tokyo 108-0022
Japan
Phone: +81 3 5730 6509
Fax: +81 3 5730 6519
Email: jpinfo@invitrogen.com

Additional international offices are listed at
www.lifetechnologies.com

These high-quality reagents and materials must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Read the Safety Data Sheet provided for each product; other regulatory considerations may apply.

Obtaining Support

For the latest services and support information for all locations, go to www.lifetechnologies.com.

At the website, you can:

- Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities
- Search through frequently asked questions (FAQs)
- Submit a question directly to Technical Support (techsupport@lifetech.com)
- Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents
- Obtain information about customer training
- Download software updates and patches

SDS

Safety Data Sheets (SDSs) are available at www.lifetechnologies.com/sds.

Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to www.lifetechnologies.com/support and search for the Certificate of Analysis by product lot number, which is printed on the product packaging (tube, pouch, or box).

Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

Disclaimer

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

Important Licensing Information

These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

The trademarks mentioned herein are the property of Life Technologies Corporation and/or its affiliate(s) or their respective owners.

©2013 Life Technologies Corporation. All rights reserved.

