

Premo™ Autophagy Sensors (LC3B-FP) *BacMam 2.0*

Catalog nos. P36235, P36236

Table 1. Contents and storage information.

Material	P36235	P36236	Concentration	Storage*
LC3B-FP (Component A)	1 mL of LC3B-GFP	1 mL of LC3B-RFP	$\sim 1 \times 10^8$ viral particles/mL	<ul style="list-style-type: none"> • 2–6°C • Protect from light • DO NOT FREEZE
Control LC3B(G120A)-FP (Component B)	200 μ L of LC3B(G120A)-GFP	200 μ L of LC3B(G120A)-RFP	$\sim 1 \times 10^8$ viral particles/mL	
Chloroquine diphosphate (Component C)	1 mL	1 mL	30 mM aqueous solution	<ul style="list-style-type: none"> • $\leq 25^\circ\text{C}$ • Protect from light
Approximate fluorescence excitation/emission maxima: GFP: 485/520 in nm; RFP: 555/584 in nm.				
*Stability: When stored as directed, this kit is stable for at least 6 months.				

Introduction

Autophagy describes the segregation and delivery of cytoplasmic cargo, including proteins and organelles, for degradation by hydrolytic enzymes in “autophagolysosomes”, also referred to as “autolysosomes”. Although first described in 1963, it has only been in the past decade that this pathway has been the subject of intense research to gain further insight into the role basal autophagy plays in cell homeostasis and development. Efforts are also directed to further elucidate the role of induced autophagy as a cell survival response to stress, microbial infection, and disease (e.g., neurodegeneration, cancer).^{1–3}

The LC3B protein plays a critical role in autophagy. Normally, this protein resides in the cytosol, but following cleavage and lipidation with phosphatidylethanolamine, LC3B associates with the phagophore. This localization can be used as a general marker for autophagic membranes (Figure 1).

The Premo™ Autophagy Sensor combines the selectivity of an LC3B-fluorescent protein (FP) chimera with the transduction efficiency of the BacMam 2.0 technology. BacMam reagents (insect **B**aculovirus with a **M**ammalian promoter) are safe to handle (Biosafety Level 1) because they are non-replicating in mammalian cells. They are also non-cytotoxic and ready-to-use. Unlike expression vectors, BacMam reagents enable titratable and reproducible expression and offer high co-transduction efficiency. Multiple BacMam reagents can therefore be readily used in the same cell. Recent improvements made to the BacMam system, BacMam 2.0, enable efficient transduction in a wider variety of cells, including neurons and neural stem cells (NSCs). The two-step protocol for imaging autophagy involves simply adding the BacMam LC3B-FP to the cells and incubating them overnight for protein expression.

Each Premo™ Autophagy Sensor Kit includes the BacMam LC3B-FP, a control BacMam LC3B(G120A)-FP, and chloroquine diphosphate to induce autophagosomes (Figure 2). The mutation on the control BacMam LC3B prevents its cleavage and subsequent lipidation during normal autophagy, and protein localization remains cytosolic and diffuse. Following treatment with chloroquine diphosphate, normal autophagic flux is disrupted and autophagosomes accumulate as a result of the increased lysosomal pH that inhibits lysosomal fusion with the autophagosomes.




Figure 1. Schematic depiction of the autophagy pathway in a eukaryotic cell. The first step involves the formation and elongation of the isolation membranes or phagophore. The second step entails the expansion and sequestering of the cytoplasm and formation of the double-membrane autophagosome and includes the association of the cytosolic LC3 protein. Fusion of lysosomes with autophagosome to generate the autolysosome is the penultimate step. In the fourth and final phase, the cargo is degraded.




Figure 2. Autophagy detection with Premo™ Autophagy Sensor. HeLa cells were plated at 5,000 cells per well and allowed to adhere overnight. Cells were then transduced with LC3B-GFP. 24 hours later, the cells were incubated with 50 μM chloroquine or left untreated (control) for 16 hours. Analysis was performed by quantifying the fluorescence from vesicular structures in the perinuclear region using the Thermo Scientific Celloomics® ArrayScan® VTI platform with Compartmental Analysis.




Figure 3. Fluorescence excitation and emission spectra for Premo™ LC3B-Green Fluorescent Protein (LC3B-GFP, panel A) and Premo™ LC3B-Red Fluorescent Protein (LC3B-RFP, panel B).

Before Starting

Materials Required but Not Provided

- Cell culture medium

Guidelines for Working with BacMam Reagents

- The following protocol is based on a 2 mL labeling volume and ~40,000 cells plated in a 35-mm dish or 1 well of a 6-well culture plate, and an MOI (multiplicity of infection) of 30 plaque forming units (pfu) per cell.
- For applications that require a larger number of cells, such as high-content screening (HCS), we recommend transducing the cells in a 10-cm dish or a T-75 flask and increasing the labeling volume to 10 mL with a proportionate increase in the volume of the virus. Following an overnight incubation under normal growth conditions, trypsinize and count cells for distribution to appropriate plates at cell number desired.
- For optimal results you may alter the MOI, cell density, temperature, and incubation time.
- We recommend transducing cells at a confluency of about 70% for best results.
- For first time users of BacMam reagents, we recommend exceptionally well-transduced cells like U-2 OS (ATCC CCL-11226).

Experimental Protocols

Protocol The following protocol was optimized using cells that have been plated. You can also treat the cells in suspension prior to plating.

- 1.1 Plate the cells at the desired density and allow them sufficient time to adhere. BacMam reagents work best when used on low-passage-number cells.
- 1.2 Calculate the volume of LC3B-FP (Component A) and LC3B(G120A)-FP (Component B) using the equation below:

$$\text{mL of LC3B-FP or LC3B(G120A)-FP} = \frac{(\text{number of cells}) (\text{MOI})}{(1 \times 10^8)}$$

where the number of cells is the estimated total number of cells at the time of cell labeling, MOI (multiplicity of infection) is the number of viral particles per cell, and 1×10^8 is the number of viral particles per mL of the reagent.

For example, to label 40,000 cells with an MOI of 30,

$$\text{mL of LC3B-FP or LC3B(G120A)-FP} = \frac{(40,000) (30)}{(1 \times 10^8)} = 0.012 \text{ mL (12 } \mu\text{L)}$$

- 1.3 Mix each LC3B reagent by inversion to ensure a homogenous solution.
- 1.4 Add the LC3B reagent directly to the cells in complete cell medium and mix gently.
- 1.5 Incubate the cells overnight (≥ 16 hours) for LC3B expression.
- 1.6 *Optional:* Treat the control cells with 30–100 μM chloroquine (Component B) for 12–16 hours.
- 1.7 Image and analyze the cells using the appropriate instrument filter sets. Refer to Figure 3 for spectral characteristics of GFP and RFP. Autophagosomes are typically located in the perinuclear region.

Note: Premo™ Autophagy Sensors were designed for use in live-cell imaging of autophagy. The cell-permeant nucleic acid stains Hoechst 33342 and HCS NuclearMask™ Blue stains are spectrally compatible with the Premo™ Autophagy Sensor fluorescence. Should you prefer fixed cell analysis, the fluorescence from GFP and RFP has been demonstrated to be resistant to fixation with 4% formaldehyde and permeabilization with 0.1% Triton® X-100. Fixation and permeabilization enables processing of labeled cells with antibodies to other cellular targets. We recommend fixed cell format for large sample sizes, such as for HCS.

- 1.8 Image the cells using the appropriate filters for GFP or RFP (see Figure 3).

Note: Overexpression of LC3B has been shown to promote aggregation of this protein in a non-autophagic dependent manner;⁴ for this reason, we recommend titrating the expression levels by varying the MOI.

References

1. Molecular Cell Biology 8, 931 (2007);
2. Genes & Development 21, 2861 (2007);
3. Drug Discovery 6, 304 (2007);
4. Meth Enz 452, 25 (2009).

Product List Current prices may be obtained from our website or from our Customer Service Department.

Cat. no.	Product Name	Unit Size
P36235	Premo™ Autophagy Sensor LC3B-GFP *BacMam 2.0*	1 kit
P36236	Premo™ Autophagy Sensor LC3B-RFP *BacMam 2.0*	1 kit
Related Products		
C10596	CellLight™ Lysosomes-GFP *BacMam 2.0*	1 mL
C10597	CellLight™ Lysosomes-RFP *BacMam 2.0*	1 mL
C10600	CellLight™ Mitochondria-GFP *BacMam 2.0*	1 mL
C10601	CellLight™ Mitochondria-RFP *BacMam 2.0*	1 mL
H3570	Hoechst 33342, trihydrochloride, trihydrate *10 mg/mL solution in water*	10 mL
H10325	HCS NuclearMask™ Blue	65 µL
L7528	LysoTracker® Red DND-99 *1 mM solution in DMSO* *special packaging*	10 × 50 µL
L7535	LysoTracker® Red DND-26 *1 mM solution in DMSO* *special packaging*	10 × 50 µL
L10382	LC3B Antibody Kit for Autophagy *rabbit polyclonal LC3B* *includes autophagosome inducer*	1 kit
M7510	MitoTracker® Orange CMTMRos *special packaging*	10 × 50 µL
M7512	MitoTracker® Red CMXRos *special packaging*	10 × 50 µL
M7514	MitoTracker® Green FM *special packaging*	10 × 50 µL
M22423	MitoTracker® Deep Red *special packaging*	10 × 50 µL

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probestech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

**For country-specific contact information,
visit www.invitrogen.com.**

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The manufacture, use, sale or import of this product may be subject to one or more patents or pending applications owned or licensed by Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity) in a manner consistent with the accompanying product literature. The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the most restrictive terms apply. Life Technologies Corporation will not assert a claim against the buyer of infringement of patents that are owned or controlled by Life Technologies Corporation and/or Molecular Probes, Inc. which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Limited Use Label License No. 21: Bac-to-Bac® and Bac-to-Bac® HT

This product is the subject of U.S. Patent No. 5,348,886. This product is sold under patent license from Monsanto for research purposes only and no license for commercial use is included. Requests for licenses for commercial manufacture or use should be directed to Director, Monsanto Corporate Research, 800 N. Lindbergh, St. Louis, Missouri 63167.

Limited Use Label License No. 127: GFP with Heterologous Promoter

This product and its use is the subject of one or more of U.S. Patent Nos. 5,491,084 and 6,146,826, and foreign equivalents. This product is sold under license from Columbia University. Rights to use this product are limited to research use only, and expressly exclude the right to manufacture, use, sell or lease this product for use for measuring the level of toxicity for chemical agents and environmental samples in cells and transgenic animals. No other rights are conveyed. Not for human use or use in diagnostic or therapeutic procedures. Inquiry into the availability of a license to broader rights or the use of this product for commercial purposes should be directed to Columbia Innovation Enterprise, Columbia University, Engineering Terrace-Suite 363, New York, New York 10027.

Limited Use Label License No. 198: Fluorescent Protein Products

This product and its use is the subject of one or more of U.S. Patent Nos. 5,777,079, 6,066,476, and 6,319,669 and foreign equivalents. Any use of this product by a commercial entity requires a separate license from either GE Healthcare or Invitrogen Corporation. For information on obtaining a commercial license to use this product, please refer to the contact information located at the bottom of this statement. No rights are conveyed to modify or clone the gene encoding GFP contained in this product. For information on licensing, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Email: outlicensing@invitrogen.com.

Limited Use Label License No. 267: Mutant GFP Products

This product and its use is the subject of one or more of U.S. Patent Nos. 6,090,919, 5,804,387, 5,994,077, and foreign equivalents.

Limited Use Label License No. 272: Humanized GFP

This product is the subject of one or more of U.S. Patent Numbers 5,786,464, 5,795,737, 5,874,304, and 6,114,148 and foreign equivalents licensed by Invitrogen Corporation. This product is sold for research use only. Not for therapeutic or diagnostic use in humans.

Limited Use Label License No. 306: Baculovirus Vectors

Certain methods that utilize the product associated with this limited use label license are covered by U.S. Patent Nos. 5,731,182 and 5,871,986. This product is for research use only by those researchers in laboratories of academic, government, industrial and/or clinical institutions engaged in the investigation of biological or biochemical processes, or research and development of biological products. This product is not to be used in the manufacture, use or sale of human or animal diagnostic, therapeutic or prophylactic products.

Limited Use Label License No. 308: WPRE Element

This product contains the Woodchuck Post-transcriptional Regulatory Element ("WPRE") which is the subject of intellectual property owned by The Salk Institute for Biological Studies, and licensed to Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; and/or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. In addition, any use of WPRE outside of this product or the product's authorized use requires a separate license from the Salk Institute. Life Technologies will not assert a claim against the buyer of infringement of patents owned by Life Technologies and claiming this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product or for a Commercial Purpose. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008, Phone (760) 603-7200. Fax (760) 602-6500, or The Salk Institute for Biological Studies, 10010 North Torrey Pines Road, La Jolla, CA 92037, Attn.: Office of Technology Management, Phone: (858) 453-4100 extension 1275, Fax: (858) 546-8093.

Limited Use Label License No. 321: TagRFP

This product is for the buyer's internal research use only and may not be used for commercial purposes. No rights are conveyed to modify or clone the gene encoding fluorescent protein contained in this product. The right to use this product specifically excludes the right to validate or screen compounds for commercial purposes. For information on commercial licensing, contact Licensing Department, Evrogen JSC, email: license@evrogen.com.

Limited Use Label License No. 332: Bacmam Virus Use

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer solely in accordance with the accompanying product literature or manual. Purchase of this product does not convey a license to expand, amplify, or otherwise propagate the provided viral particles or to otherwise modify or alter the virus by any means.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2010, Life Technologies Corporation. All rights reserved. This information is subject to change without notice.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.