

USER GUIDE

molecular  
probes®  
by *life* technologies™

# CAL-LYSE™ Lysing Solution

Whole blood lysing solution for flow cytometric applications

**Catalog Numbers** GAS-010, GAS-010S-100

**Document Part Number** L13010

**Publication Number** MAN0008970


**Revision** 1.0

CE **IVD**

*life*  
technologies™



7335 Executive Way  
Frederick, MD 21704  
E-mail: [techsupport@lifetech.com](mailto:techsupport@lifetech.com)



Life Technologies Limited  
European Regulatory Affairs,  
3 Fountain Drive, Inchinnan Business Park  
Paisley PA49RF, Scotland, UK  
Tel: +44 (0)141 81416305

Information in this document is subject to change without notice.

#### **DISCLAIMER**

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

#### **LIMITED WARRANTY**

Life Technologies and/or its affiliate(s) warrant their products as set forth in the Life Technologies General Terms and Conditions of Sale found on the Life Technologies web site at <http://www.lifetechnologies.com/termsandconditions>. If you have any questions, please contact Life Technologies at [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support).

#### **NOTICE TO PURCHASER: LIMITED USE LABEL LICENSE: 451: No right to resell**

Notice to Purchaser: No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com).

#### **TRADEMARKS**

The trademarks mentioned herein are the property of Life Technologies Corporation and/or its affiliate(s) or their respective owners.

© 2013 Life Technologies Corporation. All rights reserved.

# Contents

|                                    |          |
|------------------------------------|----------|
| <b>Introduction .....</b>          | <b>v</b> |
| Kit Contents and Storage.....      | v        |
| About the System .....             | vi       |
| <b>Methods .....</b>               | <b>1</b> |
| Procedure.....                     | 1        |
| <b>Appendix A .....</b>            | <b>3</b> |
| Interpretation of Results .....    | 3        |
| Limitations of the Procedure ..... | 4        |
| Bibliography.....                  | 5        |
| <b>Appendix B .....</b>            | <b>6</b> |
| Technical Support.....             | 6        |
| Explanation of Symbols.....        | 7        |



# Introduction

## Kit Contents and Storage

---

**Reagent provided** Sufficient Cal-Lyse™ Lysing Solution is provided for the lysis of 250 samples (GAS-010) or 1000 samples (GAS-010S-100) of anticoagulated whole blood.

| Cat. no.     | Amount              | Storage                                 |
|--------------|---------------------|---|
| GAS-010      | 25 mL (250 tests)   | Room temperature. <b>Do not freeze.</b> |
| GAS-010S-100 | 100 mL (1000 tests) | Room temperature. <b>Do not freeze.</b> |

The composition of the lysing solution includes the following ingredients:

- Polyvinylpyrrolidone
- Ethylene-bis (oxyrthylenenitruol) -tetra acetic acid
- Sodium phosphate
- Formaldehyde
- Deionized water

### Warning

Cal-Lyse™ Lysing Solution contains formaldehyde as a fixative. Formaldehyde is toxic, allergenic and a suspected carcinogen. Avoid ingestion, inhalation and contact with eyes, skin and clothing.

---

### Storage

Store at room temperature. **Do not freeze.**

---

### Evidence of deterioration

The reagent should not be used if any evidence of deterioration such as cloudiness or substantial loss of reactivity is observed.

- The normal appearance of Cal-Lyse™ Lysing Solution solution is a clear, particulate-free liquid. Do not use this reagent if discoloration occurs or a visible precipitate forms.
  - We recommend periodically checking the reagent against your standard.
- 

### Intended use

Cal-Lyse™ is a lysing solution formulated for lysis of erythrocytes in samples of anticoagulated human blood. Cal-Lyse™ Lysing Solution is intended to be used as an aid in the enumeration of leukocytes from various human biological sources including blood and bone marrow that have been stained with monoclonal antibodies for analysis using flow cytometric methods.

Results should be put within the context of other diagnostic tests, clinical history by a certified professional.

**Important: Cal-Lyse™ Lysing Solution is not a concentrate, and should not be diluted prior to use.** For optimal results, use only as directed.

**For *In Vitro* Diagnostic Use.**

---

*Continued on next page*

## About the System

---

### **Summary and principles of the test**

Cal-Lyse™ Lysing Solution is a premixed solution for use in whole blood staining procedures.

To perform analysis of leukocytes in the sample, cells bearing specific antigenic determinants are identified by incubating blood samples with fluorochrome-conjugated monoclonal antibodies. Monoclonal antibodies bind to the surfaces of viable blood cells that express the corresponding antigens.

Lysis is performed with Cal-Lyse™ Lysing Solution immediately following staining of the blood samples with fluorochrome-conjugated monoclonal antibodies. Treatment with this reagent simultaneously leads to lysis of erythrocytes and fixation of leukocytes while maintaining morphological scatter characteristics of the leukocytes. Cal-Lyse™ Lysing Solution contains formaldehyde as fixative, and no additional fixation is required.

Cells may subsequently be washed, to elimination of red cell debris as well as unbound antibody. The erythrocyte-free intact antibody-stained leukocytes are suitable for flow cytometric analysis.

---

# Methods

## Procedure

---

### Statement of warning

- Do not pipet by mouth.
  - Samples should be handled as if capable of transmitting infection. Appropriate disposal methods should be used.
  - Do not use reagent beyond the stated expiration date of the product.
  - Cal-Lyse™ Lysing Solution is not a concentrate, and should not be diluted prior to use. For optimal results, use only as directed.
  - Deviations from the recommended procedure enclosed within this product insert may invalidate the results of testing
  - The reagent should not be used if any evidence of deterioration such as cloudiness or substantial loss of reactivity is observed.
- 

### Required materials that are not supplied

- Centrifuge with swinging bucket rotor capable of  $1000 \times g$
  - Vortex mixer
  - 3–5 mL conical tubes
  - Micropipette capable of dispensing 100  $\mu\text{L}$  volumes
  - Blood collection tubes with anticoagulant (EDTA or heparin)
  - Phosphate Buffered Saline (PBS)
  - Flow cytometer
- 

*Continued on next page*

## Procedure, Continued

---

### Sample preparation

Two methods of sample preparation are provided. The first procedure includes a wash step to remove red cell debris and any unbound antibody, while the second procedure does not include the wash step.

**Note:** Failure to completely remove erythrocytes from blood samples may interfere with the subsequent flow cytometric enumeration of antibody-stained leukocytes.

---

### Antibody staining and wash lysis procedure

1. Collect blood into an appropriate anticoagulant (EDTA or heparin).
  2. Pipette 100  $\mu$ L of thoroughly mixed blood into a conical tube.
  3. Add antibody as indicated in the manufacturer's package insert to appropriately labeled tubes from step 2. Mix gently.
  4. Incubate all tubes for 15 minutes at room temperature ( $22 \pm 3^{\circ}\text{C}$ ) in the dark.
  5. Add 100  $\mu$ L of Cal-Lyse™ Lysing Solution to all tubes.
  6. Incubate all tubes for 10 minutes at room temperature ( $22 \pm 3^{\circ}\text{C}$ ) in the dark.
  7. Add 1.0 mL of deionized water to all tubes. Immediately vortex tubes.
  8. Incubate all tubes for 10 minutes at room temperature
  9. Add 3–4 mL of distilled water to each tube, and vortex. Incubate tubes for 5–10 minutes at room temperature.
  10. Centrifuge tubes for 5 minutes at  $300 \times g$ . Remove supernatant.
  11. Resuspend the cells in all tubes in 1 mL of PBS or Sheath Fluid.
  12. Analyze on a flow cytometer according immediately, or store samples at  $2^{\circ}\text{C}$ – $8^{\circ}\text{C}$  in the dark and analyze within 24 hours.
- 

### Antibody staining without wash lysis procedure

1. Collect blood into an appropriate anticoagulant (EDTA or heparin).
  2. Pipette 100  $\mu$ L of thoroughly mixed blood into a conical tube.
  3. Add antibody as indicated in the manufacturer's package insert to appropriately labeled tubes from step 2. Mix gently.
  4. Incubate all tubes for 15 minutes at room temperature ( $22 \pm 3^{\circ}\text{C}$ ) in the dark.
  5. Add 100  $\mu$ L of Cal-Lyse™ Lysing Solution to all tubes.
  6. Incubate all tubes for 10 minutes at room temperature ( $22 \pm 3^{\circ}\text{C}$ ) in the dark.
  7. Add 1.0 mL of deionized water to all tubes. Immediately vortex tubes.
  8. Incubate all tubes for 10 minutes at room temperature.
  9. Analyze on a flow cytometer according immediately, or store samples at  $2^{\circ}\text{C}$ – $8^{\circ}\text{C}$  in the dark and analyze within 24 hours.
- 

*Continued on next page*



# Appendix A


## Interpretation of Results

---

### Flow cytometry

Analyze antibody-stained cells on an appropriate flow cytometer according to the manufacturer's instructions. The right angle light scatter, or side scatter (SSC), versus forward angle light scatter (FSC) is collected to reveal the lymphocyte cell cluster. A gate is drawn for the lymphocyte cluster (lymphocyte bitmap). The fluorescence attributable to the fluorochrome-conjugated monoclonal antibody is collected, and the percentage of antibody-stained cells is determined. A fluorochrome-conjugated isotypic control may be used to estimate and correct for nonspecific binding to lymphocytes. The isotypic control must be of the same heavy chain immunoglobulin class as the specific antibody. The isotypic control should be used at approximately the same antibody concentration as the specific antibody. An analysis region is set to exclude background fluorescence and to include positively stained cells.

The following histograms are representative of cell staining from a normal donor, gated on the lymphocyte region and stained with a representative T cell-specific monoclonal antibody (CD3 FITC), and a B cell-specific monoclonal antibody (CD19 R-PE). Cells were processed with the "Antibody staining and wash lysis procedure". Red blood cells were lysed with Cal-Lyse™ Lysing Solution.



## Limitations of the Procedure

---

1. The values obtained from normal individuals may vary from laboratory to laboratory; therefore, it is recommended that each laboratory establish its own normal range.
2. When using the whole blood method, red blood cells found in some abnormal donors, as well as nucleated red cells found in normal and abnormal donors, may be resistant to lysis. Longer red cell lysis periods may be needed to avoid the inclusion of unlysed red cells.
3. Blood samples should not be refrigerated or retained at ambient temperature for an extensive period (longer than 24–30 hours) prior to incubation with monoclonal antibodies.

---

*Continued on next page*

## Bibliography

---

1. Shi, W., Kasdan, H.L., Fridge, A., Tai, Y-C. 2010. Four-part differential leukocyte count using  $\mu$ flow cytometer. *Micro Electro Mechanical Systems (MEMS), 2010 IEEE 23rd International Conference*. 1019–1022.
  2. O'Farrell, A-M., Abrams, T.J., Yuen, H.A., Ngai, T.J., Louie, S.G., Yee, K.W.H., Wong, L.M., Hong, W., Lee, L.B., Town, A., Smolich, B.D., Manning, W.C., Murray, L.J., Heinrich, M.C., and Cherrington, J.M. 2003. *Blood*. 101(9): 3597–3605.
  3. Little M.A., Al-Ani B., Ren S., Al-Nuaimi H., Leite M. Jr, et al. 2012. Anti-Proteinase 3 Anti-Neutrophil Cytoplasm Autoantibodies Recapitulate Systemic Vasculitis in Mice with a Humanized Immune System. *PLoS ONE* 7(1):e28626.
  4. Young, S-H., Wolfarth, M.G., Roberts, J.R., Kashon, M.L., and Antonini, J.M. 2013. Adjuvant effect of zymosan after pulmonary treatment in a mouse ovalbumin allergy model. *Exp. Lung Res.* 39(1): 48–57.
  5. Yassin, L.M., Londoño, J., Montoya, G., De Sanctis, J.B., Rojas, M., Ramírez, L.A., García, L.F., and Vásquez, G. 2011. Atherosclerosis development in SLE patients is not determined by monocytes ability to bind/endocytose Ox-LDL. *Autoimmunity*. 44(3): 201–210.
  6. Antonini, J.M., Zeidler-Erdely, P.C., Young, S-H., Roberts, J.R., and Erdely, A. 2012. Systemic immune cell response in rats after pulmonary exposure to manganese-containing particles collected from welding aerosols. *J. Immunotoxicology*. 9(2): 184– 192.
  7. Metzler, B., Gfeller, P., Bigaud, M., Li, J., Wieczorek, G., Heusser, C., Lake, P., and Katopodis A. 2004. Combinations of Anti-LFA-1, Everolimus, Anti-CD40 Ligand, and Allogeneic Bone Marrow Induce Central Transplantation Tolerance through Hemopoietic Chimerism, Including Protection from Chronic Heart Allograft Rejection. *J. Immunology*. 173(11): 7025 – 7036.
  8. Baribaud, F., Edwards, T.G., Sharron, M., BreLOT, A., Heveker, N., Price, K., Mortari, F., Alizon, M., Tsang, M., and Doms, R.W. 2001. Antigenically Distinct Conformations of CXCR4. *J. Virology*. 75(19): 8957 – 8967.
-

## Appendix B

### Technical Support

---

**Obtaining support** For the latest services and support information for all locations, go to [www.lifetechnologies.com](http://www.lifetechnologies.com)

At the website, you can:

- Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities
  - Search through frequently asked questions (FAQs)
  - Submit a question directly to Technical Support ([techsupport@lifetech.com](mailto:techsupport@lifetech.com))
  - Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents
  - Obtain information about customer training
  - Download software updates and patches
- 

#### **Safety Data Sheets (SDS)**

Safety Data Sheets (SDSs) are available at [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support).

---

#### **Certificate of Analysis**










The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support) and search for the Certificate of Analysis by product lot number, which is printed on the box.

---

*Continued on next page*

# Explanation of Symbols

---

| Symbol  | Description                                  | Symbol  | Description                  | Symbol  | Description                             |
|---|--|---|------------------------------|---|---|
|  | Manufacturer                                 |  | Catalog number               |  | Batch code                              |
|  | <i>In vitro</i> diagnostic medical device    |  | Use by                       |  | Temperature limitation                  |
|  | European Community authorized representative |  | Consult instructions for use |  | Caution, consult accompanying documents |

---

## Notes



**Headquarters**

5791 Van Allen Way | Carlsbad, CA 92008 USA | Phone +1 760 603 7200 | Toll Free in USA 800 955 6288

**For support visit** [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)

[www.lifetechnologies.com](http://www.lifetechnologies.com)

30 March 2013

