

## Lectin PNA Conjugates

### Quick Facts

#### Storage upon receipt:

- $\leq -20^{\circ}\text{C}$
- Desiccate
- Protect from light

**Abs/Em:** See Table 1

**Table 1.** Lectin PNA conjugates and spectral characteristics.

Catalog #	Conjugate	Abs*	Em*
L-21409	Alexa Fluor 488	495	519
L-32458	Alexa Fluor 568	579	603
L-32459	Alexa Fluor 594	590	617
L-32460	Alexa Fluor 647	650	668

\* Approximate absorption (Abs) and emission (Em) maxima, in nm.

### Introduction

Lectin PNA, isolated from peanuts, is a 110 kDa tetramer composed of identical subunits. This lectin is specific for terminal  $\beta$ -galactose, and will agglutinate human erythrocytes but only after neuraminidase treatment.<sup>1</sup> Molecular Probes offers several Alexa Fluor<sup>®</sup> dye conjugates of lectin PNA; Table 1 provides a summary of absorption and emission maxima for these conjugates.

### Contents and Storage

The fluorescent conjugates of lectin PNA are supplied lyophilized in unit sizes of 1 mg. When stored desiccated at  $\leq -20^{\circ}\text{C}$ , the lyophilized products are stable for at least one year. Solutions up to  $\sim 1$  mg/mL can be made by dissolving the protein in deionized water. Solutions, with the addition of sodium azide to a final concentration of 2 mM, can be stored at  $2-6^{\circ}\text{C}$  for at least four months with no loss of activity. For longer storage, di-

vide the solutions into aliquots and freeze at  $\leq -20^{\circ}\text{C}$ . PROTECT FROM LIGHT. AVOID REPEATED FREEZING AND THAWING OF SOLUTIONS.

It is a good practice to centrifuge the lectin conjugate solutions briefly in a microcentrifuge before use; only the supernatants should then be added to the experiment. This step will eliminate any protein aggregates that may have formed in solution, thereby reducing nonspecific background staining.

### Applications

PNA binding sites are widespread in human tissues, with staining patterns varying by tissue type.<sup>2</sup> Research has shown PNA to be selective for acrosomes in rat and human sperm,<sup>3,4</sup> and PNA serves as a marker for certain melanomas.<sup>5,6</sup> PNA has also been used to label the synaptic extracellular matrix in the study of developing neuromuscular junctions.<sup>7</sup> Since the applications of the PNA lectin are varied, researchers should consult the primary literature for protocol information.

### References

1. J Biol Chem 250, 8518 (1975); 2. Hum Pathol 15, 904 (1984); 3. Mol Reprod Dev 55, 289 (2000); 4. Histochem J 29, 583 (1997); 5. Hum Pathol 30, 556 (1999); 6. Oncol Res 5, 235 (1993); 7. J Neurosci 14, 796 (1994).

### Product List

Current prices may be obtained from our Web site or from our Customer Service Department.

Cat #	Product Name	Unit Size
L-21409	lectin PNA from <i>Arachis hypogaea</i> (peanut), Alexa Fluor <sup>®</sup> 488 conjugate .....	1 mg
L-32458	lectin PNA from <i>Arachis hypogaea</i> (peanut), Alexa Fluor <sup>®</sup> 568 conjugate .....	1 mg
L-32459	lectin PNA from <i>Arachis hypogaea</i> (peanut), Alexa Fluor <sup>®</sup> 594 conjugate .....	1 mg
L-32460	lectin PNA from <i>Arachis hypogaea</i> (peanut), Alexa Fluor <sup>®</sup> 647 conjugate .....	1 mg

---

## Contact Information

Further information on Molecular Probes' products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Leiden, the Netherlands. All others should contact our Technical Assistance Department in Eugene, Oregon.

Please visit our Web site — [www.probes.com](http://www.probes.com) — for the most up-to-date information

### **Molecular Probes, Inc.**

29851 Willow Creek Rd., Eugene, OR 97402  
Phone: (541) 465-8300 • Fax: (541) 344-6504

**Customer Service:** 6:00 am to 4:30 pm (Pacific Time)  
Phone: (541) 465-8338 • Fax: (541) 344-6504 • [order@probes.com](mailto:order@probes.com)

### **Toll-Free Ordering for USA and Canada:**

Order Phone: (800) 438-2209 • Order Fax: (800) 438-0228

**Technical Assistance:** 8:00 am to 4:00 pm (Pacific Time)  
Phone: (541) 465-8353 • Fax: (541) 465-4593 • [tech@probes.com](mailto:tech@probes.com)

### **Molecular Probes Europe BV**

PoortGebouw, Rijnsburgerweg 10  
2333 AA Leiden, The Netherlands  
Phone: +31-71-5233378 • Fax: +31-71-5233419

**Customer Service:** 9:00 to 16:30 (Central European Time)  
Phone: +31-71-5236850 • Fax: +31-71-5233419  
[eurorder@probes.nl](mailto:eurorder@probes.nl)

**Technical Assistance:** 9:00 to 16:30 (Central European Time)  
Phone: +31-71-5233431 • Fax: +31-71-5241883  
[eurotech@probes.nl](mailto:eurotech@probes.nl)

*Molecular Probes' products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.*

Several of Molecular Probes' products and product applications are covered by U.S. and foreign patents and patents pending. Our products are not available for resale or other commercial uses without a specific agreement from Molecular Probes, Inc. We welcome inquiries about licensing the use of our dyes, trademarks or technologies. Please submit inquiries by e-mail to [busdev@probes.com](mailto:busdev@probes.com). All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2003, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.