

## Isolectin IB<sub>4</sub> Conjugates

**Table 1.** Contents and Storage Information.

Material	Amount	Concentration	Storage	Stability
Isolectin IB <sub>4</sub> Alexa Fluor® dye conjugates	500 µg lyophilized powder *	NA	<ul style="list-style-type: none"> <li>• ≤-20°C</li> <li>• Desiccate</li> <li>• Protect from light</li> </ul>	1 year
Isolectin IB <sub>4</sub> Biotin conjugates	500 µg lyophilized powder *	NA	<ul style="list-style-type: none"> <li>• ≤-20°C</li> <li>• Desiccate</li> </ul>	

\* The vials are packed according to the protein content and not the dry weight, thus, it is best to solubilize the entire contents of a vial at one time.

**Approximate Fluorescence Excitation and Emission, in nm:** Alexa Fluor® 488 dye ~495/519 nm; Alexa Fluor® 568 dye ~579/603 nm; Alexa Fluor® 594 dye ~590/617 nm; Alexa Fluor® 647 dye ~650/668 nm.

### Introduction

Isolectin IB<sub>4</sub> is a 114 kDa glycoprotein and part of a family of five tetrameric type I isolectins (IA<sub>4</sub>, IA<sub>3</sub>B, IA<sub>2</sub>B<sub>2</sub>, IAB<sub>3</sub>, and IB<sub>4</sub>) isolated from the seeds of the tropical African legume *Griffonia simplicifolia* (formerly *Bandeiraea simplicifolia*). The A and B subunits comprising the family members are very similar, differing in amino acid sequence only at the N-terminus. However, the subunits display different binding specificities; the A subunit prefers *N*-acetyl-D-galactosamine end groups while the B subunit is selective for terminal α-D-galactosyl residues.<sup>1</sup> Molecular Probes offers several Alexa Fluor® dye conjugates of isolectin IB<sub>4</sub>, as well as a biotinylated form.

### Guidelines for Use

#### Preparing Isolectin IB<sub>4</sub> Conjugates

Solutions up to ~1 mg/mL can be made by dissolving the protein in an aqueous buffer at neutral pH containing 0.1–1.0 mM CaCl<sub>2</sub>. Solutions, with the addition of sodium azide to a final concentration of 2 mM, can be stored at 4°C for at least four months with no loss of activity. For longer storage, divide the solution into aliquots and freeze at ≤-20°C. AVOID REPEATED FREEZING AND THAWING.

### Using Conjugate Solutions

It is a good practice to centrifuge the isolectin conjugate solution briefly in a microcentrifuge before use; only the supernatant should then be added to the experiment. This step will eliminate any protein aggregates that may have formed in solution, thereby reducing nonspecific background staining.

Isolectin IB<sub>4</sub> specifically agglutinates blood group B erythrocytes and was originally employed for this purpose.<sup>2</sup> Subsequent work has shown that isolectin IB<sub>4</sub> is cytotoxic to several normal and tumor cell types<sup>3</sup> and has particularly strong affinity for brain microglial and perivascular cells.<sup>4</sup> It has also been particularly valuable as a histochemical and flow cytometric probe for specifically labeling endothelial cells in a number of species.<sup>5,6</sup> Isolectin IB<sub>4</sub> has been used effectively for tracing central and peripheral neuronal pathways following local injections,<sup>7,8</sup> as well as for labeling stimulated murine macrophages,<sup>9</sup> bovine thyroid cells,<sup>10</sup> various murine cell types,<sup>11,12</sup> laminin,<sup>13</sup> and thyroglobulin.<sup>14</sup> Since the applications of isolectin IB<sub>4</sub> are varied, researchers should consult the primary literature for protocol information.

## References

---

1. J Biol Chem 252, 4739 (1977); 2. Subcell Biochem 32, 127 (1999); 3. Cancer Res 42, 2977 (1982); 4. Histochemistry 102, 483 (1994); 5. Histochem J 19, 225 (1987); 6. Am J Pathol 134, 1227 (1989); 7. Neurosci Lett 222, 53 (1997); 8. Brain Res 811, 34 (1998); 9. Proc Natl Acad Sci USA 79, 166 (1982); 10. Arch Biochem Biophys 343, 73 (1997); 11. Exp Cell Res 120, 321 (1979); 12. Transplantation 61, 13 (1996); 13. Biochemistry 28, 6379 (1989); 14. J Biol Chem 259, 9858 (1984).

## Product List

Current prices may be obtained from our website or from our Customer Service Department.

---

Cat #	Product Name	Unit Size
I21411	isolectin GS-IB <sub>4</sub> from Griffonia simplicifolia, Alexa Fluor® 488 conjugate .....	500 µg
I21412	isolectin GS-IB <sub>4</sub> from Griffonia simplicifolia, Alexa Fluor® 568 conjugate .....	500 µg
I21413	isolectin GS-IB <sub>4</sub> from Griffonia simplicifolia, Alexa Fluor® 594 conjugate .....	500 µg
I32450	isolectin GS-IB <sub>4</sub> from Griffonia simplicifolia, Alexa Fluor® 647 conjugate .....	500 µg
I21414	isolectin GS-IB <sub>4</sub> from Griffonia simplicifolia, biotin-XX conjugate .....	500 µg

## Contact Information

---

**Molecular Probes, Inc.**  
29851 Willow Creek Road  
Eugene, OR 97402  
Phone: (541) 465-8300  
Fax: (541) 335-0504

**Customer Service:**  
6:00 am to 4:30 pm (Pacific Time)  
Phone: (541) 335-0338  
Fax: (541) 335-0305  
probesorder@invitrogen.com

**Toll-Free Ordering for USA:**  
Order Phone: (800) 438-2209  
Order Fax: (800) 438-0228

**Technical Service:**  
8:00 am to 4:00 pm (Pacific Time)  
Phone: (541) 335-0353  
Toll-Free (800) 438-2209  
Fax: (541) 335-0238  
probestech@invitrogen.com

**Invitrogen European Headquarters**  
Invitrogen, Ltd.  
3 Fountain Drive  
Inchinnan Business Park  
Paisley PA4 9RF, UK  
Phone: +44 (0) 141 814 6100  
Fax: +44 (0) 141 814 6260  
Email: euroinfo@invitrogen.com  
Technical Services: eurotech@invitrogen.com

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Assistance Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

### Limited Use Label License

For research use only. Not intended for any animal or human therapeutic or diagnostic use. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402. Tel: (541) 465-8300. Fax: (541) 335-0504.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2005, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.