

TaqMan® Fast Advanced Master Mix Protocol

For Research Use Only. Not intended for any animal or human therapeutic or diagnostic use.

Information in this document is subject to change without notice.

APPLIED BIOSYSTEMS DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL APPLIED BIOSYSTEMS BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT APPLIED BIOSYSTEMS IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

NOTICE TO PURCHASER: LIMITED LICENSE

A license to perform the patented 5' Nuclease Process for research is obtained by the purchase of (i) both Licensed Probe and Authorized 5' Nuclease Core Kit, (ii) a Licensed 5' Nuclease Kit, or (iii) license rights from Applied Biosystems.

The use of this product is covered by U.S. patent claims and patent claims outside the U.S. The purchase of this product includes a limited, non-transferable immunity from suit under these patent claims for using only this amount of product for the purchaser's own internal research. The right to use this product in the 5' Nuclease Process can be obtained through purchase of an Authorized 5' Nuclease Core Kit. Except under separate license rights available from Applied Biosystems, no right under any other patent claim, or to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, or to sublicense, repackage with other products, or resell in any form, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Human diagnostics uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained from the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, U.S.A.

TRADEMARKS

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

AmpliQ and TaqMan are registered trademarks of Roche Molecular Systems, Inc.

Microsoft and Excel are registered trademarks of Microsoft Corporation.

TRI Reagent is a registered trademark of Molecular Research Center, Inc.

© 2010, 2011 Life Technologies Corporation Corporation. All rights reserved.

Part Number 4444605 Rev. C

04/2011

Contents

	About This Guide	7
	Safety information	7
	Safety alert words	7
	SDSs	7
PROTOCOL	TaqMan® Fast Advanced Master Mix	9
	Product information	9
	Purpose of the product	9
	Starting template	9
	About this protocol	9
	About the kit	10
	Storage and stability	10
	Kit components	10
	Recommended real-time PCR systems	11
	Thermal-cycling conditions for your real-time PCR system	11
	Before you begin	12
	Prevent contamination	12
	Other assays or thermal-cycling conditions	12
	Section 1 TaqMan® and Custom TaqMan® Gene Expression Assays	13
	User-supplied materials	14
	Reagents for reverse transcription	14
	Reagents for PCR	15
	Optional reagents	15
	Reaction plates, consumables, and equipment	15
	Workflow	16
	Perform reverse transcription	16
	Recommended kits and reagents	16
	Perform reverse transcription	16
	Perform real-time PCR amplification	17
	PCR reagent handling and preparation	17
	Prepare the PCR reaction mix	18
	Prepare the PCR reaction plate	18
	Set up a plate document or experiment	19
	Run the PCR reaction plate	20
	Analyze the data	21
	Using custom TaqMan® primer and probe sets	21

Determine the target template and amplicon	21
Design the custom TaqMan® probes and primers	21
Determine the optimal primer concentration	22
Determine the optimal probe concentration	24
Perform real-time PCR	26
Section 2 TaqMan® Array Plates	29
User-supplied materials	29
Reagents for reverse transcription	29
Reagents for PCR	30
Optional reagents	30
Reaction plates, consumables, and equipment	30
Workflow	31
Perform reverse transcription	32
Recommended kits and reagents	32
Perform reverse transcription	32
Perform real-time PCR amplification	32
PCR reagent handling and preparation	32
Prepare the PCR reaction mix	33
Prepare the PCR reaction plate	33
Set up a plate document or experiment	34
Run the PCR reaction plate	35
Analyze the data	35
Section 3 TaqMan® MicroRNA Assays	37
User-supplied materials	37
Reagents for reverse-transcription	37
Reagents for PCR	38
Optional reagents	38
Reaction plates, consumables, and equipment	38
Workflow	39
Perform reverse transcription	40
Recommended kit	40
Perform reverse transcription	40
Perform real-time PCR amplification	40
Reagent preparation guidelines	40
Prepare the PCR reaction mix	41
Prepare the PCR reaction plate	42
Set up a plate document or experiment	42
Run the PCR reaction plate	43
Analyze the data	43
Section 4 TaqMan® Array Micro Fluidic Cards	45
User-supplied materials	45
Reagents for reverse transcription	45

	Reagents for PCR	46
	Optional reagents	46
	Reaction plates, consumables, and equipment	46
	Workflow	47
	Perform reverse transcription	48
	Recommended kits and reagents	48
	Perform reverse transcription	48
	Perform real-time PCR amplification	49
	Card handling and preparation	49
	Recommended amounts of cDNA	49
	Prepare the PCR reaction mix	50
	Prepare the TaqMan Array card	50
	Set up a plate document or experiment	51
	Run the TaqMan Array card	51
	Analyze the data	51
APPENDIX A	Additional User-Supplied Materials	53
	Reaction plates and accessories	53
	Optional reagents	55
	General laboratory equipment	56
APPENDIX B	PCR Good Laboratory Practices	57
	Sample preparation	57
	Preventing contamination	58
	False positives	58
	Uracil-N glycosylase (UNG)	58
	Prevention of PCR product carryover	59
	Fluorescent contaminants	59
APPENDIX C	About Two-Step RT-PCR	61
	Overview	62
	About TaqMan® MGB Probes	63
	About AmpliTaq® Fast DNA Polymerase	63
	About uracil-N glycosylase	63
	About ROX™ Passive Reference dye	63
	About the 5′ nuclease assay	64
APPENDIX D	Safety	67
	Chemical safety	68
	General chemical safety	68
	SDSs	68

Chemical waste safety	69
Biological hazard safety	70
Bibliography	71
Glossary	73
Documentation and Support	75
Product documentation	75
Related documentation	75
Chemistry/reagent documentation	75
Instrument documentation	76
Obtaining support	77

About This Guide

Safety information



Note: For general safety information, see this section and [Appendix D, “Safety” on page 67](#). When a hazard symbol and hazard type appear by an instrument hazard, see the “Safety” Appendix for the complete alert on the chemical or instrument.

Safety alert words

Four safety alert words appear in Applied Biosystems user documentation at points in the document where you need to be aware of relevant hazards. Each alert word—**IMPORTANT, CAUTION, WARNING, DANGER**—implies a particular level of observation or action, as defined below:



IMPORTANT! – Indicates information that is necessary for proper instrument operation, accurate chemistry kit use, or safe use of a chemical.



CAUTION! – Indicates a potentially hazardous situation that, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.



WARNING! – Indicates a potentially hazardous situation that, if not avoided, could result in death or serious injury.



DANGER! – Indicates an imminently hazardous situation that, if not avoided, will result in death or serious injury. This alert word is to be limited to the most extreme situations.

SDSs

The Safety Data Sheets (SDSs) for any chemicals supplied by Applied Biosystems or Ambion are available to you free 24 hours a day. For instructions on obtaining SDSs, see [“SDSs” on page 68](#).



IMPORTANT! For the SDSs of chemicals not distributed by Applied Biosystems or Ambion contact the chemical manufacturer.

TaqMan[®] Fast Advanced Master Mix

Product information

Purpose of the product

Use the TaqMan[®] Fast Advanced Master Mix to perform the PCR step in any gene expression experiment or other quantitative analysis, such as:

- Pathogen detection
- Differential gene expression analysis
- Viral load quantitation
- MicroRNA quantitation
- Microarray validation

Starting template

You can use the TaqMan Fast Advanced Master Mix with the DNA target of your choice, including cDNA or gDNA.

For RNA quantitation experiments, the TaqMan Fast Advanced Master Mix is used in the second step of a two-step reverse transcription–polymerase chain reaction (RT-PCR) protocol. You can use Applied Biosystems reverse transcription kits to generate the cDNA template from RNA.

About this protocol

TaqMan Fast Advanced Master Mix has been optimized for use with primers and TaqMan[®] probes that have been designed according to Applied Biosystems guidelines. This protocol describes how to perform 2-step quantitative RT-PCR, using the TaqMan Fast Advanced Master Mix in the PCR step with:

- TaqMan[®] Gene Expression Assays, Custom TaqMan[®] Gene Expression Assays, and custom TaqMan[®] primer and probe sets
- TaqMan[®] Array Plates (Fast and Standard plates)
- TaqMan[®] MicroRNA Assays
- TaqMan[®] Array Micro Fluidic Cards

Because analysis methods vary greatly between applications, this protocol provides general guidelines for analyzing data generated from experiments that use TaqMan Fast Advanced Master Mix and TaqMan[®] assays. For more detailed information about data analysis or the procedures outlined in this protocol, refer to the appropriate documentation for your instrument (see [“Product documentation” on page 75](#)).

About the kit

Storage and stability

Upon receipt, store the TaqMan® Fast Advanced Master Mix at -20 °C. After opening, store the TaqMan Fast Advanced Master Mix at 4 °C.

Performance of the TaqMan Fast Advanced Master Mix is guaranteed until the expiration date printed on the package and bottle labels.

Kit components

The TaqMan Fast Advanced Master Mix contains:

- AmpliTaq® Fast DNA Polymerase
- Uracil-N glycosylase (UNG)
- dNTPs with dUTP
- ROX™ dye (passive reference)
- Optimized buffer components

The TaqMan Fast Advanced Master Mix is supplied as a 2× concentration and is available in the following volumes:

Kit	Part number	Volume	No. of 20-µL reactions
TaqMan® Fast Advanced Master Mix Kit	4444556	1 × 1-mL	100
	4444557	1 × 5-mL	500
	4444558	1 × 50-mL	5000
	4444963	2 × 5-mL	1000
	4444964	5 × 5-mL	2500
	4444965	10 × 5-mL	5000

Recommended real-time PCR systems

You must perform the PCR step on a real-time PCR system. Traditional thermal cyclers cannot be used because they cannot detect and record the fluorescent signals generated by the cleavage of TaqMan® probes. Applied Biosystems recommends the following real-time PCR systems:

- ViiA™ 7 Real-Time PCR System
- StepOne™ and StepOnePlus™ Systems
- 7900HT/7900HT Fast System
- 7500/7500 Fast System
- 7300 System



Note: For a list of the reaction plates and accessories available for each real-time PCR system, see [page 53](#).

Thermal-cycling conditions for your real-time PCR system

The TaqMan Fast Advanced Master Mix can be used with Fast and Standard real-time PCR systems, provided the thermal-cycling profile and run mode are correctly set.

- **Thermal-cycling profile** – The thermal-cycling profile defines the temperature and time for each step. Be sure to use the appropriate thermal-cycling profile for your system:

Applied Biosystems Real-Time PCR System	Thermal-cycling profile				
	Parameter	UNG incubation [†]	Polymerase activation [‡]	PCR (40 cycles)	
		Hold	Hold	Denature	Anneal/extend
	Temp. (°C)	50	95	95	60
7900HT system	Time (mm:ss)	02:00	00:20	00:01	00:20
7900HT Fast system (Fast 96-Well, Standard 96-Well, or 384-Well Block Modules)					
ViiA™ 7 system					
StepOne™ system					
StepOnePlus™ system	Time (mm:ss)	02:00	00:20	00:03	00:30
7500 Fast system					
7500 system					
7300 system					

[†] Required for optimal UNG activity.

[‡] Required to activate the AmpliTaq® Fast DNA Polymerase.

- **Run mode** – The run mode defines the ramp rate that is used to heat or cool the sample block between temperature changes.

Applied Biosystems Real-Time PCR System	Default run mode
7900HT system	Standard
7900HT Fast system (384-Well and Standard 96-Well Block Modules)	
7500 system	
7300 system	
ViiA™ 7 system	Fast
StepOne™ system	
StepOnePlus™ system	
7900HT Fast system (Fast 96-Well Block Module)	
7500 Fast system	

Before you begin

Prevent contamination

Review [Appendix B, “PCR Good Laboratory Practices”](#) on page 57.

Other assays or thermal-cycling conditions

If you use assays other than the TaqMan® assays, or use thermal-cycling conditions other than those specified in this protocol, validate your assays and re-optimize your thermal-cycling conditions as needed.

Section 1 TaqMan® and Custom TaqMan® Gene Expression Assays

This section provides a general protocol for performing two-step RT-PCR. In the PCR step, you use the TaqMan® Fast Advanced Master Mix with one of the following assays:

- TaqMan® Gene Expression Assays
- Custom TaqMan® Gene Expression Assays
- Custom TaqMan® primer and probe sets

This section covers:

■ User-supplied materials	14
■ Workflow	16
■ Perform reverse transcription	16
■ Perform real-time PCR amplification	17
■ Analyze the data	21
■ Using custom TaqMan® primer and probe sets	21

For more information

For more detailed information about the procedures outlined in this protocol, refer to the appropriate documentation for your instrument. See [“Product documentation” on page 75](#).

User-supplied materials

Reagents for reverse transcription

Reagents		Applied Biosystems part number
High Capacity cDNA Reverse Transcription Kit	200 reactions	4368814
	200 reactions with RNase Inhibitor	4374966
	1000 reactions	4368813
	1000 reactions with RNase Inhibitor	4374967
High Capacity RNA-to-cDNA™ Kit	50 reactions	4387406
High Capacity RNA-to-cDNA Master Mix	500 reactions	4390779
	200 reactions	4390778
	50 reactions	4390777
	15 reactions	4390776
High Capacity RNA-to-cDNA Master Mix with No RT Control	500 reactions	4390713
	200 reactions	4390712
	50 reactions	4390711
	15 reactions	4390710
SuperScript® VILO™ cDNA Synthesis Kit	50 reactions	4453650
	250 reactions	4453651
Tris-EDTA (TE) buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0, made using DNase-free, RNase-free, sterile-filtered water)		AM9849
RNase inhibitor		N8080119

Reagents for PCR

Reagents		Applied Biosystems part number
TaqMan® Gene Expression Assays, inventoried		4331182
TaqMan® Gene Expression Assays, made-to-order		4351372
Custom TaqMan® Gene Expression Assays	Small-scale, 20X (144 × 50-µL reactions)	4331348
	Medium-scale, 20X (300 × 50-µL reactions)	4332078
	Large-scale, 60X (1160 × 50-µL reactions)	4332079
Custom TaqMan® probes and primers		Various†
Nuclease-free Water (not DEPC-treated)		AM9930

† Probes and primers synthesized by Applied Biosystems to your sequence and choice of quencher and reporter dyes. See Applied Biosystems web site: www.appliedbiosystems.com.

Optional reagents


See “Optional reagents” on page 55 for a list of optional user-supplied reagents.

Reaction plates, consumables, and equipment

See:

- “Reaction plates and accessories” on page 53 for a list of reaction plates that are compatible with your real-time PCR system.
- “General laboratory equipment” on page 56 for a list of required and optional laboratory consumables and equipment.

Workflow



Perform reverse transcription

Synthesis of single-stranded cDNA from total RNA samples is the first step in the two-step RT-PCR.

Recommended kits and reagents

To obtain cDNA from RNA samples, Applied Biosystems recommends the reverse transcription kits and reagents listed in “[Reagents for reverse transcription](#)” on [page 14](#).

Perform reverse transcription

For detailed RT guidelines and instructions, refer to the protocol for your reverse transcription kit. To download the protocols for Applied Biosystems kits, go to:

<http://docs.appliedbiosystems.com/search.taf>

Perform real-time PCR amplification

Target amplification using cDNA as the template is the second step in the two-step RT-PCR. In this step, the DNA polymerase amplifies target cDNA using sequence-specific primers and cleaves the TaqMan® probe to generate a fluorescent signal that is measured by the real-time PCR system. (For more information on TaqMan® chemistry, see [Appendix C on page 61](#).)



Note: If you choose to use custom TaqMan® probes and sequence detection primers, rather than a TaqMan® Gene Expression Assay or a Custom TaqMan® Gene Expression Assay, see [“Using custom TaqMan® primer and probe sets” on page 21](#).

PCR reagent handling and preparation

Follow these guidelines to ensure optimal PCR performance:

- Keep the TaqMan® assays in the freezer, away from light, until you are ready to use them. Excessive exposure to light may affect the fluorescent probes.
- Before use:
 - Thoroughly mix the TaqMan® Fast Advanced Master Mix by swirling the bottle.
 - Thaw frozen TaqMan® assays by placing them on ice. When thawed, resuspend the assays by vortexing, then briefly centrifuge.
 - Thaw frozen samples by placing them on ice. When thawed, resuspend the samples by briefly vortexing, then briefly centrifuge.

Prepare the PCR reaction mix

- Determine the total number of PCR reactions to perform. On each reaction plate, include:
 - A gene expression assay for each cDNA sample
 - No template controls (NTCs) for each assay on the plate
- !** **IMPORTANT!** You can run multiple assays on one reaction plate. Include controls for each assay that you run on a plate.

Applied Biosystems recommends that you perform 3 replicates of each reaction.

- Per the table below, calculate the total volume required for each reaction component.



Note: Include extra volume to compensate for the volume loss that occurs during pipetting.

Component	Volume (µL) for 1 reaction		Final concentration
	384-well plate	96-well and 48-well plates (both Standard and Fast)	
TaqMan® Fast Advanced Master Mix (2X)	5.0	10.0	1X
TaqMan® Gene Expression Assay (20X) or Custom TaqMan® Gene Expression Assay (20X) [†]	0.5	1.0	1X
Nuclease-free Water	3.5	7.0	—
Total volume per reaction	9.0	18.0	—

[†] See allgenes.com for assay information.

- Label a 1.5-mL microcentrifuge tube, add all components to the labeled tube, then cap the tube.
- Vortex the tube briefly to mix the components.
- Centrifuge the tube briefly to spin down the contents and eliminate air bubbles.

Prepare the PCR reaction plate

- Transfer the appropriate volume of PCR reaction mix to each well of an optical reaction plate:

Reaction plate format	Reaction volume (µL)
384-well plate	9.0
96-well and 48-well plates (both Standard and Fast)	18.0

- Cover the reaction plate with an optical adhesive film.
- Centrifuge the plate briefly to spin down the contents and eliminate air bubbles.

- Remove the optical adhesive film.
- Per the table below, add cDNA template or water to each well.

Component	Volume (µL) for 1 reaction	
	384-well plate	96-well and 48-well plates (both Standard and Fast)
cDNA template + Nuclease-free Water†	1.0	2.0
Nuclease-free Water (for the NTC reactions)	1.0	2.0
Total volume per reaction	10.0	20.0

† Use 100 ng to 1 pg of cDNA diluted in Nuclease-free Water.

- Cover the reaction plate with a new optical adhesive film.
- Centrifuge the plate briefly to spin down the contents and eliminate air bubbles.
- If required by your real-time PCR system, apply a compression pad to the plate.

Set up a plate document or experiment

When you set up a plate document or experiment, use the following thermal-cycling conditions:

- Thermal-cycling profile:

Applied Biosystems Real-Time PCR System	Thermal-cycling profile				
	Parameter	UNG incubation†	Polymerase activation‡	PCR (40 cycles)	
		Hold	Hold	Denature	Anneal/extend
Temp. (°C)	50	95	95	60	
7900HT system	Time (mm:ss)	02:00	00:20	00:01	00:20
7900HT Fast system (Fast 96-Well, Standard 96-Well, or 384-Well Block Modules)					
ViiA™ 7 system					
StepOne™ system					
StepOnePlus™ system					
7500 Fast system	Time (mm:ss)	02:00	00:20	00:03	00:30
7500 system					
7300 system					

† Required for optimal UNG activity.

‡ Required to activate the AmpliTaq® Fast DNA Polymerase.

- Run mode:

Applied Biosystems Real-Time PCR System	Default run mode
7900HT system	Standard
7900HT Fast system (384-Well and Standard 96-Well Block Modules)	
7500 system	
7300 system	
ViiA™ 7 system	Fast
StepOne™ system	
StepOnePlus™ system	
7900HT Fast system (Fast 96-Well Block Module)	
7500 Fast system	

- Sample volume:

Reaction plate format	Reaction volume (µL)
384-well plate	10.0
96-well and 48-well plates (both Standard and Fast)	20.0

Run the PCR reaction plate

1. In the system software, open the plate document or experiment that corresponds to the reaction plate.
2. Load the reaction plate into the real-time PCR system.
3. Start the run.

Analyze the data

Data analysis varies depending on the real-time PCR system that you use. The general process for analyzing gene expression quantitation data involves:

1. Viewing the amplification plots for the entire reaction plate.
2. Setting the baseline and threshold values to determine the threshold cycles (C_T) for the amplification curves.
3. Using the relative standard curve method or the comparative C_T method to analyze the data.

For detailed analysis information, refer to the appropriate documentation for your instrument. See [“Product documentation” on page 75](#).

Using custom TaqMan® primer and probe sets

To design custom probes and primers for a quantitation experiment:

- Determine the target template and amplicon (this page).
- Design the custom sequence detection primers and TaqMan® probe(s) (this page).
- Determine the optimal concentrations of the sequence detection primers (this page).
- Determine the optimal concentration of the custom TaqMan probe(s) (page 24).
- Perform real-time quantitative PCR (page 26).

Determine the target template and amplicon

A target template is a DNA sequence, including a cDNA, gDNA, or plasmid nucleotide sequence.

Design primers to amplify amplicons (short segments of DNA) within the target sequence. The shortest amplicons work the best. Consistent results are obtained for amplicon size ranges from 50 to 150 bp.

Design the custom TaqMan® probes and primers

You can design primers and probes with the Primer Express® Software. Refer to the *Primer Express® Software Version 3.0 Getting Started Guide*.

Determine the optimal primer concentration

The purpose of this procedure is to determine the minimum primer concentrations that yield the minimum threshold cycle (C_T) for each probe target without a decrease in the ΔR_n .



Note: Applied Biosystems real-time PCR systems can provide additional data for optimization using the minimum C_T . Refer to your instrument documentation for more information.

To determine the optimal primer concentration:

1. Per the table below, prepare a PCR reaction mix.

Reaction component	Volume (μ L) for 1 reaction		Final concentration
	384-well plate	96-well and 48-well plates (both Standard and Fast)	
TaqMan® Fast Advanced Master Mix	5.0	10.0	1×
Forward primer	1.0	2.0	50 to 900 nM
Reverse primer	1.0	2.0	50 to 900 nM
TaqMan® probe (2.5- μ M)	1.0	2.0	250 nM
Nuclease-free Water	1.0	2.0	—
Total volume per reaction	9.0	18.0	—

2. In an optical reaction plate, load at least 3 replicates of each of the following conditions:

Condition	Forward primer (nM)	Reverse primer (nM)
1	50	50
2	300	300
3	900	900

3. Run the reaction plate on your real-time PCR system using the following thermal-cycling profile:

Applied Biosystems Real-Time PCR System	Thermal-cycling profile				
	Parameter	UNG incubation [†]	Polymerase activation [‡]	PCR (40 cycles)	
		Hold	Hold	Denature	Anneal/ extend
	Temp. (°C)	50	95	95	60
7900HT system	Time (mm:ss)	02:00	00:20	00:01	00:20
7900HT Fast system (Fast 96-Well, Standard 96-Well, or 384-Well Block Modules)					
ViiA™ 7 system					
StepOne™ system					
StepOnePlus™ system					
7500 Fast system	Time (mm:ss)	02:00	00:20	00:03	00:30
7500 system					
7300 system					

† Required for optimal UNG activity.

‡ Required to activate the AmpliTaq® Fast DNA Polymerase.

4. At the end of run, tabulate the results for the C_T and ΔR_n . Choose the minimum forward- and reverse-primer concentrations that yield the maximum ΔR_n and the lowest C_T .

Determine the optimal probe concentration

The purpose of this procedure is to determine the minimum probe concentrations that yield the minimum C_T for each probe target without decreasing the ΔR_n .

Most TaqMan® assays are designed and run following Applied Biosystems assay development guidelines. A concentration of 900-nM primers and a 250-nM probe provides a highly reproducible and sensitive assay.

To determine the optimal probe concentration for single-probe assays:

1. Per the table below, prepare a PCR reaction mix.

Reaction component	Volume (µL) for 1 reaction		Final concentration
	384-well plate	96-well and 48-well plates (both Fast and Standard)	
TaqMan® Fast Advanced Master Mix	5.0	10.0	1X
Forward primer	1.0	2.0	Optimal [†]
Reverse primer	1.0	2.0	Optimal [†]
TaqMan® probe (2.5-µM)	1.0	2.0	50 to 250 nM
Nuclease-free Water	1.0	2.0	—
Total volume per reaction	9.0	18.0	—

[†] Use the forward- and reverse-primer concentrations determined in [“Determine the optimal primer concentration”](#) on page 22.

2. In an optical reaction plate, load at least 3 replicates at each 50-nM interval from 50 to 250 nM

3. Run the plate on your real-time PCR system using the following thermal-cycling profile:

Applied Biosystems Real-Time PCR System	Thermal-cycling profile				
	Parameter	UNG incubation [†]	Polymerase activation [‡]	PCR (40 cycles)	
		Hold	Hold	Denature	Anneal/ extend
	Temp. (°C)	50	95	95	60
7900HT system	Time (mm:ss)	02:00	00:20	00:01	00:20
7900HT Fast system (Fast 96-Well, Standard 96-Well, or 384-Well Block Modules)					
ViiA™ 7 system					
StepOne™ system					
StepOnePlus™ system					
7500 Fast system	Time (mm:ss)	02:00	00:20	00:03	00:30
7500 system					
7300 system					

† Required for optimal UNG activity.

‡ Required to activate the AmpliTaq® Fast DNA Polymerase.

4. Tabulate the results for the C_T and ΔR_n . Choose the minimum probe concentrations that yield the minimum C_T and the highest ΔR_n .

Perform real-time PCR

For routine assays that are optimized as described here, perform real-time quantitative PCR as follows:

1. Per the table below, prepare a PCR reaction mix, then transfer the appropriate volume to each well of an optical reaction plate.

Reaction component	Volume (µL) for 1 reaction		Final concentration
	384-well plate	96-well and 48-well plates (both Fast and Standard)	
TaqMan® Fast Advanced Master Mix	5.0	10.0	1×
Forward primer	1.0	2.0	Optimal [†]
Reverse primer	1.0	2.0	Optimal [†]
TaqMan® probe (2.5-µM)	1.0	2.0	Optimal [‡]
Nuclease-free Water	1.0	2.0	—
Total volume per reaction	9.0	18.0	—

[†] Use the forward- and reverse-primer concentrations determined in “Determine the optimal primer concentration” on page 22.

[‡] Use the probe concentration determined in “Determine the optimal probe concentration” on page 24.

2. Per the table below, add cDNA template or water to each well.

Component	Volume (µL) for 1 reaction	
	384-well plate	96-well and 48-well plates (both Standard and Fast)
cDNA template + Nuclease-free Water [†]	1.0	2.0
Nuclease-free Water (for the NTC reactions)	1.0	2.0
Total volume per reaction	10.0	20.0

[†] Use 100 ng to 1 pg of cDNA diluted in Nuclease-free Water.

3. Run the plate on your real-time PCR system using the following thermal-cycling profile:

Applied Biosystems Real-Time PCR System	Thermal-cycling profile				
	Parameter	UNG incubation [†]	Polymerase activation [‡]	PCR (40 cycles)	
		Hold	Hold	Denature	Anneal/ extend
Temp. (°C)	50	95	95	60	
7900HT system	Time (mm:ss)	02:00	00:20	00:01	00:20
7900HT Fast system (Fast 96-Well, Standard 96-Well, or 384-Well Block Modules)					
ViiA™ 7 system					
StepOne™ system					
StepOnePlus™ system					
7500 Fast system	Time (mm:ss)	02:00	00:20	00:03	00:30
7500 system					
7300 system					

[†] Required for optimal UNG activity.

[‡] Required to activate the AmpliTaq® Fast DNA Polymerase.

Section 2 TaqMan® Array Plates

This section provides a general protocol for performing two-step RT-PCR. In the PCR step, you use the TaqMan® Fast Advanced Master Mix with the TaqMan® Array Plates. The TaqMan Array Plates are TaqMan® Gene Expression Assays that are dried in MicroAmp® Optical 96-Well Reaction Plates and MicroAmp® Fast 96-Well Optical Reaction Plates.

This section covers:

■ User-supplied materials	29
■ Workflow	31
■ Perform reverse transcription	32
■ Perform real-time PCR amplification	32
■ Analyze the data	35

For more information

For more detailed information about the procedures outlined in this protocol, refer to the appropriate documentation for your instrument. See [“Product documentation” on page 75](#).

User-supplied materials

Reagents for reverse transcription

Reagents		Applied Biosystems part number
High Capacity cDNA Reverse Transcription Kit	200 reactions	4368814
	200 reactions with RNase Inhibitor	4374966
	1000 reactions	4368813
	1000 reactions with RNase Inhibitor	4374967
High Capacity RNA-to-cDNA™ Kit	50 reactions	4387406
High Capacity RNA-to-cDNA Master Mix	500 reactions	4390779
	200 reactions	4390778
	50 reactions	4390777
	15 reactions	4390776
High Capacity RNA-to-cDNA Master Mix with No RT Control	500 reactions	4390713
	200 reactions	4390712
	50 reactions	4390711
	15 reactions	4390710
SuperScript® VILO™ cDNA Synthesis Kit	50 reactions	4453650
	250 reactions	4453651

Reagents	Applied Biosystems part number
Tris-EDTA (TE) buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0, made using DNase-free, RNase-free, sterile-filtered water)	AM9849
RNase inhibitor	N8080119

Reagents for PCR

Reagents	Applied Biosystems part number
TaqMan® Array Standard 96 well Plates	Std_96
TaqMan® Array Fast 96 well Plates	Fast_96
TaqMan® Array Gene Signature Plates	Various [†]
Nuclease-free Water (not DEPC-treated)	AM9930

[†] See the Applied Biosystems web site: www.appliedbiosystems.com.

Optional reagents


See “Optional reagents” on page 55 for a list of optional user-supplied reagents.

Reaction plates, consumables, and equipment

See:

- “Reaction plates and accessories” on page 53 for a list of reaction plates that are compatible with your real-time PCR system.
- “General laboratory equipment” on page 56 for a list of required and optional laboratory consumables and equipment.

Workflow



Perform reverse transcription

Synthesis of single-stranded cDNA from total RNA samples is the first step in the two-step RT-PCR.

Recommended kits and reagents

To obtain cDNA from RNA samples, Applied Biosystems recommends the reverse transcription kits and reagents listed in “[Reagents for reverse transcription](#)” on [page 29](#).

Perform reverse transcription

For detailed RT guidelines and instructions, refer to the protocol for your reverse transcription kit. To download the protocols for Applied Biosystems kits, go to:

<http://docs.appliedbiosystems.com/search.taf>

Perform real-time PCR amplification

Target amplification using cDNA as the template is the second step in the two-step RT-PCR. In this step, the DNA polymerase amplifies target cDNA using sequence-specific primers and cleaves the TaqMan® probe to generate a fluorescent signal that is measured by the real-time PCR system. (For more information on TaqMan® chemistry, see [Appendix C on page 61](#).)

PCR reagent handling and preparation

Follow these guidelines to ensure optimal PCR performance:

- Keep the TaqMan® Array Plates away from light until you are ready to use them. Excessive exposure to light may affect the fluorescent probes.
- Before use:
 - Thoroughly mix the TaqMan® Fast Advanced Master Mix by swirling the bottle.
 - Thaw frozen samples by placing them on ice. When thawed, resuspend the samples by briefly vortexing, then briefly centrifuge.

Prepare the PCR reaction mix

1. For each cDNA sample, determine the total number of PCR reactions to perform (1 reaction = 1 well in the TaqMan Array Plate).
2. Per the table below, calculate the total volume required for each reaction component.



Note: Include extra volume to compensate for the volume loss that occurs during pipetting.

Component	Volume (µL) for 1 reaction	
	Fast 96-well plate	Standard 96-well plate
cDNA template + Nuclease-free Water [†]	10.0	10.0
TaqMan® Fast Advanced Master Mix (2X)	10.0	10.0
Total volume per reaction	20.0	20.0

[†] Use 5 to 50 ng of cDNA diluted in Nuclease-free Water.

Prepare the PCR reaction plate

1. Transfer the appropriate volume of PCR reaction mix to each well of an optical reaction plate:

Reaction plate format	Reaction volume (µL)
Fast 96-well plate	20.0
Standard 96-well plate	20.0

2. Cover the reaction plate with an optical adhesive film.
3. Centrifuge the plate briefly to spin down the contents and eliminate air bubbles.
4. If required by your real-time PCR system, apply a compression pad to the plate.

Set up a plate document or experiment

When you set up a plate document or experiment, use the following thermal-cycling conditions:

- Thermal-cycling profile:

Applied Biosystems Real-Time PCR System	Thermal-cycling profile				
	Parameter	UNG incubation [†]	Polymerase activation [‡]	PCR (40 cycles)	
		Hold	Hold	Denature	Anneal/ extend
Temp. (°C)	50	95	95	60	
7900HT system	Time (mm:ss)	02:00	00:20	00:01	00:20
7900HT Fast system (Fast 96-Well or Standard 96-Well)					
ViiA™ 7 system					
StepOnePlus™ system					
7500 Fast system	Time (mm:ss)	02:00	00:20	00:03	00:30
7500 system					
7300 system					

[†] Required for optimal UNG activity.

[‡] Required to activate the AmpliTaq® Fast DNA Polymerase.

- Run mode:

Applied Biosystems Real-Time PCR System	Default run mode
7900HT system	Standard
7900HT Fast system (Standard 96-Well Block Modules)	
7500 system	
7300 system	
ViiA™ 7 system	Fast
StepOnePlus™ system	
7900HT Fast system (Fast 96-Well Block Module)	
7500 Fast system	

- Sample volume: 20 µL

Run the PCR reaction plate

1. In the system software, open the plate document or experiment that corresponds to the reaction plate.
2. Load the reaction plate into the real-time PCR system.
3. Start the run.

Analyze the data

Data analysis varies depending on the real-time PCR system that you use. The general process for analyzing gene expression quantitation data involves:

1. Viewing the amplification plots for the entire reaction plate.
2. Setting the baseline and threshold values to determine the threshold cycles (C_T) for the amplification curves.
3. Using the relative standard curve method or the comparative C_T method to analyze the data.

For detailed analysis information, refer to the appropriate documentation for your instrument. See [“Product documentation” on page 75](#).

Section 3 TaqMan® MicroRNA Assays

This section provides a general protocol for performing two-step RT-PCR. In the PCR step, you use the TaqMan® Fast Advanced Master Mix with TaqMan® MicroRNA Assays. The TaqMan MicroRNA Assays are designed to detect and accurately quantitate mature microRNAs (miRNAs) using Applied Biosystems real-time PCR systems.

This section covers:

- User-supplied materials 37
- Workflow 39
- Perform reverse transcription 40
- Perform real-time PCR amplification 40
- Analyze the data 43

For more information

For more detailed information about the procedures outlined in this protocol, refer to the appropriate documentation for your instrument. See [“Product documentation” on page 75](#).

User-supplied materials

Reagents for reverse-transcription

Reagents		Applied Biosystems part number
TaqMan® MicroRNA Reverse Transcription Kit†	200 reactions	4366596
	1000 reactions	4366597
Tris-EDTA (TE) buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0, made using DNase-free, RNase-free, sterile-filtered water)		AM9849
RNase inhibitor		N8080119
Nuclease-free Water (not DEPC-treated)		AM9930

† TaqMan® MicroRNA Assays are specifically optimized to work with the TaqMan® MicroRNA Reverse Transcription Kit. Applied Biosystems cannot guarantee assay performance if you use other reverse transcription kits.

Reagents for PCR

The TaqMan MicroRNA Assays are available for a range of species. Because many mature miRNA sequences are identical across related species, many human assays are also valid for mouse and rat assays.

Reagents	Applied Biosystems part number
TaqMan® MicroRNA Assays	Various†
Custom Small RNA Assays (Early-Access)	4398987, 4398988, 4398989

† See the Applied Biosystems web site: miRNA.appliedbiosystems.com.

Optional reagents


See “Optional reagents” on page 55 for a list of optional user-supplied reagents.

Reaction plates, consumables, and equipment

See:

- “Reaction plates and accessories” on page 53 for a list of reaction plates and that are compatible with your real-time PCR system.
- “General laboratory equipment” on page 56 for a list of required and optional laboratory consumables and equipment.

Workflow



Perform reverse transcription

Synthesis of single-stranded cDNA from total RNA samples is the first step in the two-step RT-PCR.

Recommended kit

To obtain cDNA from RNA samples, Applied Biosystems recommends the TaqMan® MicroRNA Reverse Transcription Kit. TaqMan® MicroRNA Assays are specifically optimized to work with the TaqMan MicroRNA Reverse Transcription Kit. Applied Biosystems cannot guarantee assay performance if you use other reverse transcription kits. (See [“Reagents for reverse-transcription” on page 37](#) for the kit part numbers.)

Perform reverse transcription

For detailed RT guidelines and instructions, refer to the protocol provided with the TaqMan MicroRNA Reverse Transcription Kit.

Perform real-time PCR amplification


Target amplification using cDNA as the template is the second step in the two-step RT-PCR. In this step, the DNA polymerase amplifies target cDNA using sequence-specific primers and cleaves the TaqMan® probe to generate a fluorescent signal that is measured by the real-time PCR system. (For more information on TaqMan® chemistry, see [Appendix C on page 61](#).)

Reagent preparation guidelines

- Keep all TaqMan® MicroRNA Assays in the freezer, protected from light, until you are ready to use them. Excessive exposure to light may affect the fluorescent probes.
- Before use, mix the TaqMan Fast Advanced Master Mix thoroughly by swirling the bottle.
- Prepare the PCR reaction mix before transferring to the reaction plate for thermal cycling.

Prepare the PCR reaction mix

- Determine the total number of PCR reactions to perform. On each reaction plate include:
 - A mircoRNA assay for each cDNA sample
 - Endogenous control assays
 - No template controls (NTCs) for each assay on the plate

 **IMPORTANT!** You can run multiple assays on one reaction plate. Include controls for each assay that you run on a plate.

Applied Biosystems recommends that you perform 3 replicates of each reaction.

- Per the table below, calculate the total volume required for each reaction component.



Note: Include extra volume to compensate for the volume loss that occurs during pipetting.

Component	Volume (µL) for 1 reaction	
	384-well plate	96-well and 48-well plates (both Fast and Standard)
TaqMan® Fast Advanced Master Mix (2X)	5.00	10.00
Nuclease-free Water	3.83	7.67
TaqMan® MicroRNA Assay (20X)	0.50	1.00
Product from RT reaction [†]	0.67	1.33
Total volume per reaction	10.00	20.00

[†] The minimum final dilution of the RT reaction in the PCR reaction is 1:15.

- Combine the TaqMan Fast Advanced Master Mix and water:
 - Label a 1.5-mL microcentrifuge tube.
 - Add the TaqMan Fast Advanced Master Mix and water to the labeled tube, then cap the tube.
 - Mix gently, then centrifuge to bring the solution to the bottom of the tube.
 - Transfer the solution to a polypropylene tube (the PCR reaction tube). (Necessary to transfer to a different tube? Nathalie to check this.)
- Add the TaqMan® MicroRNA Assay to the PCR reaction tube.
- Add the RT reaction product to the PCR reaction tube.
- Mix gently, then centrifuge to bring the solution to the bottom of the plate.

Prepare the PCR reaction plate

1. Transfer the appropriate volume of PCR reaction mix to each well of an optical reaction plate:

Reaction plate format	Reaction volume (µL)
384-well plate	10.0
96-well and 48-well plates (both Fast and Standard)	20.0

2. Cover the reaction plate with an optical adhesive film.
3. Centrifuge the plate briefly to spin down the contents and eliminate air bubbles.
4. If required by your real-time PCR system, apply a compression pad to the plate.

Set up a plate document or experiment

When you set up a plate document or experiment, use the following thermal-cycling conditions:

- Thermal-cycling profile:

Applied Biosystems Real-Time PCR System	Thermal-cycling profile				
	Parameter	UNG incubation [†]	Polymerase activation [‡]	PCR (40 cycles)	
		Hold	Hold	Denature	Anneal/ extend
	Temp. (°C)	50	95	95	60
7900HT system	Time (mm:ss)	02:00	00:20	00:01	00:20
7900HT Fast system (Fast 96-Well, Standard 96-Well, or 384-Well Block Modules)					
ViiA™ 7 system					
StepOne™ system					
StepOnePlus™ system					
7500 Fast system	Time (mm:ss)	02:00	00:20	00:03	00:30
7500 system					
7300 system					

[†] Required for optimal UNG activity.

[‡] Required to activate the AmpliTaq® Fast DNA Polymerase.

- Run mode:

Applied Biosystems Real-Time PCR System	Default run mode
7900HT system	Standard
7900HT Fast system (384-Well and Standard 96-Well Block Modules)	
7500 system	
7300 system	
ViiA™ 7 system	Fast
StepOne™ system	
StepOnePlus™ system	
7900HT Fast system (Fast 96-Well Block Module)	
7500 Fast system	

- Sample volume:

Reaction plate format	Reaction volume (µL)
384-well plate	10.0
96-well and 48-well plates (both Standard and Fast)	20.0

Run the PCR reaction plate

1. In the system software, open the plate document or experiment that corresponds to the reaction plate.
2. Load the reaction plate into the real-time PCR system.
3. Start the run.

Analyze the data

Data analysis varies depending on the real-time PCR system that you use. The general process for analyzing gene expression quantitation data involves:

1. Viewing the amplification plots for the entire reaction plate.
2. Setting the baseline and threshold values to determine the threshold cycles (C_T) for the amplification curves.
3. Using the relative standard curve method or the comparative C_T method to analyze the data.

For detailed analysis information, refer to the appropriate documentation for your instrument. See [“Product documentation” on page 75](#).

Section 4 TaqMan® Array Micro Fluidic Cards

This section provides a general protocol for performing two-step RT-PCR. In the PCR step, you use the TaqMan® Fast Advanced Master Mix with TaqMan® Array Micro Fluidic Cards. The TaqMan Array cards are 384-well cards preloaded with TaqMan® Gene Expression Assays. The TaqMan Array cards allow you to measure gene expression using the comparative C_T ($\Delta\Delta C_T$) method of relative quantitation. You can run 1 to 8 samples per card, against 12 to 384 TaqMan Gene Expression Assay targets (including controls).

This section covers:

- User-supplied materials 45
- Workflow 47
- Perform reverse transcription 48
- Perform real-time PCR amplification 49
- Analyze the data 51

For more information

For more detailed information about the procedures outlined in this protocol, refer to the appropriate documentation for your instrument. See [“Product documentation” on page 75](#).

User-supplied materials

Reagents for reverse transcription

Reagents		Applied Biosystems part number
High Capacity cDNA Reverse Transcription Kit	200 reactions	4368814
	200 reactions with RNase Inhibitor	4374966
	1000 reactions	4368813
	1000 reactions with RNase Inhibitor	4374967
High Capacity RNA-to-cDNA™ Kit	50 reactions	4387406
High Capacity RNA-to-cDNA Master Mix	500 reactions	4390779
	200 reactions	4390778
	50 reactions	4390777
	15 reactions	4390776

Reagents		Applied Biosystems part number
High Capacity RNA-to-cDNA Master Mix with No RT Control	500 reactions	4390713
	200 reactions	4390712
	50 reactions	4390711
	15 reactions	4390710
SuperScript® VILO™ cDNA Synthesis Kit	50 reactions	4453650
	250 reactions	4453651
Tris-EDTA (TE) buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0, made using DNase-free, RNase-free, sterile-filtered water)		AM9849
RNase inhibitor		N8080119

Reagents for PCR

Reagents		Applied Biosystems part number
TaqMan® Array Custom Micro Fluidic Card	One TaqMan Array card; 384 wells per card	Various [†]
TaqMan® Gene Set	One TaqMan Array card; 384 wells per card	
TaqMan® Gene Signature Array	The endogenous control Gene Signature Arrays come in sets of two cards; the remaining Gene Signature Arrays come in sets of four cards. 384 wells per card.	
Nuclease-free Water (not DEPC-treated)		AM9930

[†] See the Applied Biosystems web site: taqmanarray.appliedbiosystems.com.

Optional reagents


See “Optional reagents” on page 55 for a list of optional user-supplied reagents.

Reaction plates, consumables, and equipment

See:

- “Reaction plates and accessories” on page 53 for a list of reaction plates and that are compatible with your real-time PCR system.
- “General laboratory equipment” on page 56 for a list of required and optional laboratory consumables and equipment.

Workflow



Perform reverse transcription

Synthesis of single-stranded cDNA from total RNA samples is the first step in the two-step RT-PCR.

Recommended kits and reagents

To obtain cDNA from RNA samples, Applied Biosystems recommends the reverse transcription kits and reagents listed in “[Reagents for reverse transcription](#)” on [page 45](#).

Perform reverse transcription

For detailed RT guidelines and instructions, refer to the:

- *TaqMan® Array Micro Fluidic Card User Guide*
and
- Protocol for your reverse transcription kit. To download the protocols for Applied Biosystems kits, go to:
<http://docs.appliedbiosystems.com/search.taf>

Perform real-time PCR amplification

Target amplification using cDNA as the template is the second step in the two-step RT-PCR. In this step, the DNA polymerase amplifies target cDNA using sequence-specific primers and cleaves the TaqMan® probe to generate a fluorescent signal that is measured by the real-time PCR system. (For more information on TaqMan® chemistry, see [Appendix C on page 61.](#))

Card handling and preparation

Follow these guidelines to ensure optimal PCR performance:

- Do not remove a TaqMan Array card from its packaging until the packaging has reached room temperature and you are ready to fill it with sample-specific PCR mix. Prolonged exposure to indoor lighting can photo-degrade the fluorescent probes contained within the card. Do not expose the card to sunlight.
- Fill each fill reservoir with sample-specific PCR mix made from a single cDNA sample.
- Add 100 µL of the sample-specific PCR mix per fill reservoir to ensure adequate filling. Volumes smaller than 100 µL will result in insufficiently filled cards.
- Do not add sample after centrifuging the cards. When you centrifuge the card, the sample-specific PCR mix resuspends the dried TaqMan® probes and primers within the wells of the card. Adding sample after centrifuging disrupts the resuspended assay positions.
- To ensure a high degree of reproducibility, Applied Biosystems recommends scheduling your runs so that each card is run as soon as possible. After sealing, there is no measurable well-to-well contamination for up to 64 hours.

Recommended amounts of cDNA

For the PCR reactions, Applied Biosystems recommends that you use:

- 30 to 1000 ng (0.3 to 10 ng/mL) of cDNA (converted from total RNA) per fill reservoir. The amount to use depends on the expression level of your target genes and the number of target copies per well that you need to detect. For example, you can use:
 - 1000 ng (10 ng/µL) per fill reservoir to detect genes with low expression. However, the cDNA concentration will be high (10 ng/µL); be sure to use high-quality cDNA without inhibitors.
 - 100 to 200 ng per fill reservoir to detect genes with moderate expression.
 - 30 to 50 ng per fill reservoir to detect genes with moderate to high expression.
- The same amount of cDNA sample for all reactions.

Prepare the PCR reaction mix

1. For each cDNA sample, determine the total number of fill reservoirs in the TaqMan Array card that you will use.
2. Per the table below, calculate the total volume required for each reaction component.



Note: Include extra volume to compensate for the volume loss that occurs during pipetting.

Reaction Component	Volume (µL) for 1 fill reservoir
cDNA template + Nuclease-free Water [†]	50.0
TaqMan® Fast Advanced Master Mix (2X)	50.0
Total volume	100.0

[†] Each sample-specific PCR reaction mix should contain 30 to 1000 ng (0.3 to 10 ng/µL) of total RNA converted to cDNA. The amount of cDNA to add depends upon the abundance of the specific gene transcript.

3. If frozen, thaw the cDNA samples on ice. Resuspend the cDNA samples by inverting the tube, then gently vortexing.
4. Mix the master mix thoroughly by swirling the bottle.
5. For each sample, label a 1.5-mL microcentrifuge tube, then add the required components to the labeled tube.
6. Cap the microcentrifuge tubes, then gently vortex the tubes to thoroughly mix the solution.
7. Briefly centrifuge the tubes to spin down the contents and eliminate air bubbles.

Prepare the TaqMan Array card

Fill, centrifuge, and seal the TaqMan Array card according to the *TaqMan® Array Micro Fluidic Card User Guide*.

Set up a plate document or experiment

When you set up a plate document or experiment, use the following thermal-cycling conditions:

- Thermal-cycling profile:

Applied Biosystems Real-Time PCR System	Thermal-cycling profile				
	Parameter	UNG incubation [†]	Polymerase activation [‡]	PCR (40 cycles)	
		Hold	Hold	Denature	Anneal/extend
	Temp. (°C)	50	92	97	62
7900HT Fast system (TaqMan® Array Micro Fluidic Card Thermal Cycling Block)	Time (mm:ss)	02:00	10:00	00:01	00:20

[†] Required for optimal UNG activity.

[‡] Required to activate the AmpliTaq® Fast DNA Polymerase.

- Default run mode: Standard
- Ramp rate: 100%
- Sample volume: 1.0 µL

Run the TaqMan Array card

1. In the system software, open the plate document or experiment that corresponds to the TaqMan Array card.
2. Load the card into the real-time PCR system.
3. Start the run.

Analyze the data

Data analysis varies depending on the real-time PCR system that you use. The general process for analyzing gene expression quantitation data involves:

1. Viewing the amplification plots for the entire reaction plate.
2. Setting the baseline and threshold values to determine the threshold cycles (C_T) for the amplification curves.
3. Using the relative standard curve method or the comparative C_T method to analyze the data.

For detailed analysis information, refer to the appropriate documentation for your instrument. See [“Product documentation” on page 75](#).

Additional User-Supplied Materials

This appendix covers:

- Reaction plates and accessories 53
- Optional reagents 55
- General laboratory equipment 56

Reaction plates and accessories

The table below lists the reaction plates and accessories available for Applied Biosystems real-time PCR systems.

Real-Time PCR System	Reaction plates and accessories
7300 system	<ul style="list-style-type: none"> • MicroAmp® Optical 96-Well Reaction Plate with Barcode: <ul style="list-style-type: none"> – 500 plates (PN 4326659) – 20 plates (PN 4306737) • MicroAmp® Optical Adhesive Film, 100 films (PN 4311971) • MicroAmp® Optical Film Compression Pad, 5 pads (PN 4312639) • MicroAmp® Optical 8-Cap Strip, 300 strips (PN 4323032)
7500 system	
7500 Fast system	<ul style="list-style-type: none"> • MicroAmp® Fast Optical 96-Well Reaction Plate with Barcode: <ul style="list-style-type: none"> – 200 plates (PN 4366932) – 20 plates (PN 4346906) • MicroAmp® Optical Adhesive Film, 100 films (PN 4311971) • MicroAmp® Optical 8-Cap Strip, 300 strips (PN 4323032)
ViiA™ 7 system, 384-Well Block Module	<ul style="list-style-type: none"> • MicroAmp® Optical 384-Well Reaction Plate with Barcode: <ul style="list-style-type: none"> – 1000 plates (PN 4343814) – 500 plates (PN 4326270) – 50 plates (PN 4309849) • MicroAmp® Optical 384-Well Reaction Plate, 1000 plates (PN 4343370) • MicroAmp® Optical Adhesive Film, 100 films (PN 4311971)
7900HT/7900HT Fast system, 384-Well Block Module	
7900HT/7900HT Fast system, Standard 96-Well Block Module	<ul style="list-style-type: none"> • MicroAmp® Optical 96-Well Reaction Plate with Barcode: <ul style="list-style-type: none"> – 500 plates (PN 4326659) – 20 plates (PN 4306737) • MicroAmp® Optical Adhesive Film, 100 films (PN 4311971) • MicroAmp® Optical Film Compression Pad, 5 pads (PN 4312639) • MicroAmp® Optical 8-Cap Strip, 300 strips (PN 4323032) • MicroAmp® Snap-On Optical Film Compression Pad, for use with the automation accessory (PN 4333292)

Appendix A Additional User-Supplied Materials
Reaction plates and accessories

Real-Time PCR System	Reaction plates and accessories
7900HT/7900HT Fast system, Fast 96-well Block Module	<ul style="list-style-type: none"> • MicroAmp® Fast Optical 96-Well Reaction Plate with Barcode: <ul style="list-style-type: none"> – 200 plates (PN 4366932) – 20 plates (PN 4346906) • MicroAmp® Optical Adhesive Film, 100 films (PN 4311971) • MicroAmp® Snap-On Optical Film Compression Pad, for use with the automation accessory (PN 4333292)
7900HT/7900HT Fast system, TaqMan® Array Micro Fluidic Card Thermal Cycling Block	TaqMan® Array Micro Fluidic Card
ViiA™ 7 system, TaqMan® Array Micro Fluidic Card Thermal Cycling Block	
StepOne™ system	<ul style="list-style-type: none"> • MicroAmp® Fast Optical 48-Well Reaction Plate, 20 plates (PN 4375816) • MicroAmp® 48-Well Optical Adhesive Film, 100 films (PN 4375323)
StepOnePlus™ system	<ul style="list-style-type: none"> • MicroAmp® Fast Optical 96-Well Reaction Plate with Barcode: <ul style="list-style-type: none"> – 200 plates (PN 4366932) – 20 plates (PN 4346906) • MicroAmp® Optical Adhesive Film, 100 films (PN 4311971)

Optional reagents

Item	Applied Biosystems part number
DNAZap™ Solution (2 bottles, 250 mL each)	AM9890
MagMAX™ AI/ND Viral RNA Isolation Kit, 50 purifications	AM1929
MagMAX™ Viral RNA Isolation Kit, 50 purifications	AM1939
mirVana™ miRNA Isolation Kit, 40 purifications	AM1560
RecoverAll™ Total Nucleic Acid Isolation Kit for FFPE, 40 purifications	AM1975
RiboPure™ Bacterial Kit	AM1925
RiboPure™ Blood Kit, 40 purifications	AM1928
RiboPure™ RNA Isolation Kit, 50 purifications	AM1924
RiboPure™ Yeast Kit	AM1926
RNAlater® ICE Solution, 25-mL	AM7030
RNAlater® Solution, 100 mL	AM7020
RNAqueous® -4PCR Kit, 30 purifications	AM1914
RNAqueous® Kit, 50 purifications	AM1912
RNaseZap® RNase Decontamination Solution, 250 mL	AM9780
RT-PCR Grade Water, 10, 1.75-mL bottles	AM9935
TRI Reagent®, 100-mL	AM9738
Turbo DNA-free™, 50 reactions	AM1907

General laboratory equipment

The following table lists required and optional laboratory equipment and consumables.

Item	Source	Part number	
Accessories for tubes of assay mixes: <ul style="list-style-type: none"> Decapper for single caps Decapper for eight caps TPE cap cluster for simultaneously capping 96 individual polypropylene tubes, 50 capmats/bag 	Micronic BV [‡]	---	
Centrifuge with plate adapter	Major laboratory suppliers (MLS)	---	
Disposable gloves	MLS	---	
Microcentrifuge	MLS	---	
Heat block or water bath or thermal cycler to 95 °C	MLS	---	
Microcentrifuge tubes, 1.5-mL	Applied Biosystems	AM12400	
Barrier (Filter) Tips	10 µL size - Pipetman™ (Ten 8 × 12 racks)	Applied Biosystems	AM12640
	10 µL size - Eppendorf® (Ten 8 × 12 racks)	Applied Biosystems	AM12635
	20 µL size (Ten 8 × 12 racks)	Applied Biosystems	AM12645
	1000 µL size (Ten 100 ct racks)	Applied Biosystems	AM12665
	200 µL size (Ten 8 x 12 racks)	Applied Biosystems	AM12655
Pipettes	Positive-displacement	MLS	---
	Air-displacement		
	Multichannel		
Vortexer	MLS	---	
Microsoft® Excel® software or equivalent spreadsheet and analysis software	Software suppliers	---	

[‡] See the Micronic BV web site: www.micronic.com. Other vendors supply similar products.

PCR Good Laboratory Practices

This appendix covers:

■ Sample preparation	57
■ Preventing contamination	58

Sample preparation

When preparing samples for PCR amplification:

- Use a positive-displacement pipette or aerosol-resistant pipette tips.
- Follow proper pipette-dispensing techniques to prevent aerosols.
- Wear clean gloves and a clean lab coat (not previously worn while handling amplified PCR products or used during sample preparation).
- Change gloves whenever you suspect that they are contaminated.
- Maintain separate areas and dedicated equipment and supplies for:
 - Sample preparation
 - PCR setup
 - PCR amplification
 - Analysis of PCR products
- Never bring amplified PCR products into the PCR setup area.
- Open and close all sample tubes carefully. Centrifuge tubes before opening. Try not to splash or spray PCR samples.
- Keep reactions and components capped as much as possible.
- Clean lab benches and equipment periodically with 10% bleach solution. Use DNAZap™ Solution (PN AM9890).

Preventing contamination

PCR assays require special laboratory practices to avoid false positive amplifications (Kwok and Higuchi, 1989). The high throughput and repetition of these assays can lead to amplification of a single DNA molecule (Saiki *et al.*, 1985; Mullis and Faloona, 1987).

False positives

Special laboratory practices are necessary in order to avoid false positive amplifications (Higuchi, *et al.*, 1989). This is because of the capability for single DNA molecule amplification provided by the PCR process (Saiki *et al.*, 1985; Mullis *et al.*, 1987; Saiki *et al.*, 1988). Because of the enormous amplification possible with PCR, amplicon carryover can result in sample contamination. Other sources of contamination could be from samples with high DNA levels or from positive control templates.

When dUTP replaces dTTP as a dNTP substrate in PCR and the method described below is used, UNG treatment can prevent the reamplification of carryover PCR products in subsequent experiments (Sninsky and Gelfand, *pers. comm.*) This method uses enzymatic and chemical reactions analogous to the restriction-modification and excision-repair systems of cells to degrade specifically PCR products from previous PCR amplifications or to degrade mis-primed, non-specific products produced prior to specific amplifications, but not degrade native nucleic acid templates.

The method used to make PCR products susceptible to degradation involves substituting dUTP for dTTP in the PCR mix and treating subsequent PCR mixes with the enzyme uracil-N glycosylase (UNG, EC 3.2.2-) prior to amplification (Longo *et al.*, 1990).

Although the protocol and reagents described here are capable of degrading or eliminating large numbers of carried over PCR products, we encourage users to continue using the specific devices and suggestions described in this protocol booklet and in Kwok (1990) and Higuchi (1989) to minimize cross-contamination from non-dU-containing PCR products or other samples.

Uracil-N glycosylase (UNG)

The UNG provided in the TaqMan[®] Fast Advanced Master Mix is a pure, nuclease-free, 26-kDa recombinant enzyme encoded by the *Escherichia coli* uracil-N glycosylase gene which has been inserted into an *E. coli* host to direct the expression of the native form of the enzyme (Kwok and Higuchi, 1989).

UNG acts on single- and double-stranded dU-containing DNA by hydrolyzing uracil-glycosidic bonds at dU-containing DNA sites. The enzyme causes the release of uracil, and creates an alkali-sensitive apyrimidic site in the DNA. Apyrimidic sites block replication by DNA polymerases. The enzyme has no activity on RNA or dT-containing DNA.

UNG incubation at 50 °C is necessary to cleave any dU-containing PCR carryover products. Ten-minute incubation at 95 °C is necessary to substantially reduce UNG activity, and to denature the native DNA in the experimental sample. Because UNG is not completely deactivated during the 95 °C incubation, it is important to keep the annealing temperatures greater than 55 °C and to refrigerate PCR products at 2 to 8 °C in order to prevent amplicon degradation.

Prevention of PCR product carryover

Use primers that contain dA nucleotides near the 3' ends so that any primer-dimer generated is efficiently degraded by UNG at least as well as any dU-containing PCR products. The further a dA nucleotide is from the 3' end, the more likely that partially degraded primer-dimer molecules may serve as templates for a subsequent PCR amplification.

Production of primer dimer could lower the amplification yield of the desired target region. If primers cannot be selected with dA nucleotides near the ends, the use of primers with 3' terminal dU-nucleotides should be considered. Single-stranded DNA with terminal dU nucleotides are not substrates for UNG (Delort et al., 1985) and thus the primers will not be degraded. Biotin-dUMP derivatives are not substrates for UNG.

The concentration of UNG and the time of the incubation step necessary to prevent amplification of contaminating dU-containing PCR product depends on the PCR conditions necessary to amplify your particular DNA sequence and the level of contamination expected. In most cases, using UNG at 1 U/100 mL reaction and incubation at 50 °C for two minutes is sufficient.

Do not attempt to use UNG in subsequent amplification of dU-containing PCR template, such as in nested-PCR protocols. The UNG will degrade the dU-containing PCR product, preventing further amplification.

Fluorescent contaminants

Since fluorescent contaminants may interfere with this assay and give false-positive results, it may be necessary to include a No-Amplification Control tube that contains sample, but no enzyme. If the absolute fluorescence of the No-Amplification Control is greater than that of the no template control (NTC) after PCR, fluorescent contaminants may be present in the sample or in the heat block of the thermal cycler.

About Two-Step RT-PCR

This appendix covers:

■ Overview	62
■ About TaqMan [®] MGB Probes.....	63
■ About AmpliTaq [®] Fast DNA Polymerase	63
■ About uracil-N glycosylase.....	63
■ About ROX [™] Passive Reference dye	63
■ About the 5' nuclease assay	64

Overview

Gene quantitation assays using TaqMan[®] Fast Advanced Master Mix and TaqMan[®] assays are performed in a two-step RT-PCR:

1. In the reverse transcription (RT) step, cDNA is reverse-transcribed from RNA.
2. In the PCR step, PCR products are quantitatively synthesized from cDNA samples using the TaqMan[®] Fast Advanced Master Mix.

The figure below illustrates two-step PCR.



 **Note:** Figure 1 does not show hybridization of the TaqMan[®] MGB probe. See Figure 2 on page 64 for details on how the TaqMan MGB probe is used in the PCR step.

Figure 1 Two-step RT-PCR



About TaqMan® MGB Probes

The TaqMan® MGB probes contain:

- A reporter dye (for example, FAM™ dye) linked to the 5' end of the probe.
- A minor groove binder (MGB) at the 3' end of the probe.

MGBs increase the melting temperature (T_m) without increasing probe length (Afonina *et al.*, 1997; Kutyaev *et al.*, 1997); they also allow for the design of shorter probes.

- A nonfluorescent quencher (NFQ) at the 3' end of the probe.

Because the quencher does not fluoresce, Applied Biosystems real-time PCR systems can measure reporter dye contributions more accurately.

About AmpliTaq® Fast DNA Polymerase

The AmpliTaq® Fast DNA Polymerase enzyme is purified through a proprietary process to reduce bacterial DNA introduced from the host organism. The purification process ensures that non-specific, false-positive DNA products due to bacterial DNA contamination are minimized during PCR.

When AmpliTaq® Fast DNA Polymerase is added to the reaction mixture at room temperature, the inactive enzyme is not capable of primer extension. Any low-stringency mispriming events that may have occurred will not be enzymatically extended and subsequently amplified. A thermal incubation step is required for activation to ensure that active enzyme is generated only at temperatures where the DNA is fully denatured.

About uracil-N glycosylase

Uracil-N glycosylase (UNG) treatment can prevent the reamplification of carryover-PCR products by removing any uracil incorporated into single- or double-stranded amplicons. (Longo *et al.*, 1990). UNG prevents reamplification of carryover-PCR products in an assay if all previous PCR for that assay was performed using a dUTP-containing master mix. See [“Preventing contamination” on page 58](#) for more information about UNG.





About ROX™ Passive Reference dye

The ROX™ Passive Reference dye provides an internal reference to which the reporter-dye signal can be normalized during data analysis. Normalization is necessary to correct for fluorescent fluctuations due to changes in concentration or volume.

About the 5' nuclease assay

The 5' nuclease assay process (Figure 2 through Figure 5) takes place during PCR amplification. This process occurs in every cycle and does not interfere with the exponential accumulation of product.

Legend for Figures 2 to 5:

-  **NFQ** = Nonfluorescent quencher
-  **MGB** = Minor groove binder
-  **R** = Reporter
-  **P** = Hot-start DNA polymerase

During PCR, the TaqMan[®] MGB probe anneals specifically to a complementary sequence between the forward and reverse primer sites.

When the probe is intact (Figure 2 and Figure 3), the proximity of the reporter dye to the quencher dye results in suppression of the reporter fluorescence primarily by Förster-type energy transfer (Förster, 1948; Lakowicz, 1983).

Figure 2 Polymerization






Figure 3 Strand displacement




The DNA polymerase cleaves only probes that are hybridized to the target (Figure 4 on page 65). Cleavage separates the reporter dye from the quencher dye; the separation of the reporter dye from the quencher dye results in increased fluorescence by the reporter. The increase in fluorescence occurs only if the target sequence is complementary to the probe and is amplified during PCR. Because of these requirements, nonspecific amplification is not detected.

Figure 4 Cleavage



Polymerization of the strand continues, but because the 3' end of the probe is blocked, no extension of the probe occurs during PCR (Figure 5).

Figure 5 Completion of polymerization



Appendix C About Two-Step RT-PCR
About the 5' nuclease assay

Safety

This appendix covers:

■ Chemical safety	68
General chemical safety	68
SDSs	68
Chemical waste safety	69
Biological hazard safety	70



Chemical safety

General chemical safety

Chemical hazard
warning



WARNING! CHEMICAL HAZARD. Before handling any chemicals, refer to the Safety Data Sheet (SDS) provided by the manufacturer, and observe all relevant precautions.



WARNING! CHEMICAL STORAGE HAZARD. Never collect or store waste in a glass container because of the risk of breaking or shattering. Reagent and waste bottles can crack and leak. Each waste bottle should be secured in a low-density polyethylene safety container with the cover fastened and the handles locked in the upright position. Wear appropriate eyewear, clothing, and gloves when handling reagent and waste bottles.

Chemical safety
guidelines

To minimize the hazards of chemicals:

- Read and understand the Safety Data Sheets (SDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials. (See “About SDSs” on page 68.)
- Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing). For additional safety guidelines, consult the SDS.
- Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (for example, fume hood). For additional safety guidelines, consult the SDS.
- Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer’s cleanup procedures as recommended in the SDS.
- Comply with all local, state/provincial, or national laws and regulations related to chemical storage, handling, and disposal.

SDSs

About SDSs

Chemical manufacturers supply current Safety Data Sheets (SDSs) with shipments of hazardous chemicals to new customers. They also provide SDSs with the first shipment of a hazardous chemical to a customer after an SDS has been updated. SDSs provide the safety information you need to store, handle, transport, and dispose of the chemicals safely.

Each time you receive a new SDS packaged with a hazardous chemical, be sure to replace the appropriate SDS in your files.

Obtaining
SDSs

The SDS for any chemical supplied by Applied Biosystems is available to you free 24 hours a day. To obtain SDSs:

1. Go to www.appliedbiosystems.com, click **Support**, then select **SDS**.
2. In the Keyword Search field, enter the chemical name, product name, SDS part number, or other information that appears in the SDS of interest. Select the language of your choice, then click **Search**.

- Find the document of interest, right-click the document title, then select any of the following:
 - Open** – To view the document
 - Print Target** – To print the document
 - Save Target As** – To download a PDF version of the document to a destination that you choose



Note: For the SDSs of chemicals not distributed by Applied Biosystems, contact the chemical manufacturer.

Chemical waste safety

Chemical waste hazards



CAUTION! HAZARDOUS WASTE. Refer to Safety Data Sheets and local regulations for handling and disposal.



WARNING! CHEMICAL STORAGE HAZARD. Never collect or store waste in a glass container because of the risk of breaking or shattering. Reagent and waste bottles can crack and leak. Each waste bottle should be secured in a low-density polyethylene safety container with the cover fastened and the handles locked in the upright position. Wear appropriate eyewear, clothing, and gloves when handling reagent and waste bottles.

Chemical waste safety guidelines

To minimize the hazards of chemical waste:


- Read and understand the Safety Data Sheets (SDSs) provided by the manufacturers of the chemicals in the waste container before you store, handle, or dispose of chemical waste.
- Provide primary and secondary waste containers. (A primary waste container holds the immediate waste. A secondary container contains spills or leaks from the primary container. Both containers must be compatible with the waste material and meet federal, state, and local requirements for container storage.)
- Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing). For additional safety guidelines, consult the SDS.
- Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (for example, fume hood). For additional safety guidelines, consult the SDS.
- Handle chemical wastes in a fume hood.
- After emptying a waste container, seal it with the cap provided.
- Dispose of the contents of the waste tray and waste bottle in accordance with good laboratory practices and local, state/provincial, or national environmental and health regulations.



Waste disposal

If potentially hazardous waste is generated, you must:

- Characterize (by analysis if necessary) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- Ensure the health and safety of all personnel in your laboratory.
- Ensure that the waste is stored, transferred, transported, and disposed of according to all local, state/provincial, and/or national regulations.

 **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.

Biological hazard safety

General biohazard



WARNING! BIOHAZARD. Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. Follow all applicable local, state/provincial, and/or national regulations. Wear appropriate protective equipment, which includes but is not limited to: protective eyewear, face shield, clothing/lab coat, and gloves. All work should be conducted in properly equipped facilities using the appropriate safety equipment (for example, physical containment devices). Individuals should be trained according to applicable regulatory and company/institution requirements before working with potentially infectious materials. Read and follow the applicable guidelines and/or regulatory requirements in the following:

- U.S. Department of Health and Human Services guidelines published in *Biosafety in Microbiological and Biomedical Laboratories* (www.cdc.gov/od/ohs/biosfty/bmb15/bmb15toc.htm).
- Occupational Safety and Health Standards, Bloodborne Pathogens (29 CFR§1910.1030; www.access.gpo.gov/nara/cfr/waisidx_01/29cfr1910a_01.html).
- Your company's/institution's Biosafety Program protocols for working with/handling potentially infectious materials.

Additional information about biohazard guidelines is available at:
www.cdc.gov

Bibliography

- Förster, V.T. 1948. Zwischenmolekulare Energiewanderung und Fluoreszenz. *Annals of Physics (Leipzig)*. 2:55–75.
- Kwok, S. and Higuchi, R. 1989. Avoiding false positives with PCR. *Nature*. 339:237–238.
- Lakowicz, J.R. 1983. *Principles of Fluorescence Spectroscopy*, ed. New York: Plenum Press. xiv, 496 pp.
- Longo, M.C., Berninger, M.S., and Hartley, J.L. 1990. Use of uracil DNA glycosylase to control carry-over contamination in polymerase chain reactions. *Gene*. 93:125–128.
- Mullis, K.B. and Faloona, F.A. 1987. Specific synthesis of DNA in vitro via a polymerase-catalyzed chain reaction. *Methods Enzymol*. 155:335–350.
- Saiki, R.K., Scharf, S., Faloona, F., et al. 1985. Enzymatic amplification of β -globin genomic sequences and restriction site analysis for diagnosis of sickle cell anemia. *Science*. 230:1350–1354.

ΔR_n	<p>The difference between the R_{n+} value and the R_{n-} value. It reliably indicates the magnitude of the signal generated by the given set of PCR conditions.</p> <p>The following equation expresses the relationship of these terms:</p> $\Delta R_n = (R_{n+}) - (R_{n-})$ <p>where:</p> $R_{n+} = \frac{\text{Emission Intensity of Reporter}}{\text{Emission Intensity of Passive Reference}} \quad \text{PCR with template}$ $R_{n-} = \frac{\text{Emission Intensity of Reporter}}{\text{Emission Intensity of Passive Reference}} \quad \text{PCR without template or early cycles of a real-time reaction}$
C_T	See threshold cycle (C_T) .
fold difference	The measured ratio of the quantity of template in Sample A over the quantity of template in Sample B, where quantity A > quantity B, so that the ratio is > 1.
full replicate	Repeated wells of the same sample with the same assay, where the contents of each well go through all experimental steps (sample preparation, reverse transcription, and PCR) separately.
minimum fold difference	The smallest fold difference that is statistically significant.
multicomponenting	The term used to distinguish the contribution each individual dye makes to the fluorescent spectra. The overlapping spectra from the dye components generate the composite spectrum, which represents one reading from one well.
normalization	The Passive Reference 1, a dye included in the 10X TaqMan [®] Buffer A, does not participate in the 5' nuclease PCR. The Passive Reference provides an internal reference to which the reporter-dye signal can be normalized during data analysis. Normalization is necessary to correct for fluorescent fluctuations due to changes in concentration or volume.
PCR/technical replicate	Identical reactions that contain identical components and volumes and evaluate the same sample.

R_n^+	<p>The R_n value of a reaction containing all components including the template.</p> <p>Normalization is accomplished by dividing the emission intensity of the reporter dye by the emission intensity of the Passive Reference to obtain a ratio defined as the R_n (normalized reporter) for a given reaction tube.</p>
R_n^-	<p>The R_n value of an unreacted sample. This value may be obtained from the early cycles of a real-time run (those cycles prior to a detectable increase in fluorescence). This value may also be obtained from a reaction not containing template.</p>
statistically significant (to a 99.7% confidence level)	<p>A result with a low probability (0.3%) of resulting from chance.</p>
threshold cycle (C_T)	<p>The PCR cycle number at which the fluorescence meets the threshold in the amplification plot.</p>

Documentation and Support

Product documentation

A portable document format (PDF) version of the document listed below is available at www.appliedbiosystems.com



Note: To open the PDF version, use the Adobe Acrobat Reader software available from www.adobe.com

Document	Part number
<i>TaqMan[®] Fast Advanced Master Mix Protocol</i>	4444605

Related documentation

To obtain the documents listed in this section or additional documentation, see “Obtaining support” on page 77.

Chemistry/reagent documentation

Document	Part number
<i>Applied Biosystems Real-Time PCR Systems Reagent Guide</i>	4387787
<i>Bioinformatic Evaluation of a Sequence for Custom TaqMan[®] Gene Expression Assays Tutorial</i>	–
<i>Custom TaqMan[®] Assays: Design and Ordering Guide</i>	4367671
<i>Custom TaqMan[®] Gene Expression Assays Protocol</i>	4334429
<i>Primer Express[®] Software Version 3.0 Getting Started Guide</i>	4362460
<i>Real-Time PCR Systems Chemistry Guide: Applied Biosystems 7900HT Fast Real-Time PCR Systems and 7300/7500/7500 Fast Real-Time PCR Systems</i>	4348358
<i>TaqMan[®] Array Micro Fluidic Card User Guide</i>	4400263
<i>TaqMan[®] Array Plates Protocol</i>	4391016
<i>TaqMan[®] Gene Expression Assays Protocol</i>	4333458
<i>TaqMan[®] MicroRNA Reverse Transcription Kit Protocol</i>	4367038
<i>TaqMan[®] Small RNA Assays Protocol</i> (For TaqMan [®] MicroRNA Assays and Custom TaqMan [®] Small RNA Assays)	4364031

Instrument documentation

7300/7500/7500 Fast system

Document	Part number
<i>Applied Biosystems 7300/7500/7500 Fast Real-Time PCR System Relative Quantification Getting Started Guide</i>	4347824
<i>Applied Biosystems 7300/7500/7500 Fast Real-Time PCR Systems Absolute Quantification Getting Started Guide</i>	4347825
<i>Applied Biosystems 7500/7500 Fast Real-Time PCR Systems Getting Started Guide for Standard Curve Experiments</i>	4387779
<i>Applied Biosystems 7500/7500 Fast Real-Time PCR Systems Getting Started Guide for Comparative C_T/Relative Standard Curve Experiments</i>	4387783

7900HT/7900HT Fast system (Fast 96-Well, Standard 96-Well, or 384-Well Block Modules)

Document	Part number
<i>Applied Biosystems 7900HT Fast Real-Time PCR System Quick Reference Card: Performing Fast Gene Quantification</i>	4351892
<i>Applied Biosystems 7900HT Fast Real-Time PCR System Relative Quantitation Using Comparative C_T Getting Started Guide</i>	4364016
<i>Applied Biosystems 7900HT Fast Real-Time PCR System Absolute Quantitation Using Standard Curve Getting Started Guide</i>	4364014
<i>Applied Biosystems 7900HT Fast Real-Time PCR System User Bulletin: Performing Fast Gene Quantitation with 394-Well Plates</i>	4369584
<i>Applied Biosystems 7900HT Fast Real-Time PCR System and SDS Enterprise Database User Guide</i>	4351684

StepOne™ and StepOnePlus™ systems

Document	Part number
<i>Applied Biosystems StepOne™ and StepOnePlus™ Real-Time PCR Systems Reagent Guide</i>	4379704
<i>Applied Biosystems StepOne™ and StepOnePlus™ Real-Time PCR Systems Relative Standard Curve and Comparative C_T Experiments Getting Started Guide</i>	4376785
<i>Applied Biosystems StepOne™ and StepOnePlus™ Real-Time PCR Systems Getting Started Guide for Standard Curve Experiments</i>	4376784

Obtaining support

For the latest services and support information for all locations, go to:

www.appliedbiosystems.com

At the Applied Biosystems web site, you can:

- Access worldwide telephone and fax numbers to contact Applied Biosystems Technical Support and Sales facilities.
- Search through frequently asked questions (FAQs).
- Submit a question directly to Technical Support.
- Order Applied Biosystems user documents, SDSs, certificates of analysis, and other related documents.
- Download PDF documents.
- Obtain information about customer training.
- Download software updates and patches.

Part Number 4444605 Rev. B 06/2010



Headquarters
5791 Van Allen Way | Carlsbad, CA 92008 USA
Phone 760.603.7200
www.lifetechnologies.com

Technical Resources and Support
For the latest technical resources and support information
for all locations, please refer to our Web site at
www.appliedbiosystems.com