

CD CHO Medium

Description

CD CHO Medium has been developed for the growth of Chinese Hamster Ovary (CHO) cells and expression of recombinant proteins in suspension culture. CD CHO is an animal origin-free (AOF), chemically defined medium that contains no proteins, hydrolysates, or components of undefined composition. CD CHO is formulated without hypoxanthine and thymidine for use in dihydrofolate reductase (DHFR) amplified systems, and without phenol red to minimize estrogen-like effects of phenol red.

Product	Catalog no.	Amount	Storage	Shelf life*
CD CHO Medium (1X), liquid	10743-011	500 mL	2°C to 8°C; Protect from light	18 months
	10743-029	1000 mL		
	10743-001	10 L (Bag)	2°C to 8°C; Protect from light	12 months
	10743-002	20 L (Bag)		
CD CHO AGT™**	12490-017	1 L	2°C to 8°C; Store dark and dry	24 months
	12490-025	1 × 10 L		
	12490-001	1 × 100 L		
	12490-003	10 kg		

* Shelf Life duration is determined from Date of Manufacture. ** AGT = Advanced Granulation Technology.

Product use

Caution: For manufacturing, processing, or repacking.

Safety information

Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Prepare media

CD CHO Medium and CD CHO AGT™ medium require supplementation with L-glutamine or GlutaMAX™-I prior to use.

1. Aseptically add L-glutamine or GlutaMAX™-I 8 mM final concentration (40 mL/L) to the medium before use.
2. If L-glutamine is not required, add 40 mL of sterile distilled water and adjust the osmolality to 320 mOsm using a sterile solution of NaCl.
3. CD CHO Medium is made without hypoxanthine and thymidine for use in dihydrofolate reductase (DHFR) amplified systems. For other applications, add 10 mL/L of HT Supplement prior to use.
4. Add 1 mL/L of Anti-Clumping Agent to media if cell clumping occurs. After any medium changes, passage cells for a minimum of 3X before use in other applications.

Note: Consider using lower levels of L-glutamine if you are using a fed batch or perfusion protocol or if the cell line in use is sensitive to ammonia.

Note: Addition of a surfactant (e.g., Pluronic® F-68) is not required.

Reconstitute CD CHO AGT™ Medium

1. Weigh out 24.3 g (equivalent to the entire contents of a 1-L package) of CD CHO AGT™ medium.
2. Add to 900 mL room temperature deionized or distilled water. Rinse inside of package to remove all traces of powder. Mix gently for 30 minutes or until medium dissolves completely.
3. Add deionized or distilled water to final volume of 1 L.
4. Filter sterilize by 0.2 µm pore size membrane filtration.
5. Supplement as described in **Prepare Media** at time of use.

Note: CD CHO AGT™ medium contains sodium bicarbonate.

Do not add additional sodium bicarbonate. CD CHO AGT™ medium is auto pH and osmolality adjusted, no further adjustment is required. For final lot pH and osmolality specifications please refer to Certificate of Analysis specification.

Culture conditions

Media: CD CHO Medium

Cell line: Chinese Hamster Ovary (CHO)

Culture type: Suspension

Culture vessels: shake flasks, spinner bottles, or bioreactor.

Procedures described are intended for use with 125-mL Erlenmeyer shake flasks.

Temperature range: 36°C to 38°C

Incubator atmosphere: Humidified atmosphere of 8–10% CO₂ in air. Ensure proper gas exchange and minimize exposure of cultures to light.

Recovery

1. Rapidly thaw (<1 minute) frozen cells in a 37°C water bath.
2. Transfer the entire contents of the cryovial into a 125-mL shake flask containing 28.5 mL of pre-warmed complete CD CHO Medium.
3. Incubate at 37°C in a humidified atmosphere of 8–10% CO₂ in air on an orbital shaker platform rotating at 125–135 rpm.
4. Subculture cells in mid-logarithmic phase 3–5 days post-thaw at a seeding density of 3 × 10⁵ viable cells/mL. Passage cells a minimum of 3X before use in other applications.

Note: Do not centrifuge CHO cells to remove DMSO as they are extremely fragile upon recovery from cryopreservation.

Subculture suspension cultures

1. Determine viable cell density using a Countess® Automated Cell Counter.
2. Seed cells at 2 × 10⁵–3 × 10⁵ viable cells/mL in sterile culture vessels containing pre-warmed complete CD CHO Medium. (30 mL per 125-mL shake flask).
3. Incubate at 37°C in a humidified atmosphere of 8–10% CO₂ in air on an orbital shaker platform rotating at 125–135 rpm. Loosen flask cap to allow for gas exchange.
4. Subculture cells when viable cell density reaches ≥1 × 10⁶ viable cells/mL into clean, sterile flask(s) with fresh pre-warmed complete CD CHO Medium.

Note: To reduce accumulation of cell debris and metabolic waste by-products in suspension cultures, gently centrifuge the cell suspension at 100 × g for 5–10 minutes and resuspend pellet in fresh complete CD CHO Medium once every 2–3 weeks.

Note: It is recommended to thaw a fresh low-passage vial of cells every 3 months or 30 passages.

Adapt CHO Cells to CD CHO Medium

It is critical that cells be in mid-logarithmic phase growth and exceed 90% viability prior to initiating adaptation procedures from conventional serum-supplemented or serum-free medium.

Direct adaptation

Transfer suspension cultures into CD CHO Medium as follows:

1. Centrifuge the cell suspension at $100 \times g$ for 5–10 minutes. Aspirate and discard the supernatant.
2. Resuspend the cell pellet in pre-warmed complete CD CHO Medium at a viable cell density of 3×10^5 – 5×10^5 cells/mL and transfer to appropriate culture vessels.
3. Return to incubator and monitor cell growth.

Note: If suboptimal cell growth is observed using the direct adaptation method, use the sequential adaptation method.

Sequential adaptation

1. Follow the procedures for **Subculture Suspension Cultures** with the following modifications.
2. During the adaptation procedure use a seeding density of 4×10^7 – 5×10^5 viable cells/mL.
3. Subculture cells into stepwise increasing ratios of complete CD CHO Medium to original medium with each subsequent passage (25:75, 50:50, 75:25, 90:10 followed by 100% CD CHO Medium). Multiple passages at each step may be needed.
4. After several passages in 100% CD CHO Medium, the viable cell count should exceed 1×10^6 – 2×10^6 cells/mL with a viability $\geq 85\%$ within 4–6 days of culture. At this stage the culture is considered to be adapted to CD CHO Medium. The seeding density may be reduced to 2×10^5 – 3×10^5 viable cells/mL during the final stages of adaptation.

Cryopreservation

Prepare the desired quantity of cells, harvesting in mid-log phase of growth with viability $>90\%$.

1. Obtain Synth-a-Freeze[®] cryopreservation medium and store at 2°C to 8°C until use.
2. Determine the viable cell density and calculate the required volume of Synth-a-Freeze[®] cryopreservation medium. Typical cell densities for cryopreservation with Synth-a-Freeze[®] medium are 0.5×10^7 – 1×10^7 viable cells/mL.
3. Harvest cells by centrifugation at $100 \times g$ for 5–10 minutes. Resuspend cell pellet in the pre-determined volume of 2°C to 8°C of Synth-a-Freeze[®] medium.
4. Immediately dispense aliquots of this suspension into cryovials according to the manufacturer's specifications
5. Achieve cryopreservation in an automated or manual controlled rate freezing apparatus following standard procedures (1°C decrease per minute).
6. Transfer frozen cells to liquid nitrogen; (vapor phase) storage at –200°C to –125°C is recommended.










Note: Check viability of cryopreserved cells 24 hours after storage of vials in liquid nitrogen. See **Recovery**.

Related products

Product	Catalog No.
L–Glutamine, 200mM (100X), liquid	25030
GlutaMAX™-I, 200mM (100X), liquid	35050
HT Supplement, (100X), liquid	11067
Anti-Clumping Agent	0010057
CD DG44 Medium, (1X) liquid	12610
CHO CD EfficientFeed™ Kit	A10241
CHO CD EfficientFeed™ A AGT™ Nutrient Supplement	A14420
CHO CD EfficientFeed™ B AGT™ Nutrient Supplement	A12456
CD EfficientFeed™ C AGT™ Nutrient Supplement	A13275
Water, Distilled	15230
CHO-S™ Cells (cGMP banked) and Media Kit	A11557
Countess [®] Automated Cell Counter	C10227

Explanation of symbols and warnings

The symbols present on the product label are explained below:

				
Catalog number	Manufacturer	Use by	Keep away from light	Batch code
				
Caution, consult accompanying documents	Consult instructions for use	Temperature limitation	Sterilized using aseptic processing techniques	

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

Limited Use Label License: Internal Research and Bioproduction Use

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product (a) to perform internal research for the sole benefit of the purchaser; and (b) to culture cells for the purpose of producing a product wherein the product will be used for any or all of the following: (i) internal research use by the purchaser; (ii) resale for internal research use by third parties; (iii) performance of research conducted by the purchaser on a fee for service or contract basis for or on behalf of third parties; (iv) resale for use as a human therapeutic agent or diagnostics product or component by third parties; (v) performance of manufacturing services conducted by the purchaser on a fee for service or contract basis for or on behalf of third parties.

The purchase of this product does not grant the purchaser any additional rights, including (without limitation) the right to transfer or resell the product in any form, the right to use the product as a therapeutic agent or diagnostics test component, or to use the product to perform other tests on a contract or fee per test basis for or on behalf of third parties. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

For additional technical information such as Safety Data Sheets (SDS), Certificates of Analysis, visit www.lifetechnologies.com/support.

For further assistance, email techsupport@lifetech.com.

© 2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Pluronic is a trademark of BASF Corporation.

DISCLAIMER: LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

www.lifetechnologies.com

 life
technologies