

USER GUIDE

**invitrogen™**  
by *life* technologies™

# pOptiVEC™-TOPO® TA Cloning Kit

For TOPO® Cloning of PCR products into a bicistronic vector

**Catalog Number** 12744-017

**Revision Date** 18 May 2011

**Part Number** 25-0977

**Publication Number** MAN0000589

**For Research Use Only. Not intended for any animal or human  
therapeutic or diagnostic use.**

*life*  
technologies™

# Table of Contents

Kit Contents and Storage.....	4
Overview.....	6
<b>Methods.....</b>	<b>8</b>
Designing PCR Primers.....	8
Producing PCR Products.....	10
Setting Up the TOPO® Cloning Reaction.....	12
Transforming One Shot® Competent <i>E. coli</i> .....	14
Analyzing Positive Clones.....	16
Next Steps.....	17
Troubleshooting.....	18
<b>Appendix.....</b>	<b>20</b>
Performing the Control Reactions.....	20
Addition of 3' A-Overhangs Post-Amplification.....	22
Map and Features of pOptiVEC™ -TOPO®.....	23
Recipes.....	25
Accessory Products.....	26
Technical Support.....	27
Purchaser Notification.....	28
Notes.....	30

## Kit Contents and Storage

Types of kits The pOptiVEC™-TOPO® TA Vector Kit is included with the Freedom™ DG44™ Kit, and as a separate catalog item:

Kit	Cat. Number
Freedom™ DG44™ Kit	A13737-01
pOptiVEC™-TOPO® TA Vector Kit	12744-017

Shipping and storage The pOptiVEC™-TOPO® TA Vector Kit is shipped on dry ice. Each kit contains two boxes as described below. Upon receipt, store boxes as detailed below.

Box	Item	Storage
1	pOptiVEC™-TOPO® TA Cloning Reagents	-20°C
2	One Shot® TOP10 Chemically Competent <i>E. coli</i>	-80°C

TOPO® TA cloning reagents The pOptiVEC™-TOPO® TA cloning reagents (Box 1) are listed below. **Note that the user must supply *Taq* polymerase. Store the contents of Box 1 at -20°C.**

Item	Concentration	Amount
pOptiVEC™ TOPO® vector, TOPO® adapted	5–10 ng linearized plasmid DNA in: 50% glycerol 50 mM Tris-HCl, pH 7.4 1 mM EDTA 1mM DTT 0.1% Triton® X-100 100 µg/mL BSA 30 µm phenol red	10 µL
10X PCR Buffer	100 mM Tris-HCl, pH 8.3 (at 42°C) 500 mM KCl 25 mM MgCl <sub>2</sub> 0.01% gelatin	100 µL
dNTP Mix	12.5 mM dATP, 12.5 mM dCTP, 12.5 mM dGTP, and 12.5 mM dTTP neutralized at pH 8.0 in water	10 µL
Salt Solution	1.2 M NaCl 0.06 M MgCl <sub>2</sub>	50 µL
Sterile Water	–	1 mL
Control PCR template	50 ng/µL in TE buffer, pH 8.0	10 µL
Control PCR primers	100 ng/µL each in TE buffer, pH 8.0	10 µL
CMV forward sequencing primer	100 ng/µL in TE buffer, pH 8.0	20 µL
EMCV IRES reverse sequencing primer	100 ng/µL in TE buffer, pH 8.0	20 µL

*Continued on next page*

## Kit Contents and Storage, continued

**Primers** The pOptiVEC™-TOPO® TA Vector Kit contains the following primers to sequence your insert.

Primer	Sequence
CMV forward	5'-CGCAAATGGGCGGTAGGCGTG-3'
EMCV IRES reverse	5'- CCTTATTCCAAGCGGCTTCG-3'

**One Shot® TOP10 reagents** The following reagents are included in the One Shot® TOP10 Chemically Competent *E. coli* kit (Box 2). Transformation efficiency is  $\geq 1 \times 10^9$  cfu/ $\mu$ g plasmid DNA. **Store the contents of Box 2 at -80°C.**

Item	Concentration	Amount
TOP10 <i>E. coli</i>	–	11 $\times$ 50 $\mu$ L
pUC19 Control DNA	10 pg/ $\mu$ L in 5 mM Tris-HCl, 0.5 mM EDTA, pH 8	50 $\mu$ L
S.O.C. Medium	2% Tryptone 0.5% Yeast extract 10 mM NaCl 2.5 mM KCl 10 mM MgCl <sub>2</sub> 10 mM MgSO <sub>4</sub> 20 mM glucose	6 mL

**Genotype of TOP10 strain** F<sup>-</sup> *mcrA*  $\Delta$ (*mrr-hsdRMS-mcrBC*)  $\Phi$ 80*lacZ* $\Delta$ M15  $\Delta$ *lac* $\chi$ 74 *recA1* *araD139*  $\Delta$ (*ara-leu*)7697 *galU galK rpsL* (Str<sup>R</sup>) *endA1 nupG*

# Overview

---

## Description

The pOptiVEC™-TOPO® vector is a TOPO®-adapted bicistronic plasmid that allows rapid cloning of a PCR product containing a mammalian secretion signal and the gene of interest downstream of the CMV promoter. In the pOptiVEC™-TOPO® vector, the transcription of the gene of interest is separated from the dihydrofolate reductase (DHFR) auxotrophic selection marker by an internal ribosome entry site (IRES), allowing transcription of the gene of interest and the selection marker on the same mRNA.

The pOptiVEC™-TOPO® vector is also used with the Freedom™ DG44™ Kit.

---

## Features of the vector

The pOptiVEC™-TOPO® vector contains the following elements:

- Human cytomegalovirus (CMV) immediate-early promoter/enhancer for high-level gene expression in a wide range of mammalian cells
- TOPO® Cloning site for rapid and efficient cloning of *Taq*-amplified PCR products
- Internal Ribosome Entry Site (IRES) from the encephalomyocarditis virus (EMCV) for cap-independent translation of DHFR
- Dihydrofolate reductase (DHFR) gene for auxotrophic selection of transfected DG44 cells and for genomic amplification of stable cell lines using methotrexate (MTX)
- The Herpes Simplex Virus thymidine kinase polyadenylation signal for proper termination and processing of the recombinant transcript
- pUC origin for high copy replication and maintenance of the plasmid in *E. coli*
- Ampicillin (*bla*) resistance gene for selection in *E. coli*

For a map and features of the pOptiVEC™-TOPO® vector, see pages 23–24.

---

## CMV promoter

The human cytomegalovirus immediate-early (HCMV IE1) gene promoter in the pOptiVEC™-TOPO® vector is 680 bp and contains the native transcriptional start site (Hennighausen & Fleckenstein, 1986). This sequence results in high levels of transgene expression.

---

## EMCV IRES

The internal ribosome entry site (IRES) from encephalomyocarditis virus (EMCV) allows cap-independent translation initiation of the DHFR selection gene. The EMCV IRES allows expression of the gene of interest and the selection marker from a single bicistronic mRNA (Gurtu *et al.*, 1996; Rees *et al.*, 1996).

---

## DHFR

Dihydrofolate reductase (DHFR) catalyzes the reduction of 5, 6-dihydrofolate to 5, 6, 7, 8-tetrahydrofolate, which is essential for DNA synthesis. CHO-derived DG44 cells lack DHFR activity and must be propagated in medium containing the purine precursors hypoxanthine and thymidine (HT) unless the cells are stably transfected with a vector that expresses DHFR.

DHFR can also function as a genomic amplification marker for your gene of interest using methotrexate (MTX) selection (Kaufman *et al.*, 1985). See the Freedom™ DG44™ Kit manual for more details on this procedure.

---

*Continued on next page*

## Overview, continued

---


### How TOPO<sup>®</sup> works

The pOptiVEC<sup>™</sup>-TOPO<sup>®</sup> vector is supplied linearized with:

- Single 3' thymidine (T) overhangs for TA Cloning<sup>®</sup>
- Topoisomerase covalently bound to the vector (this is referred to as “activated” vector)

*Taq* polymerase has a nontemplate-dependent terminal transferase activity that adds a single deoxyadenosine (A) to the 3' ends of PCR products. The linearized vector supplied in this kit has single, overhanging 3' deoxythymidine (T) residues. This allows PCR inserts to ligate efficiently with the vector.

Topoisomerase I from Vaccinia virus binds to duplex DNA at specific sites and cleaves the phosphodiester backbone after 5'-CCCTT in one strand (Shuman, 1991). The energy from the broken phosphodiester backbone is conserved by formation of a covalent bond between the 3' phosphate of the cleaved strand and a tyrosyl residue (Tyr-274) of topoisomerase I. The phospho-tyrosyl bond between the DNA and enzyme can subsequently be attacked by the 5' hydroxyl of the original cleaved strand, reversing the reaction and releasing topoisomerase (Shuman, 1994). TOPO<sup>®</sup> Cloning exploits this reaction to efficiently clone PCR products.



Once the PCR product is cloned into the pOptiVEC<sup>™</sup>-TOPO<sup>®</sup> vector and the transformants are analyzed for correct orientation and reading frame, the expression plasmid may be transfected into the CHO derived DG44 cells using the Freedom<sup>™</sup> DG44<sup>™</sup> Kit.

---

### Experimental outline

To TOPO<sup>®</sup> Clone your gene of interest into pOptiVEC<sup>™</sup>-TOPO<sup>®</sup>, you will perform the following steps:

1. Generate a PCR product containing a mammalian secretion signal and your gene of interest with *Taq* polymerase.
2. TOPO<sup>®</sup> Clone your PCR product into the pOptiVEC<sup>™</sup>-TOPO<sup>®</sup> vector and use the reaction to transform One Shot<sup>®</sup> TOP10 Chemically Competent *E. coli*
3. Pick colonies, isolate plasmid DNA, and screen for insert directionality by sequencing expression clones with the primers provided in the kit.

The following sections of this manual provide instructions and guidelines for these steps.

---

# Methods

## Designing PCR Primers

---

### Introduction

TOPO<sup>®</sup> Cloning provides a highly efficient, 5-minute, one-step cloning strategy (“TOPO<sup>®</sup> Cloning”) for the direct insertion of *Taq* polymerase-amplified PCR products into a plasmid vector. No ligase, post-PCR procedures, or PCR primers containing specific sequences are required. The section below will help you design primers to produce your PCR product for cloning into the pOptiVEC<sup>™</sup>-TOPO<sup>®</sup> vector.

---

### Points to consider when designing PCR primers

To obtain the pOptiVEC<sup>™</sup>-TOPO expression construct containing your gene of interest, your PCR primer design must include:

- A Kozak consensus sequence
- A mammalian secretion signal upstream of your gene of interest (if you wish to produce secreted protein)
- A stop codon at the end of your gene of interest

See below for more information.

---

### Kozak consensus sequence

Your gene of interest must contain a Kozak translation initiation sequence with an ATG initiation codon for proper initiation of translation (Kozak, 1987; Kozak, 1990; Kozak, 1991). An example of a Kozak consensus sequence is provided below. The ATG initiation codon is shown underlined.

(G/A)NNATGN

Other sequences are possible, but the G or A at position –3 (shown in bold) is critical for a functional Kozak sequence. At position +4 any of the four nucleotides can be present to form part of the Kozak sequence.

---

### Secretion signal

If you wish for your protein of interest to be secreted from mammalian DG44 cells, your gene of interest must include a mammalian secretion signal. To direct secretion of your protein of interest, you can include the endogenous secretion signal of your protein of interest, or add one such as the murine Ig κ-chain leader sequence (Coloma *et al.*, 1992) using PCR.

---



### Note

Do not add 5' phosphates to your primers for PCR, because the synthesized PCR product will not ligate into the vector.

Cloning efficiencies may vary depending on the primer nucleotide sequences.

---

*Continued on next page*

## Designing PCR Primers, continued

---

### TOPO<sup>®</sup> Cloning Site

Use the diagram below to help you design your PCR product for TOPO<sup>®</sup> Cloning into pOptiVEC<sup>™</sup>-TOPO<sup>®</sup>. The complete vector sequence is available from [www.invitrogen.com](http://www.invitrogen.com) or by contacting Technical Support (page 27).

```

      CAAT   CMV Forward Primer Binding Site   TATA
      |-----|-----|-----|-----|
571  GACGCAATG GCGGGTAGGC GTGTACGGTG GGAGGCTAT ATAAGCAGAG CTCGTTTAGT GAACCGTCAG
      CTGCGTTTAC CCGCCATCCG CACATGCCAC CCTCCAGATA TATTCGTCTC GAGCAAAATCA CTTGGCAGTC

641  ATCGCCTGGA GACGCCATCC ACGCTGTTTT GACCTCCATA GAAGACACCG GGACCGATCC AGCCTCCGGA
      TAGCGGACCT CTGCGGTAGG TGCGACAAAA CTGGAGGTAT CTTCGTGTGC CCTGGCTAGG TCGGAGGCCT

711  CTCTAGAGGA TCCAACCCTT PCR PRODUCT AAGGGTTGGA TCCCTACCGG TGCTGCGGCC
      GAGATCTCCT AGGTTGGGAA TTCCAACCT AGGGATGGCC ACGACCCCGG
                                     EMCV IRES Reverse Primer Binding Site
761  GCGCAGTTAA CGCGCCCCT CTCCTCCCC CCCCCCTAAC GTTACTGGCC GAAGCCGCTT GGAATAAGGC
      CGCGTCAATT GCGGCGGGGA GAGGGAGGGG GGGGGATTG CAATGACCGG CTCGGCGGAA CCTTATTCCG
```

---



# Producing PCR Products

---

## Introduction

After you have designed primers to amplify your DNA of interest, you are ready to produce your PCR product(s) for TOPO<sup>®</sup> Cloning into the pOptiVEC<sup>™</sup>-TOPO<sup>®</sup> vector.

---

## Materials supplied by user

You will need the following reagents and equipment:

- *Taq* polymerase, such as Platinum<sup>®</sup> *Taq* (see page 26)
  - Thermocycler
  - DNA template
  - Primers for PCR product
- 

## Polymerase mixtures

You may use an enzyme mixture containing *Taq* polymerase and a proofreading polymerase; however *Taq* must be used in excess of 10:1 to ensure the presence of 3' A-overhangs on the PCR product. We recommend using Platinum<sup>®</sup> *Taq* DNA Polymerase High Fidelity available from Invitrogen (see page 26 for ordering information).

If you use polymerase mixtures that do not have enough *Taq* polymerase or use a proofreading polymerase only, you can add 3' A-overhangs after amplification using the method on page 22.

---

## Producing PCR products

1. Set up the following reaction in a 50  $\mu$ L volume. Use less DNA if you are using plasmid DNA as a template and more DNA if you are using genomic DNA as a template.


Reagent	Amount
DNA template	10–100 ng
10X PCR Buffer	5 $\mu$ L
50 mM dNTPs	0.5 $\mu$ L
PCR Primers	100–200 ng each
Sterile Water	to final volume of 49 $\mu$ L
<i>Taq</i> polymerase (1 unit/ $\mu$ L)	1 $\mu$ L
<b>Total Volume</b>	<b>50 <math>\mu</math>L</b>

2. Perform amplification using the cycling parameters suitable for your primers and template. Be sure to include a 7–30 minute extension at 72° after the last cycle to ensure that all PCR products are full-length and 3' adenylated.
  3. Use agarose gel electrophoresis to verify the quality of your PCR product. You should see a single, discrete band of the correct size. If you do not see a single band, refer to the **Note** on the next page.
- 

*Continued on the next page*

## Producing PCR Products, continued

---



### Note

If you do not obtain a single, discrete band from your PCR reaction, try the following:

- The PCR Optimizer™ Kit (Catalog no. K1220-01) can help you optimize your PCR to eliminate multiple bands and smearing.
- Gel-purify your fragment before performing the TOPO® Cloning reaction using the E-Gel® CloneWell system or PureLink® Gel Extraction Kit, available separately from Invitrogen. See page 26 for more information.

Alternatively, refer to Current Protocols in Molecular Biology, Unit 2.6 (Ausubel *et al.*, 1994) for other common protocols for isolating DNA fragments.

---

# Setting Up the TOPO<sup>®</sup> Cloning Reaction

---

## Introduction

Once you have produced the desired PCR product(s), you are ready to TOPO<sup>®</sup> Clone into the pOptiVEC<sup>™</sup>-TOPO<sup>®</sup> vector and use this plasmid or plasmids for transformation of competent *E. coli*. It is important to have everything you need to set up the reaction so that you can obtain the best results. We suggest that you read this entire section and the next section about transformation before beginning. If this is the first time you have TOPO<sup>®</sup> Cloned, perform the control reactions detailed on pages 20–21 in parallel with your samples.

---



We have found that including salt (200 mM NaCl, 10 mM MgCl<sub>2</sub>) in the TOPO<sup>®</sup> Cloning reaction increases the number of transformants 2- to 3-fold. In addition, incubating the reaction mixture for greater than 5 minutes in the presence of salt can also increase the number of transformants. This is in contrast to experiments **without salt** where the number of transformants decreases as the incubation time increases beyond 5 minutes.

Including salt in the TOPO<sup>®</sup> Cloning reaction allows for longer incubation times because it prevents topoisomerase I from re-binding and potentially nicking the DNA after ligating the PCR product and dissociating from the DNA. The result is more intact molecules leading to higher transformation efficiencies.

---

## Using salt in the TOPO<sup>®</sup> cloning reaction

You will perform TOPO<sup>®</sup> Cloning in a reaction buffer containing salt (*i.e.* using the stock salt solution provided in the kit). Note that the amount of salt added to the TOPO<sup>®</sup> Cloning reaction varies depending on whether you plan to transform chemically competent cells (provided) or electrocompetent cells (see page 26 for ordering information).

- If you are transforming chemically competent *E. coli* (included with the kit), use the stock Salt Solution as supplied, and set up the TOPO<sup>®</sup> Cloning reaction as directed on the next page.
  - If you are transforming electrocompetent *E. coli* (available separately from Invitrogen; see page 26), the amount of salt in the TOPO<sup>®</sup> Cloning reaction must be reduced to 50 mM NaCl, 2.5 mM MgCl<sub>2</sub> to prevent arcing during electroporation. Dilute the stock Salt Solution 4-fold with water to prepare a 300 mM NaCl, 15 mM MgCl<sub>2</sub> Dilute Salt Solution. Use the Dilute Salt Solution to set up the TOPO<sup>®</sup> Cloning reaction as directed on the next page.
- 

*Continued on next page*

## Setting Up the TOPO<sup>®</sup> Cloning Reaction, continued

---

### Materials needed

You should have the following materials on hand before beginning:

- Your PCR product(s) (freshly prepared)
  - pOptiVEC<sup>™</sup>-TOPO<sup>®</sup> vector
  - Salt Solution or Dilute Salt Solution (see previous page)
  - Sterile Water
- 

### Performing the TOPO<sup>®</sup> cloning reaction

The table below describes how to set up your TOPO<sup>®</sup> Cloning reaction (6  $\mu$ L) to use for transformation of either chemically competent or electrocompetent *E. coli*.

**Note:** The red color of the TOPO<sup>®</sup> vector solution is normal and is used to visualize the solution.

Reagent	Chemically Competent <i>E. coli</i>	Electrocompetent <i>E. coli</i>
PCR Product	0.5 to 4 $\mu$ L	0.5 to 4 $\mu$ L
Salt Solution	1 $\mu$ L	–
Dilute Salt Solution	–	1 $\mu$ L
Sterile Water	Add to total volume of 5 $\mu$ L	Add to total volume of 5 $\mu$ L
TOPO <sup>®</sup> Vector	1 $\mu$ L	1 $\mu$ L
<b>Final Volume</b>	<b>6 <math>\mu</math>L</b>	<b>6 <math>\mu</math>L</b>

\*Store all reagents at  $-20^{\circ}\text{C}$  when finished. Salt solution and water can be stored at room temperature or  $4^{\circ}\text{C}$ .

1. Mix reaction gently and incubate for 5 minutes at room temperature ( $22^{\circ}$ – $23^{\circ}\text{C}$ ).  
**Note:** For most applications, 5 minutes will yield a sufficient number of colonies for analysis. The length of the TOPO<sup>®</sup> Cloning reaction can be varied from 30 seconds to 30 minutes. For routine subcloning of PCR products, 30 seconds may be sufficient. For larger PCR products ( $>1$  kb), increasing the reaction time may yield more colonies.
  2. Place the reaction on ice and proceed to Transforming One Shot<sup>®</sup> Competent *E. coli*, next page.  
**Note:** You may store the TOPO<sup>®</sup> Cloning reaction overnight at  $-20^{\circ}\text{C}$ .
-

# Transforming One Shot<sup>®</sup> Competent *E. coli*

---

## Introduction

Once you have performed the TOPO<sup>®</sup> Cloning reaction you are ready to use your construct(s) to transform competent *E. coli*. One Shot<sup>®</sup> TOP10 Chemically Competent *E. coli* are included with the kit (Box 2) to facilitate transformation. You may also transform One Shot<sup>®</sup> Electrocompetent cells if desired (see page 26 for ordering information). Protocols for transforming chemically competent and electrocompetent *E. coli* are provided in this section.

---

## Selecting a One Shot<sup>®</sup> chemical transformation protocol

Two protocols are provided to transform One Shot<sup>®</sup> TOP10 chemically competent *E. coli*. Consider the following factors and choose the protocol that best suits your needs.

If you wish to...	Then use the...
Maximize the number of transformants	Regular chemical transformation protocol, page 15.
Clone large PCR products (>1000 bp)	
Obtain transformants as quickly as possible	Rapid chemical transformation protocol, page 15. <b>Note:</b> This procedure is less efficient; the total number of transformants obtained may be lower than that obtained with the regular chemical transformation protocol.

---

## Materials needed

In addition to general microbiological supplies (*i.e.* plates, spreaders), you will need the following:

- TOPO<sup>®</sup> Cloning reaction(s) (From Step 2, previous page)
  - One Shot<sup>®</sup> TOP10 *E. coli*, either chemically competent (supplied with the kit, Box 2) or electrocompetent (purchased separately, see page 26)
  - S.O.C. Medium (supplied with the kit, Box 2)
  - pUC19 positive control (supplied with the kit, Box 2)
  - 42°C water bath (chemically competent cells only)
  - Electroporator with cuvettes (electrocompetent cells only)
  - 15 mL sterile, snap cap plastic culture tubes (electrocompetent cells only)
  - For each transformation reaction, 2 selective LB plates containing 100 µg/mL ampicillin. See page 25 for a recipe to prepare selective LB
  - 37°C shaking and non-shaking incubators
- 

*Continued on next page*

# Transforming One Shot<sup>®</sup> Competent *E. coli*, continued

---

## Preparing for transformation

For each transformation, you will need one vial of One Shot<sup>®</sup> competent cells and two selective LB plates.

- Equilibrate a water bath to 42°C if using chemically competent *E. coli*, or set up your electroporator if using electrocompetent *E. coli*
  - Warm the vial of S.O.C. Medium to room temperature
  - Warm selective LB plates at 37°C for 30 minutes
  - Thaw one vial of One Shot<sup>®</sup> cells **on ice** for each transformation
- 

## One Shot<sup>®</sup> chemical transformation protocol

Use the following protocol to transform One Shot<sup>®</sup> TOP10 chemically competent *E. coli*.

1. Add 2 µL of each TOPO<sup>®</sup> Cloning reaction into a vial of One Shot<sup>®</sup> Chemically Competent *E. coli* with a sterile pipette tip and mix gently. Do not mix by pipetting up and down.  
**Note:** If you are using the pUC19 control plasmid for transformation, use 1 µL (10 pg).
  2. Incubate cells/plasmid mix on ice for 5–30 minutes.  
**Note:** Longer incubations on ice seem to have a minimal effect on transformation efficiency.
  3. Heat-shock the cells for 30 seconds at 42°C without shaking.
  4. Immediately transfer the tubes to ice.
  5. Add 250 µL of room temperature S.O.C. Medium.
  6. Cap the tube tightly and shake the tube horizontally at 200 rpm in a 37°C shaking incubator for 1 hour.
  7. Spread 10–50 µL from each transformation on a prewarmed selective LB plate. To ensure even spreading of small volumes, you may add 20 µL of S.O.C. Medium to the transformation mixture. We recommend that you plate two different volumes to ensure that at least one plate contains well-spaced colonies. Incubate plates overnight at 37°C.
- 

## Rapid One Shot<sup>®</sup> chemical transformation protocol

Use the alternative protocol below to rapidly transform One Shot<sup>®</sup> TOP10 chemically competent *E. coli*. Before beginning, prewarm LB plates containing 100 µg/mL ampicillin at 37°C for 30 minutes.

1. Add 4 µL of each TOPO<sup>®</sup> Cloning reaction into a vial of One Shot<sup>®</sup> TOP10 chemically competent *E. coli* and mix gently. Do not mix by pipetting up and down.
  2. Incubate reaction on ice 5 minutes.
  3. Spread 50 µL of cells on a prewarmed selective LB plate and incubate overnight at 37°C.
-

# Analyzing Positive Clones

---

## Introduction

After transformation of your pOptiVEC™ construct(s) into *E. coli*, you will select and analyze several colonies from each transformation by sequencing using the specific primers included in the kit to determine the orientation of the insert.

---

## Analyzing positive clones

1. Pick 10 colonies from each transformation and culture them overnight in LB medium containing 100 µg/mL ampicillin.
  2. Isolate plasmid DNA using your method of choice. If you need ultra-pure plasmid DNA for automated or manual sequencing, we recommend using the PureLink® HQ Mini Plasmid Purification Kit. See page 26 for ordering information.
  3. Analyze plasmid DNA by sequencing (see below).
- 

## Sequencing

To confirm that your DNA insert is in the correct orientation, you may sequence your expression construct using the CMV forward and EMCV IRES reverse primers included with the kit. Refer to page 5 for the sequences of the primers and the diagram on page 9 for the location of the primer binding sites.

---

## Long-term storage

Once you have identified the correct clone(s), be sure to purify the colony and make a glycerol stock for long-term storage.

1. Streak the original colony out for single colonies on an LB plate containing 100 µg/mL ampicillin.
2. Isolate a single colony and inoculate into 1–2 mL of LB containing 100 µg/mL ampicillin.
3. Grow at 37°C with shaking until culture reaches stationary phase.
4. Mix 0.85 mL of culture with 0.15 mL of sterile glycerol.
5. Transfer to a cryovial and store at –80°C.

We also recommend that you store a stock of plasmid DNA at –20°C.

---

## Next Steps

---

### Introduction

Once you obtain the correct pOptiVEC™-TOPO plasmid construct(s), you will linearize and purify the plasmid prior to transfection. General guidelines for purifying the plasmid, performing a positive expression control, transfecting cells and generating stable cell lines are given below.

**Note:** If you are using the Freedom™ DG44™ Kit, please refer to the product manual supplied with the system for specific instructions on preparing the plasmid.

---

### Plasmid preparation

The pOptiVEC™-TOPO® plasmid construct must be clean, sterile, and free from contamination with phenol and sodium chloride for transfection into DG44 cells. Contaminants may kill the cells, and salt will interfere with lipid complexing, decreasing transfection efficiency. We recommend isolating plasmid DNA using the PureLink® HiPure DNA Midiprep Kit (see page 26 for ordering information).

---

### Linearizing the plasmid

Prior to transfecting DG44 cells with your pOptiVEC™ construct, you may linearize the plasmid. While linearizing your vector may not improve transfection efficiency, it increases the chances that the vector integrates into the host cell genome without disrupting the gene of interest or other elements required for expression in mammalian cells.

We suggest using *Pvu* I, which cuts once in the ampicillin resistance gene. Other unique restriction sites are possible. A complete restriction map of pOptiVEC™ TOPO® is available at [www.invitrogen.com](http://www.invitrogen.com). **Be sure that your insert does not contain the restriction enzyme site you wish to use to linearize your vector.**

After digestion, precipitate the DNA and resuspend pellet in sterile water and re-quantify using your method of choice.

---

### Transfection reagents

A wide range of transfection reagents for plasmid DNA that provide high efficiency gene expression with minimal cytotoxicity across a broad range of adherent and suspension cell lines are available from Life Technologies. Go to [www.invitrogen.com/transfection](http://www.invitrogen.com/transfection) to learn more about transfection reagents for your particular application.

---



# Troubleshooting

## Introduction

The table below lists some potential problems solutions that may help you troubleshoot your TOPO<sup>®</sup> Cloning.

Problem	Possible Cause	Solution
Few or no colonies obtained from sample reaction, but transformation control yielded colonies	Incomplete extension during PCR	Include a final extension step of 7–30 minutes during PCR. Longer PCR products will need a longer extension time.
	Excess or dilute PCR product used in the TOPO <sup>®</sup> Cloning reaction	Reduce or concentrate the amount of PCR product.
	PCR primers contain 5' phosphates	Do not add 5' phosphates to your PCR primers.
	Used a proofreading polymerase or a <i>Taq</i> /proofreading polymerase mixture for PCR	Use <i>Taq</i> polymerase to add 3' A-overhangs to your PCR product by following the method on page 22.
	Large PCR product	<ul style="list-style-type: none"> <li>• Increase the amount of PCR product used in the TOPO<sup>®</sup> Cloning reaction.</li> <li>• Increase the incubation time of TOPO<sup>®</sup> Cloning reaction from 5 minutes to 30 minutes.</li> <li>• Gel-purify the PCR product to remove primer-dimers or other artifacts.</li> </ul>
	PCR reaction contains artifacts (i.e. not a single band on an agarose gel)	<ul style="list-style-type: none"> <li>• Optimize your PCR conditions.</li> <li>• Gel-purify your PCR product.</li> </ul>
PCR product does not contain sufficient 3' A-overhangs even though you used <i>Taq</i> polymerase	<ul style="list-style-type: none"> <li>• Increase the final extension time to ensure that all 3' ends are adenylated.</li> <li>• You may redesign your primers so that they contain a 5' G instead of a 5' T.</li> </ul> <p><b>Note:</b> <i>Taq</i> polymerase is most efficient at adding a non-template 3' A next to a C, and less efficient at adding a nontemplate 3' A next to another A (Brownstein <i>et al.</i>, 1996).</p>	

*Continued on next page*

## Troubleshooting, continued

Problem	Possible Cause	Solution
Large number of incorrect inserts cloned	PCR cloning artifacts	<ul style="list-style-type: none"> <li>• Gel-purify your PCR product to remove primer-dimers and other artifacts.</li> <li>• Optimize your PCR conditions.</li> <li>• Include a final extension step of 7–30 minutes during PCR.</li> </ul>
Few or no colonies obtained from sample reaction and the transformation control gave no colonies	One Shot <sup>®</sup> competent <i>E. coli</i> stored incorrectly	<ul style="list-style-type: none"> <li>• Store One Shot<sup>®</sup> competent <i>E. coli</i> at -80°C.</li> <li>• If you are using another <i>E. coli</i> strain, follow the manufacturer's instructions.</li> </ul>
	Insufficient amount of <i>E. coli</i> plated	Increase the amount of <i>E. coli</i> plated.
	Transformants plated on selective plates with the wrong antibiotic	Use LB ampicillin plates for selection.
No protein expression in mammalian cells after transfection	PCR primer does not contain Kozak translation initiation sequence	Add a Kozak consensus site to the forward PCR primer (see page 8), resynthesize your DNA and re-clone.
	Premature stop codons	Remove stop codons by your method of choice.
	Poor secretion leader (for secreted proteins)	Include the endogenous secretion leader, if possible.
	Sequence not optimized	Optimize the codon sequence of the gene of interest.

# Appendix

## Performing the Control Reactions

---

### Introduction

We recommend performing the following control TOPO<sup>®</sup> Cloning reactions the first time you use the kit to help you evaluate your results. Performing the control reactions involves producing a control PCR product expressing the LacZ $\alpha$  fragment using the reagents included in the kit. Successful TOPO<sup>®</sup> Cloning of the control PCR product in either direction will yield > 85% blue colonies on LB plates containing ampicillin and X-gal.

---

### Before starting

The following reagents should be prepared before performing the control reaction:

- Prepare stock X-gal solution (See page 25 for recipe)
  - For each transformation, you will need two LB plates containing 100  $\mu\text{g}/\text{mL}$  ampicillin and X-gal
- 

### Producing the control PCR product

1. In a 0.5 mL microcentrifuge tube, set up the following reaction in a 50  $\mu\text{L}$  volume.

Reagent	Amount
Control DNA Template	1 $\mu\text{L}$
10X PCR Buffer	5 $\mu\text{L}$
50 mM dNTPs	0.5 $\mu\text{L}$
Control PCR Primers (0.1 $\mu\text{g}/\mu\text{L}$ each)	1 $\mu\text{L}$
Sterile Water	41.5 $\mu\text{L}$
<i>Taq</i> polymerase (1 unit/ $\mu\text{L}$ )	1 $\mu\text{L}$
<b>Total Volume</b>	50 $\mu\text{L}$

2. Amplify the control PCR product using the following cycling parameters:

Step	Time	Temp.	Cycles
Initial Denaturation	2 min.	94°C	1X
Denaturation	1 min.	94°C	25X
Annealing	1 min.	60°C	
Extension	1 min.	72°C	
Final Extension	7 min.	72°C	1X

3. Remove 10  $\mu\text{L}$  from the reaction and analyze by agarose gel electrophoresis. A discrete 500 bp band should be visible.
- 

*Continued on next page*

## Performing the Control Reactions, continued

### Control TOPO<sup>®</sup> cloning reactions

Using the control PCR product generated in the steps above and the control vector, set up two 6  $\mu\text{L}$  TOPO<sup>®</sup> Cloning reactions as described below:

Reagent	“Vector Only”	“Vector + PCR Insert”
Control PCR Product	–	1 $\mu\text{L}$
Sterile Water	4 $\mu\text{L}$	3 $\mu\text{L}$
Salt Solution or Dilute Salt Solution	1 $\mu\text{L}$	1 $\mu\text{L}$
pOptiVEC <sup>™</sup> -TOPO <sup>®</sup> vector	1 $\mu\text{L}$	1 $\mu\text{L}$

1. Incubate at room temperature for 5 minutes and place on ice.
2. Use 2  $\mu\text{L}$  of the reaction to transform two separate vials of One Shot<sup>®</sup> competent cells using the procedure on page 14.
3. Spread 10–50  $\mu\text{L}$  of each transformation mix onto LB plates containing 100  $\mu\text{g}/\text{mL}$  ampicillin and X-gal. When plating small volumes, add 20  $\mu\text{L}$  of S.O.C. Medium to ensure even spreading. Be sure to plate two different volumes to ensure that at least one plate has well-spaced colonies.
4. Incubate plates overnight at 37°C.

### What you should see

The “vector + PCR insert” reaction should produce hundreds of colonies. Greater than 85% of these will be blue.

The “vector only” reaction should yield very few colonies (<15% of the vector + PCR insert plate) and these should be white.

### Transformation control

pUC19 plasmid is included to check the transformation efficiency of the One Shot<sup>®</sup> competent cells. Transform one vial of One Shot<sup>®</sup> TOP10 cells with 10  $\mu\text{g}$  of pUC19 using the protocol on page 14. Plate 10  $\mu\text{L}$  of the transformation reaction plus 20  $\mu\text{L}$  of S.O.C. on LB plates containing 100  $\mu\text{g}/\text{mL}$  ampicillin. The transformation efficiency should be  $1 \times 10^9$  cfu/ $\mu\text{g}$  DNA.

# Addition of 3' A-Overhangs Post-Amplification

---

## Introduction

TOPO<sup>®</sup> Cloning DNA amplified by proofreading polymerases into TOPO<sup>®</sup> Cloning vectors often results in very low cloning efficiencies. Proofreading polymerases remove the 3' A-overhangs necessary for TOPO<sup>®</sup> Cloning. A method for adding 3'As post-amplification is provided below.

---

## Before starting

You will need the following items:

- *Taq* polymerase
  - A heat block equilibrated to 72°C
  - Phenol-chloroform (optional)
  - 3 M sodium acetate (optional)
  - 100% ethanol (optional)
  - 80% ethanol (optional)
  - TE buffer (optional)
- 

## Procedure

This is just one method for adding 3' A-overhangs. Other protocols may be suitable.

1. After amplification with a proofreading polymerase, place vials on ice and add 0.7–1 unit of *Taq* polymerase per tube. Mix well. It is not necessary to change the buffer.
2. Incubate at 72°C for 8–10 minutes (do not cycle).
3. Place the vials on ice. The DNA amplification product is now ready for ligation into pOptiVEC<sup>™</sup>-TOPO<sup>®</sup>.


**Note:** If you plan to store your sample overnight before proceeding with TOPO<sup>®</sup> Cloning, extract your sample with an equal volume of phenol-chloroform to remove the polymerases. Ethanol-precipitate the DNA and resuspend in TE buffer using the starting volume of the PCR.

---

# Map and Features of pOptiVEC™-TOPO®

## Map

The map below shows the elements of the pOptiVEC™-TOPO® vector. **The complete sequence is available for downloading from [www.invitrogen.com](http://www.invitrogen.com) or by contacting Technical Support (page Error! Bookmark not defined.).**



### Comments for pOptiVEC™-TOPO® 4402 nucleotides

CMV promoter:	36-715
CMV forward primer binding site:	573-593
TOPO® cloning site:	730
EMCV IRES:	772-1359
EMCV IRES reverse primer binding site:	810-829
DHFR:	1372-1935
TK polyadenylation signal:	1975-2247
pUC origin (c):	2609-3282
Ampicillin ( <i>bla</i> ) resistance gene (c):	3424-4284
<i>bla</i> promoter (c):	4279-4383

(c) = complementary strand

*Continued on next page*

## Map and Features of pOptiVEC™-TOPO®, continued

### Features

The pOptiVEC™-TOPO® vector contains the following elements. Features have been functionally tested, and the vectors have been fully sequenced.

Feature	Benefit
Full length human cytomegalovirus (CMV) immediate-early promoter/enhancer	Allows efficient, high-level expression of your recombinant protein (Andersson <i>et al.</i> , 1989; Boshart <i>et al.</i> , 1985; Hennighausen & Fleckenstein, 1986; Nelson <i>et al.</i> , 1987)
CMV forward primer	Allows sequencing of the insert
TOPO® Cloning site	Allows insertion of your PCR product
EMCV IRES reverse primer	Allows sequencing of the insert
Internal Ribosome Entry Site (IRES) from the Encephalomyocarditis virus (EMCV)	Allows cap-independent translation of DHFR (Gurtu <i>et al.</i> , 1996; Rees <i>et al.</i> , 1996)
Dihydrofolate reductase (DHFR) gene	Allows auxotrophic selection of transfected DG44 cells and for genomic amplification of stable cell lines using methotrexate (MTX) (Kaufman <i>et al.</i> , 1985)
Herpes Simplex Virus Thymidine Kinase (TK) polyadenylation signal	Allows efficient transcription termination and polyadenylation of mRNA (Cole & Stacy, 1985)
pUC origin	Allows high-copy number replication and growth in <i>E. coli</i>
Ampicillin ( <i>bla</i> ) resistance gene (β-lactamase)	Allows selection of transformants in <i>E. coli</i>

# Recipes

---

## LB (Luria-Bertani) medium and plates

### **Composition:**

1.0% Tryptone  
0.5% Yeast Extract  
1.0% NaCl  
pH 7.0

1. For 1 liter, dissolve 10 g tryptone, 5 g yeast extract, and 10 g NaCl in 950 mL deionized water.
2. Adjust the pH of the solution to 7.0 with NaOH and bring the volume up to 1 liter.
3. Autoclave on liquid cycle for 20 minutes at 15 psi. Allow solution to cool to 55°C and add antibiotic if needed.
4. Store at room temperature or at 4°C.

### **LB agar plates**

1. Prepare LB medium as above, but add 15 g/L agar before autoclaving.
  2. Autoclave on liquid cycle for 20 minutes at 15 psi.
  3. After autoclaving, cool to ~55°C, add antibiotic if needed, and pour into 10 cm plates.
  4. Let harden, then invert and store at 4°C.
- 

## X-Gal stock solution

1. To make a 40 mg/mL stock solution, dissolve 400 mg X-Gal in 10 mL dimethylformamide.
  2. Protect from light by storing in a brown bottle at -20°C.
  3. To add to previously made agar plates, warm the plate to 37°C. Pipette 40 µL of the 40 mg/mL stock solution onto the plate, spread evenly, and let dry 15 minutes.
  4. Protect plates from light.
-



## Accessory Products

**Additional Products** The products listed in this section are available separately from Life Technologies Corporation and may be used with the pOptiVEC™-TOPO® TA Cloning Kit. Ordering information is provided below.

Product	Amount	Catalog no.
Platinum® <i>Taq</i> DNA Polymerase	100 reactions	10966-018
	250 reactions	10966-026
	500 reactions	10966-034
Platinum® <i>Taq</i> DNA Polymerase High Fidelity	100 reactions	11304-011
	500 reactions	11304-029
PCR Optimizer™ Kit	100 reactions	K1220-01
One Shot® TOP10 Chemically Competent <i>E. coli</i>	10 reactions	C4040-10
	20 reactions	C4040-03
One Shot® TOP10 Electrocompetent <i>E. coli</i>	10 reactions	C4040-50
	20 reactions	C4040-52
LB Broth (1X), liquid	500 mL	10855-021
Ampicillin Sodium Salt, irradiated	200 mg	11593-027
Carbenicillin, Disodium Salt	5 g	10177-012
PureLink® HQ Mini Plasmid DNA Purification Kit	100 preps	K2100-01
PureLink® HiPure Plasmid Midiprep Kit	25 preps	K2100-04
PureLink® Quick Gel Extraction Kit	50 preps	K2100-12
E-Gel® 1.2% Starter Pak (6 gels + Powerbase™)	1 kit	G6000-01
E-Gel® 1.2% 18 Pak	18 gels	G5018-01
CHO DG44 Cells (cGMP banked) and Media Kit	1 kit	A11000-01
CD DG44 Medium	1000 mL	12610-010
OptiPRO™ SFM	100 mL	12309-050
CD OptiCHO™ Medium	1000 mL	12681-011
L-glutamine, 200 mM	100 mL	25030-081
FreeStyle™ MAX Transfection Reagent	1 mL	16447-100
Pluronic® F-68, 10%	100 mL	24040-032
Freedom™ DG44™ Kit	1 kit	A13737-01

# Technical Support

---

## Obtaining support

For the latest services and support information for all locations, go to [www.invitrogen.com](http://www.invitrogen.com). At the website, you can:

- Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities
  - Search through frequently asked questions (FAQs)
  - Submit a question directly to Technical Support ([techsupport@invitrogen.com](mailto:techsupport@invitrogen.com))
  - Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents
  - Obtain information about customer training
  - Download software updates and patches
- 

## Safety data sheets (SDS)

Safety Data Sheets (SDSs) are available at [www.invitrogen.com/sds](http://www.invitrogen.com/sds).

---

## Certificate of analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to [www.invitrogen.com/support](http://www.invitrogen.com/support) and search for the Certificate of Analysis by product lot number, which is printed on the box.

---

## Limited warranty

Life Technologies Corporation is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about a Life Technologies product or service, contact our Technical Support Representatives. All Life Technologies products are warranted to perform according to specifications stated on the certificate of analysis. The Company will replace, free of charge, any product that does not meet those specifications. This warranty limits the Company's liability to only the price of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order.

Life Technologies makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, report it to our Technical Support Representatives.

**Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.**

---

## Purchaser Notification

---

Limited Use Label  
License: Research  
Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use this product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

---

Limited Use Label  
License No. 298:  
IRES

This product is licensed under U.S. Patent No. 4,937,190 for use in research by laboratories of academic, government, industrial and/or clinical institutions engaged in the investigation of biological or biochemical processes, including applications for the enhancement of translation in living cells, but specifically excluding applications in cell-free systems or use in diagnostic or therapeutic applications or for which FDA approval is required. For information on licensing for uses other than those permitted above, contact EMD Biosciences, 441 Charmany Drive, Madison, WI, 53719 USA. Tel. (608) 442-4006.

---

## References

---

- Andersson, S., Davis, D. L., Dahlbäck, H., Jörnvall, H., and Russell, D. W. (1989) Cloning, Structure, and Expression of the Mitochondrial Cytochrome P-450 Sterol 26-Hydroxylase, a Bile Acid Biosynthetic Enzyme. *J. Biol. Chem.* *264*, 8222-8229
- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994) *Current Protocols in Molecular Biology*, Greene Publishing Associates and Wiley-Interscience, New York
- Boshart, M., Weber, F., Jahn, G., Dorsch-Häsler, K., Fleckenstein, B., and Schaffner, W. (1985) A Very Strong Enhancer is Located Upstream of an Immediate Early Gene of Human Cytomegalovirus. *Cell* *41*, 521-530
- Brownstein, M. J., Carpten, J. D., and Smith, J. R. (1996) Modulation of Non-Templated Nucleotide Addition by *Taq* DNA Polymerase: Primer Modifications that Facilitate Genotyping. *BioTechniques* *20*, 1004-1010
- Cole, C. N., and Stacy, T. P. (1985) Identification of Sequences in the Herpes Simplex Virus Thymidine Kinase Gene Required for Efficient Processing and Polyadenylation. *Mol. Cell. Biol.* *5*, 2104-2113
- Coloma, M. J., Hastings, A., Wims, L. A., and Morrison, S. L. (1992) Novel Vectors for the Expression of Antibody Molecules Using Variable Regions Generated by Polymerase Chain Reaction. *J. Imm. Methods* *152*, 89-104
- Gurtu, V., Yan, G., and Zhang, G. (1996) IRES bicistronic expression vectors for efficient creation of stable mammalian cell lines. *Biochem. Biophys. Res. Comm.* *229*, 295-298
- Hennighausen, L., and Fleckenstein, B. (1986) Nuclear factor 1 interacts with five DNA elements in the promoter region of the human cytomegalovirus major immediate early gen. *Embo J* *5*, 1367-1371
- Kaufman, R., Wasley, L., Spiliotes, A., Gossels, S., Latt, S., Larsen, G., and Kay, R. (1985) Coamplification and coexpression of human tissue-type plasminogen activator and murine dihydrofolate reductase in Chinese hamster ovary cells. *Mol Cell Biol* *5*, 1750-1759
- Kozak, M. (1987) An Analysis of 5'-Noncoding Sequences from 699 Vertebrate Messenger RNAs. *Nucleic Acids Res.* *15*, 8125-8148
- Kozak, M. (1990) Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. *Proc. Natl. Acad. Sci. USA* *87*, 8301-8305
- Kozak, M. (1991) An Analysis of Vertebrate mRNA Sequences: Intimations of Translational Control. *J. Cell Biology* *115*, 887-903
- Nelson, J. A., Reynolds-Kohler, C., and Smith, B. A. (1987) Negative and Positive Regulation by a Short Segment in the 5'-Flanking Region of the Human Cytomegalovirus Major Immediate-Early Gene. *Molec. Cell. Biol.* *7*, 4125-4129
- Rees, S., Coote, J., Stables, J., Goodson, S., Harris, S., and Lee, M. (1996) Bicistronic vector for the creation of stable cell lines that predisposes all antibiotic-resistant cells to express recombinant protein. *Biotechniques* *20*, 102-110
- Shuman, S. (1991) Recombination Mediated by Vaccinia Virus DNA Topoisomerase I in *Escherichia coli* is Sequence Specific. *Proc. Natl. Acad. Sci. USA* *88*, 10104-10108
- Shuman, S. (1994) Novel Approach to Molecular Cloning and Polynucleotide Synthesis Using Vaccinia DNA Topoisomerase. *J. Biol. Chem.* *269*, 32678-32684
- Werner, R. (2005) The development and production of biopharmaceuticals: technology and economic success factors. *Bioprocess International* *September*, 6-15

---

©2011 Life Technologies Corporation. All rights reserved.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

Pluronic is a registered trademark of BASF Corporation.

Triton is a registered trademark of Union Carbide Corporation

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

**Headquarters**

5791 Van Allen Way | Carlsbad, CA 92008 USA | Phone +1 760 603 7200 | Toll Free in USA 800 955 6288

**For support visit** [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)

[www.lifetechnologies.com](http://www.lifetechnologies.com)

