

Dynabeads[®] Human Treg Expander

Catalog no. 11129D

Store at 2 to 8°C

Rev. Date: October 2011 (Rev. 004)

Product Contents

Dynabeads[®] Human Treg Expander at 2×10^7 beads/mL concentration are supplied in phosphate buffered-saline (PBS), pH 7.4, with 0.1% human serum albumin (HSA) in a total volume of 2 mL.

Product Description

The Dynabeads[®] Human Treg Expander is intended for activation and expansion of human Treg cells isolated with the Dynabeads[®] Regulatory CD4⁺CD25⁺ T Cell Kit. The expanded Treg cells retain their regulatory capacity.

Dynabeads[®] Human Treg Expander offers a simple expansion method for Treg cells that does not require antigen-presenting cells or antigen. To expand the Treg cells, add Dynabeads[®] Human Treg Expander and recombinant IL-2 (rIL-2) to the Treg culture.

For polyclonal activation of naive T cells, use Dynabeads[®] Human T-Activator CD3/CD28. For expansion of antigen-specific T cells, use Dynabeads[®] Human T-Activator CD3/CD28/CD137 (see "Related Products")

Required Materials:

- Magnet (DynaMag[™]): See www.lifetechnologies.com/magnets for magnet recommendations.
- Buffer: PBS with 0.1% BSA, pH 7.4.
- Culture medium: 1X OpTmizer[™] T-Cell Expansion SFM.
- Round bottom tissue culture plates or tissue culture flasks.
- Humidified CO₂ incubator.
- rIL-2.

Product Use: For research use only. **Caution:** Not for human or animal therapeutic or diagnostic use.

Table 1: Volume recommendations for bead-to-cell ratio = 4:1

Specifications	For 1×10^5 Treg cells	For 1×10^6 Treg cells
Culture type	Cells/well in 96-well plate	Cells/well in 24-well plate
Dynabeads [®] Treg Expander	20 μ L	200 μ L
rIL-2	500 U/mL	500 U/mL
Vol. seeding medium	100 μ L	1 mL

General Guidelines

- Use a bead-to-cell ratio of 4:1 and the cell concentrations listed in Table 1. Other bead-to-cell ratios and cell concentrations are suitable following optimization for your application.
- Use at least 95% pure Treg cells for optimal expansion results. Low purity allows CD8⁺ and CD4⁺CD25⁺ T cells to grow and reduces the number of Treg cells in relative terms.
- If the purity is lower than recommended above, increase the amount of IL-2 up to 1000 U/mL rIL-2 to improve the expansion.

Wash Dynabeads[®]

Wash Dynabeads[®] before use.

1. Resuspend the Dynabeads[®] in the vial (i.e. vortex for >30 sec, or tilt and rotate for 5 min).
2. Transfer the desired volume of Dynabeads[®] to a tube.
3. Add an equal volume of Buffer, or at least 1 mL, and mix.
4. Place the tube on a magnet for 1 min and discard the supernatant.
5. Remove the tube from the magnet and resuspend the washed Dynabeads[®] in the same volume of culture medium as the initial volume of Dynabeads[®] taken from the vial (step 2).

Isolate Human Treg Cells

Isolate Human Treg cells using the Dynabeads[®] Regulatory CD4⁺CD25⁺ T Cell Kit.

Expand Human Treg Cells

Day 0

1. Start with 1×10^5 Treg cells in 100 μ L culture medium in a round bottom 96-well tissue culture plate.
2. Add 20 μ L pre-washed and resuspended Dynabeads[®] to the cells.
3. Add 500 U/ml rIL-2 to the cells.

Day 1

4. Add 100 μ L culture medium containing 500 U/mL rIL-2.

Day 3

5. Resuspend and split the culture in half and add 100 μ L culture medium containing 500 U/mL rIL-2 per well.

Days 5–7

6. Split the culture when needed. Examine cultures daily, noting cell size and shape. Cell shrinking and reduced proliferation rate is typically observed in exhausted cell cultures.
7. Count the cells at least twice weekly after thorough resuspension.
8. When the cell density exceeds 2×10^6 cells/mL or when the medium becomes yellow, split cultures to a density of $0.5\text{--}1 \times 10^6$ cells/mL in culture medium containing 500 U/mL rIL-2.
9. Grow the cells until the well is half-full (~500,000 cells) and transfer the cells from a 96-well plate to a 24-well plate.

Day 8

10. Remove the Dynabeads[®] by resuspending the cells and transferring the cells to a suitable tube.
11. Place the tube in a magnet for 1–2 min until the Dynabeads[®] are separated.
12. Centrifuge the supernatant and resuspend the cell pellet in fresh culture medium containing 100 U/mL rIL-2.
13. Split the cultures when needed and rest the cells until days 21–24 in culture medium with 100 U/mL rIL-2.
14. Incubate the cells in a humidified CO₂ incubator at 37°C.

Restimulation

You can restimulate the cells 15–18 days after the first stimulation or when cell shrinking and reduced rate of proliferation is observed. Restimulation guidelines are listed in Table 2.

Table 2: Volume recommendations for bead-to-cell ratio = 1:1

Specifications	For 1×10^6 Treg cells	For 1×10^6 Treg cells
Culture type	Cells/well in 96-well plate	Cells/well in 24-well plate
Dynabeads® Treg Expander	5 μ L	50 μ L
rIL-2	100 U/mL	100 U/mL
Vol. seeding medium	100 μ L	1 mL

Before restimulation

- Count the cells and split the cultures to a density of 1×10^6 cells/mL in culture medium.
- Use the volumes given in Table 2 and follow the protocol below.

Day 0

- Start with 1×10^6 Treg cells in 1 mL culture medium in a 24-well tissue culture plate.
- Add 50 μ L Dynabeads®.
- Add 100 U/mL rIL-2.

Day 1

- Add 100 μ L culture medium containing 100 U/mL rIL-2 per well.

Day 2 or 3

- Resuspend and split the culture in half and add 100 μ L culture medium containing 100 U/mL rIL-2 per well.

Days 5–7

- Split the culture when needed. Examine cultures daily, noting cell size and shape. Cell shrinking and reduced proliferation rate is typically observed in exhausted cell cultures.
- Count the cells at least twice weekly after thorough resuspension.
- When the cell density exceeds 2×10^6 cells/mL or when the medium becomes yellow, split cultures to a density of $0.5\text{--}1 \times 10^6$ cells/mL in culture medium containing 100 U/mL rIL-2.

Day 8

- Remove the Dynabeads® by resuspending the cells and transferring the cells to a suitable tube
- Place the tube in a magnet for 1–2 min until the Dynabeads® are separated.
- Centrifuge the supernatant and resuspend the cell pellet in fresh culture medium containing 100 U/mL rIL-2.
- Split the cultures when needed and rest the cells until days 21–24 in culture medium with 100 U/mL rIL-2.

Mixed Lymphocyte Reaction (MLR) Assay for Identification of Suppressive Capacity of Expanded Treg Cells

- Isolate dendritic cells. Culture with TNF- α and LPS for 2 days, and irradiate before use as stimulators in MLR.
- Establish co-cultures of 5×10^4 responder CD4⁺CD25⁻ T cells and 1×10^4 irradiated allogenic dendritic cells as stimulators in 96-well U-bottom plates.
- Add freshly purified or *in vitro* expanded CD4⁺CD25⁺ Treg cells in ratios of 1:1, 1:4, 1:8, 1:16 and 1:32 (CD25⁺:CD25⁻ ratio).
- Pulse the wells on day 6 with ³H-thymidine, and culture for an additional 18 hours before harvest.

Description of Materials

Dynabeads® Human Treg Expander are uniform 4.5 μ m, superparamagnetic polystyrene beads coated with an optimized mixture of monoclonal antibodies against the CD3 and CD28 cell surface molecules of human T cells. The CD3 antibody coated on Dynabeads® Human Treg Expander is specific for the epsilon chain of human CD3, a subunit of the TCR complex. The CD28 antibody is specific for the human CD28 co-stimulatory molecule, which is the receptor for CD80 (B7-1) and CD86 (B7-2). Both antibodies are coupled to the same bead, mimicking *in vivo* stimulation by antigen presenting cells.

Related Products

Product	Cat. No.
DynaMag™-5	12303D
DynaMag™-15	12301D
Dynabeads® Human T-Activator CD3/CD28/CD137	11162D
Dynabeads® Regulatory CD4 ⁺ CD25 ⁺ T Cell Kit	11323D
Dynabeads® Human T-Activator CD3/CD28	11131D
OpTmizer™ T-Cell Expansion SFM	0080022SA

REF on labels is the symbol for catalog number.

Limited Use Label License

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

Manufactured by Life Technologies AS, Norway. Life Technologies AS complies with the Quality System Standards ISO 9001:2008 and ISO 13485:2003.

SPEC-06044

©2011 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners, except where otherwise stated. Life Technologies and/or its affiliate(s) disclaim all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall Life Technologies and/or its affiliate(s) be liable, whether in contract, tort, warranty, or under any statute or on any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including but not limited to the use thereof.

For support visit www.lifetechnologies.com/support or email techsupport@lifetech.com

www.lifetechnologies.com

