

## Magnetic 96-Well Separator

Publication Part no. MAN0005195

Rev. Date: 17 October 2011

Catalog Number: A14179


Store at room temperature

### Product Description


The Magnetic 96-Well Separator is designed for fast and simple removal of supernatant from samples bound to magnetic particles, and may be used manually or as an integrated part of an automated workstation. The Magnetic 96-Well Separator is ideal for use in singleplex and multiplex biological assays based on the Luminex<sup>®</sup> xMAP<sup>®</sup> Technology. These assays make use of Luminex<sup>®</sup> MagPlex, or magnetic microspheres, and the MAGPIX<sup>®</sup>, Luminex 200<sup>®</sup>, or FLEXMAP 3D<sup>®</sup> instruments.

### Caution


The Magnetic 96-Well Separator is a powerful magnet and should be handled with care to avoid personal injury. Hands and skin can be pinched if caution is not taken. Keep sensitive electronic devices away from the Magnetic 96-Well Separator.



### Recommended manual washing procedure for use of the Magnetic 96-Well Separator and magnetic microspheres




1. Place a 96-well flat-bottom plate with beads onto the plate magnet.  
**Note:** A standard Greiner 96-well mylar flat-bottom plate (Cat. no. 655096), or equivalent is recommended for use with the Magnetic 96-Well Separator.
2. Allow the plate to soak for 60–90 seconds.



3. Lift the plate magnet while holding the plate securely in place on top, and invert the magnet and plate (held securely together) over an appropriate disposal container.





4. Gently shake out excess fluid.  
**Important:** When discarding fluid over a sink or waste container, use caution. Splashing may lead to contamination of plate wells and negatively impact results. Avoid adding bleach to the disposal container until washing is complete.



5. Keep the magnet and plate (held securely together) upside down, and blot the plate on a short stack of paper towels.
6. Tap the magnet and plate on dry paper towels several times to assure equal evacuation from all wells.

**For research use only. Not for human or animal therapeutic or diagnostic use.**

Manufacturing Site • 7335 Executive Way • Frederick • MD 21704 • E-mail: [techsupport@invitrogen.com](mailto:techsupport@invitrogen.com)



7. Return the magnet and plate to the right-side up position on the benchtop.
8. Remove the plate from the magnet to add additional wash solution, or for incubation.
9. Repeat, for any additional washes as needed.

### Recommended Pipetting Technique

We recommend using a modified reverse pipetting technique for all pipetting during assay setup to increase precision and limit reagent loss due to residual fluid volumes within pipette tips. A multi-channel pipette is also recommended for adding common reagents to the plate.

1. Set the pipette to the appropriate volume needed.
2. Press the plunger button slowly to the first stop and then press slightly past it (just enough past to aspirate a very small extra volume of the reagent). Alternatively, to mix the solution or sample before pipetting, pipet up and down slowly several times going to the first stop only (this allows residual fluid to build up in the tip).
3. Immerse the tip into the liquid, just below the meniscus.
4. Release the plunger button slowly and smoothly to the top resting position to aspirate the set volume of liquid. Drag the tip up the side of the tube or reservoir to remove excess volume from the outside of the tip(s).
5. Place the end of the tip against the inside wall of the recipient vessel at an angle above the fluid level.
6. Press the plunger button slowly and smoothly to the first stop. **Do not go past the first stop.** A small volume of residual liquid remains in the tip, and should not be dispensed.
7. Drag the tip out of the well while keeping the pipette pressed to the first stop and return to step 3 above if reusing tips and contamination is not an issue or change tips when contamination is possible.

### Safety Data Sheets (SDS)

Safety Data Sheets (SDSs) are available at [www.invitrogen.com/sds](http://www.invitrogen.com/sds).

### Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to [www.invitrogen.com/support](http://www.invitrogen.com/support) and search for the Certificate of Analysis by product lot number, which is printed on the box.

#### Explanation of symbols

Symbol	Description	Symbol	Description
	Catalogue Number		Batch code
	Research Use Only		<i>In vitro</i> diagnostic medical device
	Use by		Temperature limitation
	Manufacturer		European Community authorised representative
	Without, does not contain		With, contains
	Protect from light		Consult accompanying documents
	Directs the user to consult instructions for use (IFU), accompanying the product.		

### Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

©2011 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

For support visit [www.invitrogen.com/support](http://www.invitrogen.com/support) or email [techsupport@invitrogen.com](mailto:techsupport@invitrogen.com)

[www.lifetechnologies.com](http://www.lifetechnologies.com)

**life**  
technologies™