

USER GUIDE

ambion[®]
by *life* technologies[™]

MessageAmp[™] II-Biotin *Enhanced* Kit

Publication Part Number 1791M Rev. E

Revision Date November 2011

life
technologies[™]

For Research Use Only. Not intended for any animal or human therapeutic or diagnostic use.

Information in this document is subject to change without notice.

LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

NOTICE TO PURCHASER: LIMITED USE LABEL LICENSE: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

TRADEMARKS

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. LabChip is a trademark of Caliper Life Sciences Inc. GeneChip is a trademark of Affymetrix Inc.

© 2011 Life Technologies Corporation. All rights reserved.

Part Number 1791M Rev. E

November 2011

Contents

About This Guide	5
Purpose	5
User attention words	5
■ MessageAmp™ II-Biotin <i>Enhanced</i> Kit	7
Introduction	7
Product Description and Background	7
Procedure Overview	9
Materials Provided with the Kit and Storage Conditions	10
aRNA Amplification Procedure	12
Important Parameters for Successful Amplification	12
Prepare the Wash Buffer	16
Reverse Transcription to Synthesize First Strand cDNA	16
Second Strand cDNA Synthesis	18
cDNA Purification	19
In Vitro Transcription to Synthesize Biotin-labeled aRNA	20
aRNA Purification	21
Evaluation and Fragmentation of aRNA	23
aRNA Quantitation and Expected Yield	23
Analysis of aRNA Size	24
Fragmentation of Biotinylated aRNA for GeneChip® Arrays	24
Troubleshooting	26
Positive Control Reaction	26
Factors that Affect Both the Positive Control and Experimental Samples	26
Troubleshooting Low Yield and Small Average aRNA Size	27
■ APPENDIX A Materials Not Provided with the Kit	29
Materials Not Provided with the Kit	29
Related Products Available from Life Technologies	30
■ APPENDIX B Safety	31
General Safety	31
Chemical safety	31
Biological hazard safety	32

Bibliography	33
Documentation and Support	35
Obtaining SDSs	35
Obtaining support	35

About This Guide



WARNING! ABBREVIATED SAFETY ALERTS. Hazard symbols and hazard types specified in procedures may be abbreviated in this document. For the complete safety information, see the “Safety” appendix in this document.

IMPORTANT! Before using this product, read and understand the information the “Safety” appendix in this document.

Purpose

The *MessageAmp™ II-Biotin Enhanced Kit User Guide* provides detailed procedures, reference information and troubleshooting for the kit.

User attention words

Five user attention words may appear in this document. Each word implies a particular level of observation or action as described below:

Note: Provides information that may be of interest or help but is not critical to the use of the product.

IMPORTANT! Provides information that is necessary for proper instrument operation or accurate chemistry kit use.



CAUTION! Indicates a potentially hazardous situation that, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.



WARNING! Indicates a potentially hazardous situation that, if not avoided, could result in death or serious injury.



DANGER! Indicates an imminently hazardous situation that, if not avoided, will result in death or serious injury.

Except for IMPORTANTs, the safety alert words in user documentation appear with an open triangle figure that contains a hazard symbol. These hazard symbols are identical to the hazard symbols that are affixed to the instrument. See the “Safety” appendix for descriptions of the symbols.



MessageAmp™ II-Biotin *Enhanced* Kit

Introduction

Product Description and Background


The MessageAmp™ II-Biotin *Enhanced* Single Round aRNA Amplification Kit is based on the RNA amplification protocol developed in the Eberwine laboratory (Van Gelder et al. 1990). The procedure consists of reverse transcription with an oligo(dT) primer bearing a T7 promoter using ArrayScript™ reverse transcriptase (RT), engineered to produce higher yields of first-strand cDNA than wild-type enzymes. ArrayScript RT catalyzes the synthesis of virtually full-length cDNA, which is the best way to ensure production of reproducible microarray samples. The cDNA then undergoes second-strand synthesis and cleanup to become a template for in vitro transcription (IVT) in a reaction containing biotin-modified UTP and T7 RNA polymerase. To maximize biotin-labeled aRNA yield, an optimized mixture of biotin-labeled and unlabeled NTPs are supplied with the kit, and Ambion MEGAscript® IVT technology is used to generate hundreds to thousands of antisense RNA copies of each mRNA in a sample. (In this protocol the antisense amplified RNA is referred to as aRNA; it is also commonly called cRNA.) Once purified, the biotin-labeled aRNA is suitable for use on microarray gene expression systems designed for biotin-labeled antisense RNA samples.

Systematically optimized for production of biotin-labeled aRNA for array analysis

In order to develop the MessageAmp II-Biotin *Enhanced* Kit, we systematically optimized the MEGAscript T7 IVT amplification reaction for appropriate incorporation of biotin-modified UTP. We found that both the length of the linker arm connecting the biotin to the uridine base and the concentration of the modified nucleotide in the IVT reaction were important for maximal aRNA synthesis. Previous widely accepted labeling protocols included two biotin-modified nucleotides (CTP and UTP). We and others (Dorris et al.), however, found that the contribution of signal from biotin-CTP is minimal, and that the majority of the signal on arrays results from incorporated biotin-UTP. To further optimize the reaction, Life Technologies evaluated several biotin-modified UTP compounds. In these studies, biotin-11-UTP proved to have good incorporation, minimal effect on aRNA recovery during purification, and high signal on most commercial microarrays.

Figure 1 MessageAmp II-Biotin Enhanced aRNA Amplification Kit Procedure


Reverse Transcription to Synthesize First Strand cDNA



1. "Adjust RNA sample volume to 10 µL, if necessary" on page 16
2. "Add 10 µL of Reverse Transcription Master Mix and place at 42°C" on page 17
3. "Incubate for 2 hr at 42°C" on page 17

(Potential stopping point)


Second Strand cDNA Synthesis



1. "Add 80 µL Second Strand Master Mix to each sample" on page 18
2. "Incubate for 2 hr at 16°C" on page 18

(Potential stopping point)


cDNA Purification



1. "Preheat Nuclease-free Water to 55°C" on page 19
2. "Add 250 µL cDNA Binding Buffer to each sample" on page 19
3. "Pass the mixture through a cDNA Filter Cartridge" on page 19
4. "Wash with 500 µL Wash Buffer" on page 19
5. "Elute cDNA with 22 µL 55°C Nuclease-free Water" on page 19

(Potential stopping point)


In Vitro Transcription to Synthesize Biotin-labeled aRNA



1. "Add 20 µL of IVT Master Mix to each sample, and mix" on page 20
2. "Incubate 4–14 hr at 37°C" on page 21
3. "Add 60 µL Nuclease-free Water to each sample" on page 21

(Potential stopping point)

aRNA Purification



1. "Preheat Nuclease-free Water and assemble Filter Cartridges and Tubes" on page 21
2. "Add 350 µL aRNA Binding Buffer to each sample" on page 21
3. "Add 250 µL 100% ethanol and pipet 3 times to mix" on page 21
4. "Pass samples through an aRNA Filter Cartridge" on page 22
5. "Wash with 650 µL Wash Buffer" on page 22
6. "Elute aRNA with 200 µL 55°C Nuclease-free Water" on page 22
7. "Store aRNA at -80°C" on page 22

(Potential stopping point)

8. "[Optional] Concentrate the purified aRNA" on page 22



Benefits of RNA amplification

RNA amplification was originally developed as a method to expand very small RNA samples to produce enough material for array hybridization (Yue et al. 2001). Several groups have conducted studies to determine whether amplification of RNA introduces bias and they report that any bias is minimal (Li et al. 2004, Feldman et al. 2002 and Polacek et al. 2003). Additionally, among the benefits of amplification is a more reproducible expression profile from a wide range of RNA inputs. Some researchers conclude that amplification actually improves the reliability of array results regardless of whether it is needed for sample expansion (Feldman et al. 2002 and Polacek et al. 2003). As a result, RNA amplification has become the standard method for preparing RNA samples for array analysis (Kacharina et al. 1999, Pabon et al. 2001).

Procedure Overview

The MessageAmp II-Biotin aRNA amplification procedure is depicted in Figure 1.

- **Reverse Transcription to Synthesize First Strand cDNA** is primed with the T7 Oligo(dT) Primer to synthesize cDNA containing a T7 promoter sequence.
- **Second Strand cDNA Synthesis** converts the single-stranded cDNA into a double-stranded DNA (dsDNA) template for transcription. The reaction employs DNA Polymerase and RNase H to simultaneously degrade the RNA and synthesize second-strand cDNA.
- **cDNA Purification** removes RNA, primers, enzymes, and salts that would inhibit in vitro transcription.
- **In Vitro Transcription to Synthesize aRNA** with Biotin-NTP Mix generates multiple copies of biotin-modified aRNA from the double-stranded cDNA templates; this is the amplification step.
- **aRNA Purification** removes unincorporated NTPs, salts, enzymes, and inorganic phosphate to improve the stability of the biotin-modified aRNA.

To include two rounds of amplification, purchase the MessageAmp II Kit

If your experiments require two rounds of amplification, you will need to purchase the Ambion® MessageAmp II aRNA Amplification Kit (Part no. AM1751) in addition to the MessageAmp II-Biotin *Enhanced* Kit. For the first round amplification, use the MessageAmp II Kit to make unmodified aRNA.

For the second round of amplification, continue to use the MessageAmp II aRNA Amplification Kit (AM1751) until you have completed the second round cDNA synthesis step. At this point, follow the cDNA purification protocol for the MessageAmp II-Biotin *Enhanced* Kit (AM1791); this will provide the correct cDNA elution volume for the IVT reaction. Continue to follow the MessageAmp II-Biotin *Enhanced* Kit protocol for the second round amplification IVT reaction to generate biotin-labeled aRNA.

The MessageAmp II-Biotin advantage

Each step in the MessageAmp II-Biotin *Enhanced* Kit amplification procedure has been streamlined and optimized. The first-strand cDNA synthesis reaction employs ArrayScript reverse transcriptase to ensure that every cDNA bears a T7 promoter at its 5' end and that even very limited amounts of mRNA are fully converted to full-length cDNA. The second-strand cDNA synthesis reaction is designed for the efficient synthesis of full-length, double-stranded cDNAs and the complete conversion of single-stranded cDNA into double-stranded transcription templates. The cDNA purification procedure not only removes enzymes, salts, and unincorporated dNTPs,

but also efficiently removes RNA from the cDNA sample. This eliminates the heating or enzymatic digestion step commonly used in other procedures to degrade RNA (especially ribosomal RNA). The IVT reaction features MEGAscript technology to maximize transcriptional amplification and yield of aRNA. It is optimized to ensure efficient transcription of limited amounts of input DNA and synthesis of long transcripts.

Materials Provided with the Kit and Storage Conditions

The MessageAmp II-Biotin *Enhanced* Kit includes reagents for single-round amplification and biotin-labeling of 20 samples.

The MessageAmp II-Biotin Kit is configured for single-round amplification. To perform two rounds of amplification, the Ambion MessageAmp II aRNA Amplification Kit (Part no. AM1751) must be used in conjunction with the MessageAmp II-Biotin *Enhanced* Kit. See section “[Procedure Overview](#)” on page 9 for more information.

cDNA synthesis and IVT reagents

Do *not* store reagents in a frost-free freezer.

Amount		Component	Storage
20	μL	T7 Oligo(dT) Primer [†]	-20°C
22	μL	ArrayScript™ Reverse Transcriptase	-20°C
22	μL	RNase Inhibitor	-20°C
42	μL	10X First Strand Buffer	-20°C
170	μL	dNTP Mix	-20°C
210	μL	10X Second Strand Buffer	-20°C
42	μL	DNA Polymerase	-20°C
22	μL	RNase H	-20°C
84	μL	T7 Enzyme Mix	-20°C
84	μL	T7 10X Reaction Buffer	-20°C
255	μL	Biotin-NTP Mix	-20°C
10	μL	Control RNA (1 mg/mL HeLa total RNA)	-20°C
1.75	mL	Nuclease-free Water	any temp [‡]

[†] The T7 Oligo(dT) Primer is available separately from Life Technologies (Part no. AM5710).

[‡] Store the Nuclease-free Water at -20°C, 4°C, or room temperature.

Some reagents may form a precipitate when stored at -20°C. If a precipitate is visible, redissolve it by warming the solution to room temperature with gentle mixing.



cDNA and aRNA purification, and fragmentation reagents

Do *not* store reagents in a frost-free freezer.

Amount		Component	Storage
30	mL	Wash Buffer (Add 24 mL 100% ethanol before use)	4°C or room temp
7	mL	cDNA Binding Buffer	room temp [†]
9	mL	aRNA Binding Buffer	room temp
1	mL	5X Array Fragmentation Buffer	room temp
20		aRNA Filter Cartridges	room temp
40		aRNA Collection Tubes	room temp
20		cDNA Filter Cartridges + Tubes	room temp
20		cDNA Elution Tubes	room temp
10	mL	Nuclease-free Water	any temp [‡]

[†] The cDNA Binding Buffer may form a precipitate if stored colder than room temperature. If a precipitate is visible, redissolve it by warming the solution to 37°C for up to 10 min and vortexing vigorously. Cool to room temperature before use.

[‡] Store the Nuclease-free Water at -20°C, 4°C, or room temperature.

aRNA Amplification Procedure

Important Parameters for Successful Amplification

Input RNA quantity and IVT reaction incubation time

Consider both the amount of sample RNA you have and the amount of aRNA needed for your analysis when planning MessageAmp II-Biotin kit experiments. These factors will influence how much input RNA to use, whether one or two rounds of amplification should be done, and how long to incubate the IVT reaction.

Accurate quantitation

For experiments where the aRNA yield from different samples will be compared, it is *essential* to accurately quantify the input RNA used in the MessageAmp II-Biotin kit procedure. The NanoDrop® 1000A Spectrophotometer is recommended for rapid, accurate quantitation of nucleic acids; however, any reliable RNA quantitation method, such as traditional spectrophotometry or RiboGreen assay, can be used.

Recommended minimum and maximum amounts of input RNA

The following table shows the mass of total RNA that can be used in the MessageAmp II-Biotin *Enhanced* procedure. Alternatively, 10–100 ng of poly(A) selected RNA can be used in the procedure. The RNA volume must be ≤10 µL.


Recommended		Minimum		Maximum	
1000	ng	50	ng	5000	ng

Determining input RNA amount and IVT reaction incubation time

The procedure can accommodate a wide range of input RNA amounts, but for reproducible and comparable results, use a fixed amount of input RNA for all experiments. Tailor both the amount of input RNA and the amplification procedure to produce the amount of aRNA needed for your microarray hybridizations. For instance, Affymetrix GeneChip arrays require 10–15 µg of aRNA for each hybridization, but other commercial and core facility arrays may require slightly more or less aRNA.

Figure 2 shows aRNA yield data from amplification of increasing amounts of input RNA from different sources. The experiment included either a 4 hr or a 14 hr IVT reaction incubation. With most samples, amplification of 50 ng of total RNA for 14 hr produced enough aRNA for a microarray experiment. When amplifying small RNA samples (e.g., ~250 ng or less), incubating the IVT reaction for 14 hr will maximize the amount of aRNA produced.

Figure 2 aRNA Yield from Nine Different Tissue/Cell Types. Four different amounts of total RNA from nine different sources were amplified using the MessageAmp™ II-Biotin *Enhanced* Kit. Reactions were performed in triplicate using either a 4 h or a 14 h IVT reaction time. The average aRNA yields from the triplicate reactions are shown as bar graphs and in tabular format. This data is useful for determining both the amount of total RNA needed to obtain enough labeled aRNA for an array hybridization (typically ~10 µg) and the length of the IVT incubation that should be used. Note that there is a ~3-fold difference in aRNA yield between some of the samples. With most RNA sources, 50–100 ng of input total RNA amplified with the MessageAmp II-Biotin *Enhanced* Kit using a 14 hour IVT incubation will yield enough labeled aRNA for a microarray hybridization.



Input RNA	HeLa		Human Brain		Human Liver		U937 (Human)		K562 (Human)	
	4 h	14 h	4 h	14 h	4 h	14 h	4 h	14 h	4 h	14 h
50 ng	4.8	17.9	3.4	13.3	1.5	15.0	5	25.6	2.6	20.1
100 ng	9.4	45.1	7.4	31.0	2.6	37.3	12.5	47.1	4.8	45.4
500 ng	57.1	147.3	29.8	108.1	23.3	102.3	58.7	139.5	34.1	130.6
1 µg	105.8	165.6	66.0	136.5	54.7	147.3	109.1	163.2	68.8	158.9

Input RNA	Mouse Brain		Mouse Liver		Rat Thymus		Rat Liver		Overall Average	
	4 h	14 h	4 h	14 h	4 h	14 h	4 h	14 h	4 h	14 h
50 ng	1.8	7.7	1.2	8.3	1.8	16.5	2.1	18.4	2.7	15.9
100 ng	3.6	15.7	3.6	16.6	5.1	26.0	3.9	32.1	5.9	32.9
500 ng	19.1	57.7	24.6	75.0	21.1	112.2	26.9	108.6	32.7	109.0
1 µg	45.4	115.1	47.1	122.2	43.2	140.6	53.8	131.9	66.0	142.4

RNA purity


The quality of the RNA is the single most important factor affecting how efficiently an RNA sample will be amplified using the MessageAmp II-Biotin *Enhanced* Kit. RNA samples should be free of contaminating proteins, DNA, and other cellular material as well as phenol, ethanol, and salts associated with RNA isolation procedures. Impurities can lower the efficiency of reverse transcription and subsequently reduce the level of amplification. An effective measure of RNA purity is the ratio of absorbance readings at 260 and 280 nm. The ratio of A_{260} to A_{280} values should fall in the range of 1.7–2.1. RNA must be suspended in high quality water or TE (10 mM Tris-HCl, 1 mM EDTA) or THE RNA Storage Solution (Part no. AM7000, AM7001).

RNA integrity

The integrity of the RNA sample, or the proportion that is full-length, is another important component of RNA quality. Reverse transcription of partially degraded mRNAs will typically generate relatively short cDNAs that potentially lack portions of the coding region. RNA integrity can be evaluated by microfluidic analysis using the Agilent 2100 bioanalyzer and Caliper RNA LabChip® Kits. Primarily full-length RNA will exhibit a ratio of 28S to 18S rRNA bands that approaches 2:1. Using a bioanalyzer, the RIN (RNA Integrity Number) can be calculated to further evaluate RNA integrity.

Denaturing agarose gel electrophoresis and nucleic acid staining can also be used to separate and visualize the major rRNA species. When the RNA resolves into discrete rRNA bands (i.e., no significant smearing below each band), with the 28S rRNA band appearing approximately twice as intense as the 18S rRNA band, then the mRNA in the sample is likely to be mostly full-length. The primary drawback to gel electrophoresis is that it requires microgram amounts of RNA.

Figure 3 Bioanalyzer Images of Total RNA Preparations. These electropherograms (from the Agilent 2100 bioanalyzer) show RNA samples with decreasing integrity that are all of sufficient quality to use as input for the MessageAmp II-Biotin *Enhanced* Kit. The trace labeled “Intact Total RNA” represents the ideal for a bioanalyzer trace of total RNA. Notice how the ribosomal RNA peaks are at a ratio of about 2 (28S:18S) in this sample; this represents ultrahigh quality RNA in terms of integrity. Most RNA samples will more closely resemble the center trace where there are nearly equal amounts of 28S and 18S rRNA. The trace on the right shows a fairly typical human RNA prep with rRNA peaks that are lower than in the other two traces and where lower molecular weight degradation products become apparent. RNA samples with suboptimal integrity can yield meaningful array analysis results if the data are subjected to rigorous statistical analysis [Schoor et al. 2003].



Reaction incubation times should be precise and consistent

The incubation times for most of the enzymatic reactions in the procedure were optimized in conjunction with the kit reagents to ensure the maximum yield of nucleic acid product in each step—adhere to them closely. An exception is the IVT reaction, where a range of 4–14 hr incubation time is acceptable (step 2, on page 21). Refer to the table in step 2, on page 21 to help determine what incubation time to use. Although differences in IVT incubation time among samples has had very little, if any, effect on array results in our hands, we recommend using uniform IVT incubation times if aRNA yield from different samples will be compared or if you want to have equal amplification of different samples—this will provide the most reproducible amplification and array analysis.

Master mixes

We strongly recommend preparing master mixes for the MessageAmp II-Biotin *Enhanced* procedure. This approach reduces the effects of pipetting error, saves time, and improves reproducibility. Using master mixes is especially important when aRNA yield from different samples will be compared.



Thorough mixing is very important for reproducibility

Below are specific instructions for mixing kit reagents, master mixes, and individual reactions. For maximum reproducibility and aRNA yield, follow these instructions closely.

Mix each kit component after thawing.

Mix enzyme solutions by *gently* flicking the tube a few times before adding them to reactions. Thaw frozen reagents completely *at room temperature* (i.e., primers, nucleotides, and 10X buffers), then mix thoroughly by vortexing, and keep on ice before use.

Mix master mixes by gentle vortexing.

After assembling master mixes, *gently* vortex to make a homogenous mixture without inactivating the enzyme(s).

Mix individual reactions by pipetting and flicking the tube.

After adding master mixes or other reagents to individual reactions, pipet up and down 2–3 times to rinse reagents from the pipet tip. Then flick the tube with your finger 3–4 times to mix thoroughly, and finish by centrifuging briefly to collect the reaction at the bottom of the tube.

Thermal cycler recommended

The MessageAmp II-Biotin *Enhanced* procedure is very sensitive to temperature; variable or inaccurate incubation temperatures can limit aRNA synthesis. It is also very important that condensation does not form in the reaction tubes during any of the incubations. Condensation changes the composition of reaction mixtures, which can greatly reduce yield.

- ***A thermal cycler with a temperature adjustable heated lid is recommended.***

A calibrated thermal cycler, with a temperature-adjustable heated lid, is recommended, for the greatest temperature control and stability during MessageAmp II-Biotin *Enhanced* reaction incubations. Allow the thermal cycler to equilibrate to the required temperature before placing the tubes in the block for incubation. Follow the recommended settings for the lid temperatures. Too high a lid setting may inhibit the reaction; too low a setting may cause condensation.

Note: Even if you use a hybridization oven or incubator for most of the MessageAmp II-Biotin *Enhanced* reactions, a thermal cycler is strongly recommended for the 16°C second-strand synthesis reaction incubation (step 2, on page 18). Turn off the heated lid if it cannot be adjusted to match the 16°C block temperature.

If your thermal cycler does not have a temperature-adjustable lid, or a thermal cycler is unavailable, calibrated hybridization ovens or incubators (at constant temperature) may also be used. Preheat incubators so that the correct temperature has stabilized before reactions are placed in the incubator. To avoid any potential influence on the reaction temperature from the tube holder, let tube holders equilibrate in the incubator for sufficient time, or use a tube holder that doesn't touch the sides and bottoms of the tubes—for example a floating tube support.

- Heat blocks or water baths are not recommended for MessageAmp II-Biotin *Enhanced* reaction incubations.

Maintaining consistency

Procedural consistency is very important for amplification experiments. Consider implementing a detailed procedural plan that will be used by everyone in the lab to maintain consistency. This type of plan will minimize variation due to subtle procedural differences that can influence RNA amplification and may complicate gene expression studies. The plan should include basic information such as the method of RNA isolation, the amount of RNA to use in the procedure, and how long to incubate the IVT reaction. It should also address specifics that are not often included in protocols such as which tubes, tube racks, and incubators to use for each step in the process. Finally, develop a consistent work flow. For example standardize stopping points in the method. The idea is to standardize all of the variables discussed in this section of the Protocol and carefully follow all the steps in order to maximize amplification consistency among samples.

Tubes: use 0.5 mL RNase-free nonstick tubes

If a 60-well thermal cycler with temperature-adjustable lid is available, it is most convenient to conduct the MessageAmp II-Biotin *Enhanced* procedure in 0.5 mL nonstick tubes (for example, Part no. AM12350). These can be thin-wall (PCR) tubes or ordinary-weight nonstick tubes. 0.5 mL tubes are large enough to accommodate the cDNA Binding Buffer without having to transfer reactions to a larger tube. Their small size and nonstick properties also keep the reaction components at the bottom of the tube.

If your thermal cycler is equipped with a standard 96-well block, 0.2 mL non-stick tubes can be used.

Prepare the Wash Buffer

Add 24 mL 100% ethanol (ACS grade or better) to the bottle labeled Wash Buffer. Mix well and mark the label to indicate that the ethanol was added.

Reverse Transcription to Synthesize First Strand cDNA

Incubator needed

Thermal cycler programmed as follows:

Temp	Time	Cycles
42°C (50°C lid)	2 hr	1
4°C	hold	

It is important to follow the lid temperature recommendation; higher lid temperatures may inhibit reverse transcription, while lower temperatures may cause condensation, resulting in changes to the reaction composition.

1. Adjust RNA sample volume to 10 µL, if necessary
 - a. Place a maximum volume of 10 µL of total RNA (1000 ng recommended) or poly(A) selected RNA (typically 10–100 ng) into a nonstick, sterile, RNase-free, 0.5 mL tube. RNA must be in high quality water or TE. (See [“Recommended minimum and maximum amounts of input RNA”](#) on page 12.)



Note: If your experiment will include RNA spikes (e.g., Ambion ArrayControl RNA Spikes, Part no. AM1780, or Affymetrix GeneChip Poly-A Control Kit, Cat #900433), add them to samples at this step.

- b. If necessary, add Nuclease-free Water to a final volume of 10 µL, vortex briefly to mix, then centrifuge to collect the mixture at the bottom of the tube.

2. Add 10 µL of *Reverse Transcription Master Mix* and place at 42°C

- a. At room temperature, prepare *Reverse Transcription Master Mix* in a nuclease-free tube. Assemble enough to synthesize first strand cDNA from all the RNA samples in the experiment, including ≤5% overage to cover pipetting error.

At room temperature, assemble the Reverse Transcription Master Mix in the order shown:

Reverse Transcription Master Mix (for a single 20 µL reaction)		
Amount		Component
1	µL	Nuclease-free Water
1	µL	T7 Oligo(dT) Primer
2	µL	10X First Strand Buffer
4	µL	dNTP Mix
1	µL	RNase Inhibitor
1	µL	ArrayScript

- b. *Gently* vortex the tube to make a homogenous mixture without inactivating the enzyme, then centrifuge briefly (~5 sec) to collect the Reverse Transcription Master Mix at the bottom of the tube and place on ice.
- c. Transfer 10 µL of Reverse Transcription Master Mix to each RNA sample. Mix thoroughly by pipetting up and down 2–3 times, then flicking the tube 3–4 times, and centrifuge briefly to collect the reaction in the bottom of the tube.
- d. Place the samples in the thermal cycler, and start the run.

3. Incubate for 2 hr at 42°C

Incubate the reactions for 2 hr at 42°C, then centrifuge the tubes briefly (~5 sec) to collect the contents at the bottom of the tubes.

Place the tubes on ice and immediately proceed to second strand cDNA synthesis (below).

IMPORTANT! Proceed immediately to the next step.

Second Strand cDNA Synthesis

Incubator needed

Thermal cycler programmed as follows:

Temp	Time	Cycles
16°C (heat-disabled lid, or no lid)	2 hr	1
4°C	hold	

Disable the heated lid of the thermocycler if the lid temperature cannot be set in the range 16°C to room temperature. If the lid is too hot, inhibition of the reaction may occur.

1. Add 80 µL *Second Strand Master Mix* to each sample
 - a. On ice, prepare a *Second Strand Master Mix* in a nuclease-free tube in the order listed below. Assemble enough to synthesize second strand cDNA from all the samples in the experiment, including ≤5% overage to cover pipetting error. Assemble the Second Strand Master Mix on ice in the order shown:

Second Strand Master Mix (for a single 100 µL reaction)		
Amount		Component
63	µL	Nuclease-free Water
10	µL	10X Second Strand Buffer
4	µL	dNTP Mix
2	µL	DNA Polymerase
1	µL	RNase H

- b. Mix well by gently vortexing. Centrifuge briefly (~5 sec) to collect the Second Strand Master Mix at the bottom of the tube and place on ice.
 - c. Transfer 80 µL of Second Strand Master Mix to each sample. Mix thoroughly by pipetting up and down 2–3 times, then flicking the tube 3–4 times, and centrifuge briefly to collect the reaction in the bottom of the tube.
 - d. Place the tubes in a 16°C thermal cycler. It is important to cool the thermal cycler block to 16°C before adding the reaction tubes because subjecting the reactions to temperatures >16°C will compromise aRNA yield.
2. Incubate for 2 hr at 16°C

Incubate 2 hr in a 16°C thermal cycler. If the lid temperature cannot be adjusted to match the 16°C block temperature, cover the reactions with the heated lid turned off, or if the lid cannot be turned off—do not cover the tubes with it. (Do not use a water bath or a heat block in a 4°C refrigerator for this incubation because the temperature will fluctuate too much.)
3. Place reactions on ice briefly or freeze immediately

After the 2 hr incubation at 16°C, place the reactions on ice and proceed to section [“cDNA Purification” on page 19](#) (below), or immediately freeze reactions at –20°C. Do not leave the reactions on ice for more than 1 hr.

STOPPING POINT. This is a potential overnight stopping point (at –20°C), but it is better to complete the cDNA purification (next section) before stopping.



cDNA Purification

IMPORTANT! All centrifugations in this purification procedure should be done at 10,000 x g (typically ~10,000 rpm) at room temperature. cDNA Filter Cartridges should not be subjected to RCFs over 16,000 x g because the force could cause mechanical damage and/or may deposit glass filter fiber in the eluate.

1. Preheat Nuclease-free Water to 55°C

Before beginning the cDNA purification, preheat at least 24 µL per sample of Nuclease-free Water to 55°C.

IMPORTANT! Preheat the Nuclease-free Water to a maximum of 55°C; temperatures above 58°C can partially denature the cDNA, compromising final aRNA yield.

2. Add 250 µL cDNA Binding Buffer to each sample

IMPORTANT! Check the cDNA Binding Buffer for precipitation before using it. If a precipitate is visible, redissolve it by warming the solution to 37°C for up to 10 min and vortexing vigorously. Cool to room temperature before use.

Add 250 µL of cDNA Binding Buffer to each sample, and mix thoroughly by pipetting up and down 2–3 times, then flicking the tube 3–4 times. Follow up with a quick spin to collect the reaction in the bottom of the tube. Proceed quickly to the next step.

3. Pass the mixture through a cDNA Filter Cartridge

Check that the cDNA Filter Cartridge is firmly seated in its wash tube (supplied).

- Pipet the cDNA sample/cDNA Binding Buffer (from previous step) onto the center of the cDNA Filter Cartridge.
- Centrifuge for ~1 min at 10,000 x g, or until the mixture is through the filter.
- Discard the flow-through and replace the cDNA Filter Cartridge in the wash tube.

IMPORTANT! Make sure that the ethanol has been added to the bottle of Wash Buffer before using it in this step.

4. Wash with 500 µL Wash Buffer

- Apply 500 µL Wash Buffer to each cDNA Filter Cartridge.
- Centrifuge for ~1 min at 10,000 x g, or until all the Wash Buffer is through the filter.
- Discard the flow-through and spin the cDNA Filter Cartridge for an additional minute to remove trace amounts of Wash Buffer.
- Transfer cDNA Filter Cartridge to a cDNA Elution Tube.

5. Elute cDNA with 22 µL 55°C Nuclease-free Water

It is important to use Nuclease-free Water that is at 50–55°C for the cDNA elution. Colder water will be less efficient at eluting the cDNA, and hotter water (≥58°C) may result in reduced aRNA yield.

- a. Apply 22 µL of preheated Nuclease-free Water (55°C) to the center of the filter in the cDNA Filter Cartridge.
- b. Leave at room temperature for 2 min and then centrifuge for 1 min at 10,000 x g, or until all the Nuclease-free Water is through the filter. The double-stranded cDNA will now be in the eluate (~20 µL).

STOPPING POINT. The purified cDNA can be stored overnight at –20°C at this point if desired. Transfer the cDNA to a lidded, non-stick, nuclease-free tube for storage.

In Vitro Transcription to Synthesize Biotin- labeled aRNA

Incubator needed

Thermal cycler programmed as follows:

Temp	Time	Cycles
37°C (default lid; 100–105°C)	4–14 hr; see step 2.	1
4°C	hold	

It is important to use a thermal cycler with a heated lid (100–105°C), to maintain uniform temperature, and because it is extremely important that condensation does not form inside the tubes; this would change the reagent concentrations and reduce yield.

1. Add 20 µL of *IVT Master Mix* to each sample, and mix
 - a. At room temperature, prepare an *IVT Master Mix* by adding the following reagents to a nuclease-free microcentrifuge tube in the order listed below. Assemble enough for all the samples in the experiment, including ≤5% overage to cover pipetting error.

IMPORTANT! If two rounds of amplification will be done, the MessageAmp II Kit must be used for this first round transcription to make unmodified aRNA, (not biotin-labeled aRNA).

Assemble the IVT Master Mix at room temperature in the order shown:

Amount		Component
20	µL	double-stranded cDNA (from step 5. on page 19)
IVT Master Mix for a single reaction		
12	µL	Biotin-NTP Mix
4	µL	T7 10X Reaction Buffer
4	µL	T7 Enzyme Mix

- b. Mix well by gently vortexing. Centrifuge briefly (~5 sec) to collect the IVT Master Mix at the bottom of the tube and place on ice.
- c. Transfer 20 µL of IVT Master Mix to each sample. Mix thoroughly by pipetting up and down 2–3 times, then flicking the tube 3–4 times, and centrifuge briefly to collect the reaction in the bottom of the tube.
- d. Once assembled, place the tubes in the thermal cycler and start the run.



2. Incubate 4–14 hr at 37°C

The minimum recommended incubation time is 4 hr; the maximum is 14 hr. (The reactions can be held post-IVT at 4°C for up to 48 hr, for convenience.)

Use the following table as a guide to determine how long to continue your IVT reaction. There are more data and a detailed discussion of the length of the IVT incubation in section “Input RNA quantity and IVT reaction incubation time” on page 12.

aRNA Needed		Input Total RNA		IVT Incubation	
10–100	µg	1–5	µg	4	hr
1–10	µg	50 ng–1	µg	8	hr
0.1–1	µg	≤50	ng	14	hr

3. Add 60 µL Nuclease-free Water to each sample

Stop the reaction by adding 60 µL Nuclease-free Water to each aRNA sample to bring the final volume to 100 µL. Mix thoroughly by gentle vortexing.

Proceed to the aRNA purification step (below) or store at –20°C.

STOPPING POINT. The aRNA can be stored overnight at –20°C at this point if desired.

aRNA Purification

Incubator needed: heat block set at 55°C.

This purification removes enzymes, salts, and unincorporated nucleotides from the aRNA. At the end of the purification the aRNA is eluted from the filter with Nuclease-free Water.

IMPORTANT! All centrifugations in this purification procedure should be done at 10,000 × g (typically ~10,000 rpm) at room temperature. aRNA Filter Cartridges should not be subjected to RCFs over 16,000 × g because the force could cause mechanical damage and/or may deposit glass filter fiber in the eluate.

1. Preheat Nuclease-free Water and assemble Filter Cartridges and Tubes
 - Preheat a minimum of 200 µL per sample of Nuclease-free Water to 55°C.
 - For each sample, place an aRNA Filter Cartridge into an aRNA Collection Tube and set aside for use in step 4.
2. Add 350 µL aRNA Binding Buffer to each sample
 - a. Check to make sure that each IVT reaction was brought to 100 µL with Nuclease-free Water.
 - b. Add 350 µL of aRNA Binding Buffer to each aRNA sample. Proceed to the next step immediately.
3. Add 250 µL 100% ethanol and pipet 3 times to mix

Add 250 µL of ACS grade 100% ethanol to each aRNA sample, and mix by pipetting the mixture up and down 3 times. *Do NOT vortex to mix and do NOT centrifuge.*

IMPORTANT! It is crucial to follow these mixing instructions exactly, and to proceed quickly to the next step.

Proceed *immediately* to the next step as soon as you have mixed the ethanol into each sample. Any delay in proceeding could result in loss of aRNA because once the ethanol is added, the aRNA will be in a semiprecipitated state.

4. Pass samples through an aRNA Filter Cartridge
 - a. Pipet each sample mixture from the previous step onto the center of the filter in the aRNA Filter Cartridge/Collection Tube assembly.
 - b. Centrifuge for ~1 min at 10,000 x g. Continue until the mixture has passed through the filter.
 - c. Discard the flow-through and replace the aRNA Filter Cartridge back into the aRNA Collection Tube.
5. Wash with 650 µL Wash Buffer
 - a. Apply 650 µL Wash Buffer to each aRNA Filter Cartridge.
 - b. Centrifuge for ~1 min at 10,000 x g, or until all the Wash Buffer is through the filter.
 - c. Discard the flow-through and spin the aRNA Filter Cartridge for an additional ~1 min to remove trace amounts of Wash Buffer.
 - d. Transfer Filter Cartridge(s) to a fresh aRNA Collection Tube.
6. Elute aRNA with 200 µL 55°C Nuclease-free Water
 - a. To the center of the filter, add 200 µL Nuclease-free Water (preheated to 55°C).
 - b. Incubate the samples in the 55°C heat block for 10 min (recommended). Alternatively, incubate at room temperature for 2 min. This results in ~80% recovery of the aRNA.
 - c. Centrifuge for ~1.5 min at 10,000 x g, or until the Nuclease-free Water is through the filter.
 - d. The aRNA will now be in the aRNA Collection Tube in ~200 µL of Nuclease-free Water.
7. Store aRNA at -80°C

Store aRNA at -80°C for up to 1 year, and minimize repeated freeze-thawing. To prevent multiple freeze-thaw events, split samples into 5–20 µg aliquots for microarray labeling and hybridizations.
8. (Optional) Concentrate the purified aRNA

If necessary, concentrate the aRNA by vacuum centrifugation or by precipitation with ammonium acetate (NH₄OAc)/ethanol.

(Optional) Concentrate by vacuum centrifugation

If the heater on the vacuum centrifuge has different settings, use medium or low. Check the progress of drying every 5–10 min, and remove the sample from the concentrator when it reaches the desired volume.



(Optional) Precipitate with 5 M NH₄OAc and ethanol

- a. Add 0.1 volume of 5 M NH₄OAc to the purified aRNA (20 µL if the aRNA was eluted in 200 µL Nuclease-free water).
- b. Add 2.5 volumes of 100% ethanol (550 µL if the aRNA was eluted in 200 µL). Mix well and incubate at -20°C for 30 min.
- c. Microcentrifuge at top speed for 15 min at 4°C or room temperature. Carefully remove and discard the supernatant.
- d. Wash the pellet with 500 µL 70% cold ethanol, centrifuge again, and remove the 70% ethanol.
- e. To remove the last traces of ethanol, quickly respin the tube, and aspirate any residual fluid with a fine-tipped pipette or syringe needle.
- f. Air dry the pellet.
- g. Resuspend the aRNA pellet using the desired solution and volume.

Evaluation and Fragmentation of aRNA

aRNA Quantitation and Expected Yield

1. Assessing aRNA yield by UV absorbance

The concentration of an aRNA solution can be determined by measuring its absorbance at 260 nm using a spectrophotometer.

- Use a NanoDrop® spectrophotometer and measure 1.5 µL of the RNA sample directly.
- With a traditional spectrophotometer, dilute an aliquot of the RNA 1:50 to 1:100 in TE (10 mM Tris-HCl pH 8, 1 mM EDTA), and read the absorbance. (Be sure to zero the spectrophotometer with the TE used for sample dilution.) The buffer used for dilution need not be RNase-free, since slight degradation of the RNA will not significantly affect its absorbance.

Find the concentration in µg/mL by multiplying the A₂₆₀ by the dilution factor and the extinction coefficient. (1 A₂₆₀ = 40 µg RNA/mL):

$$A_{260} \times \text{dilution factor} \times 40 = \mu\text{g RNA/mL}$$

2. Assessing aRNA yield with the RiboGreen® assay

If a fluorometer or a fluorescence microplate reader is available, Molecular Probes' RiboGreen fluorescence-based assay for RNA quantitation is a convenient and sensitive way to measure RNA concentration. Follow the manufacturer's instructions for using RiboGreen assay.

3. Expected yield

The aRNA yield will depend on the amount and quality of poly(A) RNA in the input total RNA. Since the proportion of poly(A) RNA in total RNA is affected by influences such as health of the organism and the organ from which it is isolated, aRNA yield from equal amounts of total RNA may vary considerably (see Figure 2 on page 13 for empirical aRNA yield data obtained using this kit).


Analysis of aRNA Size

The size distribution of aRNA can be evaluated using an Agilent 2100 bioanalyzer with Caliper's LabChip technology, or by conventional denaturing agarose gel analysis. The bioanalyzer can provide a fast and accurate size distribution profile of aRNA samples, but aRNA yield should be determined by UV absorbance or RiboGreen assay analysis. To analyze aRNA size using a bioanalyzer, follow the manufacturer's instructions for running the assay using purified aRNA (from step 6. on page 22).

Expected aRNA size

- Agilent bioanalyzer analysis
The expected aRNA profile is a distribution of sizes 250–5500 nt with most of the aRNA 1000–1500 nt (Figure 4). To compare bioanalyzer profiles of different aRNA samples, be sure to load equal mass amounts to get an accurate comparison.
- Denaturing agarose gel analysis
Amplified aRNA should appear as a smear from 250 to 5000 nt. The average size of aRNA should be approximately 1400 nt; the average size of unmodified aRNA should be ~1150 nt.

Figure 4 Biotin-labeled and Unlabeled aRNA Made with the MessageAmp II-Biotin *Enhanced* Kit and the MessageAmp II Kits. HeLa RNA samples (1 µg) were amplified with either the MessageAmp II-Biotin *Enhanced* Kit to produce biotin-labeled aRNA, or the MessageAmp II Kit to produce unlabeled aRNA. The IVT reactions were carried out for 4 hr. The data were generated by running a sample of the reactions on an Agilent bioanalyzer using an RNA LabChip Kit.



Fragmentation of Biotinylated aRNA for GeneChip® Arrays

Most procedures for array hybridization begin with a sample fragmentation step prior to hybridization. The 5X Array Fragmentation Buffer supplied with the MessageAmp II-Biotin *Enhanced* Kit is designed for perfect compatibility with the Affymetrix GeneChip® array platform. You can use the 5X Array Fragmentation Buffer following either the procedure in the Affymetrix GeneChip Expression Analysis Technical Manual, or the equivalent procedure included here. The composition of the 5X Array Fragmentation Buffer is shown below.

Concentration		Component
200	mM	Tris Acetate, pH 8.2
500	mM	Potassium Acetate
150	mM	Magnesium Acetate

1. Assemble the aRNA fragmentation mixture



The aRNA fragmentation reaction employs metal-induced hydrolysis to fragment input aRNA.

aRNA quantity and reaction volume

Refer to the Affymetrix GeneChip Expression Analysis Technical Manual for:


- the amount of aRNA needed for hybridization with your GeneChip array format, and
- the recommended fragmentation reaction volume. This will be based on the volume of the hybridization mixture for your GeneChip array format.

40 µL rxn†		30 µL rxn†		Component
1–32	µL	1–24	µL	5–20 µg aRNA (depending on GeneChip array format)
8	µL	6	µL	5X Array Fragmentation Buffer [1X final]
to 40	µL	to 30	µL	Nuclease-free Water

† Use the reaction volume recommended for your GeneChip array platform.

- Incubate at 94°C for 35 min, then place in ice
 - Incubate the fragmentation reaction at 94°C for 35 min.
 - Place the reaction on ice immediately after the incubation.
- Optional: Evaluate a sample of the reaction on a bioanalyzer
 Analyze the size of the fragmentation reaction products by running a sample of the reaction on an Agilent bioanalyzer using an RNA LabChip Kit. Figure 5 shows a typical result of such analysis. (Follow the manufacturer’s instructions for this analysis.)
 The reaction should produce a distribution of 35–200 nt aRNA fragments with a peak at approximately 105 nt.
- Use fragmented aRNA immediately or store frozen
 Use the fragmented aRNA immediately in a GeneChip array hybridization following the instructions in the Affymetrix GeneChip Expression Analysis Technical Manual, or store undiluted, fragmented aRNA at –20°C for 1–3 days or at –80°C for long-term storage.

Figure 5 Fragmented aRNA. Agilent bioanalyzer analysis of a 1 µL sample of a 30 µL fragmentation reaction containing 10 µg of aRNA.



Troubleshooting

Positive Control Reaction

Control RNA amplification instructions

To establish if the kit is working properly, Control RNA consisting of 1 mg/mL HeLa cell total RNA is provided.

1. Use 1 μ L of the Control RNA in a MessageAmp II-Biotin *Enhanced* kit reaction; follow the procedure starting at step “Reverse Transcription to Synthesize First Strand cDNA” on page 16.
2. At step 2. on page 21, use a **14 hr** incubation for the IVT reaction.
3. Continue with the procedure through section “aRNA Purification” on page 21.

Analysis of the positive control amplification reaction

- After completing the aRNA purification, measure the A_{260} of the reaction product as described in section “Assessing aRNA yield by UV absorbance” on page 23.

The positive control reaction should produce $\geq 80 \mu\text{g}$ of aRNA.

Be aware that often the positive control reaction cannot be compared to experimental reactions, because many experimental amplification experiments will use less than the 1 μg of input RNA used in the positive control reaction, and the aRNA yield will be proportionately lower. Also the Control RNA is of exceptional quality and purity, ensuring that it will amplify with extremely high efficiency.

- Also run a 2 μg aliquot of the reaction products on a denaturing agarose gel or analyze 100–200 ng on a bioanalyzer; *the average size of the aRNA should be $\geq 1 \text{ kb}$.*

Factors that Affect Both the Positive Control and Experimental Samples

If the positive control reaction yield or amplification product size does not meet expectations, consider the following possible causes and troubleshooting suggestions. These suggestions also apply to problems with amplification of experimental RNA.

Incubation temperature(s) were incorrect

The incubation temperatures are critical for effective RNA amplification.

- Check the temperatures of all incubators used in the procedure with a calibrated thermometer.
- If a thermal cycler is used for incubation, check the accuracy of the adjustable temperature lid. If the lid temperature cannot be adjusted to match the desired reaction temperature, use the lid with the heat turned off, or do not use the lid to cover the reaction vessel(s).

Condensation formed in the tube during the reaction incubation(s)

Condensation occurs when the cap of the reaction vessel is cooler (e.g., room temperature) than the bottom of the tube. As little as 1–2 μL of condensate in an IVT reaction tube throws off the concentrations of the nucleotides and magnesium, which are crucial for good yield.

If you see condensation, spin the tube briefly and mix the reaction gently. Move the tube(s) to an incubator where condensation does not occur or is minimized.



**Troubleshooting
 Low Yield and
 Small Average
 aRNA Size**

Nuclease-contaminated tubes, tips, or equipment

Using pipettes, tubes, or equipment that are contaminated with nucleases can cleave the RNA or DNA being generated at each step in the procedure. This will reduce the size of the aRNA products and decrease aRNA yield. Both RNases and DNases can be removed from surfaces using Ambion RNaseZap® RNase Decontamination Solution (Part no. AM9780, AM9786).

Absorbance readings were inaccurate

Confirm that your spectrophotometer is accurate by measuring the absorbance of an RNA or DNA sample of known concentration. Alternatively, assess the aRNA concentration by fractionating on an agarose gel adjacent to an RNA sample whose concentration is known. Comparing the ethidium bromide staining of the aRNA and control samples can approximate the concentration of the aRNA.

Consider the following troubleshooting suggestions if the positive control reaction produced the expected results, but amplification of your experimental samples results in less or smaller (average <500 nt) aRNA than expected.

Impure RNA samples

RNA samples with significant amounts of contaminating DNA, protein, phenol, ethanol, or salts are reverse transcribed poorly and subsequently generate less aRNA than pure RNA samples. Phenol extract and ethanol precipitate your RNA, or use the Ambion MEGAclean™ Kit (Part no. AM1908) to further purify your RNA before reverse transcription.

Lower than expected input RNA concentration

Take another A₂₆₀ reading of your RNA sample or try using more RNA in the aRNA amplification procedure.

RNA integrity is compromised

RNA that is partially degraded generates cDNA that is relatively short. This will reduce the average size of the aRNA population and subsequently reduce the yield of aRNA. You can assess the integrity of an RNA sample by determining the size of the 18S and 28S rRNA bands and the relative abundance of 28S to 18S rRNA (See “[RNA integrity](#)” on page 14 for more information).

The mRNA content of your total RNA sample is lower than expected

Different RNA samples contain different amounts of mRNA. In healthy cells, mRNA constitutes 1–3% of total cellular RNA. The actual amount of mRNA depends on the cell type and the physiological state of the sample. When calculating the amount of amplification, the starting mass of mRNA in a total RNA prep should always be considered within a range of 10–30 ng per µg of total RNA (assuming good RNA quality). Most total RNA samples can be amplified up to 1000 fold producing 10–30 µg of aRNA from 1 µg of total RNA.





Materials Not Provided with the Kit

Materials Not Provided with the Kit

Lab equipment and supplies

- 100% Ethanol (to prepare the Wash Buffer)
- Thermal cycler with adjustable-temperature heated-lid (recommended), hybridization oven, or constant temperature incubators set at 70°C, 42°C, 37°C, and 16°C (See “[Thermal cycler recommended](#)” on page 15 for more information.)
- Heat block set at 55°C, for preheating the water for cDNA and aRNA purification
- Vacuum centrifuge concentrator
- Vortex mixer
- Microcentrifuge
- Non-stick RNase-free 0.5 mL microcentrifuge tubes (Part no. AM12350)
- RNase-free pipettors and tips, positive-displacement type recommended to increase the accuracy and precision of reaction inputs
- (Optional) RNA controls for microarrays analysis, such as Array Control™ RNA Spikes from Life Technologies (Part no. AM1780) or the GeneChip® Eukaryotic Poly-A RNA Control Kit from Affymetrix® (Cat #900433)
- (Optional) Non-stick RNase-free tubes for storage of cDNA (e.g., AM 12450)

Optional materials and equipment for RNA analysis

- Spectrophotometer—such as the NanoDrop® ND-1000 or ND-8000 UV-Vis Spectrophotometer. Follow the manufacturer’s instructions.
- (Optional) Agilent bioanalyzer and RNA LabChip Kits
- (Optional) Reagents and apparatus for preparation and electrophoresis of agarose gels
- (Optional) RiboGreen® RNA Quantitation Assay and Kit (Molecular Probes Inc.)

Related Products Available from Life Technologies

MessageAmp™ aRNA Amplification Kits	A full line of Ambion® MessageAmp Kits is available, tailored for different array analysis applications. The MessageAmp II Kit offers maximum flexibility; samples can be amplified using either single- or double-round amplification, and the reagent cocktails are configured to accommodate modification. For arrays requiring biotin-labeled samples, the MessageAmp Premier and MessageAmp III RNA Amplification Kits are available. For preparation of fluorescently-labeled samples, we recommend the Amino Allyl MessageAmp II Kits, which are available with and without Cy®3 and Cy®5 dyes. Bacterial RNA can be amplified using the MessageAmp II Bacteria RNA Amplification Kit. The MessageAmp II-96 and Amino Allyl MessageAmp II-96 aRNA Amplification Kits are offered for high-throughput applications.
RNA Isolation Kits	Family of kits for isolation of total or poly(A) RNA. Included in the product line are kits using classical GITC and acidic phenol, one-step disruption/denaturation, phenol-free glass fiber filter or magnetic bead binding, and combination kits. See the catalog or website (www4.appliedbiosystems.com).
GLOBINclear™ Whole Blood Globin Reduction Kits AM1980, AM1981	The GLOBINclear Whole Blood Globin Reduction Kits employ a novel, non-enzymatic technology to remove >95% of the globin mRNA from whole blood total RNA samples. The resulting mRNA is a superior template for RNA amplification and synthesis of labeled cDNA for array analysis. Kits are available for treatment of human or mouse/rat whole blood total RNA.
ArrayControl™ RNA Spikes AM1780	The ArrayControl RNA Spikes are a set of eight control RNA transcripts designed for the normalization and validation of glass microarray experiments. The RNA Spikes range in size from 750 to 2000 bases, and each transcript has a 30-base 3' poly(A) tail. The precisely quantitated RNA Spikes are designed to be added to your RNA sample before labeling, to serve as internal controls for sample labeling and hybridization efficiency.
Biotin-11-UTP and Biotin-16-UTP AM8450, AM8451, AM8452, AM8453	Biotinylated UTPs are ideal for use as substrates in vitro transcription reactions, and can be utilized by a variety of RNA polymerases, including T7, T3, and SP6 RNA polymerases. Biotinylated RNA can be used in place of radioactively labeled RNA in many applications with detection via one of a variety of streptavidin-based methods.
RNA Fragmentation Reagents AM8740	Amplified RNA is commonly fragmented prior to hybridization on oligonucleotide microarrays to improve the hybridization kinetics and signal produced on oligonucleotide microarrays. Ambion® RNA Fragmentation Reagents include a 10X Fragmentation Reagent and a Stop Solution.
Amino Allyl cDNA Labeling Kit AM1705	The Amino Allyl cDNA Labeling Kit generates cDNA for secondary fluorescent dye labeling to be used for glass array analysis. It includes all the reagents, except the amine-reactive labeling moiety (e.g. cyanine dyes) for 2-step labeling of cDNA. The reaction produces more labeled cDNA, more efficiently than direct dye incorporation.
RETROscript® Kit AM1710	First-strand cDNA synthesis kit. The RETROscript® Kit can be used to incorporate dye-modified nucleotides into cDNA using aRNA prepared with the MessageAmp™ II Kit as a template.
5-(3-aminoallyl)-dUTP AM8439	This 50 µM solution of amino allyl-modified dUTP can be used with the RETROscript® Kit (Part no. AM1710) to synthesize amine-reactive cDNA from aRNA. The amine-reactive cDNA can then be postlabeled with any amine-reactive label moiety.

General Safety



WARNING! GENERAL SAFETY. Using this product in a manner not specified in the user documentation may result in personal injury or damage to the instrument or device. Ensure that anyone using this product has received instructions in general safety practices for laboratories and the safety information provided in this document.

- Before using an instrument or device, read and understand the safety information provided in the user documentation provided by the manufacturer of the instrument or device.
- Before handling chemicals, read and understand all applicable Safety Data Sheets (SDSs) and use appropriate personal protective equipment (gloves, gowns, eye protection, etc). To obtain SDSs, see the “Documentation and Support” section in this document.

Chemical safety



WARNING! GENERAL CHEMICAL HANDLING. To minimize hazards, ensure laboratory personnel read and practice the general safety guidelines for chemical usage, storage, and waste provided below, and consult the relevant SDS for specific precautions and instructions:

Read and understand the Safety Data Sheets (SDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials. To obtain SDSs, see the “Documentation and Support” section in this document.

- Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing).
- Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (for example, fume hood).
- Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer's cleanup procedures as recommended in the SDS.
- Handle chemical wastes in a fume hood.
- Ensure use of primary and secondary waste containers. (A primary waste container holds the immediate waste. A secondary container contains spills or leaks from the primary container. Both containers must be compatible with the waste material and meet federal, state, and local requirements for container storage.).
- After emptying a waste container, seal it with the cap provided.
- Characterize (by analysis if necessary) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- Ensure that the waste is stored, transferred, transported, and disposed of



- according to all local, state/provincial, and/or national regulations.
- **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.
-

Biological hazard safety



WARNING! Potential Biohazard. Depending on the samples used on this instrument, the surface may be considered a biohazard. Use appropriate decontamination methods when working with biohazards.



WARNING! BIOHAZARD. Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. Follow all applicable local, state/provincial, and/or national regulations. Wear appropriate protective equipment, which includes but is not limited to: protective eyewear, face shield, clothing/lab coat, and gloves. All work should be conducted in properly equipped facilities using the appropriate safety equipment (for example, physical containment devices). Individuals should be trained according to applicable regulatory and company/institution requirements before working with potentially infectious materials. Read and follow the applicable guidelines and/or regulatory requirements in the following:

In the U.S.:

- U.S. Department of Health and Human Services guidelines published in Biosafety in Microbiological and Biomedical Laboratories found at: www.cdc.gov/biosafety
- Occupational Safety and Health Standards, Bloodborne Pathogens (29 CFR§1910.1030), found at: www.access.gpo.gov/nara/cfr/waisidx_01/29cfr1910a_01.html
- Your company's/institution's Biosafety Program protocols for working with/handling potentially infectious materials.
- Additional information about biohazard guidelines is available at: www.cdc.gov

In the EU:

- Check local guidelines and legislation on biohazard and biosafety precaution and refer to the best practices published in the World Health Organization (WHO) Laboratory Biosafety Manual, third edition, found at: www.who.int/csr/resources/publications/biosafety/WHO_CDS_CSR_LYO_2004_11/en/
-



Bibliography

- Dorris DR, Ramakrishnan R, Trakas D, Dudzik F, Belval R, Zhao C, Nguyen A, Domanus M, Mazumder A (2002) A highly reproducible, linear, and automated sample preparation method for DNA microarrays. *Genome Res* 12(6): 976–984.
- Feldman AJ, Costouros NG, Wang E, Qian M, Marincola FM, Alexander HR, and Libutti SK (2002) Advantages of mRNA amplification for microarray analysis. *Biotechniques* 33(4):906–914.
- Kacharmina JE, Crino PB, Eberwine J (1999) Preparation of cDNA from single cells and subcellular regions. *Methods Enzymol*, 303:3–18.
- Li Y, Li T, Liu S, Qiu M, Han Z, Jiang Z, Li R, Ying K, Xie Y, Mao Y (2004) Systematic comparison of the fidelity of aRNA, mRNA and T-RNA on gene expression profiling using cDNA microarray. *J Biotechnology* 107(1):19–28.
- Pabon C, Modrusan Z, Ruvolo MV, Coleman IM, Daniel S, Yue H, Arnold LJ Jr. (2001) Optimized T7 amplification system for microarray analysis. *Biotechniques* 31(4):874–879.
- Polacek DC, Passerini AG, Shi C, Francesco NM, Manduchi E, Grant GR, Powell S, Bischof H, Winkler H, Stoeckert CJ Jr, Davies PF (2003) Fidelity of enhanced sensitivity of differential transcription profiles following linear amplification of nanogram amounts of endothelial mRNA. *Physiol Genomics* 13:147–156.
- Schoor O, Weinschenk T, Hennenlotter J, Corvin S, Stenzl A, Rammensee H-G, and Stevanovic S (2003) Moderate degradation does not preclude microarray analysis of small amounts of RNA. *Biotechniques* 35:1192–1201.
- Van Gelder RN, von Xastrow ME, Yool A, Dement DC, Barchas JD, Eberwine JH (1990) Amplified RNA synthesized from limited quantities of heterogeneous cDNA. *Proc Natl Acad Sci USA*, 87:1663–1667.
- Yue H, Eastman PS, Wang BB, Minor J, Doctolero MH, Nuttall RL, Stack R, Becker JW, Montgomery JR, Vainer M, Johnston R (2001) An evaluation of the performance of cDNA microarrays for detecting changes in global mRNA expression. *Nucleic Acids Res* 29(8):E41–1.

Documentation and Support

Obtaining SDSs

Safety Data Sheets (SDSs) are available from www.lifetechnologies.com/sds

Note: For the SDSs of chemicals not distributed by Life Technologies, contact the chemical manufacturer.

Obtaining support

For the latest services and support information for all locations, go to:

www.lifetechnologies.com

At the website, you can:

- Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities
- Search through frequently asked questions (FAQs)
- Submit a question directly to Technical Support
- Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents
- Obtain information about customer training
- Download software updates and patches



Headquarters

5791 Van Allen Way | Carlsbad, CA 92008 USA | Phone +1 760 603 7200 | Toll Free in USA 800 955 6288

For support visit www.lifetechnologies.com/support

www.lifetechnologies.com

