

WesternDot™ 625 Western Blot Kits

Catalog nos. W10132 and W10142

Table 1. Contents and storage information.

| Material | Amount | Concentration | Storage* | Stability |
|---|--------|---|--|--|
| WesternDot™ blocking buffer (Component A) | 320 mL | 1X solution | <ul style="list-style-type: none"> • 2–6°C • DO NOT FREEZE | When stored as directed these kits are stable for 6 months |
| Wash buffer (Component B, 0.5 M Tris-HCl, 1.5 M NaCl, 0.5% Tween 20, pH 7.4) | 320 mL | 10X concentrate | | |
| Biotin-XX goat anti-mouse IgG (H+L) (Component C, for W10132, Goat Anti-Mouse Kit) OR Biotin-XX goat anti-rabbit IgG (H+L) (Component C, for W10142, Goat Anti-Rabbit Kit) | 80 µL | 2 mg/mL | | |
| Qdot® 625 streptavidin conjugate, (Component D) | 80 µL | 1 µM solution in 1 M betaine, 50 mM borate pH 8.3 with 0.05% sodium azide | | |
| WesternDot™ staining dish (Component E) | 2 each | Not applicable | | |

Number of reactions: Sufficient material is supplied for 20 mini blots (8 cm × 8 cm) using the protocol described below.

Approximate fluorescence excitation/emission maxima: Excitation: 254 to 488 nm; emission maximum: 625 nm, see Figure 2.

Introduction

The WesternDot™ 625 Western Blot Kits combine the unique properties of Qdot® 625 nanocrystals with the high affinity streptavidin-biotin interaction to allow simple, highly sensitive detection of proteins on western blots. The WesternDot™ Kits contain optimized, ready-to-use or ready-to-dilute reagents for sensitive immunodetection of proteins immobilized on nitrocellulose (NC) or polyvinylidene difluoride (PVDF) membranes.

Incorporating a standard western blotting protocol, the detection step relies on a biotinylated secondary antibody, goat anti-mouse (Cat. no. W10132) or goat anti-rabbit (Cat. no. W10142) followed by the key component, a Qdot® 625 streptavidin conjugate (Figure 1). The extremely high extinction coefficient of the Qdot® 625 nanocrystal in the UV and blue wavelengths combined with high quantum yield and an orange/red emission (Figure 2) allow for sub-nanogram sensitivity of protein detection using standard UV or blue-light based detection systems. The fluorescent signal is compatible with the commonly used modes of fluorescent detection of DNA or protein gels and does not require specialized emission filters (Table 2).

The WesternDot™ kits include a specially formulated WesternDot™ blocking buffer, Wash buffer, and two staining dishes that hold sufficient solution volumes recommended for a standard mini-blot (8 cm × 8 cm). The kits include sufficient reagents for 20 mini-blots.

Qdot® 625 Streptavidin Conjugate

The Qdot® 625 streptavidin conjugate is made from a nanometer-scale crystal of a semiconductor material (CdSe), which is coated with an additional semiconductor shell (ZnS) to improve the optical properties of the material. The core-shell material is further coated with a polymer shell that allows the materials to be conjugated to biological molecules and retain their optical properties. This polymer shell is directly coupled to streptavidin through an active ester coupling reaction yielding a material with streptavidin covalently attached on the surface (typically 5–10 streptavidins/Qdot® conjugate), which results in Qdot® streptavidin conjugates with high specific biological activity. The Qdot® 625 streptavidin conjugate is the size of a large macromolecule or protein (~15–20 nm). The Qdot® 625 streptavidin conjugate has a narrow, symmetric emission spectrum with the emission maximum near 625 nm. For details on the optical properties, spectral characteristics, and biological activity of Qdot® 625 streptavidin conjugate, refer to the manual supplied with Qdot® 625 streptavidin conjugate available separately from Invitrogen (Cat. no. A10196).




Figure 1. Schematic for the WesternDot™ kit.




Figure 2. Spectral properties of the Qdot® 625 streptavidin conjugate.

Before You Begin

Materials Required but Not Provided

- Blotted PVDF or nitrocellulose transfer membranes containing applied antigen (**Note:** We have obtained best results with low background fluorescence PVDF membranes such as Immobilon®-FL PVDF transfer membrane from Millipore.)
- Primary antibody to detect applied antigen (visit www.invitrogen.com/antibodies for a list of antibodies available from Invitrogen)
- Ultra pure water
- Orbital shaker or rocking platform
- Forceps for manipulating blotted membranes
- UV transilluminator/ethidium bromide filter/Polaroid camera or Imaging system (see Table 2 for a list of compatible instruments)

Preparing 1X Wash Buffer

Prepare 1X Wash Buffer for each blot by combining 16 mL 10X Wash Buffer (Component B) with 144 mL ultra pure water. The resulting 1X Wash Buffer contains 50 mM Tris-HCl, pH 7.4, 150 mM NaCl, 0.05% Tween®-20.

Disposal of Qdot® Conjugate

The Qdot® conjugate contains cadmium and selenium in an inorganic, crystalline form. Dispose of the material in compliance with all applicable local, state, and federal regulations. For more information on the composition of these materials, consult the Material Safety Data Sheet.

Caution

The WesternDot™ blocking buffer (Component A) contains 0.1% sodium azide. Wear gloves when handling this solution. Harmful if swallowed. In case of accident or if you feel unwell, seek medical advice immediately. Dispose of this material and its container as hazardous waste.

Experimental Protocols

General Guidelines

To obtain the best results with WesternDot™ Kits:

- Perform all steps on an orbital (rotary) shaker or rocking platform (rotating at 1 revolution/second). Be sure that the solutions cover the membrane and move freely over and around it.
- Process the blot face up (but still submerged) for good fluid flow.
- Avoid touching the working surface of the membrane, even with gloves.
- Work quickly when changing solutions as PVDF membranes dry quickly. If the membrane dries, re-wet the membrane with methanol and rinse with water before proceeding.
- Add solutions to the trays slowly, at the membrane edge, to avoid bubbles forming under the membrane. Decant from the same corner of the dish to ensure complete removal of previous solutions.

Preparing the Membrane

Immediately after transferring the proteins to a nitrocellulose or PVDF membrane, wash the membrane twice for 5 minutes with 20 mL water to remove gel and transfer buffer components. Proceed to **Immunodetection Protocol**.

If you are using water-washed and dried nitrocellulose membranes, proceed to **Immuno-detection Protocol**. If you are using water-washed and dried PVDF membranes, re-wet the membrane in 100% methanol briefly for 30 seconds, followed by two 20 mL water washes for 5 minutes, and then proceed to **Immunodetection Protocol**.

Immunodetection Protocol Perform all steps at room temperature with continuous shaking using an orbital shaker (set at ~ 1 revolution/second) or with continuous rocking for the indicated times.

- 1.1 Place the membrane in WesternDot™ staining dish (Component E) containing 8 mL WesternDot™ blocking buffer (Component A). Incubate for 60 minutes on a shaker set at 1 revolution/second.

Note: Blocking step may be performed overnight at 4°C.

- 1.2 Prepare 8 mL of the diluted primary antibody in 1X wash buffer at the appropriate concentration.

Note: Primary antibody concentration depends on the manufacturer's recommendation, and is typically 0.1 to 2.0 µg/mL.

- 1.3 Decant the blocking buffer and add diluted primary antibody solution from Step 1.2 to the blot. Incubate for 60 minutes on a shaker. Decant solution.

Note: Primary antibody incubation may be as short as 30 minutes or overnight at 4°C.

- 1.4 Wash the membrane for 5 minutes with 15 mL 1X Wash buffer, then decant. Repeat wash step 2 more times.

- 1.5 During the last wash, prepare the secondary antibody solution by diluting 4 µL Biotin-XX-Goat anti-rabbit or Biotin-XX-Goat anti-mouse (Component C) into 8 mL 1X Wash buffer.

- 1.6 Decant wash buffer and incubate the membrane in 8 mL secondary antibody solution from Step 1.5 for 30–60 minutes, then decant.

- 1.7 Wash the membrane for 5 minutes with 15 mL 1X Wash buffer, then decant. Repeat wash step 2 more times.

- 1.8 During the last wash, prepare Qdot® 625 streptavidin conjugate by diluting 4 µL Qdot® 625 streptavidin conjugate (Component D) into 8 mL WesternDot™ blocking buffer (Component A).

- 1.9 Decant wash buffer and incubate the membrane in 8 mL Qdot® 625 streptavidin conjugate solution from Step 1.8 for 30–60 minutes, then decant.

Note: Incubation with Qdot® 625 streptavidin conjugate may be as short as 15 minutes or as long as several hours to overnight. For a strong antibody and abundant antigen, signal may be visible within 5 to 10 minutes when the blot is examined *in situ* with a hand-held UV lamp.

- 1.10 Wash the membrane for 5 minutes with 15 mL 1X Wash buffer, then decant. Repeat wash step 2 more times. Perform a final wash with 15 to 20 mL ultra pure water.

Note: The blot may be stored in ultra pure water, 1X Wash buffer, TBS (Tris-buffered saline), or PBS (phosphate buffered saline) overnight with minimal signal loss.

Imaging the Blot PVDF or Immobilon®-FL membranes

For best results with PVDF membranes, dry the membrane and image with epi-illumination.

- With epi-illumination, image the membrane **wet or dry** using exposures ranging from 2 seconds to 2 minutes.
- With trans-illumination using a UV trans-illuminator, image the membrane **wet** using exposures ranging from a few milliseconds to several seconds.

Nitrocellulose membranes

- With UV epi-illumination, image the membrane **wet or dry** using exposures ranging from 2 seconds to 2 minutes.
- With trans-illumination, membranes can be imaged **wet or dry**. There will be 3- to 5-fold lower sensitivity and increased background with trans-illumination as compared to epi-illumination.

Instrumentation

See Table 2 for a list of compatible instruments.

Excitation

Deep blue or UV light in trans- or epi-illumination.

Lasers with 473 nm or 488 nm excitation are acceptable.

Emission filters

Filters that are qualified for ethidium bromide, SYPRO® Ruby, SYPRO® Red, Qdot® 605, Qdot® 625, or Qdot® 655 are suitable. Filters centered around 625 nm provide greater signal. Most orange-to-red emission filters are also acceptable.

Signal stability

Under the reaction conditions, Qdot® nanocrystals are chemically stable and photophysically stable. Repeated images may be taken over the course of several hours and the signal lasts for several days to months, though there may be a loss in sensitivity at the lower end. If a membrane begins to dry, re-wet the membrane as described in **Preparing the Membrane**. If using dried membranes, image with epi-illumination for better results.

Examples of Expected Results

Examples of expected results using the WesternDot™ 625 Western Blot Kits are shown in Figures 3 and 4.




Figure 3. Example of results using the WesternDot™ 625 Kit with a PVDF membrane: Total proteins (2-fold dilution series ranging from 10 µg to ~10 ng) from Jurkat cell extract were analyzed on a NuPAGE® Novex® 4–12% Bis-Tris gel and then transferred to an Immobilon®-FL PVDF membrane. Immunodetection of glyceraldehyde-3-phosphate dehydrogenase (GAPDH), an endogenous “housekeeping” protein in the Jurkat cell extract was performed with WesternDot™ 625 Goat Anti-Mouse Western Blot Kit (Cat. no. W10132) using a mouse monoclonal anti-GAPDH antibody (Invitrogen Cat. no. 39-8600) at 1 µg/mL as described in this manual. The wet membrane was imaged using an Alpha Innotech HD2 instrument with a SYPRO® Red emission filter (620 nm ±40 nm) and excitation at 302 nm trans-illumination with an exposure time of 300 milliseconds.




Figure 4. Example of results using the WesternDot™ 625 Kit with a nitrocellulose membrane: Jurkat cell extract was spiked with biotin-XX- transferrin (Invitrogen Cat. no. T23363) and rabbit IgG, and was analyzed on a NuPAGE® Novex® 4–12% Bis-Tris gel. Proteins were transferred to nitrocellulose membrane using the iBlot™ Gel Transfer System. Immunodetection of biotin-XX- transferrin and rabbit IgG in the Jurkat cell extract was performed with WesternDot™ 625 Goat Anti-Rabbit Western Blot Kit (Cat. no. W10142) as described in this manual. The membrane was imaged using an Alpha Innotech HD2 instrument with a SYPRO® Red emission filter (620 nm ±40 nm) and excitation at 254 nm epi-illumination with an exposure time of 45 seconds.

The biotin-XX-transferrin is directly detected using the Qdot® 625 streptavidin conjugate (Component D) while biotinylated secondary antibody (Component C) serves as a primary antibody for rabbit IgG detection and illustrates the utility of the Qdot® 625 streptavidin conjugate for detection via any biotinylated primary antibody. Biotinylated proteins that are intrinsic to the extract (see lane 3) will be detected with the WesternDot™ Kits.

Lane 1: 1 ng biotin-XX- transferrin; Lane 2: 10 ng rabbit IgG; Lane 3: 5 mg Jurkat whole-cell extract; Lane 4: 5 mg Jurkat whole-cell extract spiked with 10 ng rabbit IgG and 1 ng biotin-XX- transferrin; Lanes 5-11: 2-fold dilution series of the sample applied to Lane 4.

Table 2. Imaging platforms recommended for detecting signal using WesternDot™ 625 Western Blot Kits.

| Manufacturer | Instrument |
|-------------------------|---|
| Alpha Innotech | AlphaDigiDoc PRO |
| | Alphamager® EP |
| | Alphamager® HP |
| | FluorChem® HD2 |
| BioRad | ChemiDoc XRS |
| | GelDoc XR |
| | PharosFX Systems |
| | VersaDoc MP 4000 |
| Fuji Film Life Sciences | FLA-7000 |
| | FLA-8000 |
| | FLA-9000 |
| | LAS-3000 |
| | LAS-4000 |
| GE Healthcare | Typhoon™ |
| | Storm™ with Blue LED |
| Invitrogen | Safe Imager™ |
| Kodak | Gel Logic |
| | Image Station 4000 MM and PRO |
| | Image Station 4000R and PRO |
| Syngene | DigiGenius |
| | G: Box fluorescence and chemiluminescence |
| | InGenius |
| | U: Genius |
| UVP | Benchtop UV Transilluminators |
| | BioDocIt™ Imaging System |
| | DigiDocIt |
| | EC3™ Imaging System |
| | FirstLight® UV illuminator |
| | Visi-Blue™ Transilluminator |

Product List Current prices may be obtained from our website or from our Customer Service Department.

| Cat. no. | Product Name | Unit Size |
|-------------------------|--|-------------|
| W10132 | WesternDot™ 625 Goat Anti-Mouse Western Blot Kit *20 minigel blots* | 1 kit |
| W10142 | WesternDot™ 625 Goat Anti-Rabbit Western Blot Kit *20 minigel blots* | 1 kit |
| <i>Related Products</i> | | |
| A10196 | Qdot® 625 streptavidin conjugate *1 µM solution* | 200 µL |
| B2763 | Biotin-XX goat anti-mouse IgG (H+L) *2 mg/mL* | 0.5 mL |
| B2770 | Biotin-XX goat anti-rabbit IgG (H+L) *2 mg/mL* | 0.5 mL |
| IB1001 | iBlot™ Gel Transfer Device | 1 each |
| IB3010-02 | iBlot™ Transfer Stack, Mini (Nitrocellulose) | 10 sets/box |
| LC5800 | Novex® Sharp Pre-stained Protein Standard | 2 × 250 µL |
| LC5925 | SeeBlue® Plus2 Pre-Stained Standard | 500 µL |
| S37102 | Safe Imager™ Blue Light Transilluminator | 1 each |

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probetech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

For country-specific contact information,
visit www.invitrogen.com.

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2008, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.