


RNAsecure™ Reagent

Catalog Number AM7005, AM7006

Pub. No. 4386507 Rev. B

Cat. no.	Concentration	Volume	Storage conditions
AM7005	25X	1 mL	Store at -20°C.
AM7006	25X	10 mL	

 **WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from www.lifetechnologies.com/support.

Product description

RNAsecure™ Reagent is designed to use instead of DEPC in molecular biology reagents, in solutions that are incompatible with DEPC (such as Tris or MOPS buffers), in solutions that cannot be autoclaved, and in enzymatic reactions prior to adding the enzyme. No inhibition of the following enzymatic reactions was seen when treated with RNAsecure™ Reagent prior to enzyme addition: T7, T3, and SP6 RNA Polymerase, in vitro transcription (37°C); MMLV and AMV, reverse transcription (42°C); and SuperTaq™ Polymerase, PCR (94°C).

Treat solutions with RNAsecure™ Reagent

1. Dilute RNAsecure™ Reagent to 1X final concentration in the buffer or solution to be treated, and mix well.
2. Incubate the mixture at 60°C for 10 minutes, then cool to room temperature.
The buffer/solution is then ready to use for RNA isolation and analysis. Alternatively, the treated solution can be stored at room temperature, 4°C or -20°C until use.
3. To eliminate any RNase contamination that may occur after initial treatment, reheat the treated solution to 60°C for 10–20 minutes.

Note: RNAsecure™ Reagent is active above 45°C, and may inactivate some enzymes in reactions incubated above this temperature. Therefore it may be critical to treat reaction buffers with RNAsecure™ Reagent prior to addition of some enzymes and then to perform the enzymatic incubations below 45°C. This does not apply to SuperTaq™ Polymerase.

Quality control

Nonspecific endonuclease activity: A sample is incubated with supercoiled plasmid DNA and analyzed by agarose gel electrophoresis.

Exonuclease activity: A sample is incubated with labeled double-stranded DNA, followed by PAGE analysis.

Functional testing: RNAsecure™ Reagent is tested functionally for its ability to inactivate RNase in solution.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

For Research Use Only. Not for use in diagnostic procedures.

The information in this guide is subject to change without notice.

Disclaimer: **DISCLAIMER:** LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

© 2013 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation and/or its affiliate(s) or their respective owners.

For support visit lifetechnologies.com/support or email techsupport@lifetech.com

lifetechnologies.com

16 May 2013

