

PROTOCOL

ambion® | RNA
by *life* technologies™

MessageAmp™ II aRNA Amplification Kit

life
technologies™

For Research Use Only. Not intended for any animal or human therapeutic or diagnostic use.

Information in this document is subject to change without notice.

APPLIED BIOSYSTEMS DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL APPLIED BIOSYSTEMS BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT APPLIED BIOSYSTEMS IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

LIMITED USE LABEL LICENSE: RESEARCH USE ONLY

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

TRADEMARKS

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

Cy and CyDye are trademarks of GE Healthcare.

© 2011 Life Technologies Corporation. All rights reserved.

Part Number 1751M Rev. E
02/2011

MessageAmp™ II aRNA Amplification Kit

(Part Number AM1751)

Protocol

I.	Introduction	1
	A. Product Description and Background	
	B. Procedure Overview	
	C. Materials Provided with the Kit and Storage Conditions	
	D. Materials Not Provided with the Kit	
	E. Related Products	
II.	aRNA Amplification Procedure	8
	A. Important Parameters for Successful Amplification	
	B. Prepare the Wash Buffer	
	C. Reverse Transcription to Synthesize First Strand cDNA	
	D. Second Strand cDNA Synthesis	
	E. cDNA Purification	
	F. In Vitro Transcription to Synthesize aRNA	
	G. aRNA Purification	
III.	Assessing aRNA Yield and Quality	22
	A. aRNA Quantitation and Expected Yield	
	B. Analysis of aRNA Size	
IV.	(Optional) Second Round Amplification	24
	A. Synthesis of First Strand cDNA (Second Round)	
	B. Synthesis of Second Strand cDNA (Second Round)	
V.	Troubleshooting	28
	A. Positive Control Reaction	
	B. Factors that Affect Both the Positive Control and Experimental Samples	
	C. Troubleshooting Low Yield and Small Average aRNA Size	
	D. aRNA is Not Efficiently Reverse Transcribed	
VI.	Appendix	31
	A. References	
	B. Quality Control	
	C. Safety Information	
VII.	Documentation and Support	34
	A. Obtaining SDSs	
	B. Obtaining support	

I. Introduction

**IMPORTANT**

Before using this product, read and understand the “Safety Information” in the appendix in this document.

A. Product Description and Background

The MessageAmp™ II aRNA Amplification Kit is based on the RNA amplification protocol developed in the Eberwine laboratory (Van Gelder et al. 1990). The procedure consists of reverse transcription with an oligo(dT) primer bearing a T7 promoter using ArrayScript™ reverse transcriptase (RT), engineered to produce higher yields of first-strand cDNA than wild-type enzymes. ArrayScript RT catalyzes the synthesis of virtually full-length cDNA, which is the best way to ensure production of reproducible microarray samples. The cDNA then undergoes second-strand synthesis and cleanup to become a template for in vitro transcription (IVT) with T7 RNA polymerase. To maximize aRNA yield, Ambion MEGAscript® IVT technology is used to generate hundreds to thousands of antisense RNA copies of each mRNA in a sample. (In this protocol the antisense amplified RNA is referred to as aRNA; it is also commonly called cRNA.) The IVT can be configured to synthesize either biotin-labeled RNA, or unlabeled aRNA that can subsequently be labeled by reverse transcription (for example with fluorescently labeled dNTPs).


The resulting aRNA is suitable for use on most commercially available microarray gene expression systems. RNA samples of limited amounts (0.1–100 ng) can be put through two rounds of amplification if desired. This strategy makes the production of microarray samples from picogram amounts of total RNA entirely possible (Luo et al. 1999).

Benefits of RNA amplification

RNA amplification was originally developed as a method to expand very small RNA samples to produce enough material for array hybridization (Yue et al. 2001). Several groups have conducted studies to determine whether amplification of RNA introduces bias and they report that any bias is minimal (Li et al. 2004, Feldman et al. 2002 and Polacek et al. 2003). Additionally, among the benefits of amplification is a more reproducible expression profile from a wide range of RNA inputs. Some researchers conclude that amplification actually improves the reliability of array results regardless of whether it is needed for sample expansion (Feldman et al. 2002 and Polacek et al. 2003). As a result, RNA amplification has become the standard method for preparing RNA samples for array analysis (Kacharmina et al. 1999, Pabon et al. 2001).

MessageAmp™ II aRNA Amplification Kit

Figure 1. MessageAmp™ II aRNA Amplification Procedure



B. Procedure Overview


The MessageAmp II aRNA amplification procedure is depicted in Figure 1.

- **Reverse Transcription to Synthesize First Strand cDNA** is primed with the T7 Oligo(dT) Primer to synthesize cDNA containing a T7 promoter sequence.
- **Second Strand cDNA Synthesis** converts the single-stranded cDNA into a double-stranded DNA (dsDNA) template for transcription. The reaction employs DNA Polymerase and RNase H to simultaneously degrade the RNA and synthesize second-strand cDNA.
- **cDNA Purification** removes RNA, primers, enzymes, and salts that would inhibit in vitro transcription.
- **In Vitro Transcription to Synthesize aRNA** generates multiple copies of aRNA from the double-stranded cDNA templates; this is the amplification step.
- **aRNA Purification** removes unincorporated NTPs, salts, enzymes, and inorganic phosphate to improve the stability of the aRNA and to facilitate subsequent enzymatic reactions.

Optional second round of amplification

Additional amplification of an RNA sample can be achieved by subjecting the aRNA to a second round of amplification (see Figure 2). The reagents and methodology used in the first and second rounds of amplification are only slightly different (see section IV starting on page 24.)

Figure 2. Second Round Amplification



The MessageAmp II advantage

Each step in the MessageAmp II aRNA Amplification Kit amplification procedure has been streamlined and optimized. The first-strand cDNA synthesis reaction employs ArrayScript reverse transcriptase to ensure that every cDNA bears a T7 promoter at its 5' end and that even very limited amounts of mRNA are fully converted to full-length cDNA. The second-strand cDNA synthesis reaction is designed for the efficient synthesis of full-length, double-stranded cDNAs and the complete conversion of single-stranded cDNA into double-stranded transcription templates. The cDNA purification procedure not only removes enzymes, salts, and unincorporated dNTPs, but also efficiently removes RNA from the cDNA sample. This eliminates the heating or enzymatic digestion step commonly used in other procedures to degrade RNA (especially ribosomal RNA). The IVT reaction features MEGAscript technology to maximize transcriptional amplification and yield of aRNA. It is optimized to ensure efficient transcription of limited amounts of input DNA and synthesis of long transcripts.

The NTPs for IVT are provided separately so that modified nucleotides (e.g., biotinylated UTP, or cyanine 3/cyanine 5 CTP and UTP) can be readily incorporated into aRNA. The simple, rapid aRNA purification procedure prepares the aRNA for downstream applications (reverse transcription or post-labeling reactions).

C. Materials Provided with the Kit and Storage Conditions

The MessageAmp II aRNA Amplification Kit includes reagents for single-round amplification of 20 samples or two-round amplification of 10 samples.

cDNA synthesis and IVT reagents

Do *not* store reagents in a frost-free freezer.

Amount	Component	Storage
60 µL	T7 Oligo(dT) Primer	-20°C
22 µL	ArrayScript™ Reverse Transcriptase	-20°C
22 µL	RNase Inhibitor	-20°C
42 µL	10X First Strand Buffer	-20°C
170 µL	dNTP Mix	-20°C
210 µL	10X Second Strand Buffer	-20°C
42 µL	DNA Polymerase	-20°C
22 µL	RNase H	-20°C
84 µL	T7 Enzyme Mix	-20°C
84 µL	T7 10X Reaction Buffer	-20°C
84 µL	T7 ATP	-20°C

Amount	Component	Storage
84 µL	T7 CTP	-20°C
84 µL	T7 GTP	-20°C
84 µL	T7 UTP	-20°C
40 µL	Second Round Primers	-20°C
10 µL	Control RNA (1 mg/mL HeLa total RNA)	-20°C
1.75 mL	Nuclease-free Water	any temp*

* Store the Nuclease-free Water at -20°C, 4°C, or room temperature.

Some reagents may form a precipitate when stored at -20°C. If a precipitate is visible, redissolve it by warming the solution to room temperature with gentle mixing.

cDNA and aRNA purification

Do **not** store reagents in a frost-free freezer.

Amount	Component	Storage
30 mL	Wash Buffer (Add 24 mL 100% ethanol before use)	4°C or room temp
7 mL	cDNA Binding Buffer	room temp*
9 mL	aRNA Binding Buffer	room temp
20	aRNA Filter Cartridges	room temp
40	aRNA Collection Tubes	room temp
20	cDNA Filter Cartridges + Tubes	room temp
20	cDNA Elution Tubes	room temp
10 mL	Nuclease-free Water	any temp†

* The cDNA Binding Buffer may form a precipitate if stored colder than room temperature. If a precipitate is visible, redissolve it by warming the solution to 37°C for up to 10 min and vortexing vigorously. Cool to room temperature before use.

† Store the Nuclease-free Water at -20°C, 4°C, or room temperature.

D. Materials Not Provided with the Kit

Lab equipment and supplies

- 100% Ethanol (to prepare the Wash Buffer)
- Thermal cycler with adjustable-temperature heated-lid (recommended), hybridization oven, or constant temperature incubators set at 70°C, 42°C, 37°C, and 16°C (See [Thermal cycler recommended](#) on page 11 for more information.)
- Heat block set at 55°C, for preheating the water for cDNA and aRNA purification
- Vacuum centrifuge concentrator
- Vortex mixer
- Microcentrifuge
- Non-stick RNase-free 0.5 mL microcentrifuge tubes (P/N AM12350)

- RNase-free pipettors and tips, positive-displacement type recommended to increase the accuracy and precision of reaction inputs
- (Optional) RNA controls for microarrays analysis, such as Array Control™ RNA Spikes from Life Technologies (P/N AM1780) or the GeneChip® Eukaryotic Poly-A RNA Control Kit from Affymetrix® (Cat #900433)
- (Optional) Non-stick RNase-free tubes for storage of cDNA (e.g., AM12450)
- Spectrophotometer—such as the NanoDrop ND-1000 or ND-8000 UV-Vis Spectrophotometer. Follow the manufacturer's instructions.
- (Optional) Agilent bioanalyzer and RNA LabChip Kits
- (Optional) Reagents and apparatus for preparation and electrophoresis of agarose gels
- (Optional) Quant-iT™ RiboGreen® RNA Assay Kit from Invitrogen (R11490) for use with a fluorescence microplate reader, standard spectrofluorometer, or filter fluorometer

Optional materials and equipment for RNA analysis

(Optional) Biotin-labeled UTP

Biotin-labeled UTP can be added to the in vitro transcription reaction to synthesize biotin-labeled aRNA. Biotin-11-UTP (P/N AM8451, 75 mM) is recommended because it gives good incorporation, has minimal effect on aRNA recovery during purification, and results in high signal on most commercial microarrays.

E. Related Products

MessageAmp™ aRNA Amplification Kits

See web or print catalog for P/Ns

A full line of Ambion® MessageAmp Kits is available, tailored for different array analysis applications. The MessageAmp II Kit offers maximum flexibility; samples can be amplified using either single- or double-round amplification, and the reagent cocktails are configured to accommodate modification. For arrays requiring biotin-labeled samples, the MessageAmp Premier and MessageAmp III RNA Amplification Kits are available. For preparation of fluorescently-labeled samples, we recommend the Amino Allyl MessageAmp II Kits, which are available with and without Cy™3 and Cy5 dyes. Bacterial RNA can be amplified using the MessageAmp II Bacteria RNA Amplification Kit. The MessageAmp II-96 and Amino Allyl MessageAmp II-96 aRNA Amplification Kits are offered for high-throughput applications.

FirstChoice® Total RNA

High-quality total RNA from a variety of human, mouse and rat tissues and from human cell lines. DNA is removed with a stringent DNase treatment, and the purity and integrity of these RNAs are verified by Agilent bioanalyzer evaluation, denaturing agarose gel electrophoresis, or Northern analysis. FirstChoice Total RNA is prepared by methods that quantitatively recover small RNAs (miRNA, siRNA, and snRNA). FirstChoice Total RNAs are ready for use in any application that requires highly purified, intact RNA. See the catalog or website (www.invitrogen.com/ambion) for a complete listing of available FirstChoice RNAs.

RNA Isolation Kits See web or print catalog for P/Ns	Family of kits for isolation of total RNA. Included in the product line are kits using classical GITC and acidic phenol, one-step disruption/denaturation, phenol-free glass fiber filter or magnetic bead binding, and combination kits. See the catalog or website (www.invitrogen.com/ambion).
GLOBINclear™ Whole Blood Globin Reduction Kits P/N AM1980, AM1981	The GLOBINclear Whole Blood Globin Reduction Kits employ a novel, non-enzymatic technology to remove >95% of the globin mRNA from whole blood total RNA samples. The resulting mRNA is a superior template for RNA amplification and synthesis of labeled cDNA for array analysis. Kits are available for treatment of human or mouse/rat whole blood total RNA.
ArrayControl™ RNA Spikes P/N AM1780	The ArrayControl RNA Spikes are a set of eight control RNA transcripts designed for the normalization and validation of glass microarray experiments. The RNA Spikes range in size from 750 to 2000 bases, and each transcript has a 30-base 3' poly(A) tail. The precisely quantitated RNA Spikes are designed to be added to your RNA sample before labeling, to serve as internal controls for sample labeling and hybridization efficiency.
Biotin-11-UTP and Biotin-16-UTP P/N AM8450, AM8451, AM8452	Biotinylated UTPs are ideal for use as substrates in vitro transcription reactions, and can be utilized by a variety of RNA polymerases, including T7, T3, and SP6 RNA polymerases. Biotinylated RNA can be used in place of radioactively labeled RNA in many applications with detection via one of a variety of streptavidin-based methods.
RNA Fragmentation Reagents P/N AM8740	Amplified RNA is commonly fragmented prior to hybridization on oligonucleotide microarrays to improve the hybridization kinetics and signal produced on oligonucleotide microarrays. Ambion® RNA Fragmentation Reagents include a 10X Fragmentation Reagent and a Stop Solution.
Amino Alkyl cDNA Labeling Kit P/N AM1705	The Amino Alkyl cDNA Labeling Kit generates cDNA for secondary fluorescent dye labeling to be used for glass array analysis. It includes all the reagents, except the amine-reactive labeling moiety (e.g. cyanine dyes) for 2-step labeling of cDNA. The reaction produces more labeled cDNA, more efficiently than direct dye incorporation.
RETROscript® Kit P/N AM1710	First-strand cDNA synthesis kit. The RETROscript® Kit can be used to incorporate dye-modified nucleotides into cDNA using aRNA prepared with the MessageAmp™ II Kit as a template.
5-(3-aminoalkyl)-dUTP P/N AM8439	This 50 µM solution of amino alkyl-modified dUTP can be used with the RETROscript® Kit (P/N AM1710) to synthesize amine-reactive cDNA from aRNA. The amine-reactive cDNA can then be postlabeled with any amine-reactive label moiety.

II. aRNA Amplification Procedure

A. Important Parameters for Successful Amplification

Input RNA quantity and IVT reaction incubation time

Consider both the amount of sample RNA you have and the amount of aRNA needed for your analysis when planning MessageAmp II kit experiments. These factors will influence how much input RNA to use, whether one or two rounds of amplification should be done, and how long to incubate the IVT reaction.

Accurate quantitation

For experiments where the aRNA yield from different samples will be compared, it is *essential* to accurately quantify the input RNA used in the MessageAmp II kit procedure. The NanoDrop 1000A Spectrophotometer is recommended for rapid, accurate quantitation of nucleic acids; however, any reliable RNA quantitation method, such as traditional spectrophotometry or RiboGreen, can be used.

Recommended minimum and maximum amounts of input RNA

Table 1 shows the mass of total RNA that can be used in the MessageAmp II aRNA Amplification procedure. The RNA volume must be $\leq 10 \mu\text{L}$.

Table 1. Total RNA Input for MessageAmp II Procedure

Amplification	Recommended	Minimum	Maximum
single round	1000 ng	100 ng	5000 ng
two rounds	100 ng	0.1 ng	100 ng

Determining input RNA amount and IVT reaction incubation time

The procedure can accommodate a wide range of input RNA amounts, but for reproducible and comparable results, use a fixed amount of input RNA for all experiments. Tailor both the amount of input RNA and the amplification procedure to produce the amount of aRNA needed for your microarray hybridizations. For instance, Affymetrix GeneChip arrays require 10–15 μg of aRNA for each hybridization, but other commercial and core facility arrays may require slightly more or less aRNA.

Figure 3 shows aRNA yields from different amounts of Control RNA amplified with increasing IVT incubation times. aRNA yields from additional RNA sources are tabulated in Table 3 on page 22. These data show that when using 1000 ng or more RNA, a 4 hr IVT incubation produces enough aRNA for several microarray hybridizations. When amplifying samples with limited amounts of RNA however, (e.g. ~250 ng or less) incubating the IVT reaction for 14 hr will maximize the amount of aRNA produced.

Typically 100 ng of total RNA input is the lower limit for synthesizing ~10 µg of aRNA using a 4 hr IVT reaction. If your total RNA input will be ~100 ng or less, we advise conducting a preliminary amplification from a representative sample to determine how much aRNA you can expect from your experimental samples. If the preliminary experiment does not produce enough aRNA for two array hybridizations, repeat the experiment using an IVT incubation as long as 14 hr and/or consider using two rounds of amplification.




Figure 3. aRNA Yield vs. IVT Incubation Time and Total RNA Input.

The indicated amounts of Control RNA (HeLa total RNA) were amplified for the IVT incubation times shown. Although longer IVT incubation times produced more aRNA, the amount of aRNA needed for most microarrays was obtained using a 4-hour incubation time with ≥1000 ng input RNA. Note that this is empirical data obtained using the Control RNA provided with the kit; aRNA yield from experimental samples may be considerably different.

RNA purity

The quality of the RNA is the single most important factor affecting how efficiently an RNA sample will be amplified using the MessageAmp II aRNA Amplification Kit. RNA samples should be free of contaminating proteins, DNA, and other cellular material as well as phenol, ethanol, and salts associated with RNA isolation procedures. Impurities can lower the efficiency of reverse transcription and subsequently reduce the level of amplification. An effective measure of RNA purity is the ratio of absorbance readings at 260 and 280 nm. The ratio of A₂₆₀ to A₂₈₀ values should

fall in the range of 1.7–2.1. RNA must be suspended in high quality water or TE (10 mM Tris-HCl, 1 mM EDTA) or THE RNA Storage Solution (P/N AM7000, AM7001).

RNA integrity

The integrity of the RNA sample, or the proportion that is full-length, is another important component of RNA quality. Reverse transcription of partially degraded mRNAs will typically generate relatively short cDNAs that potentially lack portions of the coding region. RNA integrity can be evaluated by microfluidic analysis using the Agilent 2100 bioanalyzer and RNA LabChip® Kits. Primarily full-length RNA will exhibit a ratio of 28S to 18S rRNA bands that approaches 2:1. Using a bioanalyzer, the RIN (RNA Integrity Number) can be calculated to further evaluate RNA integrity.

Denaturing agarose gel electrophoresis and nucleic acid staining can also be used to separate and visualize the major rRNA species. When the RNA resolves into discrete rRNA bands (i.e., no significant smearing below each band), with the 28S rRNA band appearing approximately twice as intense as the 18S rRNA band, then the mRNA in the sample is likely to be mostly full-length. The primary drawback to gel electrophoresis is that it requires microgram amounts of RNA.




Figure 4. Bioanalyzer Images of Total RNA Preparations

These electropherograms (from the Agilent 2100 bioanalyzer) show RNA samples with decreasing integrity that are all of sufficient quality to use as input for the MessageAmp II aRNA Amplification Kit. The trace labeled “Intact Total RNA” represents the ideal for a bioanalyzer trace of total RNA. Notice how the ribosomal RNA peaks are at a ratio of about 2 (28S:18S) in this sample; this represents ultrahigh quality RNA in terms of integrity. Most RNA samples will more closely resemble the center trace where there are nearly equal amounts of 28S and 18S rRNA. The trace on the right shows a fairly typical human RNA prep with rRNA peaks that are lower than in the other two traces and where lower molecular weight degradation products become apparent. RNA samples with suboptimal integrity can yield meaningful array analysis results if the data are subjected to rigorous statistical analysis (Schoor et al. 2003).

Reaction incubation times should be precise and consistent

The incubation times for most of the enzymatic reactions in the procedure were optimized in conjunction with the kit reagents to ensure the maximum yield of nucleic acid product in each step—adhere to them closely. An exception is the IVT reaction, where a range of 4–14 hr incubation time is acceptable (step [II.F.3](#) on page 19). Refer to Figure 3 on page 9 and Table 2 on page 19 to help determine what incubation time to use. Although differences in IVT incubation time among samples has

had very little, if any, effect on array results in our hands, we recommend using uniform IVT incubation times if aRNA yield from different samples will be compared or if you want to have equal amplification of different samples—this will provide the most reproducible amplification and array analysis.

Master mixes

We strongly recommend preparing master mixes for the MessageAmp II aRNA Amplification procedure. This approach reduces the effects of pipetting error, saves time, and improves reproducibility. Using master mixes is especially important when aRNA yield from different samples will be compared.

Thorough mixing is very important for reproducibility

Below are specific instructions for mixing kit reagents, master mixes, and individual reactions. For maximum reproducibility and aRNA yield, follow these instructions closely.

Mix each kit component after thawing.

Mix enzyme solutions by *gently* flicking the tube a few times before adding them to reactions. Thaw frozen reagents completely *at room temperature* (i.e., primers, nucleotides, and 10X buffers), then mix thoroughly by vortexing, and keep on ice before use.

Mix master mixes by gentle vortexing.

After assembling master mixes, *gently* vortex to make a homogenous mixture without inactivating the enzyme(s).

Mix individual reactions by pipetting and flicking the tube.

After adding master mixes or other reagents to individual reactions, pipet up and down 2–3 times to rinse reagents from the pipet tip. Then flick the tube with your finger 3–4 times to mix thoroughly, and finish by centrifuging briefly to collect the reaction at the bottom of the tube.

Thermal cycler recommended

The MessageAmp II aRNA Amplification procedure is very sensitive to temperature; variable or inaccurate incubation temperatures can limit aRNA synthesis. It is also very important that condensation does not form in the reaction tubes during any of the incubations. Condensation changes the composition of reaction mixtures, which can greatly reduce yield.

- *A thermal cycler with a temperature adjustable heated lid is recommended.*

A calibrated thermal cycler, with a temperature-adjustable heated lid, is recommended, for the greatest temperature control and stability during MessageAmp II aRNA Amplification reaction incubations. Allow the thermal cycler to equilibrate to the required temperature before placing the tubes in the block for incubation.



NOTE

Even if you use a hybridization oven or incubator for most of the MessageAmp II aRNA Amplification reactions, a thermal cycler is strongly recommended for the 16°C second-strand synthesis reaction incubation (step II.D.2 on page 15). Turn off the heated lid if it cannot be adjusted to match the 16°C block temperature.

Maintaining consistency

Tubes: use 0.5 mL RNase-free nonstick tubes

Follow the recommended settings for the lid temperatures. Too high a lid setting may inhibit the reaction; too low a setting may cause condensation.

If your thermal cycler does not have a temperature-adjustable lid, or a thermal cycler is unavailable, calibrated hybridization ovens or incubators (at constant temperature) may also be used. Preheat incubators so that the correct temperature has stabilized before reactions are placed in the incubator. To avoid any potential influence on the reaction temperature from the tube holder, let tube holders equilibrate in the incubator for sufficient time, or use a tube holder that doesn't touch the sides and bottoms of the tubes—for example a floating tube support.

- *Heat blocks or water baths are not recommended for MessageAmp II aRNA Amplification reaction incubations.*

Procedural consistency is very important for amplification experiments. Consider implementing a detailed procedural plan that will be used by everyone in the lab to maintain consistency. This type of plan will minimize variation due to subtle procedural differences that can influence RNA amplification and may complicate gene expression studies. The plan should include basic information such as the method of RNA isolation, the amount of RNA to use in the procedure, and how long to incubate the IVT reaction. It should also address specifics that are not often included in protocols such as which tubes, tube racks, and incubators to use for each step in the process. Finally, develop a consistent work flow. For example standardize stopping points in the method. The idea is to standardize all of the variables discussed in this section of the Protocol and carefully follow all the steps in order to maximize amplification consistency among samples.

If a 60-well thermal cycler with temperature-adjustable lid is available, it is most convenient to conduct the MessageAmp II aRNA Amplification procedure in 0.5 mL nonstick tubes (for example, P/N AM12350). These can be thin-wall (PCR) tubes or ordinary-weight nonstick tubes. 0.5 mL tubes are large enough to accommodate the cDNA Binding Buffer without having to transfer reactions to a larger tube. Their small size and nonstick properties also keep the reaction components at the bottom of the tube.

If your thermal cycler is equipped with a standard 96-well block, 0.2 mL non-stick tubes can be used.

B. Prepare the Wash Buffer

Add 24 mL 100% ethanol (ACS grade or better) to the bottle labeled Wash Buffer. Mix well and mark the label to indicate that the ethanol was added.

C. Reverse Transcription to Synthesize First Strand cDNA

Incubator needed

Thermal cycler programmed as follows:

Temp	Time	Cycles
42°C (50°C lid)	2 hr	1
4°C	hold	

It is important to follow the lid temperature recommendation; higher lid temperatures may inhibit reverse transcription, while lower temperatures may cause condensation, resulting in changes to the reaction composition.

1. Adjust RNA sample volume to 10 μ L, if necessary

- a. Place a maximum volume of 10 μ L of total RNA (1000 ng recommended) into a nonstick, sterile, RNase-free, 0.5 mL tube. RNA must be in high quality water or TE. (See Table 1 on page 8 for minimum and maximum RNA input amounts.)



NOTE

If your experiment will include RNA spikes (e.g., Ambion ArrayControl RNA Spikes, P/N AM1780, or Affymetrix GeneChip Poly-A Control Kit, Cat #900433), add them to samples at this step.

2. Add 10 μ L of **Reverse Transcription Master Mix** and place at 42°C

- b. If necessary, add Nuclease-free Water to a final volume of 10 μ L, vortex briefly to mix, then centrifuge to collect the mixture at the bottom of the tube.
- a. At room temperature, prepare **Reverse Transcription Master Mix** in a nuclease-free tube. Assemble enough to synthesize first strand cDNA from all the RNA samples in the experiment, including $\leq 5\%$ overage to cover pipetting error. Assemble the Reverse Transcription Master Mix in the order shown:

Reverse Transcription Master Mix (for a single 20 μ L reaction)

Amount	Component
1 μ L	Nuclease-free Water
1 μ L	T7 Oligo(dT) Primer
2 μ L	10X First Strand Buffer
4 μ L	dNTP Mix
1 μ L	RNase Inhibitor
1 μ L	ArrayScript

- b. *Gently* vortex the tube to make a homogenous mixture without inactivating the enzyme, then centrifuge briefly (~5 sec) to collect the Reverse Transcription Master Mix at the bottom of the tube and place on ice.
- c. Transfer 10 µL of Reverse Transcription Master Mix to each RNA sample. Mix thoroughly by pipetting up and down 2–3 times, then flicking the tube 3–4 times, and centrifuge briefly to collect the reaction in the bottom of the tube.
- d. Place the samples in the thermal cycler, and start the run.

3. Incubate for 2 hr at 42°C

Incubate the reactions for 2 hr at 42°C, then centrifuge the tubes briefly (~5 sec) to collect the contents at the bottom of the tubes.

Place the tubes on ice and immediately proceed to second strand cDNA synthesis (below).



IMPORTANT

Proceed immediately to the next step.

D. Second Strand cDNA Synthesis

Incubator needed

Thermal cycler programmed as follows:

Temp	Time	Cycles
16°C (heat-disabled lid, or no lid)	2 hr	1
4°C	hold	

Disable the heated lid of the thermocycler if the lid temperature cannot be set in the range 16°C to room temperature. If the lid is too hot, inhibition of the reaction may occur.

1. Add 80 µL *Second Strand Master Mix* to each sample

- a. On ice, prepare a *Second Strand Master Mix* in a nuclease-free tube in the order listed below. Assemble enough to synthesize second strand cDNA from all the samples in the experiment, including ≤5% overage to cover pipetting error. Assemble the Second Strand Master Mix on ice in the order shown:

Second Strand Master Mix (for a single 100 µL reaction)

Amount	Component
63 µL	Nuclease-free Water
10 µL	10X Second Strand Buffer
4 µL	dNTP Mix
2 µL	DNA Polymerase
1 µL	RNase H

- b. Mix well by gently vortexing. Centrifuge briefly (~5 sec) to collect the Second Strand Master Mix at the bottom of the tube and place on ice.
- c. Transfer 80 μ L of Second Strand Master Mix to each sample. Mix thoroughly by pipetting up and down 2–3 times, then flicking the tube 3–4 times, and centrifuge briefly to collect the reaction in the bottom of the tube.
- d. Place the tubes in a 16°C thermal cycler. It is important to cool the thermal cycler block to 16°C before adding the reaction tubes because subjecting the reactions to temperatures >16°C will compromise aRNA yield.

2. Incubate for 2 hr at 16°C

Incubate 2 hr in a 16°C thermal cycler. If the lid temperature cannot be adjusted to match the 16°C block temperature, cover the reactions with the heated lid turned off, or if the lid cannot be turned off—do not cover the tubes with it. (Do not use a water bath or a heat block in a 4°C refrigerator for this incubation because the temperature will fluctuate too much.)

3. Place reactions on ice briefly or freeze immediately

After the 2 hr incubation at 16°C, place the reactions on ice and proceed to section [E. cDNA Purification](#) (below), or immediately freeze reactions at –20°C. Do not leave the reactions on ice for more than 1 hr.



STOPPING POINT

This is a potential overnight stopping point (at –20°C), but it is better to complete the cDNA purification (next section) before stopping.

E. cDNA Purification



IMPORTANT

All centrifugations in this purification procedure should be done at 10,000 x g (typically ~10,000 rpm) at room temperature. cDNA Filter Cartridges should not be subjected to RCFs over 16,000 x g because the force could cause mechanical damage and/or may deposit glass filter fiber in the eluate.

1. Preheat Nuclease-free Water to 55°C

Before beginning the cDNA purification, preheat at least 20 µL per sample of Nuclease-free Water to 55°C.



IMPORTANT

Preheat the Nuclease-free Water to a maximum of 55°C; temperatures above 58°C can partially denature the cDNA, compromising final aRNA yield.

2. Add 250 µL cDNA Binding Buffer to each sample



IMPORTANT

Check the cDNA Binding Buffer for precipitation before using it. If a precipitate is visible, redissolve it by warming the solution to 37°C for up to 10 min and vortexing vigorously. Cool to room temperature before use.

Add 250 µL of cDNA Binding Buffer to each sample, and mix thoroughly by pipetting up and down 2–3 times, then flicking the tube 3–4 times. Follow up with a quick spin to collect the reaction in the bottom of the tube. Proceed quickly to the next step.

3. Pass the mixture through a cDNA Filter Cartridge

Check that the cDNA Filter Cartridge is firmly seated in its wash tube (supplied).

- Pipet the cDNA sample/cDNA Binding Buffer (from step [2](#)) onto the center of the cDNA Filter Cartridge.
- Centrifuge for ~1 min at 10,000 x g, or until the mixture is through the filter.
- Discard the flow-through and replace the cDNA Filter Cartridge in the wash tube.



IMPORTANT

Make sure that the ethanol has been added to the bottle of Wash Buffer before using it in this step.

4. Wash with 500 μ L Wash Buffer

- Apply 500 μ L Wash Buffer to each cDNA Filter Cartridge.
- Centrifuge for \sim 1 min at 10,000 \times g, or until all the Wash Buffer is through the filter.
- Discard the flow-through and spin the cDNA Filter Cartridge for an additional minute to remove trace amounts of Wash Buffer.
- Transfer cDNA Filter Cartridge to a cDNA Elution Tube.

5. Elute cDNA with 18 μ L 55°C Nuclease-free Water

It is important to use Nuclease-free Water that is at 50–55°C for the cDNA elution. Colder water will be less efficient at eluting the cDNA, and hotter water (\geq 58°C) may result in reduced aRNA yield.

- Apply 18 μ L of preheated Nuclease-free Water (55°C) to the center of the filter in the cDNA Filter Cartridge.
- Leave at room temperature for 2 min and then centrifuge for 1 min at 10,000 \times g, or until all the Nuclease-free Water is through the filter. The double-stranded cDNA will now be in the eluate (\sim 16 μ L).



STOPPING POINT

The purified cDNA can be stored overnight at -20°C at this point if desired. Transfer the cDNA to a lidded, non-stick, nuclease-free tube for storage.

F. In Vitro Transcription to Synthesize aRNA

Incubator needed

Thermal cycler programmed as follows:

Temp	Time	Cycles
37°C (default lid; 100–105°C)	4–14 hr; see step 3	1
4°C	hold	

It is important to use a thermal cycler with a heated lid (100–105°C), to maintain uniform temperature, and because it is extremely important that condensation does not form inside the tubes; this would change the reagent concentrations and reduce yield.

The IVT reaction can be set up to synthesize unmodified aRNA, or biotin-labeled UTP (or other modified nucleotides) can be incorporated into the aRNA during the IVT. Samples that will undergo two rounds of amplification *cannot* be labeled with biotin at this step because it would interfere with the second round of amplification.

For the highest aRNA yield, conduct the IVT in a 40 μ L final reaction volume. We provide instructions for a 20 μ L biotin-labeled reaction for users who want to reduce the reaction cost by using only half the quantity of biotin-labeled nucleotide.

1a. For 20 µL IVT reactions, concentrate the cDNA

If you will be conducting a 20 µL IVT reaction, concentrate the eluted cDNA from step [II.E.5](#) on page 17 in a vacuum centrifuge concentrator until the volume is reduced to 8 µL. This should take only a few minutes. Avoid drying the mixture to completion.

1b. For 40 µL IVT reactions, proceed directly to step 2

If you will be conducting a 40 µL IVT reaction, using unmodified or biotin-labeled UTP, transfer the eluted cDNA from step [II.E.5](#) on page 17 to a PCR tube and proceed directly to step [2](#) below.

2. Add IVT Master Mix to each sample, and mix

a. At room temperature, prepare an **IVT Master Mix** by adding the following reagents to a nuclease-free microcentrifuge tube in the order listed below. Assemble enough for all the samples in the experiment, including ≤5% overage to cover pipetting error.



IMPORTANT

*If two rounds of amplification will be done, this first round transcription should be **unmodified**, containing **only unmodified UTP**.*

Assemble the IVT Master Mix at room temperature in the order shown:

Biotin labeled		Unmodified		Component
40 µL rxn	20 µL rxn	40 µL rxn	40 µL rxn	
--	(8 µL)	--	--	concentrated cDNA from step 1a on page 18
(16 µL)	--	(16 µL)	--	double-stranded cDNA (from step II.E.5 on page 17)
IVT Master Mix for a single reaction				
4 µL	2 µL	4 µL	4 µL	T7 ATP
4 µL	2 µL	4 µL	4 µL	T7 CTP
4 µL	2 µL	4 µL	4 µL	T7 GTP
2.6 µL	1.3 µL	4 µL	4 µL	T7 UTP
1.4 µL	0.7 µL	--	--	Biotin-11-UTP*, 75 mM
4 µL	2 µL	4 µL	4 µL	T7 10X Reaction Buffer
4 µL	2 µL	4 µL	4 µL	T7 Enzyme Mix

* Biotin-16-UTP can be used instead of Biotin-11-UTP, if needed.

- b. Mix well by gently vortexing. Centrifuge briefly (-5 sec) to collect the IVT Master Mix at the bottom of the tube and place on ice.
- c. Transfer IVT Master Mix to each sample following the volume guidelines below. Mix thoroughly by pipetting up and down 2–3 times, then flicking the tube 3–4 times, and centrifuge briefly to collect the reaction in the bottom of the tube.

**IMPORTANT**

If the cDNA was vacuum dried, resuspend the template by washing the sides of the tube with the IVT Master Mix as you add it to the tube. It is important to make sure that all of the material is resuspended in the IVT reaction mix before the mixture is transferred to a PCR tube for IVT.

Biotin labeled		Unmodified	
40 μ L rxn	20 μ L rxn	40 μ L rxn	
24 μ L	12 μ L	24 μ L	volume IVT Master Mix

d. Once assembled, place the tubes in the thermal cycler and start the run.

3. Incubate 4–14 hr at 37°C

The minimum recommended incubation time is 4 hr; the maximum is 14 hr. (The reactions can be held post-IVT at 4°C for up to 48 hr, for convenience.)

Use the following table as a guide to determine how long to continue your IVT reaction. There are more data and a detailed discussion of the length of the IVT incubation in section [II.A Input RNA quantity and IVT reaction incubation time](#) starting on page 8.

Table 2. Recommended IVT Incubation Times

aRNA Needed	Input Total RNA	IVT Incubation
10–100 μ g	1–5 μ g	4 hr
1–10 μ g	0.1–1 μ g	8 hr
0.1–1 μ g	\leq 100 ng	14 hr

4. Add Nuclease-free Water to bring each sample to 100 μ L

Stop the reaction by adding Nuclease-free Water to each aRNA sample to bring the final volume to 100 μ L. Mix thoroughly by gentle vortexing.

Proceed to the aRNA purification step (below) or store at -20°C .

**STOPPING POINT**

The aRNA can be stored overnight at -20°C at this point if desired.

G. aRNA Purification

Incubator needed: heat block set at 55°C.

This purification removes enzymes, salts, and unincorporated nucleotides from the aRNA. At the end of the purification the aRNA is eluted from the filter with Nuclease-free Water.



IMPORTANT

All centrifugations in this purification procedure should be done at 10,000 x g (typically ~10,000 rpm) at room temperature.

aRNA Filter Cartridges should not be subjected to RCFs over 16,000 x g because the force could cause mechanical damage and/or may deposit glass filter fiber in the eluate.

1. Preheat Nuclease-free Water and assemble Filter Cartridges and Tubes

- Preheat a minimum of 200 µL per sample of Nuclease-free Water to 55°C.
- For each sample, place an aRNA Filter Cartridge into an aRNA Collection Tube and set aside for use in step 4.

2. Add 350 µL aRNA Binding Buffer to each sample

- a. Check to make sure that each IVT reaction was brought to 100 µL with Nuclease-free Water.
- b. Add 350 µL of aRNA Binding Buffer to each aRNA sample. Proceed to the next step immediately.

3. Add 250 µL 100% ethanol and pipet 3 times to mix



IMPORTANT

It is crucial to follow these mixing instructions exactly, and to proceed quickly to the next step.

Add 250 µL of ACS grade 100% ethanol to each aRNA sample, and mix by pipetting the mixture up and down 3 times. **Do NOT vortex to mix and do NOT centrifuge.**

Proceed **immediately** to the next step as soon as you have mixed the ethanol into each sample. Any delay in proceeding could result in loss of aRNA because once the ethanol is added, the aRNA will be in a semi-precipitated state.

4. Pass samples through an aRNA Filter Cartridge

- a. Pipet each sample mixture from step 3 onto the center of the filter in the aRNA Filter Cartridge/Collection Tube assembly.
- b. Centrifuge for ~1 min at 10,000 x g. Continue until the mixture has passed through the filter.
- c. Discard the flow-through and replace the aRNA Filter Cartridge back into the aRNA Collection Tube.

5. Wash with 650 µL Wash Buffer

- a. Apply 650 µL Wash Buffer to each aRNA Filter Cartridge.
- b. Centrifuge for ~1 min at 10,000 x g, or until all the Wash Buffer is through the filter.
- c. Discard the flow-through and spin the aRNA Filter Cartridge for an additional ~1 min to remove trace amounts of Wash Buffer.
- d. Transfer Filter Cartridge(s) to a fresh aRNA Collection Tube.

6. Elute aRNA with 200 µL 55°C Nuclease-free Water

- a. To the center of the filter, add 200 µL Nuclease-free Water (preheated to 55°C).
- b. Incubate the samples in the 55°C heat block for 10 min (recommended).

Alternatively, incubate at room temperature for 2 min. This results in ~80% recovery of the aRNA.

- c. Centrifuge for ~1.5 min at 10,000 × g, or until the Nuclease-free Water is through the filter.
- d. The aRNA will now be in the aRNA Collection Tube in ~200 μL of Nuclease-free Water.

7. Store aRNA at –80°C

Store aRNA at –80°C for up to 1 year, and minimize repeated freeze-thawing. To prevent multiple freeze-thaw events, split samples into 5–20 μg aliquots for microarray labeling and hybridizations.

8. (Optional) Concentrate the purified aRNA

If necessary, concentrate the aRNA by vacuum centrifugation or by precipitation with ammonium acetate (NH₄OAc)/ethanol.

(Optional) Concentrate by vacuum centrifugation

If the heater on the vacuum centrifuge has different settings, use medium or low. Check the progress of drying every 5–10 min, and remove the sample from the concentrator when it reaches the desired volume.

(Optional) Precipitate with 5 M NH₄OAc and ethanol

- a. Add 0.1 volume of 5 M NH₄OAc to the purified aRNA (20 μL if the aRNA was eluted in 200 μL Nuclease-free water).
- b. Add 2.5 volumes of 100% ethanol (550 μL if the aRNA was eluted in 200 μL). Mix well and incubate at –20°C for 30 min.
- c. Microcentrifuge at top speed for 15 min at 4°C or room temperature. Carefully remove and discard the supernatant.
- d. Wash the pellet with 500 μL 70% cold ethanol, centrifuge again, and remove the 70% ethanol.
- e. To remove the last traces of ethanol, quickly respin the tube, and aspirate any residual fluid with a fine-tipped pipette or syringe needle.
- f. Air dry the pellet.
- g. Resuspend the aRNA pellet using the desired solution and volume.

III. Assessing aRNA Yield and Quality

A. aRNA Quantitation and Expected Yield

1. Assessing aRNA yield by UV absorbance

The concentration of an aRNA solution can be determined by measuring its absorbance at 260 nm using a spectrophotometer.

- Use a NanoDrop spectrophotometer and measure 1.5 µL of the RNA sample directly.
- With a traditional spectrophotometer, dilute an aliquot of the RNA 1:50 to 1:100 in TE (10 mM Tris-HCl pH 8, 1 mM EDTA), and read the absorbance. (Be sure to zero the spectrophotometer with the TE used for sample dilution.) The buffer used for dilution need not be RNase-free, since slight degradation of the RNA will not significantly affect its absorbance.

Find the concentration in µg/mL by multiplying the A_{260} by the dilution factor and the extinction coefficient. ($1 A_{260} = 40 \mu\text{g RNA/mL}$):

$$A_{260} \times \text{dilution factor} \times 40 = \mu\text{g RNA/mL}$$

2. Assessing aRNA yield with the RiboGreen® assay

If a fluorometer or a fluorescence microplate reader is available, Invitrogen's RiboGreen fluorescence-based assay for RNA quantitation is a convenient and sensitive way to measure RNA concentration. Follow the manufacturer's instructions for using RiboGreen.

3. Expected yield

The aRNA yield will depend on the amount and quality of poly(A) RNA in the input total RNA. Since the proportion of poly(A) RNA in total RNA is affected by influences such as health of the organism and the organ from which it is isolated, aRNA yield from equal amounts of total RNA may vary considerably. The following table shows empirical aRNA yield data obtained using this kit.

Table 3. aRNA Yields from Different Amounts of Selected RNAs Amplified Using a 4 hr IVT Incubation.

RNA source	Input total RNA used in MessageAmp II					Amount of aRNA synthesized (µg)
	50 ng	100 ng	500 ng	1 µg	3 µg	
Human brain	4.2	12.1	47.1	127.0	204.6	
Human heart	2.3	6.7	23.3	57.1	126.8	
Human kidney	2.0	4.0	21.3	39.0	93.9	
HeLa S3 cells	5.3	11.5	59.1	87.3	89.7	
Rat thymus	2.4	4.0	23.1	46.5	100.5	
Universal Human Reference RNA*	—	—	47.3	75	105.9	

* Stratagene

B. Analysis of aRNA Size

The size distribution of aRNA can be evaluated using an Agilent 2100 bioanalyzer with LabChip technology, or by conventional denaturing agarose gel analysis. The bioanalyzer can provide a fast and accurate size distribution profile of aRNA samples, but aRNA yield should be determined by UV absorbance or RiboGreen analysis. To analyze aRNA size using a bioanalyzer, follow the manufacturer's instructions for running the assay using purified aRNA (from step [II.G.6](#) on page 20). For more information, go to www.invitrogen.com/site/us/en/home/support/technical-support.html

Expected aRNA size

Agilent bioanalyzer analysis

The expected aRNA profile is a distribution of sizes 250–5500 nt with most of the aRNA 1000–1500 nt (Figure 5). To compare bioanalyzer profiles of different aRNA samples, be sure to load equal mass amounts to get an accurate comparison.

Denaturing agarose gel analysis

Amplified aRNA should appear as a smear from 250 to 5000 nt. The average size of biotin labeled aRNA should be approximately 1400 nt; the average size of unmodified aRNA should be ~1150 nt.




Figure 5. Bioanalyzer Analysis of a Biotin-Labeled MessageAmp™ II aRNA Amplification Kit Reaction.

This electropherogram displays the aRNA size distribution from 1 round of amplification of 1 µg of the Control RNA using biotin labeled NTPs in a 40 µL IVT reaction that was incubated for 14 hr. Using the Agilent bioanalyzer mRNA smear assay, 50% of the aRNA is calculated to be larger than 1470 nt. (The vertical grey lines mark 65 and 1470 bases. The area outside these lines represents 50% of the area under the curve produced by the product of the positive control reaction.)

IV. (Optional) Second Round Amplification

If one cycle of amplification does not yield the amount of aRNA necessary for your experiments, a second round of amplification can be conducted to generate additional aRNA (see Figure 2 on page 3 for an overview of second round amplification). In order to conduct two rounds of amplification, the first round of amplification must contain *only* unmodified CTP and UTP; ***modified aRNA cannot undergo a second round of amplification.***

The procedure is similar to the first round of amplification, and the reaction products are equivalent, but different primers are used for first and second strand synthesis, and the reaction setup is slightly different.

Second round amplification products are typically shorter than first round amplification products, but we have not seen that this has adverse effects on array hybridization results if all samples are prepared using two rounds of amplification.

A. Synthesis of First Strand cDNA (Second Round)

Incubator needed

Thermal cycler programmed as follows:

Temp	Time	Cycles
37°C (default lid; 100–105°C)	30 min	1
4°C	hold	

It is important to use a thermal cycler with a heated lid (100–105°C), to maintain uniform temperature, and because it is extremely important that condensation does not form inside the tubes; this would change the reagent concentrations and reduce yield.

1. Mix ≤ 2 μ g aRNA with 2 μ L Second Round Primers, and adjust the volume to 12 μ L

- a. Place up to 2 μ g of purified aRNA from the first round of amplification into a sterile RNase-free microcentrifuge tube. With very small RNA samples, we have dried the entire first round amplification reaction, and used it as starting material for the second round amplification.



IMPORTANT

The volume of the aRNA must be ≤ 10 μ L. If necessary, concentrate the aRNA in a vacuum centrifuge. Do not dry the aRNA completely, as this could impede reverse transcription.

- b. Add 2 μ L of Second Round Primers.
- c. Add Nuclease-free Water to bring the volume to 12 μ L, vortex briefly to mix, and centrifuge briefly to collect the reaction in the bottom of the tube.

2. Incubate 10 min at 70°C

- a. Program a thermal cycler for the annealing step:

Temp	Time	Cycles
70°C (default lid; 100–105°C)	10 min	1
4°C	hold	

- b. Place the samples in the equilibrated thermal cycler, start the run, and incubate for 10 minutes at 70°C.

It is important to follow the lid temperature recommendation; higher lid temperatures may inhibit reverse transcription, while lower temperatures may cause condensation, resulting in changes to the reaction composition.

- c. Remove the RNA samples from the 70°C incubation and centrifuge briefly (~5 sec) to collect the reaction at the bottom of the tube. Place the reaction on ice briefly before starting step 3.

3. Add 8 µL of Reverse Transcription Master Mix and place at 42°C

- a. While the samples are incubating at 70°C, prepare
- Reverse Transcription Master Mix*
- in a nuclease-free tube at room temperature. Assemble enough master mix for all of the samples in the experiment, including ≤5% overage to cover pipetting error. Assemble the Reverse Transcription Master Mix in the order shown:

Reverse Transcription Master Mix (for a single 20 µL reaction)

Amount	Component
2 µL	10X First Strand Buffer
4 µL	dNTP Mix
1 µL	RNase Inhibitor
1 µL	ArrayScript

- b. Mix well by gently vortexing. Centrifuge briefly (~5 sec) to collect the Reverse Transcription Master Mix at the bottom of the tube and place on ice.

- c. Transfer 8 µL of Reverse Transcription Master Mix to each sample. Mix thoroughly by pipetting up and down 2–3 times, then flicking the tube 3–4 times, and centrifuge briefly to collect the reaction in the bottom of the tube.

- d. Program the thermal cycler for the reverse transcription reaction:

Temp	Time	Cycles
42°C (50°C lid)	2 hr	1
4°C	hold	

- e. Place the tubes in the thermal cycler which has equilibrated to 42°C and start the run.

4. Incubate 2 hr at 42°C

Incubate reactions for 2 hr at 42°C, then centrifuge briefly (~5 sec) to collect the reaction at the bottom of the tube.

5. Add 1 µL RNase H and incubate 30 min at 37°C

RNase H specifically degrades the aRNA leaving only the cDNA as template for second strand synthesis. This helps assure that the second strand synthesis reaction will be primed exclusively by the T7 Oligo(dT) Primer.

- a. Add 1 µL RNase H to each reaction. Mix thoroughly by pipetting up and down 2–3 times, then flicking the tube 3–4 times, and centrifuge briefly to collect the reaction in the bottom of the tube.
- b. Incubate for 30 min at 37°C in a thermal cycler. After the incubation, proceed directly to Second Strand Synthesis (below).

B. Synthesis of Second Strand cDNA (Second Round)

Incubator needed:

Thermal cycler with a temperature-adjustable lid

1. Add 5 µL T7 Oligo(dT) Primer to each sample

- a. Add 5 µL T7 Oligo(dT) Primer to cDNA sample.
- b. Mix well by gently vortexing, then centrifuge briefly (~5 sec) to collect the sample at the bottom of the tube.

2. Incubate for 10 min at 70°C, then place on ice

- a. Program a thermal cycler for the annealing step:

Temp	Time	Cycles
70°C (default lid; 100–105°C)	10 min	1
4°C	hold	

- b. Place the samples in the equilibrated thermal cycler, start the run, and incubate 10 min at 70°C in a thermal cycler.
- c. Place the reaction on ice briefly before adding the remaining second strand cDNA synthesis reagents.

3. Add 74 µL *Second Strand Master Mix* to each sample

- a. On ice, Assemble *Second Strand Master Mix* for all the samples in the experiment, including ≤5% overage to cover pipetting error. Assemble the Second Strand Master Mix on ice in the order shown:

Second Strand Master Mix (for a single 100 µL reaction)	
Amount	Component
58 µL	Nuclease-Free Water
10 µL	10X Second Strand Buffer
4 µL	dNTP Mix
2 µL	DNA Polymerase

- b. Mix well by gently vortexing. Centrifuge briefly (~5 sec) to collect the Second Strand Master Mix at the bottom of the tube and place on ice.

- c. Transfer 74 μL of Second Strand Master Mix to each sample. Mix thoroughly by pipetting up and down 2–3 times, then flicking the tube 3–4 times, and centrifuge briefly to collect the reaction in the bottom of the tube.
- d. Program a thermal cycler for the reverse transcription reaction:

Temp	Time	Cycles
16°C (heat-disabled lid, or no lid)	2 hr	1
4°C	hold	

Disable the heated lid of the thermocycler if the lid temperature cannot be set in the range 16°C to room temperature. If the lid is too hot, inhibition of the reaction may occur.

- e. Place the tubes in the pre-cooled 16°C thermal cycler. It is important to cool the thermal cycler block to 16°C before placing the reaction tubes; subjecting the reactions to temperatures >16°C could compromise aRNA yield.

4. Incubate 2 hr at 16°C

Incubate 2 hr in a 16°C thermal cycler. If the lid temperature cannot be adjusted to match the 16°C block temperature, cover the reactions with the heated lid turned off, or if the lid cannot be turned off—do not cover the reactions. (Do not use a water bath or a heat block in a 4°C refrigerator for this incubation because the temperature will fluctuate too much.)



STOPPING POINT

This is a potential overnight stopping point, but it is better to complete the cDNA purification (section [II.E](#) on page 16) before stopping.

5. Continue with the procedure at section [II.E](#) starting on page 16

After the 2 hr incubation at 16°C, place the reactions on ice and proceed to section [II.E. cDNA Purification](#) on page 16, or immediately freeze reactions at -20°C . Do not leave the reactions on ice for long periods of time.

Complete the rest of the second round amplification; start at section [II.E. cDNA Purification](#) on page 16, and continue through the remainder of section [II](#).

V. Troubleshooting

A. Positive Control Reaction

Control RNA amplification instructions

To establish if the kit is working properly, Control RNA consisting of 1 mg/mL HeLa cell total RNA is provided.

1. Use 1 µL of the Control RNA in a single-round MessageAmp II kit reaction; follow the procedure starting at step [II.C.1](#) on page 13.
2. At step [II.F.3](#) on page 19, use a **14 hr** incubation for the IVT reaction.
3. Continue with the procedure for making unmodified aRNA through section [II.G.](#)

Analysis of the positive control amplification reaction

- After completing the aRNA purification, measure the A_{260} of the reaction product as described in section [III.A.1](#) on page 22. ***The positive control reaction should produce ≥ 50 µg of aRNA.***
Be aware that often the positive control reaction cannot be compared to experimental reactions, because many experimental amplification experiments will use less than the 1 µg of input RNA used in the positive control reaction, and the aRNA yield will be proportionately lower. Also the Control RNA is of exceptional quality and purity, ensuring that it will amplify with extremely high efficiency.
- Also run a 2 µg aliquot of the reaction products on a denaturing agarose gel or analyze 100–200 ng on a bioanalyzer; ***the average size of the aRNA should be ≥ 1 kb.***

B. Factors that Affect Both the Positive Control and Experimental Samples

If the positive control reaction yield or amplification product size does not meet expectations, consider the following possible causes and troubleshooting suggestions. These suggestions also apply to problems with amplification of experimental RNA.

Incubation temperature(s) were incorrect

The incubation temperatures are critical for effective RNA amplification.

- Check the temperatures of all incubators used in the procedure with a calibrated thermometer.
- If a thermal cycler is used for incubation, check the accuracy of the adjustable temperature lid. If the lid temperature cannot be adjusted to match the desired reaction temperature, use the lid with the heat turned off, or do not use the lid to cover the reaction vessel(s).

Condensation formed in the tube during the reaction incubation(s)

Condensation occurs when the cap of the reaction vessel is cooler (e.g., room temperature) than the bottom of the tube. As little as 1–2 μL of condensate in an IVT reaction tube throws off the concentrations of the nucleotides and magnesium, which are crucial for good yield.

If you see condensation, spin the tube briefly and mix the reaction gently. Move the tube(s) to an incubator where condensation does not occur or is minimized.

Nuclease-contaminated tubes, tips, or equipment

Using pipettes, tubes, or equipment that are contaminated with nucleases can cleave the RNA or DNA being generated at each step in the procedure. This will reduce the size of the aRNA products and decrease aRNA yield. Both RNases and DNases can be removed from surfaces using Ambion RNaseZap[®] RNase Decontamination Solution (P/N AM9780, AM9786).

Absorbance readings were inaccurate

Confirm that your spectrophotometer is accurate by measuring the absorbance of an RNA or DNA sample of known concentration. Alternatively, assess the aRNA concentration by fractionating on an agarose gel adjacent to an RNA sample whose concentration is known. Comparing the ethidium bromide staining of the aRNA and control samples can approximate the concentration of the aRNA.

C. Troubleshooting Low Yield and Small Average aRNA Size

Consider the following troubleshooting suggestions if the positive control reaction produced the expected results, but amplification of your experimental samples results in less or smaller (average <500 nt) aRNA than expected.

Impure RNA samples

RNA samples with significant amounts of contaminating DNA, protein, phenol, ethanol, or salts are reverse transcribed poorly and subsequently generate less aRNA than pure RNA samples. Phenol extract and ethanol precipitate your RNA, or use the Ambion MEGAClear[™] Kit (P/N AM1908) to further purify your RNA before reverse transcription.

Lower than expected input RNA concentration

Take another A_{260} reading of your RNA sample or try using more RNA in the aRNA amplification procedure.

RNA integrity is compromised

RNA that is partially degraded generates cDNA that is relatively short. This will reduce the average size of the aRNA population and subsequently reduce the yield of aRNA. You can assess the integrity of an RNA sample by determining the size of the 18S and 28S rRNA bands and the relative abundance of 28S to 18S rRNA (See [RNA integrity](#) on page 10 for more information).

The mRNA content of your total RNA sample is lower than expected

Different RNA samples contain different amounts of mRNA. In healthy cells, mRNA constitutes 1–3% of total cellular RNA. The actual amount of mRNA depends on the cell type and the physiological state of the sample. When calculating the amount of amplification, the starting mass of mRNA in a total RNA prep should always be considered within a range of 10–30 ng per μg of total RNA (assuming good RNA quality). Most total RNA samples can be amplified up to 1000 fold producing 10–30 μg of aRNA from 1 μg of total RNA.

D. aRNA is Not Efficiently Reverse Transcribed

The cDNA procedure relies on oligo(dT) priming

The aRNA has a poly(U) tract near the 5' end but lacks a poly(A) tract at its 3' end. Thus any reverse transcription procedures that rely on oligo(dT) primers will not effectively convert aRNA to cDNA. Try using gene specific or random primers.

The filter in the aRNA Filter Cartridge was not completely dried after the wash steps

If the aRNA contains ethanol carried over from the Wash Buffer, it can inhibit reverse transcription. Make sure that the filter is completely dry at step [II.G.5.c](#) on page 20 just before eluting the aRNA.

To remove ethanol from an aRNA sample, ethanol precipitate it following the instructions in step [II.G.8](#) on page 21. Be sure to include the double spin described in step [8.e](#) to remove the last traces of ethanol before air drying the aRNA pellet.

Absorbance readings are inaccurate

Confirm that your spectrophotometer is accurate by measuring the absorbance of an RNA or DNA sample of known concentration. Alternatively, assess the quantity of aRNA by a different method such as fractionating on an agarose gel adjacent to an RNA sample whose concentration is known and comparing the ethidium bromide staining or using a sensitive RNA dye such as RiboGreen.

VI. Appendix

A. References

- Feldman AJ, Costouros NG, Wang E, Qian M, Marincola FM, Alexander HR, and Libutti SK (2002) Advantages of mRNA amplification for microarray analysis. *Biotechniques* 33(4):906–914.
- Kacharina JE, Crino PB, Eberwine J (1999) Preparation of cDNA from single cells and subcellular regions. *Methods Enzymol*, 303:3–18.
- Li Y, Li T, Liu S, Qiu M, Han Z, Jiang Z, Li R, Ying K, Xie Y, Mao Y (2004) Systematic comparison of the fidelity of aRNA, mRNA and T-RNA on gene expression profiling using cDNA microarray. *J Biotechnology* 107(1):19–28.
- Luo L, Salunga RC, Guo H, Bittner A, Joy KC, Galindo JE, Xiao H, Rogers KE, Wan JS, Jackson MR, Erlander MG (1999) Gene expression profiles of laser-captured adjacent neuronal subtypes. *Nat Med* 5(1):117–122.
- Pabon C, Modrusan Z, Ruvolo MV, Coleman IM, Daniel S, Yue H, Arnold LJ Jr. (2001) Optimized T7 amplification system for microarray analysis. *Biotechniques* 31(4):874–879.
- Polacek DC, Passerini AG, Shi C, Francesco NM, Manduchi E, Grant GR, Powell S, Bischof H, Winkler H, Stoeckert CJ Jr, Davies PF (2003) Fidelity of enhanced sensitivity of differential transcription profiles following linear amplification of nanogram amounts of endothelial mRNA. *Physiol Genomics* 13:147–156.
- Schoor O, Weinschenk T, Hennenlotter J, Corvin S, Stenzl A, Rammensee H-G, and Stevanovic S (2003) Moderate degradation does not preclude microarray analysis of small amounts of RNA. *Biotechniques* 35:1192–1201.
- Van Gelder RN, von Xastrow ME, Yool A, Dement DC, Barchas JD, Eberwine JH (1990) Amplified RNA synthesized from limited quantities of heterogeneous cDNA. *Proc Natl Acad Sci USA*, 87:1663–1667.
- Yue H, Eastman PS, Wang BB, Minor J, Doctolero MH, Nuttall RL, Stack R, Becker JW, Montgomery JR, Vainer M, Johnston R (2001) An evaluation of the performance of cDNA microarrays for detecting changes in global mRNA expression. *Nucleic Acids Res* 29(8):E41–1.

B. Quality Control

Functional testing

The Control RNA is used in a MessageAmp II kit reaction following the instructions in section [V.A](#) on page 28. The aRNA yield is assessed by measuring the A_{260} on the NanoDrop ND1000A spectrophotometer. The median size of the aRNA is assessed using the mRNA smear assay on the Agilent 2100 bioanalyzer.

Nuclease testing

Relevant kit components are tested in the following nuclease assays:

RNase activity

Meets or exceeds specification when a sample is incubated with labeled RNA and analyzed by PAGE.

Nonspecific endonuclease activity

Meets or exceeds specification when a sample is incubated with super-coiled plasmid DNA and analyzed by agarose gel electrophoresis.

Exonuclease activity

Meets or exceeds specification when a sample is incubated with labeled double-stranded DNA, followed by PAGE analysis.

Protease testing

Meets or exceeds specification when a sample is incubated with protease substrate and analyzed by fluorescence.

C. Safety Information



WARNING

GENERAL SAFETY. Using this product in a manner not specified in the user documentation may result in personal injury or damage to the instrument or device. Ensure that anyone using this product has received instructions in general safety practices for laboratories and the safety information provided in this document.

- Before using an instrument or device, read and understand the safety information provided in the user documentation provided by the manufacturer of the instrument or device.
- Before handling chemicals, read and understand all applicable Safety Data Sheets (SDSs) and use appropriate personal protective equipment (gloves, gowns, eye protection, etc). To obtain SDSs, see the “Documentation and Support” section in this document.

1. Chemical safety



WARNING

GENERAL CHEMICAL HANDLING. To minimize hazards, ensure laboratory personnel read and practice the general safety guidelines for chemical usage, storage, and waste provided below, and consult the relevant SDS for specific precautions and instructions:

- Read and understand the Safety Data Sheets (SDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials. To obtain SDSs, see the “Documentation and Support” section in this document.
- Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing).
- Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (for example, fume hood).
- Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer's cleanup procedures as recommended in the SDS.
- Handle chemical wastes in a fume hood.

- Ensure use of primary and secondary waste containers. (A primary waste container holds the immediate waste. A secondary container contains spills or leaks from the primary container. Both containers must be compatible with the waste material and meet federal, state, and local requirements for container storage.)
- After emptying a waste container, seal it with the cap provided.
- Characterize (by analysis if necessary) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- Ensure that the waste is stored, transferred, transported, and disposed of according to all local, state/provincial, and/or national regulations.
- **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.

2. Biological hazard safety



WARNING

Potential Biohazard. Depending on the samples used on the instrument, the surface may be considered a biohazard. Use appropriate decontamination methods when working with biohazards.



WARNING

BIOHAZARD. Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. Follow all applicable local, state/provincial, and/or national regulations. Wear appropriate protective equipment, which includes but is not limited to: protective eyewear, face shield, clothing/lab coat, and gloves. All work should be conducted in properly equipped facilities using the appropriate safety equipment (for example, physical containment devices). Individuals should be trained according to applicable regulatory and company/institution requirements before working with potentially infectious materials. Read and follow the applicable guidelines and/or regulatory requirements in the following:

In the U.S.:

- U.S. Department of Health and Human Services guidelines published in Biosafety in Microbiological and Biomedical Laboratories found at: www.cdc.gov/biosafety
- Occupational Safety and Health Standards, Bloodborne Pathogens (29 CFR§1910.1030), found at: www.access.gpo.gov/nara/cfr/waisidx_01/29cfr1910a_01.html
- Your company's/institution's Biosafety Program protocols for working with/handling potentially infectious materials.
- Additional information about biohazard guidelines is available at: www.cdc.gov

In the EU:

Check local guidelines and legislation on biohazard and biosafety precaution and refer to the best practices published in the World Health Organization (WHO) Laboratory Biosafety Manual, third edition, found at: www.who.int/csr/resources/publications/biosafety/WHO_CDS_CSR_LYO_2004_11/en/

VII. Documentation and Support

A. Obtaining SDSs

Safety Data Sheets (SDSs) are available from:

www.invitrogen.com/sds

or

www.appliedbiosystems.com/sds

B. Obtaining support

For the latest services and support information for all locations, go to:

www.invitrogen.com

or

www.appliedbiosystems.com

At the website, you can:

- Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities
- Search through frequently asked questions (FAQs)
- Submit a question directly to Technical Support
- Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents
- Obtain information about customer training
- Download software updates and patches



1751ME

Headquarters

5791 Van Allen Way | Carlsbad, CA 92008 USA

Phone +1 760 603 7200 | Toll Free in USA 800 955 6288

For support visit www.appliedbiosystems.com/support

www.lifetechnologies.com

