

Human Resistin ELISA Kit

Catalog. no. KHP0051

Pub. No. MAN0009568 Rev 2.0

Description

The Human Resistin ELISA Kit is a solid-phase sandwich Enzyme-Linked Immunosorbent Assay (ELISA) designed to detect and quantify the level of human resistin in serum, plasma, and cell culture supernatant. The assay will recognize both natural and recombinant human resistin.

Resistin is an adipocyte-derived peptide first identified during a search for targets of thiazolidinediones. It has been shown to affect various physiological processes, including insulin response and inflammatory processes within the human body.

Contents and storage

The components included in the ELISA kit are listed below. Upon receipt, store the kit at 2°C to 8°C.

Components	Cat. no. KHP0051 96 tests
Hu Resistin Antibody Coated Wells. 96 well strip plate.	1 plate
Hu Resistin Standard. Lyophilized	1 vial
Wash Buffer 10X.	2 × 30 mL
ELISA Buffer 10X.	2 × 30 mL
Detection Antibody.	30 µL
HRP Labeled Streptavidin. Lyophilized.	30 µL
TMB Substrate Solution.	12 mL
Stop Solution.	12 mL
Adhesive Plate Covers.	2



CAUTION! This kit contains materials with small quantities of Proclin™ 300. Proclin™ 300 is toxic, corrosive, and a skin irritant, so avoid ingestion and contact with eyes, skin and mucous membranes. Observe all federal, state, and local regulations for disposal.

Materials required but not provided

- Deionized water
- Microtiter plate reader with software capable of measurement at or near 450 nm
- Plate washer – automated or manual (squirt bottle, manifold dispenser, or equivalent)
- Calibrated adjustable precision pipettes and glass or plastic tubes for diluting solutions

General Guidelines

- Review the **Procedural guidelines** and **Plate washing directions** in the *ELISA Technical Guide* available at thermofisher.com for details prior to starting the procedure.
- Reagents are lot-specific. Do not mix or interchange different reagent lots from various kit lots.
- **Allow reagents to reach room temperature before use.** Mix to redissolve any precipitated salts.

For Research Use Only. Not for use in diagnostic procedures.

Manufacturing Site • 7335 Executive Way • Frederick • MD 21704 • E-mail: techsupport@lifetech.com

Prepare 1X Wash Buffer

1. Dilute 10X Wash Buffer 1:10 (e.g., 50 mL of 10X Wash Buffer with 450 mL of deionized water). Label as 1X Wash Buffer.
2. Store the concentrate and 1X Wash Buffer in the refrigerator. Use the diluted buffer within 14 days.

Sample preparation guidelines

- Collect samples in pyrogen/endotoxin-free tubes. Thirty minutes after collection, centrifuge for 15 minutes at $1,000 \times g$.
- Freeze samples in aliquots after collection if samples will not be tested immediately. Avoid multiple freeze-thaw cycles of frozen samples. Thaw completely and mix well (do not vortex) prior to analysis.
- Avoid the use of hemolyzed or lipemic sera. If large amounts of particulate matter are present in the sample, centrifuge or filter sample prior to analysis.

Dilute samples

Dilute **serum** and **plasma** samples at least 20-fold in 1X ELISA Buffer.

Because conditions may vary, it is recommended that each investigator determine the optimal dilution to be used for each application.


Prepare 1X ELISA Buffer

Dilute 10X ELISA Buffer 1:10 (e.g., 20 mL 10X ELISA Buffer with 180 mL of deionized water). Label as 1X ELISA Buffer.

Dilute standards

Note: Use glass or plastic tubes for diluting standards.

1. Reconstitute Hu Resistin Standard with 1 mL of deionized water. Swirl or mix gently and allow the contents to sit for a minimum of 15 minutes to ensure complete reconstitution. Label as 16 ng/mL Hu resistin.
2. Add 300 μ L 1X ELISA Buffer to each of 8 tubes labeled as follows: 8, 4, 2, 1, 0.5, 0.25, 0.125 and 0 ng/mL Hu resistin.
3. Mix the reconstituted standard well and make serial dilutions of the standard as shown in the dilution diagram below. Mix thoroughly between steps.
4. Discard any remaining reconstituted standard.



Prepare 1X Detection Antibody solution

Dilute Detection Antibody 1:2,000 in 1X ELISA Buffer (e.g., 5 μ L Detection Antibody with 10 mL 1X ELISA Buffer). Label as 1X Detector Antibody.

Diluted Detection Antibody is not stable and cannot be stored.

Prepare 1X HRP solution

1. Reconstitute HRP Labeled Streptavidin with 100 μ L 1X ELISA Buffer. After reconstitution, aliquots can be stored at -20°C . Avoid freeze/thaw cycles.
2. Dilute reconstituted HRP Labeled Streptavidin to the working concentration by adding 25 μ L in 10 mL (1:400) of 1X ELISA Buffer.


Prepare within 15 minutes of use. Diluted HRP-labeled Streptavidin is not stable and cannot be stored.

ELISA procedure

Allow all reagents to reach room temperature before use. Mix all liquid reagents prior to use. **Total assay time is 3.2 hours.**

IMPORTANT! Perform a standard curve with each assay.

Determine the number of 16-well strips required for the assay. Insert the strips in the frames for use. Re-bag any unused strips and frames, and store at 4°C for up to 1 month.




Bind antigen

1. Add 100 µL of the different standards in duplicate to the wells. Leave wells for chromogen blanks empty.
2. Add 100 µL of diluted serum, plasma, or cell culture supernatant samples (see page 2) to the appropriate wells.
3. Cover the plate with plate cover and incubate for 1 hour at 37°C.
4. Thoroughly aspirate the solution and wash wells 3 times with 1X Wash Buffer (300 µL/wash).




Add detector antibody

5. Add 100 µL Detection Antibody solution into each well except the chromogen blanks.
6. Cover the plate with plate cover and incubate for 1 hour at 37°C.
7. Thoroughly aspirate the solution and wash wells 3 times with 1X Wash Buffer (300 µL/wash).




Add 1X HRP Labeled Streptavidin

8. Add 100 µL 1X HRP labeled Streptavidin solution into each well.
9. Cover the plate with plate cover and incubate for 1 hour at 37°C.
10. Thoroughly aspirate the solution and wash wells 5 times with 1X Wash Buffer (300 µL/wash).



Add TMB substrate

11. Add 100 µL TMB Substrate Solution to each well. The substrate solution will begin to turn blue.
12. Incubate for 10 minutes at room temperature **in the dark**.
Note: TMB should not touch aluminum foil or other metals.



Add stop solution

13. Add 100 µL Stop Solution to each well. Tap side of the plate gently to mix. The solution in the wells changes from blue to yellow.



Read the plate and generate the standard curve

1. Read the absorbance at 450 nm. Read the plate within 30 minutes after adding the Stop Solution.
2. Use curve-fitting software to generate the standard curve. A four parameter algorithm provides the best standard curve fit. Optimally, the background absorbance may be subtracted from all data points, including standards, unknowns and controls, prior to plotting.
3. Read the concentrations for unknown samples and controls from the standard curve. Multiply value(s) obtained for sample(s) by the appropriate factor to correct for the sample dilution.


Note: Dilute samples producing signals greater than that of the highest standard in 1X ELISA Buffer and reanalyze. Multiply the concentration by the appropriate dilution factor.

Performance characteristics

Standard curve (example)

Typical standard curve over the range of 0–8 ng/mL Hu resistin.

Standard Hu resistin (ng/mL)	Optical Density (450 nm)
8	2.82
4	1.91
2	1.27
1	0.70
0.5	0.34
0.25	0.23
0.125	0.13
0	0



Specificity

This ELISA is specific for the measurement of natural and recombinant human resistin. It does not cross-react with **Human** leptin, RELM-β or adiponectin; **Mouse** resistin, RELM-α or RELM-β; **Rat** resistin or RELM-α.

Analytical Sensitivity

The minimum detectable dose of Hu resistin is 100 pg/mL

Performance characteristics, continued

Intra-assay precision

Five samples of known Hu resistin concentration were assayed in replicates of 10 to determine precision within an assay.

Sample	Average (µg/mL)	SD	%CV
1	10.78	0.31	2.86
2	19.23	0.99	5.17
3	21.49	0.67	3.12
4	5.19	0.20	3.77
5	12.57	0.47	3.73

SD = Standard Deviation; CV = Coefficient of Variation

Recovery

When samples (serum or plasma) are spiked with known concentrations of Hu resistin, the recovery averages 96% (range 93-108%).

Type	Sample	Average % Recovery	% Range
Serum	1	94.4	93-96
	2	96.6	95-97
	3	98.9	92-108

Expected values

Resistin levels in plasma and serum range from 1 to >20 ng/mL (from healthy donors).

Inter-assay precision

Five samples of known Hu resistin concentration were assayed in replicates of 10 to determine precision between assays.

Sample	Average (µg/mL)	SD	%CV
1	6.80	0.49	7.20
2	22.96	1.24	5.40
3	6.49	0.27	4.20
4	15.32	1.07	6.97
5	25.66	1.11	4.35

SD = Standard Deviation; CV = Coefficient of Variation

Linearity of Dilution








Different human serum samples containing resistin were diluted several fold (1:10 to 1:40).

Samples	Sample Dilution	Expected (ng/mL)	Observed (ng/mL)	% of Expected
1	1:10	5.59	5.59	100
	1:20	2.80	2.63	94.15
	1:40	1.40	1.23	88.29
2	1:10	6.13	6.13	100
	1:20	3.06	3.71	121.11
	1:40	1.53	1.80	117.13
3	1:10	6.33	6.33	100
	1:20	3.16	3.66	115.66
	1:40	1.58	1.63	102.90

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

Product label explanation of symbols and warnings

 Catalog Number	 Batch code	 Temperature limitation	 Use by	 Manufacturer	 Consult instructions for use	 Caution, consult accompanying documents
--------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------

DISCLAIMER: LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Corporate entity: Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For support visit thermofisher.com/support or contact techsupport@lifetech.com.

thermofisher.com

26 April 2016

ThermoFisher
SCIENTIFIC