



ELISA Kit
Catalog #KHO0741

Cleaved PARP
[214/215]

www.invitrogen.com

Invitrogen Corporation

542 Flynn Road, Camarillo, CA 93012

Tel: 800-955-6288

E-mail: techsupport@invitrogen.com

TABLE OF CONTENTS

Introduction	4
Principle of the Method.....	5
Reagents Provided.....	6
Supplies Required But Not Provided	7
Procedural Notes/Lab Quality Control.....	7
Safety	9
Directions for Washing.....	9
Procedure for Extraction of Proteins from Cells.....	10
Reagent Preparation and Storage.....	12
Reconstitution and Dilution of Cleaved PARP [214/215] Standard	12
Dilution of Cleaved PARP [214/215] Standard.....	13
Storage and Final Dilution of Anti-Rabbit IgG HRP (100X).....	14
Dilution of Wash Buffer	15
Assay Method.....	15
Typical Data	18
Limitations of the Procedure.....	19
Performance Characteristics	20
Sensitivity	20
Precision	21
Recovery.....	22
Parallelism	22
Linearity of Dilution	23
Specificity.....	24
References	26

INTRODUCTION

Poly (ADP-Ribose) Polymerase (PARP, MW: 116 kDa) is one of the most abundant proteins in the nucleus. It catalyzes the polymerization of ADP-ribose units from donor NAD^+ molecules on target nuclear proteins, resulting in the attachment of linear or branched polymers. This 116 kDa protein consists of three main domains: an amino N-terminal DNA-binding domain (DBD), an automodification domain and a carboxy C-terminal catalytic domain.

PARP plays diverse roles in many molecular and cellular processes, including DNA damage detection and repair, chromatin modification, transcription, and cell death pathways. These processes are critical for many physiological and pathological outcomes, including genome maintenance, carcinogenesis, aging, inflammation, and neuronal function.

In particular, PARP has been implicated in both apoptosis and necrosis. PARP is recognized by cell biologists as the ‘death substrate’, because PARP was one of the first identified substrates of caspases. During apoptosis, caspase 7 and caspase 3 cleave PARP between Asp214 and Gly215, which results in two fragments: p85 and p25. The cleavage of PARP separates its DBD from its catalytic domain and inactivates the enzyme. This process eliminates PARP activation in response to DNA fragmentation during apoptosis, which prevents its futile attempts at DNA repair and subsequent ATP depletion and necrotic cell death. Therefore, PARP cleavage helps to commit cells to the apoptotic pathway and is generally viewed as a hallmark for apoptosis.

The Invitrogen Cleaved PARP [214/215] ELISA is designed to detect and quantify the levels of PARP p85 fragment, which is produced from

full-length PARP by cleavage between Asp214 and Gly215 during apoptosis. This assay is intended to detect cleaved PARP [214/215] from human cell lysates or tissue homogenates. It doesn't recognize mouse or rat cleaved PARP [214/215].

For Research Use Only. CAUTION: Not for human or animal therapeutic or diagnostic use.

READ ENTIRE PROTOCOL BEFORE USE

PRINCIPLE OF THE METHOD

The Invitrogen Cleaved PARP [214/215] kit is a solid phase sandwich Enzyme Linked-Immuno-Sorbent Assay (ELISA). A monoclonal antibody specific for PARP (116 kDa and 85 kDa) has been coated onto the wells of the microtiter strips provided. Samples, including a standard containing cleaved PARP [214/215] (p85 fragment), control specimens, and unknowns, are pipetted into these wells together with a rabbit (detection) antibody specific for cleaved PARP [214/215]. During the first incubation, the cleaved PARP [214/215] antigen binds simultaneously to the immobilized (capture) antibody and to the solution phase rabbit (detection) antibody. After washing, a horseradish peroxidase-labeled anti-rabbit IgG (anti-rabbit IgG-HRP) is added. This binds to the detection antibody to complete the four-member sandwich. After the second incubation and washing to remove all the excess anti-rabbit IgG-HRP, a substrate solution is added, which is acted upon by the bound enzyme to produce color. The intensity of this colored product is directly proportional to the concentration of cleaved PARP [214/215] present in the original specimen.

REAGENTS PROVIDED

Note: Store all reagents at 2 to 8°C.

<i>Reagent</i>	<i>96 Test Kit</i>
<i>Cleaved PARP [214/215] Standard.</i> Contains 0.1% sodium azide. Refer to vial label for quantity and reconstitution volume.	2 vials
<i>Standard Diluent Buffer.</i> Contains 0.1% sodium azide; red dye*; 25 mL per bottle.	1 bottle
<i>Cleaved PARP Antibody Coated Wells, 12 x 8 Well Strips.</i>	1 plate
<i>Cleaved PARP [214/215] Detection Antibody.</i> Contains 0.1% sodium azide; blue dye*; 6.0 mL per bottle.	1 bottle
<i>Anti-Rabbit IgG HRP (100X).</i> Contains 3.3 mM thymol; 0.125 mL per vial.	1 vial
<i>HRP Diluent.</i> Contains 3.3 mM thymol; yellow dye*; 25 mL per bottle.	1 bottle
<i>Wash Buffer Concentrate (25X);</i> 100 mL per bottle.	1 bottle
<i>Stabilized Chromogen, Tetramethylbenzidine (TMB);</i> 25 mL per bottle.	1 bottle
<i>Stop Solution;</i> 25 mL per bottle.	1 bottle
<i>Plate Covers, adhesive strips.</i>	3
* In order to help our customers avoid any mistakes in pipetting the ELISAs, we provide colored <i>Standard Diluent Buffer</i> , <i>Detection Antibody</i> , and <i>HRP Diluent</i> to help monitor the addition of solutions to the reaction wells. This does not in any way interfere with the test results.	

Disposal Note: This kit contains materials with small quantities of sodium azide. Sodium azide reacts with lead and copper plumbing to form explosive metal azides. Upon disposal, flush drains with a large volume of water to prevent azide accumulation. Avoid ingestion and contact with eyes, skin and mucous membranes. In case of contact, rinse affected area with plenty of water. Observe all federal, state and local regulations for disposal.

SUPPLIES REQUIRED BUT NOT PROVIDED

1. Microtiter plate reader capable of measurement at or near 450 nm.
2. Calibrated adjustable precision pipettes, preferably with disposable plastic tips. (A manifold multi-channel pipette is desirable for large assays.)
3. Cell extraction buffer (see Recommended Formulation, p. 10).
4. Distilled or deionized water.
5. Plate washer: automated or manual (squirt bottle, manifold dispenser, etc.).
6. Data analysis and graphing software.
7. Glass or plastic tubes for diluting and aliquoting standard.
8. Absorbent paper towels.
9. Calibrated beakers and graduated cylinders in various sizes.
10. Orbital shaker with maximum speed up to 600 rpm.

PROCEDURAL NOTES/LAB QUALITY CONTROL

1. When not in use, kit components should be refrigerated. All reagents should be warmed to room temperature before use.
2. **Microtiter plates should be allowed to come to room temperature before opening the foil bags.** Once the desired number of strips has been removed, immediately reseal the bag and store at 2 to 8°C to maintain plate integrity.

3. Samples should be frozen if not analyzed shortly after collection. Avoid multiple freeze-thaw cycles of frozen samples. Thaw completely and mix well prior to analysis.
4. If particulate matter is present, centrifuge or filter prior to analysis.
5. All standards, controls and samples should be run in duplicate.
6. Samples that are greater than the highest standard point should be diluted with *Standard Diluent Buffer* and retested.
7. When pipetting reagents, maintain a consistent order of addition from well-to-well. This ensures equal incubation times for all wells.
8. Cover or cap all reagents when not in use.
9. **Do not mix or interchange different reagent lots from various kit lots.**
10. Do not use reagents after the kit expiration date.
11. Read absorbances within 2 hours of assay completion.
12. In-house controls should be run with every assay. If control values fall outside pre-established ranges, the accuracy of the assay is suspect.
13. All residual wash liquid must be drained from the wells by efficient aspiration or by decantation followed by tapping the plate forcefully on absorbent paper. *Never* insert absorbent paper directly into the wells.
14. Because *Stabilized Chromogen* is light sensitive, avoid prolonged exposure to light. Also avoid contact between *Stabilized Chromogen* and metal, or color may develop.

SAFETY

All blood components and biological materials should be handled as potentially hazardous. Follow universal precautions as established by the Centers for Disease Control and Prevention and by the Occupational Safety and Health Administration when handling and disposing of infectious agents.

DIRECTIONS FOR WASHING

Incomplete washing will adversely affect the test outcome. All washing must be performed with *Wash Buffer Concentrate (25X)* provided.

Washing can be performed manually as follows: completely aspirate the liquid from all wells by gently lowering an aspiration tip (aspiration device) into the bottom of each well. Take care not to scratch the inside of the well.

After aspiration, fill the wells with at least 0.4 mL of diluted wash solution. Let soak for 15 to 30 seconds, then aspirate the liquid. Repeat as directed under **ASSAY METHOD**. After the washing procedure, the plate is inverted and tapped dry on absorbent tissue.

Alternatively, the wash solution may be put into a squirt bottle. If a squirt bottle is used, flood the plate with wash buffer, completely filling all wells. After the washing procedure, the plate is inverted and tapped dry on absorbent tissue.

If using an automated washer, the operating instructions for washing equipment should be carefully followed. If your automated washer allows, 30 second soak cycles should be programmed into the wash cycle.

PROCEDURE FOR EXTRACTION OF PROTEINS FROM CELLS

A. Recommended Formulation of Cell Extraction Buffer:

10 mM Tris, pH 7.4

100 mM NaCl

1 mM EDTA

1 mM EGTA

1 mM NaF

20 mM $\text{Na}_4\text{P}_2\text{O}_7$

2 mM Na_3VO_4

1% Triton X-100

10% glycerol

0.1% SDS

0.5% deoxycholate

1 mM PMSF (stock is 0.3 M in DMSO)

Protease inhibitor cocktail (e.g., Sigma Cat. # P-2714; reconstituted according to manufacturer's guideline). Add 100 μL per 1 mL Cell Extraction Buffer.

This buffer is stable for 2 to 3 weeks at 4°C or for up to 6 months when aliquoted (without protease inhibitors and PMSF added) and stored at -20°C. When stored frozen, the Cell Extraction Buffer should be thawed on ice. This buffer (minus protease inhibitor cocktail and PMSF) can be obtained from Invitrogen Cat. # FNN0011. **Important:** add the protease inhibitors just before using. The stability of protease inhibitor supplemented Cell Extraction Buffer is 24 hours at 4°C. PMSF is very unstable and must be added prior to use, even if added previously.

B. Protocol for Cell Extraction

This protocol has been successfully applied to several cell lines. Researchers should optimize the cell extraction procedures for their own applications.

1. Collect cells in PBS by centrifugation (non-adherent) or scraping from culture flasks (adherent).
2. Wash cells twice with cold PBS.
3. Remove and discard the supernatant and collect the cell pellet. (At this point the cell pellet can be frozen at -80°C and lysed at a later date).
4. Lyse the cell pellet in Cell Extraction Buffer for 30 minutes, on ice, with vortexing at 10 minute intervals. The volume of Cell Extraction Buffer depends on the cell number in cell pellet and expression of PARP. For example, 1×10^6 Jurkat cells grown in RPMI supplemented with 10% FBS and treated with $1 \mu\text{M}$ staurosporine can be extracted in 1 mL of Extraction Buffer. Under these conditions, use of 1-10 μL of the clarified cell extract diluted to a volume of 50 μL /well in *Standard Diluent Buffer* (See **Assay Method**) is sufficient for the detection of cleaved PARP [214/215].
5. Transfer extract to microcentrifuge tubes and centrifuge at 13,000 rpm for 10 minutes at 4°C .
6. Aliquot the clear lysate to clean microfuge tubes. These samples are ready for assay. Lysates can be stored at -80°C . Avoid multiple freeze/thaw cycles.

REAGENT PREPARATION AND STORAGE

A. Reconstitution and Dilution of Cleaved PARP [214/215] Standard

Note: This *Cleaved PARP [214/215] Standard* is prepared from *in vitro* caspase 3-cleaved PARP recombinant protein.

1. Reconstitute *Cleaved PARP [214/215] Standard* with *Standard Diluent Buffer*. Refer to standard vial label for instructions. Swirl or mix gently and allow to sit for 10 minutes to ensure complete reconstitution. Label as 10 ng/mL cleaved PARP [214/215]. Use the standard within 1 hour of reconstitution.
2. Add 0.15 mL of *Standard Diluent Buffer* to each of 6 tubes labeled 5, 2.5, 1.25, 0.625, 0.312 and 0.156 ng/mL cleaved PARP [214/215].
3. Make serial dilutions of the standard as described in the following dilution table. Mix thoroughly between steps.

B. Dilution of Cleaved PARP [214/215] Standard

Standard:	Add:	Into:
10 ng/mL	Prepare as described in step 1	
5 ng/mL	0.15 mL of the 10 ng/mL std.	0.15 mL of the Diluent Buffer
2.5 ng/mL	0.15 mL of the 5 ng/mL std.	0.15 mL of the Diluent Buffer
1.25 ng/mL	0.15 mL of the 2.5 ng/mL std.	0.15 mL of the Diluent Buffer
0.625 ng/mL	0.15 mL of the 1.25 ng/mL std.	0.15 mL of the Diluent Buffer
0.312 ng/mL	0.15 mL of the 0.625 ng/mL std.	0.15 mL of the Diluent Buffer
0.156 ng/mL	0.15 mL of the 0.312 ng/mL std.	0.15 mL of the Diluent Buffer
0 ng/mL	0.15 mL of the Diluent Buffer	An empty tube

Remaining reconstituted standard should be discarded or frozen in aliquots at -80°C for further use. Standard can be frozen and thawed one time only without loss of immunoreactivity.

C. Storage and Final Dilution of Anti-Rabbit IgG HRP (100X)

Please Note: The *Anti-Rabbit IgG HRP (100X)* is in 50% glycerol. This solution is viscous. To ensure accurate dilution, allow *Anti-Rabbit IgG HRP (100X)* to reach room temperature. Gently mix. Pipette *Anti-Rabbit IgG HRP (100X)* slowly. Remove excess concentrate solution from pipette tip by gently wiping with clean absorbent paper.

1. Dilute 10 μL of this 100x concentrated solution with 1 mL of *HRP Diluent* for each 8-well strip used in the assay. Label as Anti-Rabbit IgG HRP Working Solution.

For Example:

# of 8-Well Strips	Volume of Anti-Rabbit IgG HRP (100X)	Volume of Diluent
2	20 μL solution	2 mL
4	40 μL solution	4 mL
6	60 μL solution	6 mL
8	80 μL solution	8 mL
10	100 μL solution	10 mL
12	120 μL solution	12 mL

2. Return the unused *Anti-Rabbit IgG HRP (100X)* to the refrigerator.

D. Dilution of Wash Buffer

Allow the *Wash Buffer Concentrate (25X)* to reach room temperature and mix to ensure that any precipitated salts have redissolved. Dilute 1 volume of the *Wash Buffer Concentrate (25X)* with 24 volumes of deionized water (e.g., 50 mL may be diluted up to 1.25 liters, 100 mL may be diluted up to 2.5 liters). Label as Working Wash Buffer.

Store both the concentrate and the Working Wash Buffer in the refrigerator. The diluted buffer should be used within 14 days.

ASSAY METHOD: PROCEDURE AND CALCULATIONS

Be sure to read the *Procedural Notes/Lab Quality Control* section before carrying out the assay.

Allow all reagents to reach room temperature before use. Gently mix all liquid reagents prior to use.

Note: A standard curve must be run with each assay.

1. Determine the number of 8-well strips needed for the assay. Insert these in the frame(s) for current use. (Re-bag extra strips and frame. Store these in the refrigerator for future use.)
2. Add 50 μL of the *Standard Diluent Buffer* to zero wells. Well(s) reserved for chromogen blank should be left empty.
3. Add 50 μL of standards, samples or controls to the appropriate microtiter wells. Tap gently on side of plate to mix. Samples prepared in Cell Extraction Buffer must be diluted 1:5 or greater in *Standard Diluent Buffer* (for example, 10 μL sample into 40 μL buffer). While a 1:5 sample dilution has been found to be satisfactory, higher dilutions such as 1:10 or 1:20 may be optimal.

The dilution chosen should be optimized for each experimental system. (See **REAGENT PREPARATION AND STORAGE**, Section B.)

4. Pipette 50 μL *Cleaved PARP [214/215] Detection Antibody* solution into each well except the chromogen blank(s). Tap gently on the side of the plate to thoroughly mix.
5. Cover wells with *plate cover* and incubate for **3 hours at room temperature on an orbital shaker, shaking at 400-600 rpm.**
6. Thoroughly aspirate or decant solution from wells and discard the liquid. Wash wells 4 times. See **DIRECTIONS FOR WASHING.**
7. Add 100 μL *Anti-Rabbit IgG HRP Working Solution* to each well except the chromogen blank(s). (Prepare the working dilution as described in **REAGENT PREPARATION AND STORAGE**, Section C.)
8. Cover wells with the *plate cover* and incubate for **30 minutes at room temperature.**
9. Thoroughly aspirate or decant solution from wells and discard the liquid. Wash wells 4 times. See **DIRECTIONS FOR WASHING.**
10. Add 100 μL of *Stabilized Chromogen* to each well. The liquid in the wells will begin to turn blue.
11. Incubate for **30 minutes at room temperature and in the dark.** **Please Note: Do not cover the plate with aluminum foil or metalized mylar.** The incubation time for chromogen substrate is often determined by the microtiter plate reader used. Many plate readers have the capacity to record a maximum optical density (O.D.) of 2.0. The O.D. values should be monitored and the substrate reaction stopped before the O.D. of the positive wells

exceed the limits of the instrument. The O.D. values at 450 nm can only be read after the *Stop Solution* has been added to each well. If using a reader that records only to 2.0 O.D., stopping the assay after 20 to 25 minutes is suggested.

12. Add 100 μL of *Stop Solution* to each well. Tap side of plate gently to mix. The solution in the wells should change from blue to yellow.
13. Read the absorbance of each well at 450 nm having blanked the plate reader against a chromogen blank composed of 100 μL each of *Stabilized Chromogen* and *Stop Solution*. Read the plate within 2 hours after adding the *Stop Solution*.
14. Plot the absorbances of the standards against the standard concentrations. (Optimally, the background absorbance may be subtracted from *all* data points, including standards, unknowns and controls, prior to plotting.) Draw the best smooth curve through these points to construct the standard curve. If using curve fitting software, the four parameter algorithm provides the best curve fit.
15. Read the cleaved PARP [214/215] concentrations for unknown samples and controls from the standard curve plotted in step 14. **Multiply value(s) obtained for sample(s) by the appropriate dilution factor to correct for the dilution with *Standard Diluent Buffer*.** (Samples producing signals higher than the highest standard [10 ng/mL] should be further diluted in *Standard Diluent Buffer* and reanalyzed, multiplying the concentration by the appropriate dilution factor.)

TYPICAL DATA

The following data were obtained for the various standards over the range of 0 to 10 ng/mL cleaved PARP [214/215].

Standard Cleaved PARP [214/215] (ng/mL)	Optical Density (450 nm)
10	3.01
5	2.19
2.5	1.43
1.25	0.85
0.625	0.49
0.312	0.30
0.156	0.20
0	0.09

LIMITATIONS OF THE PROCEDURE

Do not extrapolate the standard curve beyond the 10 ng/mL standard point; the dose-response is non-linear in this region and accuracy is difficult to obtain. Dilute samples >10 ng/mL with *Standard Diluent Buffer*; reanalyze these and multiply results by the appropriate dilution factor.

The influence of various extraction buffers has not been thoroughly investigated. The rate of degradation of native PARP in various matrices has not been investigated.

For Research Use Only. CAUTION: Not for human or animal therapeutic or diagnostic use.


PERFORMANCE CHARACTERISTICS

SENSITIVITY

The analytical sensitivity of this assay is <0.062 ng/mL of cleaved PARP [214/215]. This was determined by adding two standard deviations to the mean O.D. obtained when the zero standard was assayed 30 times. Using Jurkat cells treated with 1 μ M staurosporine for 3 hours, this level of sensitivity was equivalent to the detection of cleaved PARP [214/215] in 50 cells.

The sensitivity of this ELISA was compared to Western blotting using known quantities of cleaved PARP [214/215]. The data presented in Figure 1 show that the sensitivity of the ELISA is about 10x greater than that of Western blotting. The bands shown in the Western blotting data were developed using rabbit anti-cleaved PARP [214/215], and an alkaline phosphatase conjugated anti-rabbit IgG followed by chemiluminescent substrate and autoradiography. A different amount of cell lysate was used in Western blotting and ELISA, due to the much higher sensitivity of ELISA.

Figure 1: Detection of Cleaved PARP [214/215] by ELISA vs Western Blot: Jurkat Cells Treated with Staurosporine



PRECISION

1. Intra-Assay Precision

Samples of known cleaved PARP [214/215] concentrations were assayed in replicates of 16 to determine precision within an assay.

	Sample 1	Sample 2	Sample 3
Mean (ng/mL)	5.09	2.14	0.92
SD	0.26	0.078	0.046
%CV	5.11	3.64	5.00

SD = Standard Deviation
CV = Coefficient of Variation

2. Inter-Assay Precision

Samples were assayed 48 times in multiple assays to determine precision between assays.

	Sample 1	Sample 2	Sample 3
Mean (ng/mL)	5.09	2.12	0.91
SD	0.25	0.09	0.062
%CV	4.91	4.25	6.81

SD = Standard Deviation
CV = Coefficient of Variation


RECOVERY

To evaluate recovery, *Cleaved PARP [214/215] Standard* was spiked at 3 different concentrations into 20% Cell Extraction Buffer. The average recovery was 105%.

PARALLELISM

Natural cleaved PARP [214/215] from staurosporine-induced apoptotic Jurkat cell lysate was serially diluted in *Standard Diluent Buffer*. The optical density of each dilution was plotted against the cleaved PARP [214/215] standard curve. Parallelism demonstrated by the figure below indicates that the standard accurately reflects cleaved PARP [214/215] content in samples.

Figure 2



LINEARITY OF DILUTION

Jurkat cells were grown in RPMI with 10% FBS, treated with 1 μ M staurosporine for 3 hours and lysed with Cell Extraction Buffer. This lysate was diluted in *Standard Diluent Buffer* over the range of the assay and measured for cleaved PARP [214/215]. Linear regression analysis of sample values versus the expected concentration yielded a correlation coefficient of 0.99.

Dilution	Cell Lysate		
	Measured (ng/mL)	Expected (ng/mL)	% Expected
Neat	2.83	2.83	100
1/2	1.31	1.41	92.9
1/4	0.65	0.71	91.5

SPECIFICITY

The Invitrogen Cleaved PARP [214/215] ELISA specifically recognizes human cleaved PARP p85 fragment. It does not recognize full-length or p25 fragment of PARP protein. It does not detect mouse or rat cleaved PARP. In Figure 3, Jurkat cells were treated with staurosporine (Invitrogen Cat. # PHZ1271) at 1 μ M for 3 hours. Untreated cells were used as controls. Cell extracts were prepared and 5 μ g/mL of total cell lysates were analyzed with the Cleaved PARP [214/215] ELISA. The results show that the level of cleaved PARP [214/215] is highly elevated in staurosporine-induced apoptotic cells.

Figure 3





Figure 4 shows that PARP cleavage is the hallmark of apoptosis, which is dependent on the dose of staurosporine treatment. In this experiment, Jurkat cells were treated with staurosporine at varying concentrations (0 - 1,000 nM) for 3 hours. Cell extracts were prepared and 5 $\mu\text{g/mL}$ of total cell lysates were assayed for the levels of cleaved PARP [214/215]. PARP cleavage between Asp214 and Gly215 increases with increasing dose of staurosporine.

Figure 4
















REFERENCES

1. Lazebnik, Y.A., et al. (1994) Cleavage of poly(ADP-ribose) polymerase by a proteinase with properties like ICE. *Nature* 371 (6495):346-347.
2. Tewari, M., et al. (1995) Yama/CPP32 beta, a mammalian homolog of CED-3, is a CrmA-inhibitable protease that cleaves the death substrate poly(ADP-ribose) polymerase. *Cell* 81(5):801-809.
3. Chiarugi, A., et al. (2002) Cell biology. PARP-1--a perpetrator of apoptotic cell death? *Science* 297(5579):200-201.
4. Schreiber, V., et al. (2006) Poly(ADP-ribose): novel functions for an old molecule. *Nat. Rev. Mol. Cell Biol.* 7(7):517-528.
5. Soldani, C., et al. (2001) Poly(ADP-ribose) polymerase cleavage during apoptosis: when and where? *Exp. Cell Res.* 269(2):193-201.
6. Soldani, C., et al. (2002) Poly(ADP-ribose) polymerase-1 cleavage during apoptosis: an update. *Apoptosis* 7(4):321-328.
7. Farmer, H., et al. (2006) Targeting the DNA repair defect in BRCA mutant cells as a therapeutic strategy. *Nature* 434(7035):917-921.
8. Bryant, H.E., et al. (2006) Specific killing of BRCA2-deficient tumors with inhibitors of poly(ADP-ribose) polymerase. *Nature* 434(7035):913-917.
9. Kim, M.Y., et al. (2004) NAD⁺-dependent modulation of chromatin structure and transcription by nucleosome binding properties of PARP-1. *Cell* 119(6):803-814.
10. Lis, J.T., et al. (2006) Promoter cleavage: a topoIIbeta and PARP-1 collaboration. *Cell* 125(7):1225-1227.

11. Midorikawa, R., et al. (2006) KIF4 motor regulates activity-dependent neuronal survival by suppressing PARP-1 enzymatic activity. *Cell* 125(2):371-383.

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

Explanation of symbols

Symbol	Description	Symbol	Description
	Catalogue Number		Batch code
	Research Use Only		<i>In vitro</i> diagnostic medical device
	Use by		Temperature limitation
	Manufacturer		European Community authorised representative
	Without, does not contain		With, contains
	Protect from light		Consult accompanying documents
	Directs the user to consult instructions for use (IFU), accompanying the product.		

Copyright © Invitrogen Corporation. 25 May 2010

Cleaved PARP [214/215] Assay Summary

