

## Contents

1. Description
  - 1.1 Background information
  - 1.2 Applications
  - 1.3 Reagent requirements
2. General protocol for immunofluorescent staining
3. Examples of immunofluorescent staining with Annexin V fluorochrome and conjugates
4. References

## Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

## 1. Description

This product is for research use only.

**Components** Annexin V conjugated to:

Conjugate	Order no. 1 mL (100 tests)	Order no. 300 µL (30 tests)
FITC	130-093-060	–
PE	130-108-077	130-108-112
Biotin	130-092-773	–

**Capacity** 1 mL: 100 tests or up to  $10^8$  total cells  
300 µL: 30 tests or up to  $3 \times 10^7$  total cells.

**Product format** Annexin V conjugates are supplied in buffer containing stabilizer and 0.05% sodium azide.

**Storage** Store protected from light at 2–8 °C. Do not freeze. The expiration date is indicated on the vial label.


### 1.1 Background information

In most normal, viable eukaryotic cells, the negatively charged phospholipid phosphatidylserine (PS) is located in the cytosolic leaflet of the plasma membrane lipid bilayer.<sup>1</sup> PS redistribution from the inner to the outer leaflet is an early and widespread event during apoptosis.<sup>1,2</sup> However, in necrosis, PS becomes accessible due to the disruption of membrane integrity.<sup>2</sup> Apart from necrosis and apoptosis, PS also becomes accessible in activated platelets<sup>3</sup>, in certain cell anomalies such as sickle cell anemia<sup>4</sup>, in erythrocyte senescence<sup>5</sup>, upon degranulation of mast cells<sup>6</sup>, and in certain stages of B cell differentiation<sup>7</sup>. PS exposure also serves as a trigger for the recognition and removal of apoptotic cells by macrophages.<sup>8,9</sup> Annexin V is a 36 kDa phospholipid-binding protein and has a high

affinity to PS in the presence of physiological concentrations of calcium ( $\text{Ca}^{2+}$ ).<sup>10</sup>

MACS® Annexin V fluorochrome and biotin conjugates have been developed for the detection and discrimination of apoptotic and dead cells. Apoptotic cells, which are otherwise undetectable by staining with propidium iodide (PI) or 4',6-diamidino-2-phenylindole (DAPI), can be directly detected through their staining with fluorochrome-conjugated Annexin V. Dead cells are stained with both Annexin V and PI or DAPI, whereas viable cells cannot be stained with either.

## Staining procedure



### 1.2 Applications

- Studies on cell death (apoptosis and/or necrosis).
- Evaluation of MACS Separations with the Annexin V MicroBead Kit (# 130-090-201) and the Dead Cell Removal Kit (# 130-090-101).

### 1.3 Reagent requirements

- **Buffer:** Prepare 1× Annexin V Binding Buffer from the Annexin V Binding Buffer 20× Stock Solution (# 130-092-820): For  $10^6$  total cells, dilute 500 µL of the Annexin V Binding Buffer (20× Stock Solution) with 9.5 mL of sterile, distilled water. Alternatively, prepare 1× Annexin V Binding Buffer by diluting 25 mL of the 20× Stock Solution with 475 mL of sterile, distilled water. Store at 2–8 °C.
- ▲ **Note:** Handle under sterile conditions.
- (Optional) Conjugated anti-biotin antibodies, e.g., Anti-Biotin-PE (# 130-090-756) as secondary antibody reagent in combination with Annexin V-Biotin.
- (Optional) Propidium Iodide Solution (# 130-093-233), DAPI, or 7-AAD for flow cytometric exclusion of dead cells without cell fixation. For cell fixation and flow cytometric exclusion of dead cells, the Fixation and Dead Cell Discrimination Kit (# 130-091-163) is recommended.

## 2. General protocol for immunofluorescent staining

▲ Volumes for fluorescent labeling given below are for **up to**  $10^6$  nucleated cells. When working with fewer than  $10^6$  cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes, accordingly (e.g. for  $2 \times 10^6$  nucleated cells, use twice the volume of all indicated reagent volumes and total volumes).


▲ When working with cell samples containing platelets (e.g. blood samples), wash samples carefully at a low centrifugation speed ( $200 \times g$ ) in order to remove platelets. Use buffer containing the ion chelator EDTA for these washing steps. Finally, wash with Annexin V Binding Buffer to avoid chelation of  $Ca^{2+}$  ions. Activated platelets expose PS and therefore bind Annexin V.<sup>5</sup>

1. Wash  $10^6$  cells in 1 mL of  $1 \times$  Annexin V Binding Buffer and centrifuge at  $300 \times g$  for 10 minutes. Aspirate supernatant completely.
2. (Optional) Repeat washing step.
3. Resuspend  $10^6$  cells in 100  $\mu$ L of  $1 \times$  Annexin V Binding Buffer.
4. Add 10  $\mu$ L of Annexin V conjugate.
5. Mix well and incubate for 15 minutes in the dark at room temperature.  
▲ **Note:** Lower temperatures may require increased incubation times. Higher temperatures and/or longer incubation times may lead to non-specific cell labeling.
6. Wash cells by adding 1 mL of  $1 \times$  Annexin V Binding Buffer per  $10^6$  cells and centrifuge at  $300 \times g$  for 10 minutes. Aspirate supernatant completely.
7. (Optional) Repeat washing step.
8. (Optional) If Annexin V-Biotin was used, resuspend the cell pellet in 100  $\mu$ L of  $1 \times$  Annexin V Binding Buffer per  $10^6$  cells, add 10  $\mu$ L of anti-biotin antibody, incubate for 10 minutes in the dark in the refrigerator ( $2-8^\circ C$ ), and continue as described in step 6.
9. Resuspend cell pellet in 500  $\mu$ L of  $1 \times$  Annexin V Binding Buffer per  $10^6$  **total** cells.
10. (Optional) Add 1  $\mu$ g/mL of Propidium Iodide Solution or 0.01  $\mu$ g of DAPI immediately prior to analysis by flow cytometry or fluorescence microscopy.


## 3. Examples of immunofluorescent staining with Annexin V fluorochromes and conjugates

Jurkat cells, cultured with staurosporine (50 nM) for 15 hours, were stained with Annexin V conjugated to FITC (A) or PE (B) followed by staining with DAPI and analyzed by flow cytometry using the MACSQuant® Analyzer.


A)



B)



Jurkat cells were cultured with and without different concentrations of staurosporine (200 nM up to 2  $\mu$ M) for different times (from 4 up to 20.5 hours). Cells were then pooled and stained with Annexin V-Biotin followed by Anti-Biotin-APC and PI and analyzed by flow cytometry.



## 4. References

1. Koopman, G. *et al.* (1994) Annexin V for flow cytometric detection of phosphatidylserine expression on B cells undergoing apoptosis. *Blood* 84: 1415–1420.
2. Martin, S.J. *et al.* (1995) Early redistribution of plasma membrane phosphatidylserine is a general feature of apoptosis regardless of the initiating stimulus: inhibition by overexpression of Bcl-2 and Abl. *J. Exp. Med.* 182: 1545–1556.
3. Thiagarajan, P. and Tait, J.F. (1990) Binding of annexin V/placental anticoagulant protein I to platelets. Evidence for phosphatidylserine exposure in the procoagulant response of activated platelets. *J. Biol. Chem.* 265: 17420–17423.
4. Kuypers, F.A. *et al.* (1996) Detection of altered membrane phospholipid asymmetry in subpopulations of human red blood cells using fluorescently labeled Annexin V. *Blood* 87: 1179–1187.
5. Schroit, A.J. and Zwaal, R.F.A. (1991) Transbilayer movement of phospholipids in red cells and platelet membranes. *Acta Biochim. Biophys.* 1071: 313–329.
6. Demo, S.D. *et al.* (1999) Quantitative measurement of mast cell degranulation using a novel flow cytometric Annexin-V binding assay. *Cytometry* 36: 340–348.
7. Dillon, S.R. *et al.* (2001) Annexin V binds to positively selected B cells. *J. Immunol.* 166: 58–71.
8. Fadok, V.A. *et al.* (1992) Exposure of phosphatidylserine on the surface of apoptotic lymphocytes triggers specific recognition and removal by macrophages. *J. Immunol.* 148: 2207–2216.
9. Fadok, V.A. *et al.* (2000) A receptor for phosphatidylserine-specific clearance of apoptotic cells. *Nature* 405: 85–90.
10. Moss, S.E. *et al.* (1991) Diversity in the Annexin family. In *Novel Calcium Binding Proteins*, Springer Verlag, 535–566.

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols.

## Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

MACS and MACSQuant are registered trademarks of Miltenyi Biotec GmbH.

Copyright © 2015 Miltenyi Biotec GmbH.