

# Exosome Isolation Kits

## human

Exosome Isolation Kit CD9	130-110-913
Exosome Isolation Kit CD63	130-110-918
Exosome Isolation Kit CD81	130-110-914
Exosome Isolation Kit Pan	130-110-912

### Contents

1. Description
  - 1.1 Principle of the MACS® Separation
  - 1.2 Background information
  - 1.3 Applications
  - 1.4 Reagent and instrument requirements
2. Protocol
  - 2.1 Sample preparation
  - 2.2 Magnetic labeling
  - 2.3 Magnetic separation
  - 2.4 Elution
3. Troubleshooting
4. References

### Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

### 1. Description

**This product is for research use only.**

<b>Components</b>	<b>Exosome Isolation Kit CD9, human:</b>
	1 mL Exosome Isolation MicroBeads CD9, human
	4 mL Equilibration Buffer
	25 mL Isolation Buffer
	20 µ Columns and plungers
	or
	<b>Exosome Isolation Kit CD63, human:</b>
	1 mL Exosome Isolation MicroBeads CD63, human
	4 mL Equilibration Buffer
	25 mL Isolation Buffer
	20 µ Columns and plungers
	or
	<b>Exosome Isolation Kit CD81, human:</b>
	1 mL Exosome Isolation MicroBeads CD81, human
	4 mL Equilibration Buffer
	25 mL Isolation Buffer
	20 µ Columns and plungers
	or

### Exosome Isolation Kit Pan, human:

1 mL Exosome Isolation MicroBeads Pan, human (Cocktail of MicroBeads conjugated to CD9, CD63, and CD81)  
 4 mL Equilibration Buffer  
 25 mL Isolation Buffer  
 20 µ Columns and plungers

### Size

For 20 isolations.

### Product format

Exosome Isolation MicroBeads are supplied in buffer containing stabilizer and 0.05% sodium azide.

### Storage

Store Exosome Isolation MicroBeads, Equilibration Buffer, and Isolation Buffer protected from light at 2–8 °C. Do not freeze. Store µ Columns and plungers dry at 10–35 °C. The expiration dates are indicated on the labels.

#### 1.1 Principle of the MACS® Separation

The isolation of exosomes or extracellular vesicles (EVs) is performed by positive selection using MicroBeads recognizing the tetraspanin proteins CD9, CD63, or CD81. First, EVs are magnetically labeled during a short incubation period. The labeled EVs are loaded onto a µ Column, which is placed in the magnetic field of a µMACS™ Separator. The magnetically labeled EVs are retained within the column, while the unlabeled vesicles and cell components run through the column. After removing the column from the magnetic field, the intact EVs can either be collected by elution with Isolation Buffer, or directly lysed in the column and the protein in the lysates can be analysed, e.g., by Western blotting.

#### 1.2 Background information

Exosomes are extracellular vesicles released from living cells in an energy-dependent process. Exosomes are specifically loaded with nucleic acids, lipids, and proteins from their parental cell. Therefore, the constitution of extracellular vesicles reflects the type and status of the originating cell. Exosomes are secreted by many cell types<sup>1</sup> into diverse body fluids such as blood<sup>2</sup>, semen<sup>3</sup>, urine<sup>4</sup>, saliva<sup>5</sup>, breast milk<sup>6</sup>, ascites fluid<sup>7</sup>, and cerebrospinal fluid<sup>8</sup>.

The main difference to other extracellular vesicles such as apoptotic vesicles or membrane vesicles is the endocytic origin of exosomes. Exosomes are released from intact cells after inward budding of multivesicular bodies and fusion with the plasma membrane. They have the same membrane orientation as the originating cell, i.e., displaying extracellular domains on their surface<sup>9,10</sup>.

CD9, CD63, and CD81 are three of the most-studied members of the tetraspanin protein family and can be used to isolate exosomes.

#### 1.3 Applications

- Isolation of exosomes or extracellular vesicles from cell culture supernatant.
- Isolation of exosomes or extracellular vesicles from body fluids like plasma, urine, or ascites.

## 1.4 Reagent and instrument requirements

- $\mu$ MACS™ Separator (# 130-042-602)
- MACS® MultiStand (# 130-042-303)
- Centrifuge
- 1.5 mL tubes
- Exosome lysis and elution buffer for direct loading of exosomal proteins on sodium dodecyl sulfate (SDS) gels, e.g., 50 mM Tris-Cl, pH 6.8, with 2% (w/v) SDS, 8% (v/v) glycerin, and 0.005% (w/v) bromophenol blue.
- (Optional) Serum-free cell culture medium
- (Optional) EDTA or citrate tubes
- (Optional) Phosphate-buffered saline (PBS)
- (Optional) Filter with 0.22  $\mu$ m nylon mesh

## 2. Protocol

▲ Kit components should not be substituted or mixed with components from other kits or lots.

### 2.1 Sample preparation

The Exosome Isolation Kits can be used for pre-cleared cell culture supernatant or pre-cleared body fluids like plasma, ascites, or urine as well as to isolate exosomes carrying the respective markers from EV preparations, e.g., from ultracentrifugation or density gradient centrifugations.

#### Pre-clearing cell culture supernatant


1. Incubate the cells of interest in serum-free medium for 12–72 hours depending on the cell line. Adjust the incubation conditions to an apoptosis rate of less than 5%.
2. Remove cells, cell debris, and larger vesicles by serial centrifugations at 300 $\times$ g for 10 minutes, 2000 $\times$ g for 30 minutes, and 10,000 $\times$ g for 45 minutes. Take off the respective supernatant for the next centrifugation step.
3. (Optional) Filter the supernatant through a filter with 0.22  $\mu$ m nylon mesh.

#### Pre-clearing plasma

1. Add blood to EDTA or citrate tubes.
2. Separate plasma by centrifugation at 1,000 $\times$ g for 10 minutes.
3. Dilute plasma with an equal volume of PBS.
4. Remove cell debris and larger vesicles by serial centrifugations at 2000 $\times$ g for 30 minutes and 10,000 $\times$ g for 45 minutes. Take off the respective supernatant for the next centrifugation step.

#### Pre-clearing urine or ascites

1. Remove cells, cell debris, and larger vesicles by serial centrifugations at 300 $\times$ g for 10 minutes, 2000 $\times$ g for 30 minutes, and 10,000 $\times$ g for 45 minutes. Take off the respective supernatant for the next centrifugation step.
2. (Optional) Filter the supernatant through a filter with 0.22  $\mu$ m nylon mesh.



### 2.2 Magnetic labeling

▲ Volumes for magnetic labeling given below are for up to 2 mL exosome-containing sample, such as pre-cleared cell culture supernatant, plasma, urine, ascites, etc.

1. Add 50  $\mu$ L of Exosome Isolation MicroBeads to exosome-containing sample and vortex.
2. Incubate for 1 hour at room temperature.
3. Proceed to magnetic separation (2.3).



### 2.3 Magnetic separation

▲ Always wait until the column reservoir is empty before proceeding to the next step.

▲ Pre-warm Equilibration Buffer and Isolation Buffer to room temperature before use.

1. Place a  $\mu$  Column in the magnetic field of the  $\mu$ MACS Separator that is attached to the MACS MultiStand.
2. Prepare column by applying 100  $\mu$ L Equilibration Buffer on top of the column.
3. Rinse column with 3 $\times$ 100  $\mu$ L of Isolation Buffer. Only add new buffer when the column reservoir is empty.
4. Apply magnetically labeled sample onto the column and let it run through.
5. Wash column with 4 $\times$ 200  $\mu$ L of Isolation Buffer. Only add new buffer when the column reservoir is empty.
6. (Optional) Remove any residual drop at the column tip by touching the column tip with a pipette tip.
7. Proceed to elution of exosome (2.4).

## 2.4 Elution

Perform elution of exosomes with column outside of the magnetic separator to obtain intact vesicles.

### Elution of intact exosomes

▲ Eluted exosomes remain attached to MicroBeads and might interfere with downstream analysis like nanoparticle tracking analysis or electron microscopy as both, the exosomes and MicroBeads, will be analyzed.

1. Remove the column from the magnetic separator and place column onto a 1.5 mL tube. Add 100 µL Isolation Buffer to the column and immediately flush out the magnetically labeled vesicles by firmly pushing the plunger into the column.
2. Eluted exosomes can be used for, e.g., RNA isolation or Western blot analysis as well as for analysis of exosomes using the MACSPlex Exosome Kit.

### Elution of exosome lysate for Western blot analysis

1. Remove the column from the magnetic separator and place column onto a 1.5 mL tube. Add 100 µL Exosome lysis and elution buffer to the column and immediately flush out the vesicle lysate by firmly pushing the plunger into the column.
2. The lysed exosome samples can directly be used for further Western blot analysis.

## 3. Troubleshooting

### Slow column flow

▲ Air bubble formation within the column can impair column flow. To prevent air bubble formation, use room-temperature buffers for the wash steps or, where possible, degas the buffers before use.

### No or low exosome recovery

▲ Exosome binding to the Exosome Isolation MicroBeads mainly depend on the exosome concentration. Low recoveries can be indicative for low exosome concentration. Extended culture times could improve the yield. Alternatively, exosomes can be pre-concentrated from larger volumes, e.g., by ultracentrifugation.

▲ Vesicles might have no epitopes used for isolation. Although, most exosome carry the tetraspanin markers CD9, CD63, or CD81, some exosome populations have been shown to not equally carry all of the three markers<sup>11</sup>.

▲ Columns are still magnetized. The µ Columns must be removed from the µMACS™ Separator to enable successful elution.

### Exosomes still present in the flow through

▲ The binding capacity of the Exosome Isolation MicroBeads could be overloaded in case of very high concentration of exosomes, e.g., from cell culture medium or pre-enriched exosome samples. Repeat the experiment with diluted exosome samples.

▲ Contaminations, e.g., from cell culture medium could interfere with exosome binding to the Exosome Isolation MicroBeads. Repeat the experiment with diluted exosome samples or use isolated exosomes instead of cell culture supernatant.

### Co-purification of negative exosomes

▲ High concentration of exosomes or contaminations, e.g., from cell culture medium can give rise to non-specific binding of exosomes to the MicroBeads. Repeat the experiment with diluted exosome samples or use isolated exosomes instead of cell culture supernatant.

## 4. References

1. Raposo, G. and Stoorvogel, W. (2013) Extracellular vesicles: exosomes, microvesicles, and friends. *J. Cell Biol.* 200 (4): 373–383.
2. Caby, M. P. *et al.* (2005) Exosomal-like vesicles are present in human blood plasma. *Int. Immunol.* 17 (7): 879–887.
3. Aalberts, M. *et al.* (2012) Identification of distinct populations of prostasomes that differentially express prostate stem cell antigen, annexin A1, and GLIPR2 in humans. *Biol. Reprod.* 86 (3): 82.
4. Pisitkun, T. *et al.* (2004) Identification and proteomic profiling of exosomes in human urine. *Proc. Natl. Acad. Sci. U. S. A.* 101 (36): 13368–13373.
5. Ogawa, Y. *et al.* (2011) Proteomic analysis of two types of exosomes in human whole saliva. *Biol. Pharm. Bull.* 34 (1): 13–23.
6. Admyre, C. *et al.* Exosomes with immune modulatory features are present in human breast milk. *J. Immunol.* 179 (3): 1969–1978.
7. Andre, F. *et al.* (2002) Malignant effusions and immunogenic tumour-derived exosomes. *Lancet* 360(9329): 295–305.
8. Street, J. M. *et al.* (2012) Identification and proteomic profiling of exosomes in human cerebrospinal fluid. *J. Transl. Med.* 10: 5.
9. Thery, C. *et al.* (2009) Membrane vesicles as conveyors of immune responses. *Nat. Rev. Immunol.* 9 (8): 581–593.
10. Thery, C. *et al.* Exosomes: composition, biogenesis and function. *Nat. Rev. Immunol.* 2 (8): 569–579.
11. Koliha, N. *et al.* (2016) A novel multiplex bead-based platform highlights the diversity of extracellular vesicles. *J. Extracell. Vesicles* 5: 29975.
12. Koliha, N. *et al.* (2016) Melanoma affects the composition of blood cell-derived extracellular vesicles. *Front. Immunol* 7: 282.

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols.

### Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

MACS, the MACS logo, and MultiMACS are registered trademarks or trademarks of Miltenyi Biotec GmbH and/or its affiliates in various countries worldwide.

Copyright © 2018 Miltenyi Biotec GmbH and/or its affiliates. All rights reserved.