

# GlutaMAX™ Media



# Get the most from your cell cultures with GlutaMAX™ media.

## Spontaneous Decomposition of L-glutamine




Figure 1. L-glutamine spontaneously decomposes into ammonia and pyrrolidone carboxylic acid at a decomposition rate dependent on pH and temperature.<sup>1</sup>

## GlutaMAX™ Media

- Increase media stability
- Minimize toxic ammonia build-up
- Maximize cell performance

L-glutamine is an essential nutrient in cell cultures for energy production as well as protein and nucleic acid synthesis. However, L-glutamine in cell culture media spontaneously degrades (figure 1). This generates ammonia as a byproduct, which is toxic to the cells<sup>2</sup> and can affect protein glycosylation<sup>3,4</sup> and cell viability, lowering protein production and changing glycosylation patterns.

Lower ammonia concentrations can be advantageous in attaining high cell yields, particularly for cells that are sensitive to ammonia toxicity.<sup>5</sup> Cells can be sensitive to ammonia even at non-toxic levels, creating artifacts.

## Media Stability Keeps Cells Healthier, Longer

GIBCO™ GlutaMAX™ media are standard cell culture media that contain a stabilized form of L-glutamine, the dipeptide L-alanyl-L-glutamine, that prevents degradation and ammonia build-up even during long-term cultures (figures 2 and 3). Extremely stable in aqueous solution, the dipeptide does not degrade in storage or incubation.

## Stability of GlutaMAX™-I vs. L-glutamine in D-MEM




Figure 2. D-MEM was supplemented with GlutaMAX™-I or L-glutamine, aliquoted into vials and stored at 37°C. Samples were taken daily and frozen at -20°C. Levels of GlutaMAX™-I and L-glutamine were determined by HPLC.

## Stable Dipeptide Maximizes Cell Performance

The dipeptide makes the difference.

The GlutaMAX™ dipeptide is split by aminopeptidases, releasing L-glutamine and L-alanine from the dipeptide.

The mechanism of dipeptide utilization involves the gradual release of peptidase during culture to allow the gradual hydrolysis of the dipeptide in the medium (figure 4). This can be compared to the strategy of a fed-batch culture in which L-glutamine is continuously fed into the culture but maintained at low concentration. The result is an efficient energy metabolism and a high-growth yield.

The GlutaMAX™ dipeptide (GlutaMAX™-I Supplement\*) can be used as a direct substitute for L-glutamine with minimal or no adaptation.<sup>6</sup>

Now you can increase media stability, minimize toxic ammonia build-up, and maximize cell performance. Get the most from your cell cultures with GlutaMAX™ media.

## Ammonia Levels in Supplemented Media




Figure 3. D-MEM was supplemented with GlutaMAX™-I or L-glutamine, aliquoted into vials and stored at 37°C. Samples were taken daily and frozen at -20°C. Levels of ammonia were determined by HPLC.

## Stable Dipeptide Delivery






Figure 4. Controlled delivery of L-glutamine from media to cells in culture.

## Growth of AE-1 Cells



**Figure 5. Cell Growth.** AE-1 mouse myeloma cells were seeded at  $1 \times 10^5$  cells/ml. Samples were taken daily after three days and evaluated in triplicate for cell density. Cell viability was determined by trypan blue exclusion. Cell density continued to increase in samples containing GlutaMAX™-I past the point when L-glutamine containing samples decreased in cell density.

## IgG<sub>1</sub> Production by AE-1 Cells



**Figure 6. IgG<sub>1</sub> Production.** The samples from figure 5 were also examined by ELISA for IgG<sub>1</sub> production. GlutaMAX™-I containing samples outperformed L-glutamine containing samples.

## Obtain Equivalent or Improved Production Yields

While results may vary depending upon the cell line, using GlutaMAX™-I Supplement in place of L-glutamine in cell culture can improve cell viability and growth, potentially increasing productivity levels. In an application-specific example, figures 5 and 6 show an AE-1 cell growth curve and recombinant IgG<sub>1</sub> production. The GlutaMAX™-I culture demonstrates improved cell numbers and productivity.

GlutaMAX™-I Supplement can also extend cell culture life, which may reduce the number of times the cells must be passaged. Figure 7 compares MDBK cells cultured in D-MEM with 10% FBS and L-glutamine or GlutaMAX™-I Supplement. Cells cultured in GlutaMAX™-I reach peak density two days later and viability declines less rapidly than that observed in cultures with L-glutamine supplementation. The slight increase of the lag phase is attributed to the time needed to release the peptidase and digest the dipeptide. This allows a gradual increase in availability of L-glutamine to the cells.<sup>2</sup>

## Choose From Many Formulations

We offer many widely used liquid media formulations in which the GlutaMAX™ dipeptide substitutes for L-glutamine. They include D-MEM, MEM, IMDM, RPMI, Opti-MEM®, and others. For details, see the insert in this brochure or the GIBCO™ catalog.

## Do It Yourself with GlutaMAX™-I Supplement

You can purchase the GlutaMAX™ dipeptide as a stand-alone supplement. Use this 200 mM solution, GlutaMAX™-I Supplement, as a direct substitute for L-glutamine at equimolar concentrations in your current cell culture media formulation.

**Note:** This supplement is suitable for mammalian cell cultures. It is not recommended for insect cell cultures.

## Applications for GlutaMAX™ Media

GlutaMAX™ media and media supplemented with GlutaMAX™-I are suitable for both adherent and suspension mammalian cell cultures including:

- Culture systems requiring long periods of incubations without feeding (e.g., cloning assays)
- Long-term studies requiring optimum standardization of media (e.g., cancer cell lines, long-term cultures passaged over time, toxicity testing)
- Culture systems sensitive to ammonia (e.g., high-density bioreactors)

## Growth of MDBK Cells in D-MEM Supplemented with L-glutamine or GlutaMAX™-I and 10% FBS




Figure 7. MDBK cells were seeded at approximately  $1 \times 10^5$  cells/flask in D-MEM with 10% FBS and L-glutamine or GlutaMAX™-I in 25 cm<sup>2</sup> T-flasks.

## Common Cell Lines Cultured with GlutaMAX™-I

MDBK	Bovine Kidney
MDCK	Canine Kidney
HELA	Human Ovary
Per. C6	Human Embryonic Retinoblastoma
293	Human Embryonic Kidney
AE-1	Mouse Hybridoma
3D9	Mouse Hybridoma
CHO	Hamster Ovary
BHK	Hamster Kidney

For a complete list of citations, go to [www.invitrogen.com](http://www.invitrogen.com)

## Get What You Need for Better Cell Culture

Better cell culture begins with GIBCO™ products and services.

From the world's largest manufacturer of products for cell culture, GIBCO™ media, sera, reagents, and technical support have set the global standard for over 40 years.

- Publications reference GIBCO™ media more than any other media.
- The biopharmaceutical industry relies on us more than any other supplier for its critical research and production needs.
- Researchers rate our service the best, citing quality, technical support, reliable information, on time delivery, and on-site Supply Centers.

Scientists worldwide trust our quality, rely on our service, and welcome our innovations, now even more powerful through integration with Invitrogen tools for molecular biology.


You can depend on us to help you save time and money, meet regulatory requirements, and improve experimental results.

You get more with GIBCO™ products and services from Invitrogen.

You get what you need for better cell culture. †

## References

1. Tritsch, G.L. and Moore, G.E. (1962) Spontaneous Decomposition of Glutamine in Cell Culture Media. *Experimental Research* **28**, 360-364.
2. Hassell, T., Gleave, S., and Butler, M. (1991) Growth Inhibition in Cell Culture. *Applied Biochemistry and Biotechnology* **30**, 30-41.
3. Yang, M. and Butler, M. (2002) Effects of Ammonia and Glucosamine on the Heterogeneity of Erythropoietin Glycoforms. *Biotechnology Progress* **18**, 129-138.
4. Yang, M. and Butler, M. (2000) Effects of Ammonia on the Glycosylation of Human Recombinant Erythropoietin in Culture. *Biotechnology Progress* **16**, 751-759.
5. Christie, A. and Butler, M. (1994) Growth and Metabolism of a Murine Hybridoma in Cultures Containing Glutamine-based Dipeptides. *FOCUS*® **16**, 1, 9.
6. Brand, K., Feki, W., Hintzentern, J., von Langer, K., Lupp, P., and Schroener, C. (1989) *Metabolism* **38**, 29.



For a free sample of GlutaMAX™ Media or GlutaMAX™-I Supplement, call Cell Culture Technical Services, 1-800-955-6288.

## ORDERING INFORMATION

Description	Classical Media		Size	GlutaMAX™ Media	
	With L-Glutamine Cat. No.	Without L-Glutamine Cat. No.		Cat. No.	Size
<b>Dulbecco's Modified Eagle Medium (D-MEM) (1X), liquid</b> Low glucose, contains sodium pyruvate.	11885-076		1,000 ml	<b>10567-014</b>	500 ml
	11885-084		500 ml		
	11885-092		10 × 500 ml		
<b>Dulbecco's Modified Eagle Medium (D-MEM) (1X), liquid</b> High glucose, contains sodium pyruvate.	11995-040		1,000 ml	<b>10569-010</b>	500 ml
	11995-081		6 × 1,000 ml		
	11995-065	10313-021	500 ml		
	11995-073		10 × 500 ml		
<b>Dulbecco's Modified Eagle Medium (D-MEM) (1X), liquid</b> High glucose, contains no sodium pyruvate.	11965-084	11960-051	1,000 ml	<b>10566-016</b>	500 ml
	11965-126	11960-077	6 × 1,000 ml		
	11965-092	11960-044	500 ml		
	11965-118	11960-069	10 × 500 ml		
<b>Dulbecco's Modified Eagle Medium (D-MEM) (1X), liquid</b> High glucose, contains HEPES buffer, but no sodium pyruvate.	12430-047		1,000 ml	<b>10564-011</b>	500 ml
	12430-054		500 ml		
	12430-062		10 × 500 ml		
<b>D-MEM/F-12, (1X) liquid, 1:1</b>	11320-033		500 ml	<b>10565-018</b>	500 ml
<b>F-12 Nutrient Mixture, (Ham) (1X), liquid</b>	11765-047		1,000 ml	<b>31765-035</b>	500 ml
	11765-070		6 × 1,000 ml		
	11765-054		500 ml		
	11765-062		10 × 500 ml		
<b>Iscove's Modified Dulbecco's Medium (IMDM) (1X), liquid</b>	12440-046		1,000 ml	<b>31980-030</b>	500 ml
	12440-053		500 ml		
	12440-061		10 × 500 ml		
<b>Minimum Essential Medium (MEM) alpha (1X), liquid</b> Contains no ribonucleosides or deoxyribonucleosides.	12561-049		1,000 ml	<b>32561-037</b>	500 ml
	12561-056		500 ml		
<b>Minimum Essential Medium (MEM) alpha (1X), liquid</b> Contains ribonucleosides and deoxyribonucleosides.	12571-048		1,000 ml	<b>32571-036</b>	500 ml
	12571-063		500 ml		
	12571-071		10 × 500 ml		
<b>Minimum Essential Medium (MEM), liquid</b> Contains Earle's Salts.	11095-072	11090-073	1,000 ml	<b>41090-036</b>	500 ml
	11095-080	11090-081	500 ml		
	11095-098	11090-099	10 × 500 ml		
<b>Minimum Essential Medium (MEM), liquid</b> Contains Earle's Salts and HEPES buffer.		12360-038	500 ml	<b>42360-032</b>	500 ml
<b>Opti-MEM® I Reduced-Serum Medium (1X), liquid</b>	31985-062		100 ml	<b>51985-034</b>	500 ml
	31985-070		500 ml		
	31985-088		10 × 500 ml		
<b>RPMI Medium 1640 (1X), liquid</b>	11875-085	21870-084	1,000 ml	<b>61870-036</b>	500 ml
	11875-135		6 × 1,000 ml		
	11875-093	21870-076	500 ml		
	11875-119	21870-092	10 × 500 ml		
	11875-101		100 ml		
	11875-127		20 × 100 ml		
<b>RPMI Medium 1640 (1X), liquid</b> Contains HEPES buffer.	22400-071		1,000 ml	<b>72400-047</b>	500 ml
	22400-089		500 ml		
	22400-105		10 × 500 ml		
Reagents	Cat. No.		Size	GlutaMAX™-I Supplement*	
<b>L-Glutamine-200 mM (100X), liquid</b>	25030-149		20 ml	<b>35050-061</b>	100 ml
	25030-081				



www.invitrogen.com

**Corporate Headquarters:** Invitrogen Corporation • 1600 Faraday Avenue • Carlsbad, California 92008 U.S.A.  
Tel: 1 760 603 7200 • Tel (Toll Free): 1 800 955 6288 • Toll Free Fax: 1 800 331 2286 • E-mail: tech\_service@invitrogen.com

**European Headquarters:** Invitrogen Ltd • 3 Fountain Drive • Inchinnan Business Park • Paisley PA4 9RF, UK  
Tel: +44(0) 141 814 6100 • Fax: +44(0) 141 814 6260 • E-mail: eurotech@invitrogen.com



These products are for in vitro diagnostic use and are not intended for human or animal therapeutic use. Uses other than the labeled intended use may be a violation of local law. \*This product is for research use, and where appropriate, as raw material components in further cell culture manufacturing applications. It is not intended for human or animal diagnostic, therapeutic, or other clinical uses, unless otherwise stated. © 2003 Invitrogen Corporation PAS02-074MS Part No. 332-22195





#### United States Headquarters

Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, California 92008  
Phone: 1 760 603 7200  
Toll Free: 1 800 955 6288  
Toll Free Fax: 1 800 331 2286  
E-mail: [tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)

#### European Headquarters

Invitrogen Ltd  
3 Fountain Drive  
Inchinnan Business Park  
Paisley PA4 9RF, UK  
Tel: +44 (0) 141 814 6100  
Fax: +44 (0) 141 814 6260  
E-mail: [eurotech@invitrogen.com](mailto:eurotech@invitrogen.com)

#### International Offices

Argentina 5411 4556 0844  
Australia 1 800 331 627  
Austria 0800 20 1087  
Belgium 0800 14894  
Brazil 0800 11 0575  
Canada 800 263 6236  
China 10 6849 2578  
Denmark 80 30 17 40  
France 0800 23 20 79  
Germany 0800 083 0902  
Hong Kong 2407 8450  
India 11 577 3282  
Italy 02 98 22 201  
Japan 03 3663 7974  
The Netherlands 0800 099 3310  
New Zealand 0800 600 200  
Norway 22 83 16 00  
Spain & Portugal 900 181 461  
Sweden 020 26 34 52  
Switzerland 0800 848 800  
Taiwan 2 2651 6156  
For other countries see our website.

[www.invitrogen.com](http://www.invitrogen.com)

## Maximize your cell cultures with GlutaMAX™ media.

Increase media stability, minimize toxic ammonia build-up, and maximize cell performance. Use GlutaMAX™ media to get the most from your cell cultures.



[www.invitrogen.com](http://www.invitrogen.com)



Printed on recycled paper

These products are for *in vitro* diagnostic use and are not intended for human or animal therapeutic use. Use other than the labeled intended use may be a violation of local law.

\*This product is for research use, and where appropriate, as raw material components in further cell culture manufacturing applications. It is not intended for human or animal diagnostic, therapeutic, or other clinical uses, unless otherwise stated.