

Fast and Reliable: Single-Tube Real-Time RT-PCR

High-quality, ready-to-use reagents for amplification of animal pathogen RNA

- Obtain results in approximately 1 hour
- Consistently amplify RNA targets with high sensitivity
- Versatile design allows for use with multiple PCR cyclers, primers, and probes
- Minimize cross-contamination potential with fast and simple reaction setup




Figure 1. AgPath-ID™ One-Step RT-PCR Procedure Overview.

The AgPath-ID™ One-Step RT-PCR Kit is designed for sensitive, robust amplification of RNA targets using a rapid, single-tube real-time reverse transcription PCR (RT-PCR) strategy. You supply the PCR primers and TaqMan® probe (optional); the kit includes the buffers and enzymes needed for RT-PCR.

The kit is configured for fast and simple reaction setup, and reactions are assembled in a single tube, minimizing sample handling errors and expediting setup (Figure 1). Once the reactions are assembled, results are available in approximately 1 hour. The 25X RT-PCR Enzyme Mix is composed of Ambion's highly efficient ArrayScript™ Reverse Transcriptase, an MMLV RT capable of producing high cDNA yields, and

AmpliTaQ Gold® DNA Polymerase, the preferred hot-start DNA polymerase for specific target amplification. The 2X RT-PCR Buffer contains optimized reagents for efficient, robust reverse transcription and PCR. It also contains the passive reference dye, ROX™ Dye, for quantitative fluorescent signal normalization. Detection Enhancer is provided as an optional reagent for amplification of templates with high GC content or persistent secondary structure. Finally, a tube of ultrapure Nuclease-free Water is also provided.

Rigorous Linearity, Unsurpassed Sensitivity

To illustrate the consistent performance of the AgPath-ID™ One-Step RT-PCR Kit, serial dilutions of Virus A Control RNA containing 5 to 5×10^6 copies were amplified (Figure 2). The amplification plots show the consistent set of curves expected from highly efficient PCR, and the inset graph shows the reliability and efficiency of the PCR across a wide range of input template amounts.




Figure 2. The AgPath-ID™ One-Step RT-PCR Kit Consistently Amplifies Across a Wide Dynamic Range of RNA Input. qRT-PCR targeting serially diluted Virus A Control RNA transcript (5 to 5×10^6 copies) demonstrates highly efficient and consistent performance. Reactions were performed on an Applied Biosystems 7500 Fast Real-Time PCR System.

Figure 3 shows amplification of a serial dilution of another control RNA, Virus B. Amounts of RNA were kept low (20 to 40,000 copies) in order to compare the analytical sensitivity of target amplification of the AgPath-ID™ kit and a competitor's RT-PCR reagent kit. The AgPath-ID™ One-Step RT-PCR Kit provided better analytical sensitivity than the competitor's kit across the dilution range. If your research demands consistent, reliable performance and maximum sensitivity, the AgPath-ID™ One-Step RT-PCR Kit delivers.




Figure 3. The AgPath-ID™ One-Step RT-PCR Kit Provides Higher Analytical Sensitivity Than a Competitor's Kit. Serially diluted Virus B Control RNA (20 to 40,000 copies) was amplified using the AgPath-ID™ One-Step RT-PCR Kit or a leading competitor's kit (Competitor Q). Reactions were performed on an Applied Biosystems 7500 Fast Real-Time PCR System. The AgPath-ID™ One-Step RT-PCR Kit provided earlier C_t values and better analytical sensitivity than the competitor's kit.

ORDERING INFORMATION

| Description | Size | P/N |
|--------------------------------|-----------|---------|
| AgPath-ID™ One-Step RT-PCR Kit | 100 rxns | AM1005 |
| AgPath-ID™ One-Step RT-PCR Kit | 500 rxns | 4387424 |
| AgPath-ID™ One-Step RT-PCR Kit | 1000 rxns | 4387391 |

Limited Use Label License: ENVIRONMENTAL TESTING, QUALITY CONTROL/QUALITY ASSURANCE TESTING, FOOD AND AGRICULTURAL TESTING

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product (a) to perform internal research for the sole benefit of the purchaser; and (b) for environmental testing, quality control/quality assurance testing, food and agricultural testing, including reporting results of purchaser's activities in environmental testing, quality control/quality assurance testing, food and agricultural testing for a fee or other commercial consideration. No other right is hereby granted expressly, by implication, or by estoppel. This product is for environmental testing, quality control/ quality assurance testing, food and agricultural testing and research purposes only.

The purchase of this product does not grant the purchaser any additional rights, including (without limitation) the right to transfer or resell the product in any form or the right to use the product as a therapeutic agent or diagnostics test component. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

Human diagnostic uses require a separate license from Roche.

© 2009, 2011 Life Technologies Corporation. SYBR is a registered trademark of Molecular Probes, Inc. TaqMan and AmpliTaq Gold are registered trademarks of Roche Molecular Systems, Inc. All other trademarks are the sole property of their respective owners. This product is compatible with the 5' nuclease detection and dsDNA-binding dye processes covered by patents owned or licensable by Applied Biosystems. No license under these patents is conveyed expressly, by implication, or by estoppel to the purchaser by the purchase of this product. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Printed in the USA, 06/2011 Publication B37PBII-03