

# Path-ID™ Multiplex One-Step RT-PCR Kit

## Sensitive Multiplex Amplification of RNA Targets

- **Efficient**—save time and money by amplifying up to four different RNA targets simultaneously
- **Simple**—single-tube reaction minimizes cross-contamination
- **Sensitive**—optimized for amplification of low-abundance RNA
- **Wide dynamic range**—efficient amplification over 7 logs of RNA input amount

The Path-ID™ Multiplex One-Step RT-PCR Kit is designed for the sensitive, robust amplification and multiplex quantitation of up to four RNA targets simultaneously using a simple real-time RT-PCR strategy (Figure 1). The kit includes the buffers and enzymes needed for RT-PCR; you supply only the RNA sample, PCR primers, and TaqMan® probes.

### Simple, Robust Reactions

Path-ID™ Multiplex One-Step RT-PCR Kit reactions are run using a single-tube, one-step procedure to reverse transcribe the RNA and amplify your targets (Figure 1). Single-tube reactions minimize sample handling errors, producing more consistent results. Reactions are driven with Multiplex Enzyme Mix, which includes ArrayScript™ Reverse Transcriptase, an MMLV RT capable of producing high cDNA yields. The mix also includes AmpliTaq Gold® DNA Polymerase, the preferred hot-start DNA polymerase for specific target amplification, and RT-PCR Buffer with optimized reagents for efficient, robust results from both the reverse




Figure 1. Path-ID™ Multiplex One-Step RT-PCR Kit Procedure Overview.




Figure 2. The Path-ID™ Multiplex One-Step RT-PCR Kit Consistently Amplifies Four RNA Targets in a Single Reaction. Xeno™ RNA Control and control RNAs for virus A, virus B, and virus C were amplified in a single multiplex reaction using the Path-ID™ Multiplex One-Step RT-PCR Kit and run on the Applied Biosystems® 7500 Real-Time PCR System. A sample set with fixed amounts of three of the targets and a dilution series of the virus B control RNA (red curve) were included.

transcription reaction and the PCR. The buffer also contains ROX™ passive reference dye for quantitative fluorescent signal normalization and increased data precision.

### Increase Efficiency With Multiplexing

The Path-ID™ Multiplex One-Step RT-PCR Kit is designed for the amplification of up to four RNA targets simultaneously, resulting in both time and cost savings. Figure 2 shows


the amplification of four targets by multiplex RT-PCR using the kit. The quantities of three of the targets in the experiment were held constant, but the fourth target was serially diluted to show the dynamic range of multiplex target amplification with the kit. Furthermore, Figure 3 shows that similar sensitivity is obtained in multiplexed reactions and in singleplex reactions.

### Better Sensitivity


In comparison to a competitor's product, the Path-ID™ Multiplex One-Step RT-PCR Kit shows superior sensitivity with more consistent results (Figure 4).

### Streamline Your Research

The Path-ID™ Multiplex One-Step RT-PCR Kit brings the sensitivity, specificity, and wide dynamic range of real-time PCR to amplification of nucleic acid target sequences. The kit is configured to amplify up to four RNA targets in a single RT-PCR, without compromising results. With a single-tube procedure, potential sample handling errors are minimized and reaction setup is fast and easy. The kit has been shown to amplify as few as 20 copies of target, enabling unsurpassed sensitivity in a streamlined workflow.



**Figure 3. Comparison of the Path-ID™ Multiplex One-Step RT-PCR Kit Amplifying Targets in Singleplex and Duplex RT-PCRs.** Virus E RNA was reverse transcribed and PCR amplified in a singleplex reaction, and virus D and virus E were reverse transcribed and co-amplified in a duplex reaction using the Path-ID™ Multiplex One-Step RT-PCR Kit. No difference in sensitivity was observed between the singleplex and multiplex reactions.



**Figure 4. The Path-ID™ Multiplex One-Step RT-PCR Kit Provides Higher Target Analytical Sensitivity Than a Competitor's Kit.** The quadruplex RT-PCR experiment described in Figure 2 was performed using the Path-ID™ Multiplex One-Step RT-PCR Kit and a kit from Supplier X. Only data for the virus B target are shown. The Path-ID™ Multiplex One-Step RT-PCR Kit amplified the lower amounts of target with better sensitivity (lower Ct values) than the kit from Supplier X. Note that at <1,000 copies/rxn, only the Path-ID™ Multiplex One-Step RT-PCR Kit was able to amplify virus B RNA.

### ORDERING INFORMATION

Description	Size	Part Number
Path-ID™ Multiplex One-Step RT-PCR Kit	100 rxns	4428206
Path-ID™ Multiplex One-Step RT-PCR Kit	500 rxns	4428207
Path-ID™ Multiplex One-Step RT-PCR Kit	1,000 rxns	4440022

For Research Use Only. Not intended for human or animal therapeutic or diagnostic use.

© 2010 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. AmpliTaq Gold and TaqMan are registered trademarks of Roche Molecular Systems, Inc. Printed in the U.S.A. 04/2010 Publication C012758