

FilmTracer™ SYPRO® Ruby Biofilm Matrix Stain

Catalog no. F10318

Table 1. Contents and storage information.

Material	Amount	Storage	Stability
FilmTracer™ SYPRO® Ruby biofilm matrix stain	200 mL	<ul style="list-style-type: none"> • 2–25°C • Desiccate • Protect from light • DO NOT FREEZE 	When stored as directed, the product is stable for at least 9 months.

Number of reactions: Sufficient material is supplied for 1000 reactions, based on the protocol below.

Approximate fluorescence excitation/emission maxima: FilmTracer™ SYPRO® Ruby biofilm matrix stain: 280, 450/610 nm.

Introduction

Biofilms present a unique set of challenges for fluorescent staining and subsequent imaging. A typical biofilm not only exhibits heterogeneous thickness throughout the surface, placing stringent restrictions on stain penetration, but also contains regions of widely varying environmental conditions. Evidence suggests that bacterial cells exist in various physiological states within these biofilm microenvironments. Furthermore, biofilms contain many undefined components (*e.g.*, the extracellular polymeric matrix) that differ with species and conditions.

SYPRO® Ruby stain labels most classes of proteins, including glycoproteins, phosphoproteins, lipoproteins, calcium binding proteins, fibrillar proteins, and other proteins that are difficult to stain.¹ SYPRO® Ruby stain has been tested to stain the matrix of *Pseudomonas aeruginosa* (ATCC 15442) and some strains of *E. coli* (not *E. coli* K-12, which does not produce cellulose). Staining patterns may vary depending upon the organism and the matrix composition.

Before You Begin

Materials Required but Not Provided

- Biofilm samples: Biofilms may be grown on coupons in a biofilm reactor, as colony biofilms, in flow-cell system, or in drip-flow reactors. This protocol describes staining biofilms grown on glass coupons in a CDC reactor. For more information on the CDC reactor or other reactor types, refer to the BioSurface Technologies website (www.imt.net/~mitbst/Products.html), contact Center for Biofilm Engineering, Montana State University, Bozeman, Montana, or refer to the standard protocols outlined in the following ASTM methods: ASTM E2647, ASTM E2562, ASTM E2196.

- Fluorescence microscope with appropriate excitation/emission filters
- 0.2 µm filter-sterilized water
- Staining dishes (*e.g.*, 60-mm dish, 6-well plate, etc.)

Preparing Reagents The FilmTracer™ SYPRO® Ruby biofilm matrix stain comes as a ready to use stain and does not require further dilution.

Caution SYPRO® Ruby biofilm matrix stain is characterized as an irritant due entirely to the solvent system and buffer salts in the product. The heavy metal component of the stain (ruthenium) is not a regulated hazardous heavy metal in the United States. Independent toxicity tests found SYPRO® Ruby stain to be nontoxic at 5000 mg/kg in Sprague-Dawley® rats. Nevertheless, SYPRO® Ruby stain should be handled with care, consistent with good laboratory practices.

Federal regulation in the US does not consider ruthenium or the solvent system and buffers in SYPRO® Ruby biofilm matrix stain to be hazardous waste. Nevertheless, the solvent may be subject to more stringent regulations local to the user, so we recommend that you mix used staining solution with flammable waste and dispose by destructive incineration and not by sewerage or landfill.

Experimental Protocols

- Guidelines for Staining**
- We recommend performing staining in water as the phosphates in buffers may interfere with fluorescent staining.
 - If you need to stain and image multiple samples, do not stain more than two samples at a time. Evidence suggests that, in many cases, stain might be drawn from cells over time as they sit in water. Stagger staining, so that samples are stained, rinsed, and imaged following the same schedule. Image immediately following rinsing.
 - For imaging biofilm on CDC reactor coupons, use glass coupons only. In particular, avoid polycarbonate coupons for imaging purposes as polycarbonate is autofluorescent, and the rough surface interferes with imaging.
 - If you follow the protocol below, you do not need to use fixatives on the biofilm.

Staining Procedure The protocol below describes how to stain biofilms grown on glass CDC reactor coupons. For any other growth surface, you may need to adjust the staining volumes.

- 1.1. Add 200 µL (or appropriate volume) of staining solution onto the biofilm sample. Add the stain very gently so as not to disturb the biofilm. Immediately add the stain before the biofilm dries.
- 1.2. Incubate the sample for up to 30 minutes at room temperature, **protected from light**.
- 1.3. Rinse the sample gently with filter-sterilized water. Remove all excess stain and rinse water from the base of the support material.
- 1.4. For best results with reactor coupons, place coupon in a 60-mm dish, fill the dish with filter-sterilized water to cover the coupon surface by 1–3 mm, and observe on the microscope using a 40X 0.7NA 3.3 mm WD water objective or a 63X 0.9NA 2.2 mm WD water immersion objective.

References

1. Electrophoresis 21, 2509 (2000).

Product List **Current prices may be obtained from our website or from our Customer Service Department.**

Cat. no.	Product Name	Unit Size
F10318	FilmTracer™ SYPRO® Ruby biofilm matrix stain.....	200 mL
Related Products		
F10317	FilmTracer™ FM® 1-43 green biofilm cell stain.....	1 mg
F10319	FilmTracer™ calcein red-orange biofilm stain.....	20 × 50 µg
F10320	FilmTracer™ calcein violet biofilm stain.....	20 × 25 µg
F10322	FilmTracer™ calcein green biofilm stain.....	20 × 50 µg
L10316	FilmTracer™ LIVE/DEAD™ Biofilm Viability Kit.....	1 kit
S34854	SYTO® 9 green fluorescent nucleic acid stain *5 mM solution in DMSO*.....	100 µL

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probetech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the most restrictive terms apply. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents that are owned or controlled by Invitrogen Corporation and/or Molecular Probes, Inc. which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2009, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.