


## Anti-Mouse CD8a eFluor® 605NC

**Catalog Number:** 93-0081

**Also Known As:** CD8 alpha, Ly-2, Ly-35, Ly-B, Lyt-2

**RUO: For Research Use Only. Not for use in diagnostic procedures.**



Staining of BALB/c splenocytes with Anti-Mouse CD4 APC (cat. 17-0041) and Rat IgG2a K Isotype Control eFluor® 605NC (cat. 93-4321) (left) or Anti-Mouse CD8a eFluor® 605NC (right). Cells in the lymphocyte gate were used for analysis.

### Product Information

**Contents:** Anti-Mouse CD8a eFluor® 605NC

**REF** **Catalog Number:** 93-0081

**Clone:** 53-6.7

**Concentration:** 5 µL

**Host/Isotype:** Rat IgG2a, kappa

**Formulation:** aqueous buffer, 0.09% sodium azide, may contain carrier protein/stabilizer

**Temperature Limitation:** Store at 2-8°C. Light sensitive material. This product is guaranteed for 6 months upon receipt when stored properly.

**LOT** **Batch Code:** Refer to Vial

**Use By:** Refer to Vial

**Caution, contains Azide**

### Description

The 53-6.7 monoclonal antibody reacts with the mouse CD8a molecule. CD8a is an approximately 32-34 kDa cell surface receptor expressed either as a heterodimer with the CD8 beta chain (CD8 alpha beta) or as a homodimer (CD8 alpha alpha). A majority of thymocytes and a subpopulation of mature alpha beta TCR T cells express CD8 alpha beta while gamma delta; TCR T cells, a subpopulation of intestinal intraepithelial lymphocytes (IELs) and dendritic cells express CD8 alpha alpha. CD8 binds to MHC class I and through its association with protein tyrosine kinase p56lck plays a role in T cell development and activation of mature T cells.

### Applications Reported

This 53-6.7 antibody has been reported for use in flow cytometric analysis.

### Applications Tested

This 53-6.7 antibody has been pre-titrated and tested by flow cytometric analysis of mouse splenocytes. This can be used at 5 µL per test. A test is defined as the amount of antibody that will stain a cell sample in a final volume of 100 µL. Cell number should be determined empirically but can range from 10<sup>5</sup> to 10<sup>8</sup> cells/test.

The isotype control eFluor 605NC rat IgG2a (cat. 93-4321) should be used at 5 µL/test.

**Laser/Filter Recommendation:** When using eFluor 605NC, we recommend excitation with the 405nm violet laser with an appropriate filter set, such as the 595LP dichroic mirror with the 605/40 bandpass filter. An acceptable alternative is the 610/20 bandpass filter. For instruments not equipped with a violet laser, the eFluor 605NC is also excited by the 488 nm blue laser and can be used as a PE-Texas Red alternative.

**Fixation Recommendation:** When fixing samples that have been stained with nanocrystal reagents, we recommend keeping the total volume at approximately 200 µL of IC Fixation Buffer (cat. 00-8222) and the exposure time 30-60 minutes to preserve the optimal fluorescent signal from the nanocrystal reagent.

For answers about fixation and other questions, please refer to Nanocrystal Frequently Asked Questions or contact eBioscience Technical Support.

### References

Mochimaru H, Usui T, Yaguchi T, Nagahama Y, Hasegawa G, Usui Y, Shimmura S, Tsubota K, Amano S, Kawakami Y, Ishida S. Suppression of alkali burn-induced corneal neovascularization by dendritic cell vaccination targeting VEGF receptor 2. Invest Ophthalmol Vis Sci. 2008 May;49(5):2172-7. (53-6.7, in vivo depletion, PubMed)

Yang Z, Day YJ, Toufektsian MC, Xu Y, Ramos SI, Marshall MA, French BA, Linden J. Myocardial infarct-sparing effect of adenosine A2A receptor

activation is due to its action on CD4+ T lymphocytes. *Circulation*. 2006 Nov 7;114(19):2056-64. (53-6.7, in vivo depletion, PubMed)

Taylor JL, Ordway DJ, Trout J, Gonzalez-Juarrero M, Basaraba RJ, Orme IM. Factors associated with severe granulomatous pneumonia in *Mycobacterium tuberculosis*-infected mice vaccinated therapeutically with hsp65 DNA. *Infect Immun*. 2005 Aug;73(8):5189-93. (53-6.7, IHC frozen)

Grabbe S, Varga G, Beissert S, Steinert M, Pendl G, Seeliger S, Bloch W, Peters T, Schwarz T, Sunderkötter C, Scharffetter-Kochanek K. Beta2 integrins are required for skin homing of primed T cells but not for priming naïve T cells. *J Clin Invest*. 2002 Jan;109(2):183-92. (53-6.7, IHC frozen)

Ledbetter JA, Rouse RV, Micklem HS, Herzenberg LA. T cell subsets defined by expression of Lyt-1,2,3 and Thy-1 antigens. Two-parameter immunofluorescence and cytotoxicity analysis with monoclonal antibodies modifies current views. *J Exp Med*. 1980 Aug 1;152(2):280-95.

Ledbetter, J. A. and L. A. Herzenberg. Xenogeneic monoclonal antibodies to mouse lymphoid differentiation antigens. *Immunol Rev*. 1979;47:63-90.

#### **Related Products**

00-4222 Flow Cytometry Staining Buffer

17-0041 Anti-Mouse CD4 APC (GK1.5)

93-4321 Rat IgG2a K Isotype Control eFluor® 605NC (eBR2a)

#### **Legal**

Under patent number: US 7,939,170 and additional pending patent application(s)

---

Not for further distribution without written consent.

Copyright © 2000-2012 eBioscience, Inc.

Tel: 888.999.1371 or 858.642.2058 • Fax: 858.642.2046 • [www.eBioscience.com](http://www.eBioscience.com) • [info@eBioscience.com](mailto:info@eBioscience.com)