

# NativePAGE™ Cathode Buffer Additive (20X)

(Novex®)

货号 BN2002

包装 250 ml

NativePAGE™ Cathode Buffer Additive (20X) is added to the NativePAGE™ Running Buffer and the resulting buffer is used in the cathode chamber of an XCell SureLock® Mini Cell for running NativePAGE™ Gels.

## 详细说明

## 常用规格

<b>Gel Type:</b>	IEF Gel
<b>Buffer Type:</b>	Anode & Cathode Buffers
<b>Concentrated:</b>	20 X
<b>Product Size:</b>	250 ml
<b>Gel Compatibility:</b>	NativePAGE™ Gels
<b>Shipping Condition:</b>	Room Temperature
<b>Regulatory Statement:</b>	For Research Use Only. Not for use in diagnostic procedures.

## 内容和储存

Store at room temperature