

Recombinant Rat Interleukin-10

Publication Number MAN0004425 Rev. 1.00

Catalog Number:	PRC0104	PRC0105
Quantity:	10 µg	25 μg
Lot Number:	See product label.	
Molecular Weight:	18.7 kDa, 160 amino acid residues.	
Purity:	>98% as determined by SDS-PAGE and HPLC analyses.	
Biological Activity:	$ED_{50} \le 10 \text{ ng/mL}$ ($\ge 1 \times 10^5 \text{ units/mg}$). The biological activity is determined in an assay measuring dose–dependent inhibition of antigen–specific T cell proliferation. For most in vitro applications, IL–10 is effective at a concentration range of 0.2–20 ng/mL. The optimal concentration should be determined for each specific application.	
Formulation:	Lyophilized from 5 mM Tris, pH 7.2. Sterile filtered prior to lyophilization.	
Endotoxin:	$< 0.1 \mathrm{ng}/\mathrm{\mu g}$	
Production:	Recombinant rat IL–10 is produced in <i>E. coli</i> and purified via affinity chromatography.	
Reconstitution Recommendation:	We recommend that the vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute lyophilized rat IL–10 in 5 mM Tris, pH 7.2, to a concentration of 0.1–1.0 mg/mL. Further dilution should be made in medium or buffered solution containing carrier protein, such as PBS with 0.1% BSA.	
Storage:	This lyophilized preparation is stable at 2°C to 8°C, but should be kept at –20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2°C to 8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at –20°C. Avoid repeated freeze/thaw cycles.	
Expiration Date:	Expires one year from date of receipt when stored as instructed.	
References:	Hirano, T. (1994) In: The Cytokine Handbook, 2nd Edition, (A. Thomson, editor) Academic Press Limited, San Diego, Ca, pp. 223–237. Terrone, D.A., B.K. Rinehart, J.P. Granger, P.S. Barrilleaux, J.N. Martin, and W.A. Bennett (2001) Interleukin–10 administration and bacterial endotoxin-induced preterm birth in a rat model. Obstetrics and Gynecology 98(3):476–480.	

Explanation of Symbols

The symbols present on the product label are explained below:

Symbol	Description
REF	Catalog Number
RUO	Research Use Only
	Use by
***	Manufacturer
[-]	Without, does not contain
from Light	Protect from light
<u> </u>	Directs the user to consult instructions for use (IFU), accompanying the product.

Symbol	Description	
LOT	Batch code	
IVD	In vitro diagnostic medical device	
1	Temperature limitation	
EC REP	European Community authorized representative	
[+]	With, contains	
À	Consult accompanying documents	

Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

 $\label{lem:continuous} \textbf{For Research Use Only. Caution: Not for human or animal the rapeutic or diagnostic use.}$

 $Manufacturing\ site:\ 7335\ Executive\ Way\ |\ Frederick,\ MD\ 21704\ |\ Toll\ Free\ in\ USA\ 800.955.6288$

© 2011 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

 $\label{prop:comsupport} \textbf{For support visit} \ www. life technologies.com/support or email techsupport@life tech.com/support or email techsupport@life technologies.com/support or email techsupport@life technologies.com/support@life tec$

Revision Date 24 May 2011

