

The purchase of this product conveys to the end user the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the end user (whether the end user is an academic or for-profit entity). The end user cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The end user may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components for contract manufacturing services; (2) use of the product or its components to provide a service, information, or data; (3) use of the product itself or its components as a therapeutic, diagnostic or prophylactic; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. If the purchaser is not willing to accept the limitations of this limited use statement, PeproTech is willing to accept return of the product if notified within thirty (30) days of purchase, in exchange for full refund or credit. For information on purchasing a license to this product for purposes other than research, contact

PeproTech, Inc.  
PO Box 275  
Rocky Hill, NJ 08553  
Tel.: 609-497-0253  
Email: [info@peprotech.com](mailto:info@peprotech.com)



Princeton Business Park, 5 Crescent Ave., P O Box 275  
Rocky Hill, NJ 08553 USA  
Tele: 609.497.0253 Fax: 609.497.0321  
[www.peprotech.com](http://www.peprotech.com) Email: [info@peprotech.com](mailto:info@peprotech.com)

**Human EG-VEGF  
ELISA Development Kit  
900-K244 Lot# 0908244  
Expiration one year from date of receipt**

**Description:** Human EG-VEGF ELISA development kit contains the key components required for the quantitative measurement of natural and/or recombinant hEG-VEGF in a sandwich ELISA format within the range of 16–2000pg/ml. Using the ELISA protocol described below, the recommended microplates, reagents and solutions, the components supplied in this kit are sufficient to assay hEG-VEGF in approximately 1000 ELISA plate wells.

**RECONSTITUTION & STORAGE**

**Capture Antibody:** 25µg of antigen-affinity purified goat anti-hEG-VEGF + 2.5mg D-mannitol. Centrifuge vial prior to opening. Reconstitute in 0.25ml sterile water for a concentration of 100µg/ml.

**Detection Antibody:** 25µg of biotinylated antigen-affinity purified goat anti-hEG-VEGF + 2.5mg D-mannitol. Centrifuge vial prior to opening. Reconstitute in 0.25ml sterile water for a concentration of 100µg/ml.

**Human EG-VEGF Standard:** 1µg of recombinant hEG-VEGF + 2.2mg BSA + 11.0mg D-mannitol. Centrifuge vial prior to opening. Reconstitute in 1ml sterile water for a concentration of 1µg/ml.

**Note:** The reconstituted components are stable for 2 weeks when stored at 2–8°C. Components that have been reconstituted and aliquoted can be stored at –20°C for up to 6 months.

**Avidin-HRP Conjugate:** 60µl vial. Upon receipt, avidin-HRP conjugate should be aliquoted into ten 6µl vials and stored at ≤–20°C. Aliquots stored frozen at ≤–20°C are stable for up to 2 years from date of receipt. Avoid more than one freeze-thaw cycle. **Avidin should be used in conjunction with ABTS only.**

**RECOMMENDED MATERIALS (or purchase  
PeproTech's ELISA Buffer Kit: Cat.# 900-K00)**

ELISA microplates (Nunc MaxiSorp Prod. # 439454, or Corning Prod # 3590);  
Tween-20 (Sigma Cat. # P-7949);  
BSA (Sigma Cat # A-7030);  
ABTS Liquid Substrate Solution (Sigma Cat. # A3219);  
Dulbecco's PBS [10x] (Gibco BRL Cat. # 14200-075).

**RECOMMENDED SOLUTIONS**

*All solutions should be at ambient temperature prior to use.*

**PBS:** dilute 10xPBS to 1xPBS, pH 7.20 in sterile water.

**Wash Buffer:** 0.05% Tween-20 in PBS

**Block Buffer:** 1% BSA in PBS \*

**Diluent:** 0.05% Tween-20, 0.1% BSA in PBS \*

**\* Sterile filter and store at 4°C for up to 1 week.**

**PLATE PREPARATION**

1. Dilute capture antibody with PBS to a concentration of 0.25µg/ml. Immediately, add 100µl to each ELISA plate well. Seal the plate and incubate overnight at room temperature.

2. Aspirate the wells to remove liquid and wash the plate 4 times using 300µl of wash buffer per well. After the last wash invert plate to remove residual buffer and blot on paper towel.
3. Add 300µl block buffer to each well. Incubate for at least 1 hour at room temperature.
4. Aspirate and wash plate 4 times.

**ELISA PROTOCOL**

**Standard/Sample:** dilute standard from 2ng/ml to zero in diluent. Immediately add 100µl of standard or sample to each well in triplicate. Incubate at room temperature for at least 2 hours.

**Detection:** Aspirate and wash plate 4 times. Dilute detection antibody in diluent to a concentration of 0.25µg/ml. Add 100µl per well. Incubate at room temperature for 2 hours.

**Avidin-HRP Conjugate:** Aspirate and wash plate 4 times. Dilute 5.5µl of avidin-HRP conjugate 1:2000 in diluent for a total volume of 11ml. Add 100µl per well. Incubate 30 minutes at room temperature.

**ABTS Liquid Substrate:**

*(ABTS Substrate should be at ambient temperature prior to use)*

Aspirate and wash plate 4 times. Add 100µl of substrate solution to each well. Incubate at room temperature for color development. Monitor color development with an ELISA plate reader at 405 nm with wavelength correction set at 650 nm.

**NOTE:** Reliable standard curves are obtained when either O.D. readings do not exceed 0.2 units for the zero standard concentrations, or 1.2 units for the highest standard concentration. The plate should be monitored at 5-minute intervals for approximately 20 minutes. *O.D. readings may vary.*

**CROSS REACTIVITY**

When tested at 50ng/ml the following antigens did not exhibit significant cross reactivity:

Human	EGF, FGF-16, flt3-L, GM-CSF, PDGF-AA, PDGF-AB, PDGF-BB, PIGF, Prokineticin-2, RANTES, SCF, VEGF <sub>165</sub> , VEGF <sub>121</sub>
Murine	EGF, flt3-L, GM-CSF, SCF, VEGF
Rat	EGF, SCF

