Technical Data Sheet

Purified Mouse IgE Isotype Control

Product Information

554118 **Material Number:** Size: 0.1 mg 0.5 mg/ml**Concentration:** Clone: IgE-3

Mouse (BALB/c) IgE, κ Isotype: Reactivity: QC Testing: Mouse

Aqueous buffered solution containing ≤0.09% sodium azide. Storage Buffer:

Description

This reagent is used as a mouse IgE isotype (negative) control. It is specific for TNP (trintrophenol).

Preparation and Storage

Store undiluted at 4°C.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

Application Notes

Application

Isotype control	Routinely Tested
Intracellular staining (flow cytometry)	Routinely Tested

Product Notices

- 1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
- Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.

BD Biosciences

bdbiosciences.com

United States Canada Asia Pacific Latin America/Caribbean Europe Japan 877.232.8995 888.259.0187 32.53.720.550 0120.8555.90 65.6861.0633

For country-specific contact information, visit bdbiosciences.com/how_to_order/

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2008 BD



Page 1 of 1

554118 Rev. 8