# **Technical Data Sheet**

# Alexa Fluor® 647 Mouse IgG2a, κ Isotype Control

## **Product Information**

Material Number:	558053
Size:	50 tests
Vol. per Test:	20 µl
Clone:	MOPC-173
Isotype:	Mouse (BALB/c) IgG2a, κ
Storage Buffer:	Aqueous buffered solution containing BSA and $\leq 0.09\%$ sodium azide.

#### Description

The MOPC-173 antibody has unknown specificity. The transplantable plasmacytoma MOPC-173 was induced by intraperitoneal injection of mineral oils into BALB/c mice. In the absence of antigen-specific binding, this immunoglobulin may bind non-specifically to Fc receptors.

#### **Preparation and Storage**

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated to Alexa Fluor® 647 under optimum conditions, and unreacted Alexa Fluor® 647 was removed.

#### **Application Notes**

Application				
	Intracellular staining (flow cytometry)	Routinely Tested		

#### **Recommended Assay Procedure:**

This antibody conjugate is suitable for intracellular staining of human whole blood (using BD Phosflow™ Lyse/Fix Buffer) and peripheral blood mononuclear cells (using BD Cytofix™ Fixation Buffer or BD Phosflow Fix Buffer I). Any of the three BD Phosflow permeabilization buffers may be used.

## **Suggested Companion Products**

Catalog Number	Name	Size	Clone
554655	Fixation Buffer	100 ml	(none)
557870	Fix Buffer I	250 ml	(none)
558049	Lyse/Fix Buffer 5X	250 ml	(none)
557885	Perm/Wash Buffer I	125 ml	(none)
558052	Perm Buffer II	125 ml	(none)
558050	Perm Buffer III	125 ml	(none)

#### Product Notices

- 1. This reagent has been pre-diluted for use at the recommended Volume per Test. We typically use  $1 \times 10^{6}$  cells in a 100-µl experimental sample (a test).
- 2. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- 3. Source of all serum proteins is from USDA inspected abattoirs located in the United States.
- 4. Alexa Fluor® 647 fluorochrome emission is collected at the same instrument settings as for allophycocyanin (APC).
- 5. The Alexa Fluor®, Pacific Blue<sup>™</sup>, and Cascade Blue<sup>®</sup> dye antibody conjugates in this product are sold under license from Molecular Probes, Inc. for research use only, excluding use in combination with microarrays, or as analyte specific reagents. The Alexa Fluor® dyes (except for Alexa Fluor® 430), Pacific Blue<sup>™</sup> dye, and Cascade Blue® dye are covered by pending and issued patents.
- 6. Alexa Fluor® is a registered trademark of Molecular Probes, Inc., Eugene, OR.
- 7. For fluorochrome spectra and suitable instrument settings, please refer to our Multicolor Flow Cytometry web page at www.bdbiosciences.com/colors.
- 8. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.

#### References

Sibinovic KH, Potter M, Hoostelaere, Rode B, Wax J, ed. Catalogue of plasmacytomas and other tumors of the lymphoreticular system, 3rd edition. Kensington, Maryland: Litton Bionetics, Inc; 1976:1-33. (Clone-specific)

#### **BD Biosciences**

bdbiosciences.com

bublosciences.com							
	United States	Canada	Europe	Japan	Asia Pacific	Latin America/Caribbean	
	877.232.8995	800.268.5430	32.2.400.98.95	0120.8555.90	65.6861.0633	55.11.5185.9995	
For country contact information, visit <b>bdbiosciences.com/contact</b>							
Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation							
of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the							
	use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone						



be of our product or as a component of another product. Any any fight to reserve that the permitted use without the express product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited. For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale. Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2014 BD