Technical Data Sheet

Alexa Fluor® 488 Mouse anti-PDGFRβ (CD140b) (pY1021)

Product Information

Material Number:	558419		
Size:	50 tests		
Vol. per Test:	20 µl		
Clone:	J105-412		
Immunogen:	Phosphorylated Human PDGFRβ (pY1021)		
Isotype:	Mouse (BALB/c) IgG1, ĸ		
Reactivity:	QC Testing: Mouse		
	Tested in Development: Human		
Storage Buffer:	Aqueous buffered solution containing BSA and ≤0.09% sodium azide.		

Description

Platelet-derived growth factor (PDGF) is a potent mitogen for cells of mesenchymal origin and exerts its effects by binding to the PDGF receptor (PDGFR), a transmembrane protein tyrosine kinase. PDGFR is composed of PDGFR α (CD140a) and/or PDGFR β (CD140b) polypeptides. Both PDGF and PDGFR consist of subunits that form homo- or heterodimers with varying specificities: PDGF-AA binds only to $\alpha\alpha$ PDGFR, PDGF-AB binds to both $\alpha\alpha$ and $\alpha\beta$ PDGFR, and PDGF-BB binds to all three PDGFRs. Ligand binding induces dimerization and activation of the receptor. Upon activation, CD140b is phosphorylated at multiple tyrosine sites and, in turn, an intracellular phosphorylation cascade is initiated. PDGFR localizes primarily to membrane invaginations termed caveolae, compartments that are enriched in several of its downstream effectors, including phosphatidylinositol 3'-kinase, Src, and phospholipase C- γ (PLC- γ).

The J105-412 monoclonal antibody recognizes the phosphorylated tyrosine 1021 (pY1021) in the C-terminal noncatalytic region of CD140b, which interacts primarily with PLC- γ . The orthologous phosphorylation site in mouse PDGFR β is Y1020.



Analysis of PDGFRβ (CD140b) (pY1021) in mouse embryonic fibroblasts. Serum-starved NIH/373 cells were either stimulated with PDGF-BB (Cat. No. 354051, shaded histogram) or unstimulated (open histogram). The cells were fixed (BD Phosflow™ Fix Buffer I, Cat. No. 557870) for 10 minutes at 37°C, then permeabilized (BD Phosflow Perm Buffer III, Cat. No. 558050) on ice for at least 30 minutes, and then stained with Alexa Fluor® 488 anti-PDGFRβ (CD140b) (pY1021). Flow cytometry was performed on a BD™ FACSCalibur flow cytometry system.

Preparation and Storage

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated to Alexa Fluor® 488 under optimum conditions, and unreacted Alexa Fluor® 488 was removed.

Application Notes

Application

Intracellular staining (flow cytometry)		Routinely Tested		sted			
BD Bioscie	ences						
bdbiosciences.com							
United States 877.232.8995	Canada 800.268.5430	Europe 32.2.400.98.95	Japan 0120.8555.90	Asia Pacific 65.6861.0633	Latin America/Caribbean 55.11.5185.9995		M B D
For country cor	ntact informatio	on, visit bdbiosci e	ences.com/conta	ct			
Conditions: The in of any patents. BL use of our produc product or as a co written authoriza	formation disclose D Biosciences will n ts. Purchase does n imponent of anoth tion of Becton Dic	d herein is not to be ot be held responsil ot include or carry a er product. Any use kinson and Compan	e construed as a reco ole for patent infring any right to resell or of this product othe y is stirtly prohibited	ommendation to use gement or other vio transfer this producer than the permittee d	e the above product in violation lations that may occur with the ct either as a stand-alone ed use without the express		

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale. Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2014 BD

Suggested Companion Products

Catalog Number	Name	Size	Clone
557870	Fix Buffer I	250 ml	(none)
558050	Perm Buffer III	125 ml	(none)

Product Notices

- This reagent has been pre-diluted for use at the recommended Volume per Test. We typically use 1×10^{6} cells in a 100-µl experimental 1. sample (a test).
- 2. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- Source of all serum proteins is from USDA inspected abattoirs located in the United States. 3
- 4. Alexa Fluor® 488 fluorochrome emission is collected at the same instrument settings as for fluorescein isothiocyanate (FITC).
- Species testing during development may have been performed with a different format of the same clone. Selected applications have been 5. tested for cross-reactivity.
- Alexa Fluor® is a registered trademark of Molecular Probes, Inc., Eugene, OR. 6.
- The Alexa Fluor®, Pacific Blue™, and Cascade Blue® dye antibody conjugates in this product are sold under license from Molecular 7. Probes, Inc. for research use only, excluding use in combination with microarrays, or as analyte specific reagents. The Alexa Fluor® dyes (except for Alexa Fluor® 430), Pacific Blue™ dye, and Cascade Blue® dye are covered by pending and issued patents.
- For fluorochrome spectra and suitable instrument settings, please refer to our Multicolor Flow Cytometry web page at 8.
- www.bdbiosciences.com/colors.
- 9 Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.

References

Claesson-Welsh L. Platelet-derived growth factor receptor signals. J Biol Chem. 1994; 269(51):32023-32026. (Biology) Liu J, Oh P, Horner T, Rogers RA, Schnitzer JE. Organized endothelial cell surface signal transduction in caveolae distinct from glycosylphosphatidylinositol-anchored protein microdomains. J Biol Chem. 1997; 272(11):7211-7222. (Biology)

BD Biosciences

bdbiosciences.com United States

 Canada
 Europe
 Japan

 800.268.5430
 32.2.400.98.95
 0120.8555.90
Asia Pacific 65.6861.0633 877.232.8995 For country contact information, visit bdbiosciences.com/contact

For Country Contact Information, Visit **Dabiosciences.com/contact** Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited. For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale. Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2014 BD

Latin America/Caribbean

55.11.5185.9995

