

## Technical Data Sheet

## Functional Bead A9

## Product Information

Material Number:	558579
Size:	0.9 ml
Storage Buffer:	Aqueous buffered solution containing ≤0.09% sodium azide.

## Description

The BD™ CBA Functional Beads are color-coded 7.5 µm polystyrene beads that can be classified using a flow cytometer. The Functional Beads can be used in conjunction with the BD CBA Functional Bead Conjugation Buffer Set (Cat. No. 558556) to covalently link water-soluble proteins to the surface of the beads via a sulfo-SMCC chemistry. The beads, once linked to a protein, can be used as a capture matrix for analytes of interest. This vial of Functional Beads contains enough beads for 6,000 tests (6 - 12 conjugations of 500 - 1000 tests each). It is recommended that the user carefully read the BD CBA Functional Bead Conjugation Buffer Set instruction manual before purchasing or using the BD CBA Functional Beads.

## Preparation and Storage

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

## Application Notes

## Application

Array	Routinely Tested
-------	------------------

## Recommended Assay Procedure:

For details on the specific conjugation procedure used to covalently link proteins to the BD CBA Functional Beads, guidelines for developing bead-based immunoassays and how to confirm the presence of immunoglobulins conjugated to the beads, refer to the BD CBA Functional Bead Conjugation Buffer Set instruction manual.

## Suggested Companion Products

Catalog Number	Name	Size	Clone
558556	Functional Bead Conjugation Buffer Set	15 reactions	(none)
558550	PE Goat anti-Mouse Ig Detector	25 tests	Polyclonal
558551	PE Goat anti-Rat Ig Detector	25 tests	Polyclonal
558553	PE Mouse anti-Rabbit IgG Detector	25 tests	C101-359

## Product Notices

- Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- Please refer to [www.bdbiosciences.com/pharmingen/protocols](http://www.bdbiosciences.com/pharmingen/protocols) for technical protocols.

## BD Biosciences

[bdbiosciences.com](http://bdbiosciences.com)

United States 877.232.8995 Canada 888.268.5430 Europe 32.53.720.550 Japan 0120.8555.90 Asia Pacific 65.6861.0633 Latin America/Caribbean 0800.771.7157

For country-specific contact information, visit [bdbiosciences.com/how\\_to\\_order/](http://bdbiosciences.com/how_to_order/)

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2011 BD

