

Technical Data Sheet


PE Rat Anti-Mouse CD131

Product Information

Material Number:	559920
Alternate Name:	β IL-3R, β c
Size:	0.2 mg
Concentration:	0.2 mg/ml
Clone:	JORO50
Immunogen:	Bone Marrow-Derived C4-77 Pro-T Lymphocyte Clone
Isotype:	Rat IgG1, κ
Reactivity:	QC Testing: Mouse
Storage Buffer:	Aqueous buffered solution containing $\leq 0.09\%$ sodium azide.

Description

The JORO50 antibody reacts with the extracellular regions of the mouse CD131 (β IL-3 (AIC2A) and β c (AIC2B) cytokine receptor subunits). A variety of mouse cell types, including multipotential hematopoietic stem cells, mast cells, megakaryocytes, eosinophils, erythroblasts, pre-B cells, and osteoclasts, express β IL-3 and β c subunits. Either β IL-3 or β c can combine with the IL-3 α chain to form two distinct, functional (i.e., signaling), high affinity receptors for mouse IL-3. The β IL-3 subunit by itself can bind mouse IL-3 with low affinity whereas the β c subunit can not. β c (but not β IL-3) can combine with the IL-5R α (CD125) and GM-CSFR α (CD116) subunits to form high affinity, signaling receptors for mouse IL-5 or GM-CSF, respectively. The immunogen used to generate the JORO50 hybridoma was the bone marrow-derived, C4-77 pro-T lymphocyte clone.



High level expression of β IL-3 and β c by Mac-1⁺ bone marrow cells. BALB/c bone marrow cells (BMC) were stained (20 minutes at 4°C) with phycoerythrin-conjugated JORO50 antibody (1 μ g mAb/10⁶ cells; Cat. No. 559920) and FITC-anti-CD11b (Cat. No. 553310).

Preparation and Storage

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated with R-PE under optimum conditions, and unconjugated antibody and free PE were removed.

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

Application Notes

Application

Flow cytometry	Routinely Tested
----------------	------------------

Recommended Assay Procedure:

Immunofluorescent Staining and Flow Cytometric Analysis: The phycoerythrin-conjugated form of JORO50 (Cat. No. 559920) can be used for the immunofluorescent staining (≤ 1 μ g antibody/10⁶ cells) and flow cytometric analysis of normal mouse cells or cell lines to measure their expressed levels of β IL-3 and β c. An appropriate phycoerythrin-conjugated immunoglobulin isotype control is clone R3-34 (Cat. No. 554685). PE-JORO50 can be used for the direct immunofluorescent staining and flow cytometric analysis of cells that express β IL-3 or β c as described in the image legend.

BD Biosciences

bdbiosciences.com

United States	Canada	Europe	Japan	Asia Pacific	Latin America/Caribbean
877.232.8995	888.259.0187	32.53.720.550	0120.8555.90	65.6861.0633	55.11.5185.9995

For country-specific contact information, visit bdbiosciences.com/how_to_order/

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2008 BD



Suggested Companion Products

Catalog Number	Name	Size	Clone
554685	PE Rat IgG1, κ Isotype Control	0.1 mg	R3-34
553310	FITC Rat Anti-Mouse CD11b	0.5 mg	M1/70

Product Notices

1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
2. Please refer to www.bdbiosciences.com/pharming/en/protocols for technical protocols.
3. For fluorochrome spectra and suitable instrument settings, please refer to our Fluorochrome Web Page at www.bdbiosciences.com/colors.
4. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.