# **Technical Data Sheet**

# V450 Hamster IgG2, κ Isotype Control

#### **Product Information**

560563 Material Number anti-KLH Alternate Name:  $0.1 \, \text{mg}$ 0.2 mg/ml **Concentration:** B81-3 Clone:

Immunogen: Keyhole Limpet Hemocyanin Armenian Hamster IgG2, κ Isotype:

Storage Buffer: Aqueous buffered solution containing protein stabilizer and ≤0.09% sodium

azide.

#### Description

The B81-3 monoclonal antibody is specific for keyhole limpet hemocyanin (KLH), an antigen not expressed by mammalian cells or cell lines.

The antibody is conjugated to BD Horizon<sup>TM</sup> V450, which has been developed for use in multicolor flow cytometry experiments and is available exclusively from BD Biosciences. It is excited by the Violet laser Ex max of 406 nm and has an Em Max at 450 nm. Conjugates with BD Horizon<sup>TM</sup> V450 can be used in place of Pacific Blue<sup>TM</sup> conjugates.

## **Preparation and Storage**

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated with BD Horizon<sup>TM</sup> V450 under optimum conditions, and unreacted BD Horizon<sup>TM</sup> V450 was removed.

## Application Notes

#### Application

| F | low cytometry  | Routinely Tested |
|---|----------------|------------------|
| Ι | sotype control | Routinely Tested |

#### **Product Notices**

- Since applications vary, each investigator should titrate the reagent to obtain optimal results.
- Although hamster immunoglobulin isotypes have not been well defined, BD Biosciences Pharmingen has grouped Armenian and Syrian hamster IgG monoclonal antibodies according to their reactivity with a panel of mouse anti-hamster IgG mAbs. A table of the hamster IgG groups, Reactivity of Mouse Anti-Hamster Ig mAbs, may be viewed at http://www.bdbiosciences.com/documents/hamster\_chart\_11x17.pdf.
- BD Horizon<sup>TM</sup> V450 has a maximum absorption of 406 nm and maximum emission of 450 nm. Before staining with this reagent, please confirm that your flow cytometer is capable of exciting the fluorochrome and discriminating the resulting fluorescence.
- Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- Pacific Blue<sup>TM</sup> is a trademark of Molecular Probes, Inc., Eugene, OR.
- For fluorochrome spectra and suitable instrument settings, please refer to our Fluorochrome Web Page at www.bdbiosciences.com/colors.
- Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.

## References

Mendrick DL, Kelly DM. Temporal expression of VLA-2 and modulation of its ligand specificity by rat glomerular epithelial cells in vitro. Lab Invest. 1993; 69(6):690-702. (Immunogen)

# **BD Biosciences**

bdbiosciences.com

United States 877.232.8995 888.268.5430 32.53.720.550 0120.8555.90 65.6861.0633 0800.771.7157

For country-specific contact information, visit bdbiosciences.com/how\_to\_order/

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.
For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.
BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2011 BD

