## **Technical Data Sheet**

# Alexa Fluor® 647 Rat IgM, κ Isotype Control

#### **Product Information**

Material Number: Size: Clone: Immunogen: Isotype: Reactivity: Storage Buffer:

# Description

The R4-22 antibody has unknown specificity.

This antibody is routinely tested by flow cytometric analysis. Other applications were tested at BD Biosciences Pharmingen during antibody development only or reported in the literature.

560892

0.1 mg

R4-22

Pooled Mouse Ig

QC Testing: Negative Human and Mouse

Aqueous buffered solution containing ≤0.09% sodium azide.

Rat IgM, ĸ

### **Preparation and Storage**

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated to Alexa Fluor® 647 under optimum conditions, and unreacted Alexa Fluor® 647 was removed.

### **Application Notes**

Application						
	Flow cytometry	Routinely Tested				
	Isotype control	Routinely Tested				

#### Product Notices

- 1. This reagent has been pre-diluted for use at the recommended Volume per Test. We typically use  $1 \times 10^{6}$  cells in a 100-µl experimental sample (a test).
- 2. Alexa Fluor® 647 fluorochrome emission is collected at the same instrument settings as for allophycocyanin (APC).
- 3. Alexa Fluor® is a registered trademark of Molecular Probes, Inc., Eugene, OR.
- 4. The Alexa Fluor®, Pacific Blue™, and Cascade Blue® dye antibody conjugates in this product are sold under license from Molecular Probes, Inc. for research use only, excluding use in combination with microarrays, or as analyte specific reagents. The Alexa Fluor® dyes (except for Alexa Fluor® 430), Pacific Blue™ dye, and Cascade Blue® dye are covered by pending and issued patents.
- 5. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- 6. For fluorochrome spectra and suitable instrument settings, please refer to our Fluorochrome Web Page at www.bdbiosciences.com/colors.
- 7. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.

BD Biosciences bdbiosciences.com								
877.232.8995	888.268.5430	32.53.720.550	0120.8555.90	65.6861.0633	0800.771.7157			
For country-spe	cific contact inf	ormation, visit I	odbiosciences.co	m/how_to_orde	r/			
of any patents. BL use of our produc product or as a co written authoriza For Research Use	) Biosciences will n ts. Purchase does n mponent of anoth tion of Becton Dick Dnly. Not for use in	ot be held responsil ot include or carry a er product. Any use inson and Company diagnostic or thera	ble for patent infring any right to resell or of this product oth y is strictly prohibite peutic procedures.	gement or other vio transfer this produce than the permitted.	e the above product in violation lations that may occur with the ct either as a stand-alone ed use without the express 11 BD			