Technical Data Sheet

Alexa Fluor[®] 488 Mouse anti-Human IL-17F

Product Information

Material Number:	561331
Alternate Name:	Interleukin-17F; cytokine ML-1; ML-1; IL-24
Size:	100 tests
Vol. per Test:	5 µl
Clone:	O33-782
Immunogen:	Human IL-17F
Isotype:	Mouse IgG1, ĸ
Reactivity:	QC Testing: Human
Storage Buffer:	Aqueous buffered solution containing BSA and ≤0.09% sodium azide.

Description

The O33-782 monoclonal antibody specifically binds to Interleukin-17F (IL-17F). IL-17F is a member of the IL-17 family of cytokines. IL-17F is encoded by the *IL17F* gene located in chromosome 6 (location: 6p12). IL-17F is a proinflammatory cytokine that is produced by activated T cells including differentiated CD4+ T helper 17 (Th17) cells. Activated Th17 cells can express disulfide-linked IL-17F and IL-17A homodimers as well as IL-17A/IL-17F heterodimers. These IL-17 dimers act by binding to and signaling through IL-17 receptor complexes (IL-17R). IL-17R are comprised of transmembrane IL-17RA and IL-17-RC protein subunits that are expressed by a variety of target cells including epithelial and endothelial cells, keratinocytes, fibroblasts, and granulocytes. IL-17F can induce target cells to produce proinflammatory cytokines such as IL-1 β , IL-6, G-CSF, GM-CSF, and TNF and chemokines including CXCL1/Gro- α , CXCL2/Gro- β , and CXCL8/IL-8 that attract and activate leukocytes, eg, neutrophils. Th17 and other IL-17F-producing cells play protective roles in the clearance of extracellular pathogens, including bacteria and fungi. IL-17F can also play adverse roles in inflammation associated with asthma and autoimmune diseases.



Flow cytometric analysis of IL-17F expression by resting and activated human peripheral blood CD4+ T cells and Th17 polarized cells. Human peripheral blood mononuclear cells were either unstimulated (Left Panels) or stimulated with Phorbol 12-Myristate 13-Acetate (PMA; Sigma P-8139) plus lonomycin (Sigma I-0634) in the presence of BD GolgiStop™ Protein Transport Inhibitor (Cat. No. 554724) for 5 hours (Middle Panels) or were cultured in Th17 polarization conditions and restimulated with PMA and lonomycin in the presence of BD GolgiStop™ for 5 hours (Right Panels). Cells were then fixed and permeabilized using BD Cytofix/Cytoperm™ reagents (Cat. No. 554714) followed by staining with Alexa Fluor® 488 Mouse anti-Human IL-17F (Cat. No. 561331), PerCP-Cy5.5 Mouse anti-Human CD4 (Cat. No. 341654), PE Mouse Anti-Human IFN-γ (Cat. No. 554701) and Alexa Fluor® 647 Mouse anti-Human IL-17A (Cat. No. 56490). Two-color flow cytometric dot plots showing the correlated expression patterns of IL-17F versus IL-17A or IFN-γ were derived from CD4+ gated events with the forward and side light-scatter characteristics of intact lymphocytes. Flow cytometry was performed using a BD™ LSR II Flow Cytometry System. Other compatible fixation and permeabilized in the "Recommended Assay Procedure."

BD Biosciences

bdbiosciences.com							
United States	Canada	Europe	Japan	Asia Pacific	Latin America/Caribbean		
877.232.8995	888.268.5430	32.53.720.550	0120.8555.90	65.6861.0633	0800.771.7157		
For country-specific contact information, visit bdbiosciences.com/how_to_order/							
Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation							
of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the							
use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone							
product or as a component of another product. Any use of this product other than the permitted use without the express							
written authorization of Becton Dickinson and Company is strictly prohibited.							
For Research Use C	only. Not for use in	diagnostic or thera	peutic procedures. No	ot for resale.			
BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2011 BD							



Preparation and Storage

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography. The antibody was conjugated to Alexa Fluor® 488 under optimum conditions, and unreacted Alexa Fluor® 488 was removed. Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

Application Notes

Application

11		
Intracellular staining (flow cytometry)	Routinely Tested	

Recommended Assay Procedure:

This antibody conjugate is suitable for intracellular staining of human peripheral blood mononuclear cells using BD Cytofix/Cytoperm[™] reagents, BD Pharmingen[™] Human FoxP3 Buffer Set or BD[™] Phosflow fixation and permeabilization buffers (Fix buffer I with Perm/Wash Buffer I, Perm Buffer II, or Perm Buffer III).

Suggested Companion Products

Catalog Number	Name	Size	Clone	
554724	Protein Transport Inhibitor (Containing Monensin)	0.7 ml	(none)	
554714	BD Cytofix/Cytoperm [™] Fixation/Permeablization Kit	250 tests	(none)	
341654	PerCP-Cy5.5 Mouse anti-Human CD4	50 tests	SK3	
554701	PE Mouse Anti-Human IFN-γ	0.1 mg	B27	
560490	Alexa Fluor® 647 Mouse anti-Human IL-17A	100 tests	N49-653	
557721	Alexa Fluor® 488 Mouse IgG1 κ Isotype Control	100 tests	MOPC-21	
557870	Fix Buffer I	250 ml	(none)	
557885	Perm/Wash Buffer I	125 ml	(none)	
558052	Perm Buffer II	125 ml	(none)	
558050	Perm Buffer III	125 ml	(none)	
560098	Human FoxP3 Buffer Set	100 tests	(none)	

Product Notices

- 1. This reagent has been pre-diluted for use at the recommended Volume per Test. We typically use 1×10^{6} cells in a 100-µl experimental sample (a test).
- 2. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.
- 3. The Alexa Fluor®, Pacific Blue™, and Cascade Blue® dye antibody conjugates in this product are sold under license from Molecular Probes, Inc. for research use only, excluding use in combination with microarrays, or as analyte specific reagents. The Alexa Fluor® dyes (except for Alexa Fluor® 430), Pacific Blue™ dye, and Cascade Blue® dye are covered by pending and issued patents.
- 4. Alexa Fluor® 488 fluorochrome emission is collected at the same instrument settings as for fluorescein isothiocyanate (FITC).
- 5. For fluorochrome spectra and suitable instrument settings, please refer to our Fluorochrome Web Page at www.bdbiosciences.com/colors.
- 6. Alexa Fluor® is a registered trademark of Molecular Probes, Inc., Eugene, OR.
- 7. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- 8. An isotype control should be used at the same concentration as the antibody of interest.
- 9. Source of all serum proteins is from USDA inspected abattoirs located in the United States.
- 10. Licensed for Research Purposes Only. Commercial use requires license from Boyce Thompson Institute for Plant Research.

References

Fouser LA, Wright JF, Dunussi-Joannopoulos K, Collins M. Th17 cytokines and their emerging roles in inflammation and autoimmunity. *Immunol Rev.* 2008; 226:87-102. (Biology)

Shen F, Gaffen SL. Structure-function relationships in the IL-17 receptor: implications for signal transduction and therapy. *Cytokine*. 2008; 41(2):92-104. (Biology) Starnes T, Robertson MJ, Sledge G, et al.. Cutting edge: IL-17F, a novel cytokine selectively expressed in activated T cells and monocytes, regulates angiogenesis and endothelial cell cytokine production. *J Immunol.* 2001; 167(8):4137-4140. (Biology)

Wang YH, Liu YJ. The IL-17 cytokine family and their role in allergic inflammation. Curr Opin Immunol. 2008; 20(6):697-702. (Biology)

Wright JF, Bennett F, Li B, et al. The human IL-17F/IL-17A heterodimeric cytokine signals through the IL-17RA/IL-17RC receptor complex. J Immunol. 2008; 181(4):2799-2805. (Biology)

Wright JF, Guo Y, Quazi A, et al. Identification of an interleukin 17F/17A heterodimer in activated human CD4+ T cells. J Biol Chem. 2007; 282(18):13447-13455. (Biology)