Technical Data Sheet

BV421 Mouse Anti-Human IL-2

Product Information

Material Number: 562914

Alternate Name: IL2; Interleukin-2; T-cell growth factor; TCGF

 Size:
 50 tests

 Vol. per Test:
 5 μ l


 Clone:
 5344.111


Storage Buffer: Aqueous buffered solution containing BSA and ≤0.09% sodium azide.

Description

The 5344.111 antibody specifically binds to the multifunctional cytokine, human interleukin-2 (IL-2). IL-2 is produced by activated T cells and has multiple functions that can affect the growth, proliferation and differentiation of many different target cell types including T cells, B cells, NK cells, monocytes and macrophages. The immunogen used to generate the 5344.111 hybridoma was affinity-purified, recombinant human IL-2 protein.

The antibody was conjugated to BD HorizonTM BV421 which is part of the BD HorizonTM Brilliant VioletTM family of dyes. With an Ex Max of 407-nm and Em Max at 421-nm, BD HorizonTM BV421 can be excited by the violet laser and detected in the standard Pacific BlueTM filter set (eg, 450/50-nm filter). BD HorizonTM BV421 conjugates are very bright, often exhibiting a 10 fold improvement in brightness compared to Pacific BlueTM conjugates.





Multicolor flow cytometric analysis of IL-2 expressed in stimulated human peripheral blood mononuclear cells. HiCK-1 Human Cytokine Positive Control Cells (Cat. No. 555061) were permeabilized with BD Perm/Wash™ Buffer (Cat. No. 554723). The cells were then stained with FITC Mouse Anti-Human CD3 antibody (Cat. No. 555916) and either BD Horizon™ BV421 Mouse IgG1, κ Isotype Control (Cat No. 562438, Left Panel) or with BD Horizon™ BV421 Mouse Anti-Human IL-2 antibody (Cat No. 562914, Right Panel). Two-color flow cytometric dot plots show IL-2 expression (or Ig Isotype control staining) versus CD3 derived from gated events with the forward and side light-scatter characteristics of intact lymphocytes. Flow cytometry was performed using a BD™ LSR II Flow Cytometer System.

Preparation and Storage

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated with BD Horizon™ BV421 under optimum conditions, and unconjugated antibody and free BD Horizon™ BV421 were removed.

Application Notes

Application

Intracellular staining (flow cytometry)

Routinely Tested

BD Biosciences

bdbiosciences.com

 United States
 Canada
 Europe
 Japan
 Asia Pacific
 Latin America/Caribbean

 877.232.8995
 800.979.9408
 32.53.720.550
 0120.8555.90
 65.6861.0633
 55.11.5185.9995

For country contact information, visit bdbiosciences.com/contact

Conditions: The information disclosed herein is not to be constructed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be help responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited. For Research Use Only, Not for use in diagnostic or therapeutic procedures. Not for resale.

Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2011 BD



562914 Rev. 1 Page 1 of 2

Suggested Companion Products

Catalog Number	Name Name	Size	Clone
554656	Stain Buffer (FBS)	500 ml	(none)
562438	BV421 Mouse IgG1, k Isotype Control	50 μg	X40
555061	HiCK-1 Human Cytokine Positive Control Cells	1.0 ml	(none)
554723	Perm/Wash Buffer	100 ml	(none)
555916	FITC Mouse Anti-Human CD3	100 tests	UCHT1

Product Notices

- This reagent has been pre-diluted for use at the recommended Volume per Test. We typically use 1 × 10⁶ cells in a 100-µl experimental sample (a test).
- 2. An isotype control should be used at the same concentration as the antibody of interest.
- 3. Source of all serum proteins is from USDA inspected abattoirs located in the United States.
- 4. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.
- Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- 6. For fluorochrome spectra and suitable instrument settings, please refer to our Multicolor Flow Cytometry web page at www.bdbiosciences.com/colors.
- 7. Pacific BlueTM is a trademark of Molecular Probes, Inc., Eugene, OR.
- 8. Brilliant VioletTM 421 is a trademark of Sirigen.
- 9. This product is manufactured and sold under license from Pestka Biomedical Laboratories, Inc. (d/b/a PBL InterferonSource) and may be used solely as indicated. This product may not be resold or incorporated in any manner into another product for resale. Any use for therapeutics is strictly prohibited. This product is covered by U.S. Patent No. 5,597,901 and Bulgarian Patent No. BG1895.

References

Liao W, Lin JX, Leonard WJ. IL-2 family cytokines: new insights into the complex roles of IL-2 as a broad regulator of T helper cell differentiation. Curr Protoc Immunol. 2011; 23(5):598-604. (Biology)

Malek TR. The biology of interleukin-2. Annu Rev Immunol. 2008; 26:453-479. (Biology)

Waldmann TA. The biology of interleukin-2 and interleukin-15: implications for cancer therapy and vaccine design. Nat Rev Immunol. 2006; 6(8):595-601. (Biology)

BD Biosciences

bdbiosciences.com

 United States
 Canada
 Europe
 Japan
 Asia Pacific
 Latin America/Caribbean

 877.232.8995
 800.979.9408
 32.53.720.550
 0120.8555.90
 65.6861.0633
 55.11.5185.9995

For country contact information, visit bdbiosciences.com/contact

Conditions: The information disclosed herein is not to be constructed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be help responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited. For Research Use Only, Not for use in diagnostic or therapeutic procedures. Not for resale.

Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2011 BD



562914 Rev. 1 Page 2 of 2