Technical Data Sheet

BV650 Rat IgG2a, κ Isotype Control

Product Information

Material Number:	563236
Size:	50 µg
Concentration:	0.2 mg/ml
Clone:	R35-95
Immunogen:	Mouse Pooled Immunoglobulin
Isotype:	Rat (LOU) IgG2a, κ
Storage Buffer:	Aqueous buffered solution containing BSA and ≤0.09% sodium azide.

Description

The R35-95 hybridoma was generated by hybridization of Y3 myeloma cells with spleen cells from LOU rats immunized with mouse immunoglobulins. The R35-95 hybridoma produces rat IgG2a, κ immunoglobulin that has no measurable reactivity with mouse immunoglobulins. The R35-95 immunoglobulin was selected as an isotype control following screening for low background binding on a variety of mouse and human tissues.

The antibody was conjugated to BD Horizon[™] BV650 which is part of the BD Horizon[™] Brilliant Violet[™] family of dyes. This dye is a tandem fluorochrome of BD Horizon[™] BV421 with an Ex Max of 405-nm and an acceptor dye with an Em Max at 650-nm. BD Horizon[™] BV650 can be excited by the violet laser and detected in a filter used to detect APC-like dyes (eg, 660/20-nm filter). Due to the excitation and emission characteristics of the acceptor dye, there will be spillover into the APC and Alexa Fluor® 700 detectors. However, the spillover can be corrected through compensation as with any other dye combination.

Preparation and Storage

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated with BD Horizon[™] BV650 under optimum conditions, and unconjugated antibody and free BD Horizon[™] BV650 were removed.

Application Notes

Ар	plication	
F	low cytometry	Routinely Tested
Is	sotype control	Routinely Tested

Suggested Companion Products

Catalog Number	Name	Size	Clone
554656	Stain Buffer (FBS)	500 ml	(none)
— • • • • •			

Product Notices

- 1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
- 2. An isotype control should be used at the same concentration as the antibody of interest.
- 3. Brilliant Violet[™] 650 is a trademark of Sirigen.
- 4. Alexa Fluor® is a registered trademark of Molecular Probes, Inc., Eugene, OR.
- 5. Source of all serum proteins is from USDA inspected abattoirs located in the United States.
- 6. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- 7. For fluorochrome spectra and suitable instrument settings, please refer to our Multicolor Flow Cytometry web page at www.bdbiosciences.com/colors.
- 8. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.

BD Biosciences bdbiosciences.com							
United States 877.232.8995	Canada 800.979.9408	Europe 32.53.720.550	Japan 0120.8555.90	Asia Pacific 65.6861.0633	Latin America/Caribbean 55.11.5185.9995		
For country con	ntact informatio	on, visit bdbiosci	ences.com/conta	ct			
of any patents. Bl use of our produc product or as a co written authoriza	D Biosciences will n ts. Purchase does n imponent of anoth tion of Becton, Dic	ot be help responsil ot include or carry a er product. Any use kinson and Compan	ble for patent infrin any right to resell o	gement or other vio r transfer this producer than the permitte d.	use the above product in violatio lations that may occur with the ct either as a stand-alone ed use without the express		

Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2011 BD

