Technical Data Sheet

BV711 Rat IgG1, κ Isotype Control

Product Information

Material Number: 563283

Alternate Name: IgG1, kappa; Ighg1, Igkc

 Size:
 50 μg

 Concentration:
 0.2 mg/ml

 Clone:
 R3-34

Immunogen: Mouse immunoglobulin

Isotype: Rat IgG1, κ

Storage Buffer: Aqueous buffered solution containing BSA and ≤0.09% sodium azide.

Description

Following immunization of a rat with mouse immunoglobulin (Ig), the Ig from the R3-34 hybridoma was identified as a non-reactive clone. The R3-34 immunoglobulin was selected as an Ig isotype control following screening for low background staining on a variety of mouse and human cells and tissues.

The antibody was conjugated to BD HorizonTM BV711 which is part of the BD HorizonTM Brilliant VioletTM family of dyes. This dye is a tandem fluorochrome of BD HorizonTM BV421 with an Ex Max of 405-nm and an acceptor dye with an Em Max at 711-nm. BD HorizonTM BV711 can be excited by the violet laser and detected in a filter used to detect CyTM5.5 / Alexa Fluor® 700-like dyes (eg, 712/20-nm filter). Due to the excitation and emission characteristics of the acceptor dye, there may be moderate spillover into the Alexa Fluor® 700 and PerCP-CyTM5.5 detectors. However, the spillover can be corrected through compensation as with any other dye combination.

Preparation and Storage

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated with BD Horizon™ BV711 under optimum conditions, and unconjugated antibody and free BD Horizon™ BV711 were removed.

Application Notes

Application

rr ·······		
Flow cytometry	Routinely Tested	
Isotype control	Routinely Tested	

Suggested Companion Products

Catalog Number	Name	Size	Clone
554656	Stain Buffer (FBS)	500 ml	(none)

Product Notices

- 1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
- 2. Source of all serum proteins is from USDA inspected abattoirs located in the United States.
- 3. An isotype control should be used at the same concentration as the antibody of interest.
- 4. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.
- 5. Alexa Fluor® is a registered trademark of Molecular Probes, Inc., Eugene, OR.
- 6. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- For fluorochrome spectra and suitable instrument settings, please refer to our Multicolor Flow Cytometry web page at www.bdbiosciences.com/colors.
- 8. Cy is a trademark of Amersham Biosciences Limited.
- 9. Brilliant Violet™ 711 is a trademark of Sirigen.

BD Biosciences

bdbiosciences.com

United States Canada Europe Japan Asia Pacific Latin America/Caribbean 877.232.8995 800.979.9408 32.53.720.550 0120.8555.90 65.6861.0633 55.11.5185.9995

For country contact information, visit bdbiosciences.com/contact

Conditions: The information disclosed herein is not to be constructed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be help responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited. For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale. Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2011 BD



Page 1 of 1