



61-9111 2 ml
FITC-Rabbit Anti-
Porcine IgG (H+L)

FITC-RABBIT ANTI-PORCINE IgG (H+L)

INVITROGEN AFFINITY PURIFIED ANTIGEN SPECIFIC ANTIBODY

FORM: Liquid. The product is in 10 mM phosphate buffered saline, pH 7.4 (PBS), containing 1% BSA, 40% glycerol and 0.1% sodium azide as preservative.

AFFINITY PURIFIED RABBIT ANTI-PORCINE IgG (H+L) SPECIFIC Antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and is free of unconjugated antibodies or uncoupled FITC.

APPLICATIONS: Fluorescent-antibody conjugates can be used to localize or to label antigens in cells, tissue sections or other fluorescent immunoassay applications.

SUGGESTED WORKING DILUTION: 1:20 to 1:50 with PBS containing 10% normal rabbit serum. Working dilution for specific application should be determined by the investigator to obtain the best conditions.

STORAGE: Store at 4°C.

REMARKS: For research use only.

Research Use Only

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288

E-mail: techsupport@invitrogen.com

PI619111

(Rev 12/08) DCC-08-1089

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.



61-9111 2 ml
FITC-Rabbit Anti-
Porcine IgG (H+L)

FITC-RABBIT ANTI-PORCINE IgG (H+L)

INVITROGEN AFFINITY PURIFIED ANTIGEN SPECIFIC ANTIBODY

FORM: Liquid. The product is in 10 mM phosphate buffered saline, pH 7.4 (PBS), containing 1% BSA, 40% glycerol and 0.1% sodium azide as preservative.

AFFINITY PURIFIED RABBIT ANTI-PORCINE IgG (H+L) SPECIFIC Antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and is free of unconjugated antibodies or uncoupled FITC.

APPLICATIONS: Fluorescent-antibody conjugates can be used to localize or to label antigens in cells, tissue sections or other fluorescent immunoassay applications.

SUGGESTED WORKING DILUTION: 1:20 to 1:50 with PBS containing 10% normal rabbit serum. Working dilution for specific application should be determined by the investigator to obtain the best conditions.

STORAGE: Store at 4°C.

REMARKS: For research use only.

Research Use Only

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288

E-mail: techsupport@invitrogen.com

PI619111

(Rev 12/08) DCC-08-1089

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.