7500 Fast Real-Time PCR Systems Spectral Calibration Kit I

Catalog Number 4360788

Pub. No. 4362281 Rev. C

WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from **thermofisher.com/support**.

Contents and storage

Contents	Amount	Storage
Background Plate sealed with an optical cover	1	–25°C to –15°C
Spectral Calibration Plates sealed with optical covers	7	
Region of Interest (ROI) Calibration Plate sealed with an optical cover	1	

Related documentation

For detailed information on instrument setup and the calibration process, refer to the *Applied Biosystems*[™] 7300/7500/7500 Fast Real-Time PCR System Installation and Maintenance Guide (Pub. no. 4347828).

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at **www.thermofisher.com/us/en/home/global/terms-and-conditions.html**. If you have any questions, please contact Life Technologies at **www.thermofisher.com/support**.

For support visit thermofisher.com/support or email techsupport@lifetech.com

The information in this guide is subject to change without notice.

DISCLAIMER

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Corporate entity: Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

©2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

