

PI AHF0082

## Mouse (monoclonal) Anti-Cyclin D1 Unconjugated

## **PRODUCT ANALYSIS SHEET**

| Catalog Number:                  | AHF0082   |  |
|----------------------------------|---|--|
| Lot Number:                      | See product label   |  |
| Expiration Date:                 | See product label   |  |
| Quantity/Volume:                 | 100 µg/0.5 mL   |  |
| Clone Number:                    | DCS-6   |  |
| Isotype:                         | IgG2a kappa (mouse)   |  |
| Form of Antibody:                | Purified immunoglobulin in phosphate buffered saline, pH 7.4, with 0.2% bovine serum albumin.   |  |
| Preservation:                    | 0.09% sodium azide (Caution: sodium azide is a poisonous and hazardous substance.<br>Handle with care and dispose of properly.)   |  |
| Purification:                    | Purified from ascites by Protein A affinity chromatography.   |  |
| Immunogen:                       | Human recombinant full length cyclin D1 protein.  |  |
| Myeloma/Fusion<br>Partners:      | Produced by fusion between female BALB/c splenocytes and mouse myeloma NS-2 cells.  |  |
| Specificity:                     | This monoclonal antibody recognizes a protein of 36 kDa, identified as cyclin D1 (also known as PRAD-1 or bcl-1). Cyclins are regulatory subunits of the cyclin-dependent kinases (cdk), which control transition at different specific phases of the cell cycle. In mammalian cells, cyclin D1 plays a critical role in $G_1$ to S transition. Cyclin D1 is a putative proto-oncogene overexpressed in a wide variety of human neoplasms.  |  |
| Species Reactivity:              | Human, monkey, mouse, rat, and dog cyclin D1. Other species were not tested.  |  |
| Applications:                    | This antibody is suitable for use in immunoprecipitation (both non-denatured and denatured), immunofluorescence, and is recommended for immunohistochemistry with cryostat sections, and formalin-fixed or methacarn-fixed paraffin embedded sections. Staining of formalin-fixed paraffin embedded tissue sections requires boiling in 10 mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at room temperature for 20 minutes prior to addition of antibody. This antibody is also highly recommended for Western blotting, biological blockade, and flow cytometry. |  |
| Suggested Working<br>Dilutions:  | For immunoprecipitation, use 10 $\mu$ L/mg of protein lysate; for immunohistochemistry, use 1:50-1:100; for flow cytometry, use 0.25 $\mu$ g to stain 10 <sup>6</sup> human PBMCs fixed in 100% methanol; and for Western blotting, 1:250-1:500 is recommended. The optimal antibody concentration should be determined for each specific application.  |  |
| Recommended Positive<br>Control: | T47D, ZR75, BT474, SKBR3, MCF-7, MDA-MB-453, U-2-OS, HT29, Ramos, Jurkat, KNRK, C32 or A431. About 50-70% of mantle cell lymphomas and 40% of breast carcinomas are positive.   |  |

For Research Use Only. CAUTION: Not for human or animal therapeutic or diagnostic use.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 • E-mail: techsupport@invitrogen.com

(Rev 04/10) DCC-10-0604

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, <u>www.invitrogen.com</u>). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

Store at 2-8°C.

**References:** 

Storage:

PI AHF0082

Bartkova, J., *et al.* (1996) The p16-cyclin D/Cdk4-pRb pathway as a functional unit frequently altered in melanoma pathogenesis. Cancer Research 56:5475-5483.

Bartkova, J., *et al.* (1995) Cyclin D1 oncoprotein aberrantly accumulates in malignancies of diverse histogenesis. Oncogene 10:775-778.

Bartek, J., *et al.* (1996) Abundance and subcellular localization of cyclin D3 in human tumours. Int. J. Cancer 65:323-327.

De Boer, C.J., *et al.* (1995) Cyclin D1 protein analysis in the diagnosis of mantle cell lymphoma. Blood 86(7):2715-2723.

Depoortere, F., *et al.* (1998) A requirement for cyclin D3-cyclin-dependent kinase (cdk)-4 assembly in the cyclic adenosine monophosphate-dependent proliferation of thyrocytes. The Journal of Cell Biology 140(6):1427-1439.

Jadayel, D.M., *et al.* (1997) Potential role for concurrent abnormalities of the cyclin D1, p16<sup>CDKN2</sup> and p15<sup>CDKN2B</sup> genes in certain B cell non-Hodgkin's lymphomas. Functional studies in a cell line (Granta 519). Leukemia 11:64-72.

Lukas, J., *et al.* (1996) Deregulated expression of E2F family members induces S-phase entry and overcomes p16<sup>INK4A</sup>-mediated growth suppression. Molecular and Cellular Biology 16(3):1047-1057.

Lukas, J., *et al.* (1994) Cyclin D1 protein oscillates and is essential for cell cycle progression in human tumour cell lines. Oncogene 9:707-718.

Niculescu III, A., *et al.* (1998) Effects of p21Cip1/Waf1 at both the G1/S and the G2/M cell cycle transitions: pRb is a critical determinant in blocking DNA replication and in preventing endoreduplication. Molecular and Cellular Biology 18(1):629-643.

Ott, M.M., *et al.* (1997) Cyclin D1 expression in mantle cell lymphoma is accompanied by downregulation of cyclin D3 and is not related to the proliferative activity. Blood 90(8):3154-3159.

Pokrovskaja, K., *et al.* (1996) Phenotype-related differences in the expression of D-type cyclins in human B cell-derived lines. Cell Growth and Differentiation 7:1723-1732.

Steiner, P., *et al.* (1995) Identification of a Myc-dependent step during the formation of active  $G_1$  cyclin-cdk complexes. The EMBO Journal 14(19):4814-4826.

| Explanation of symbols            |   |             |  |  |
|-----------------------------------|---|-------------|--|--|
| Symbol                            | Description   | Symbol      | Description                                  |  |
| REF                               | Catalogue Number  | LOT         | Batch code                                   |  |
| RUO                               | Research Use Only   | IVD         | In vitro diagnostic medical device           |  |
| X                                 | Use by  | ł           | Temperature limitation                       |  |
| 1                                 | Manufacturer  | EC REP      | European Community authorised representative |  |
| [-]                               | Without, does not contain   | [+]         | With, contains                               |  |
| er til konstante<br>from Light    | Protect from light  | $\triangle$ | Consult accompanying documents               |  |
| $\begin{bmatrix} 1 \end{bmatrix}$ | Directs the user to consult instructions for use (IFU), accompanying the product. |             |  |  |

## For Research Use Only. CAUTION: Not for human or animal therapeutic or diagnostic use.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 • E-mail: techsupport@invitrogen.com

(Rev 04/10) DCC-10-0604

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, <u>www.invitrogen.com</u>). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.