

# **Incubation Trays**

Catalog no. LC2102

Store at room temperature

## Description

The Incubation Tray is designed to facilitate the processing of western blots. The tray is reusable, optically clear with a sealable lid, and is suitable for most western detection and immunostaining applications.

**Important:** The Incubation Tray is **NOT** heat-resistant and cannot be heated in a microwave oven or autoclaved.

## **Specifications**

Quantity: 8/pk

Tray Dimensions: 10 cm x 14 cm

Volume Capacity:

Tray Material: APET (amorphous

 $400 \, \text{ml}$ 

polyethylene terephthalate)

#### **Chemical Resistance**

The Incubation Tray is compatible with most detection reagents including weak acids, alcohols and organic solvents but is not compatible with strong acids or alkali.

Part No. LC2102.pps

Rev. Date: 20 Dec 2004

This product is distributed for laboratory research only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.

For technical questions about this product, call the Invitrogen Tech-Line™ U.S.A. 800 955 6288

#### **Additional Products**

The following products are available separately from Invitrogen for use with E-PAGE™ Gels. For more information on these products, visit our website at <a href="https://www.invitrogen.com">www.invitrogen.com</a> or contact Technical Service.

Product	Quantity	Catalog no.
E-PAGE™ 96 6% Gels	8 gels	EP096-06
E-PAGE™ 96 Starter Kit	1 kit	EPST96-06
Mother E-Base™ Device	1 unit	EB-M03
Daughter E-Base™ Device	1 unit	EB-D03
E-PAGE™ 48 8% Gels	8 gels	EP048-08
E-PAGE™ 48 Starter Kit	1 kit	EPST48-08
E-PAGE <sup>™</sup> MagicMark <sup>™</sup> Unstained Protein Standard	250 μl	LC5701
E-PAGE <sup>™</sup> SeeBlue <sup>®</sup> Pre-stained Standard	500 μl	LC5700
Invitrolon™ Filter Paper Sandwich (0.45 µm)	16/pk	LC2007
Nitrocellulose Filter Paper Sandwich (0.45 µm)	16/pk	LC2006
Nitrocellulose Filter Paper Sandwich (0.2 µm)	16/pk	LC2009
NuPAGE® Transfer Buffer (20X)	1 L	NP0006-1
NuPAGE® Antioxidant	15 ml	NP0005
WesternBreeze® Chromogenic Immunodetection System-Anti-Mouse	1 kit	WB7103
WesternBreeze® Chemiluminescent Immunodetection System-Anti-Mouse	1 kit	WB7104

©2004 Invitrogen Corporation. All rights reserved