



# M0328S 🖌 🖼 🕅 🛤

A NEW ENGLAND BioLabs

1-800-632-7799

info@neb.com

www.neb.com

BioLabs

1-800-632-7799

info@neb.com

www.neb.com

💓 RR 65° M&

500 units 5,000 U/ml Lot: 0041209 RECOMBINANT Store at -20°C Exp: 9/14

**Description:** Bst DNA Polymerase, Full Length is the full length polymerase from Bacillus stearothermophilus. It has  $5' \rightarrow 3'$  polymerase and double-strand specific  $5' \rightarrow 3'$  exonuclease activity, but lacks  $3' \rightarrow 5'$  exonuclease activity (2).

**Source:** An *E. coli* strain that carries the DNA Polymerase gene from *Bacillus stearothermophilus*. (H. Kong)

Supplied in: 50 mM KCl, 10 mM Tris-HCl (pH 7.5), 0.1 mM EDTA, 1 mM dithiothreitol, 0.1% Triton X-100 and 50% glycerol.

Bst DNA Polymerase, Full Length



# M0328S

500 units 5,000 U/ml Lot: 0041209 RECOMBINANT Store at -20°C Exp: 9/14

**Description:** Bst DNA Polymerase, Full Length is the full length polymerase from Bacillus stearothermophilus. It has  $5' \rightarrow 3'$  polymerase and double-strand specific  $5' \rightarrow 3'$  exonuclease activity, but lacks  $3' \rightarrow 5'$  exonuclease activity (2).

**Source:** An *E. coli* strain that carries the DNA Polymerase gene from *Bacillus stearothermophilus*. (H. Kong)

Supplied in: 50 mM KCl, 10 mM Tris-HCl (pH 7.5), 0.1 mM EDTA, 1 mM dithiothreitol, 0.1% Triton X-100 and 50% glycerol. **Reagents Supplied with Enzyme:** 10X ThermoPol<sup>™</sup> Reaction Buffer

Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 65°C.

## 1X ThermoPol Reaction Buffer:

20 mM Tris-HCl 10 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> 10 mM KCl 2 mM MgSO<sub>4</sub> 0.1% Triton<sup>®</sup> X-100 pH 8.8 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme that will incorporate 10 nmol of dNTP into acid insoluble material in 30 minutes at 65°C.

Unit Assay Conditions: 50 mM KCl, 20 mM Tris-HCl (pH 8.8), 10 mM MgCl<sub>2</sub>, 30 nM M13mp18 SS DNA, 70 nM M13 sequencing primer (-47) 24-mer, 200  $\mu$ M dATP, 200  $\mu$ M dCTP, 200  $\mu$ M dGTP, 100  $\mu$ M dTTP including [<sup>3</sup>H]-dTTP, 100  $\mu$ g/ml BSA and enzyme.

Specific Activity: 80,000 units/mg

Heat Inactivation: 80°C for 20 minutes.

**Reagents Supplied with Enzyme:** 10X ThermoPol<sup>™</sup> Reaction Buffer

Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 65°C.

#### 1X ThermoPol Reaction Buffer:

20 mM Tris-HCl 10 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> 10 mM KCl 2 mM MgSO<sub>4</sub> 0.1% Triton<sup>®</sup> X-100 pH 8.8 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme that will incorporate 10 nmol of dNTP into acid insoluble material in 30 minutes at 65°C.

Unit Assay Conditions: 50 mM KCl, 20 mM Tris-HCl (pH 8.8), 10 mM MgCl<sub>2</sub>, 30 nM M13mp18 SS DNA, 70 nM M13 sequencing primer (-47) 24-mer, 200  $\mu$ M dATP, 200  $\mu$ M dCTP, 200  $\mu$ M dGTP, 100  $\mu$ M dTTP including [<sup>3</sup>H]-dTTP, 100  $\mu$ g/ml BSA and enzyme.

Specific Activity: 80,000 units/mg

Heat Inactivation: 80°C for 20 minutes.

# Quality Control Assays

**Endonuclease Activity:** Incubation of a 50 µl reaction in ThermoPol Reaction Buffer containing a minimum of 50 units of *Bst* DNA Polymerase, Full Length with 1 µg of supercoiled  $\phi$ X174 DNA for 4 hours at either 37°C or 65°C results in < 10% conversion to the nicked form as determined by agarose gel electrophoresis.

**Physical Purity:** Purified to > 95% homogeneity as determined by SDS-PAGE analysis using Coomassie Blue detection.

Notes On Use: Bst DNA Polymerase, Full Length does not exhibit  $3' \rightarrow 5'$  exonuclease activity. Reaction temperatures above 70°C are not recommended.

Cannot be used for thermal cycle sequencing or PCR.

#### References:

**Quality Control Assays** 

Coomassie Blue detection.

recommended.

References:

76-87.

185-195.

PCR.

1. Mead, D.A. et al. (1991) *BioTechniques*, 11, 76–87.

Endonuclease Activity: Incubation of a 50 µl

Full Length with 1 µg of supercoiled  $\phi$ X174

in < 10% conversion to the nicked form as

determined by agarose gel electrophoresis.

as determined by SDS-PAGE analysis using

reaction in ThermoPol Reaction Buffer containing

a minimum of 50 units of *Bst* DNA Polymerase,

DNA for 4 hours at either 37°C or 65°C results

**Physical Purity:** Purified to > 95% homogeneity

Notes On Use: Bst DNA Polymerase, Full Length

Cannot be used for thermal cycle sequencing or

1. Mead, D.A. et al. (1991) *BioTechniques*, 11,

2. Aliotta, J.M. et al. (1996) Genetic Analysis, 12,

does not exhibit  $3' \rightarrow 5'$  exonuclease activity.

Reaction temperatures above 70°C are not

 Aliotta, J.M. et al. (1996) *Genetic Analysis*, 12, 185–195.

### **Companion Products Sold Separately:**

Magnesium Sulfate (MgSO<sub>4</sub>) Solution #B1003S 6.0 ml

ThermoPol Reaction Buffer Pack #B9004S 6.0 ml

ThermoPol II (Mg-free) Reaction Buffer Pack #B9005S 6.0 ml

ThermoPol DF (Detergent-free) Reaction Buffer Pack #B9013S 6.0 ml

Deoxynucleotide Solution Set #N0446S 25 μmol each

DeoxynucleotideSolution Mix#N0447S8 μmol each#N0447L40 μmol each

U.S. Patent No. 5,814,506.

THERMOPOL<sup>™</sup> is a trademark of New England Biolabs, Inc. TRITON<sup>®</sup> is a registered trademark of Union Carbide Corporation.



CERTIFICATE OF ANALYSIS

# **Companion Products Sold Separately:**

Magnesium Sulfate (MgS0<sub>4</sub>) Solution #B1003S 6.0 ml

ThermoPol Reaction Buffer Pack #B9004S 6.0 ml

ThermoPol II (Mg-free) Reaction Buffer Pack #B9005S 6.0 ml

ThermoPol DF (Detergent-free) Reaction Buffer Pack #B9013S 6.0 ml

Deoxynucleotide Solution Set #N0446S 25 μmol each

DeoxynucleotideSolution Mix#N0447S8 μmol each#N0447L40 μmol each

U.S. Patent No. 5,814,506

THERMOPOL<sup>™</sup> is a trademark of New England Biolabs, Inc. TRITON<sup>®</sup> is a registered trademark of Union Carbide Corporation.

