HeLa Genomic DNA



1-800-632-7799 info@neb.com www.neb.com

N4006S

15 μg Lot: 0041211 Exp: 11/14 100 μg/ml Store at -20°C

Description: Human female HeLa genomic DNA.

Source: HeLa (cervix adenocarcinoma) cells were grown to confluency in DMEM plus 10% fetal bovine serum. Genomic DNA was isolated by a standard genomic purification protocol (1), phenol extracted and equilibrated to 10 mM Tris-HCI (pH 7.5) and 1 mM EDTA.

Applications:

- A control for Methylation-Specific PCR (MSP) (2), Bisulfite sequencing, Methylationsensitive Single-Nucleotide Primer Extension (Ms-SNuPE), Combined Bisulfite Restriction Analysis (COBRA), Bisulfite treatment and PCR-Single-Strand Conformation Polymorphism Analysis (Bisulfite-PCR-SSCP/BiPS).
- PCR, SNP Analysis, Southern Blotting
- Genomic DNA library construction

Supplied in: 10 mM Tris-HCl (pH 7.5) and 1 mM EDTA. Avoid multiple freeze/thaw cycles.

Quality Assurance: Purified free of contaminating proteins and RNA.

A_{260/280} Ratio: 1.87

Reference:

- Sambrook, J. and Russell, D. (2001) *Molecular Cloning: A Laboratory Manual*, (3rd ed.), (pp. 6.4–6.12). Cold Spring Harbor: Cold Spring Harbor Laboratory Press.
- Herman, J.G. and Baylin, S.B. (1996). U.S. Patent No. 5,786,146. John Hopkins University School of Medicine.

CERTIFICATE OF ANALYSIS

HeLa Genomic DNA



1-800-632-7799 info@neb.com www.neb.com

N4006S

15 μg Lot: 0041211 Exp: 11/14 100 μg/ml Store at -20°C

Description: Human female HeLa genomic DNA.

Source: HeLa (cervix adenocarcinoma) cells were grown to confluency in DMEM plus 10% fetal bovine serum. Genomic DNA was isolated by a standard genomic purification protocol (1), phenol extracted and equilibrated to 10 mM Tris-HCI (pH 7.5) and 1 mM EDTA.

Applications:

- A control for Methylation-Specific PCR (MSP) (2), Bisulfite sequencing, Methylationsensitive Single-Nucleotide Primer Extension (Ms-SNuPE), Combined Bisulfite Restriciton Analysis (COBRA), Bisulfite treatment and PCR-Single-Strand Conformation Polymorphism Analysis (Bisulfite-PCR-SSCP/BiPS).
- PCR, SNP Analysis, Southern Blotting
- Genomic DNA library construction

Supplied in: 10 mM Tris-HCl (pH 7.5) and 1 mM EDTA. Avoid multiple freeze/thaw cycles.

Quality Assurance: Purified free of contaminating proteins and RNA.

A_{260/280} Ratio: 1.87

Reference:

- Sambrook, J. and Russell, D. (2001)
 Molecular Cloning: A Laboratory Manual,
 (3rd ed.), (pp. 6.4–6.12). Cold Spring Harbor:
 Cold Spring Harbor Laboratory Press.
- 2. Herman, J.G. and Baylin, S.B. (1996). U.S. Patent No. 5,786,146. John Hopkins University School of Medicine.